

## DNeasy<sup>®</sup> 96 Plant Kit

The DNeasy 96 Plant Kit (cat. no. 69181) can be stored at room temperature (15–25°C) for up to 1 year.

### Further information

- *DNeasy Plant Handbook*: [www.qiagen.com/HB-1166](http://www.qiagen.com/HB-1166)
- *TissueLyser Handbook*: [www.qiagen.com/HB-0487](http://www.qiagen.com/HB-0487)
- Safety Data Sheets: [www.qiagen.com/safety](http://www.qiagen.com/safety)
- Technical assistance: [support.qiagen.com](http://support.qiagen.com)

### Notes before starting

- If necessary, redissolve any precipitates in Buffer AP1 and Buffer AW1 concentrates.
- Add ethanol to Buffer AW1 and Buffer AW2 concentrates.
- Preheat Buffer AP1 to 65°C.
- Prepare a fresh working lysis solution: For 2 x 96 samples, combine 90 ml Buffer AP1, 225 µl RNase A, and 225 µl Reagent DX.

### Procedure

1. Place up to 50 mg leaves into each tube in 2 collection microtube racks.
2. Add 1 tungsten carbide bead to each collection microtube.
3. Pipet 400 µl working lysis solution into each collection microtube. Tightly seal the microtubes using the caps provided.
4. Assemble each rack of collection microtubes into the TissueLyser.
5. Grind the sample for 1.5 minutes at 30 Hz.
6. Reassemble the racks so that the collection microtubes nearest the TissueLyser in steps 4 and 5 are now furthest from the TissueLyser.
7. Grind the samples for another 1.5 min at 30 Hz.



8. Centrifuge to collect any solution from the caps.
9. Add 130  $\mu$ l Buffer P3 to each collection microtube and reseal using new caps.
10. Place a clear cover over each rack and shake vigorously up and down for 15 s.  
Centrifuge to collect any solution from the caps.
11. Incubate the collection-microtube racks for 10 min at  $-20^{\circ}\text{C}$ .
12. Centrifuge the collection-microtube racks for 5 min at  $5788 \times g$  (6000 rpm).
13. Transfer 400  $\mu$ l of each supernatant to a new collection microtube.
14. Add 600  $\mu$ l of Buffer AW1 to each sample. Close microtubes with new caps.
15. Place a clear cover over each rack and shake vigorously up and down for 15 s.  
Centrifuge to collect any solution from the caps.
16. Place 2 DNeasy 96 plates on top of S-Blocks. Mark the DNeasy 96 plates for later sample identification.
17. Transfer 1 ml of each sample to each well of the DNeasy 96 plates.
18. Seal each DNeasy 96 plate with an AirPore Tape Sheet. Centrifuge for 4 min at  $5788 \times g$ .  
If lysate remains in the DNeasy 96 plates after centrifugation, centrifuge for another 4 min.
19. Remove the AirPore Tape Sheet. Add 800  $\mu$ l Buffer AW2 to each sample.
20. Centrifuge for 15 min at  $5788 \times g$  without AirPore Tape Sheet to dry the membranes.
21. Place each DNeasy 96 plate on a new Elution Microtubes RS rack.
22. Add 100  $\mu$ l Buffer AE and seal with new AirPore Tape Sheets. Incubate for 1 min at room temperature ( $15$ – $25^{\circ}\text{C}$ ). Centrifuge for 2 min at  $5788 \times g$ .
23. Repeat step 22. Seal the Elution Microtubes RS with new caps to store DNA.

## Document Revision History

Date	Changes
04/2021	Changed speed condition " $3800 \times g$ " to " $5788 \times g$ ". Revised the "Notes before starting" section.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN<sup>®</sup>, Sample to Insight<sup>®</sup>, DNeasy<sup>®</sup> (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law. 1123238 04/2021 HB-0544-003 © 2021 QIAGEN, all rights reserved.