

QIAamp[®] DSP Virus Spin Kit

The QIAamp DSP Virus Spin Kit is designed for manual or automated purification of viral nucleic acids from human plasma and serum. The kit is designed to provide rapid and reliable viral DNA and RNA purification while minimizing cross-contamination risks.

Kit performance is not guaranteed for each virus species and must be validated by the user. It is the user's responsibility to validate system performance for any procedures used in their laboratory that are not covered by the QIAGEN[®] performance evaluation studies.

Performance characteristics

The performance of the QIAcube for the QIAamp DSP Virus Spin Kit was compared using hepatitis C virus (HCV) RNA as an example virus. The tests were performed with a dilution of quantified virus panels made in HCV-negative human plasma and serum (n=15).

HCV RNA was detected using an in-house RT-PCR assay (Figure 1). Viral nucleic acids were purified from 200 or 400 μ l samples with an elution volume of 60 μ l.




Figure 1. Performance comparison of the QIAamp DSP Virus Spin Kit. The performance of the QIAamp DSP Virus Spin Kit in combination with different protocols was compared by analysis of serum and plasma samples. The QIAcube protocols used were: QIAamp DSP Virus Spin — Large Plasma Samples — Standard (400 μ l); QIAamp DSP Virus Spin — Large Plasma Samples — Manual Lysis (400 μ l); QIAamp DSP Virus Spin — Standard (200 μ l, plasma or serum samples); and QIAamp DSP Virus Spin — Manual Lysis (200 μ l, plasma or serum samples). Viral RNA was detected using viral dilution series and an in-house RT-PCR assay for HCV RNA.

Linear range

Hepatitis C virus (HCV) RNA was extracted using the QIAamp DSP Virus Spin Kit with one manual and two automated QIAcube protocols: QIAamp DSP Virus Spin (200 μ l); automated QIAamp DSP Virus Spin — Large Plasma Samples — Standard (400 μ l); and automated QIAamp DSP Virus Spin — Standard (200 μ l). Performance of the QIAamp DSP Virus Spin Kit was compared to the QIAamp DSP Virus Kit using the manual QIAamp DSP Virus protocol (500 μ l) together with the QIAvac 24 Plus system. The tests were performed with dilutions of quantified virus panels made in HCV-negative human plasma. Dilution series with 5 different virus titers were tested, with 12 replicates each. The linear range of the QIAamp DSP Virus Spin Kit procedure was determined for HCV using an in-house RT-PCR assay (Figure 2). Viral nucleic acids were purified from sample volumes of 200–500 μ l, with a 60 μ l elution volume.




Figure 2. The linear range of the QIAamp DSP Virus Spin Kit. The linear range of yields using the QIAamp DSP Virus Spin Kit (cat. no. 61704) with 1 manual and 2 different automated protocols, compared to the manual QIAamp DSP Virus Kit (cat. no. 60704) procedure using the QIAvac 24 Plus. The linear range of the protocols was determined using viral dilution series and an in-house RT-PCR assay for HCV RNA.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, QIAamp®, QIAcube® (QIAGEN Group).

Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

© 2012 QIAGEN, all rights reserved.

www.qiagen.com

Australia ■ 1-800-243-800

Austria ■ 0800/281010

Belgium ■ 0800-79612

China ■ 021-51345678

Denmark ■ 80-885945

Finland ■ 0800-914416

France ■ 01-60-920-930

Germany ■ 02103-29-12000

Hong Kong ■ 800 933 965

Ireland ■ 1800 555 049

Italy ■ 800 787980

Japan ■ 03-5547-0811

Korea (South) ■ 1544 7145

Luxembourg ■ 8002 2076

The Netherlands ■ 0800 0229592

Norway ■ 800-18859

Singapore ■ 65-67775366

Spain ■ 91-630-7050

Sweden ■ 020-790282

Switzerland ■ 055-254-22-11

UK ■ 01293-422-911



Sample & Assay Technologies