



BioMag[®] Goat Anti-Mouse IgG

BioMag[®] suspensions are a series of cell-sorting products containing superparamagnetic particles. BioMag particles are irregularly shaped, with an average diameter of 1.6 μm , and consist of an iron oxide core with a silane coating. The surface of the particle is coated with amine or carboxyl groups, facilitating the covalent attachment of proteins, glycoproteins, secondary antibodies, and other ligands, with retention of biological activity. The particles' irregular shape provides a large surface area, ensuring a high binding efficiency.

BioMag particles are superparamagnetic — they respond to magnetic fields but do not retain magnetic properties upon removal of the magnetic field. This inability to become magnetized permits magnetic extraction without magnetically induced aggregation. Rapid and efficient removal of BioMag particles from suspension is achieved by the application of an external magnetic field.

Depending on antigen availability and the total cell population, cell sorting may require between 10 and 50 magnetic particles per cell. Multiple cell sorts can also be performed. BioMag suspensions contain approximately 1×10^8 magnetic particles per milligram, and are supplied in concentrations of 1 or 5 mg/ml. The particle-to-cell ratio is based on the total cell population, as demonstrated in the following sample calculation.

Sample variables

- A system with 1×10^7 total cells
- 10–50 magnetic particles per cell (total cell population)
- BioMag suspension with 1 mg BioMag particles/ml

Therefore, 1×10^7 cells \times 10–50 particles per cell = $1\text{--}5 \times 10^8$ magnetic particles are required.

1 mg/ml BioMag suspension contains 1×10^8 magnetic particles per milligram, which is equivalent to 1×10^8 magnetic particles per milliliter. Therefore, the volume of BioMag suspension required is 1.0–5.0 ml per 10^7 cells.

Applications

BioMag Goat Anti-Mouse IgG is a standard BioMag particle coated with polyclonal goat anti-mouse IgG antibodies and is highly suited for use in cell sorting methods where a mouse IgG antibody is used as a primary antibody. BioMag Goat Anti-Mouse IgG can be used to separate the cells of interest from a heterogeneous cell population using negative selection. BioMag Goat Anti-Mouse IgG can also be used as a secondary antibody in enzyme immunoassays and radioassays that utilize a mouse IgG primary monoclonal antibody.

Specifications

Form	Magnetic particle suspension (1 mg/ml) in phosphate-buffered saline (PBS) containing EDTA and 0.08% sodium azide.
BioMag concentration	1 mg/ml
Binding capacity	1 ml (1 mg) BioMag Goat Anti-Mouse IgG binds >0.20 mg mouse IgG. The specificity of goat anti-mouse for IgG subclasses may vary.
Storage and stability	BioMag Goat Anti-Mouse IgG can be stored at 2–8°C until the expiration date. Do not freeze. Do not dry. Centrifugation should only be used if it is the last step of a procedure, i.e., if resuspension of the BioMag particles is not required. Freezing, drying, and centrifugation result in extensive aggregation of the BioMag particles and a loss of binding activity.
Safety information	When working with chemicals always wear a suitable lab coat, disposable gloves, and protective goggles. For more information please consult the appropriate material safety data sheet, available online at www.qiagen.com/ts/msds.asp

Preparation of cells

Lymphocytes can be partially purified using a Ficoll® gradient, nylon wool, or other similar method, although it may be possible to sort cells directly from whole blood or other sources. Since any given cell source will have unique purification requirements, we recommend determining the optimal conditions. For optimal results, we recommend using the most dilute cellular suspension possible.

General protocol considerations

The protocol for the indirect method of cell separation is a general guide to cell separation using the BioMag System, and is applicable in most cases. We recommend optimizing the conditions for your cells. A linear scaleup or -down of volumes and cell numbers is possible. Guidelines for the direct method of cell separation are included. We recommend using the indirect method, since less primary antibody is required.

Protocol for indirect method of cell separation

1. **Wash required amount of BioMag particles 2–3 times in appropriate sterile culture medium or buffer. Resuspend in 500 μ l sterile medium. Use a magnetic separator (>20 megaoersted) to pull the magnetic particles to the side of the tube.**

Washing of BioMag particles removes the sodium azide preservative. Low protein (5–10% FCS or 0.5% BSA) media and buffers are recommended to reduce nonspecific binding. If coagulation of the cells is a problem, include an anticoagulant, such as 2 mM EDTA in the sterile culture medium or buffer. Separation should be performed for 2–3 min using a magnetic separator. Appropriate magnets are available from QIAGEN; see Ordering Information.

Important: Do not centrifuge the BioMag suspension during wash steps. Centrifugation results in extensive aggregation and loss of binding activity.

2. **Prepare 10^7 total cells in 100–500 μ l sterile culture medium.**

3. **Add the appropriate amount of primary antibody to the cells and incubate for 20–30 min at 4°C.**

For optimal results, we recommend optimizing the amount of primary antibody. As little as 0.1–1.0 μ g primary antibody per 10^7 target cells may be sufficient for effective selection of target cells using this protocol.

4. **Remove unbound primary antibody from the cells by centrifugation at 200–300 \times g for 5–10 min. Resuspend the cells in 1 ml sterile medium and centrifuge at 200–300 \times g for 5–10 min to wash the cells. Repeat wash step 1–2 times, using 1 ml sterile medium for each wash.**

5. **Resuspend cells in 1 ml sterile medium and add the washed BioMag particles prepared in step 1.**

Note: We recommend sorting cells in total volumes \geq 1 ml (including BioMag suspension and cell volume). Where volumes < 1 ml are used, additional medium or buffer should be added to a final volume of 1 ml.

6. **Incubate at 4°C for 15–30 min. Swirl reaction vessel occasionally, or place eppendorf tubes on a rotating wheel during incubation.**

Longer incubation is not recommended as magnetic particles may detach from the target cells as a result of cell surface changes over time.

Room temperature (15–25°C) or 37°C are optimal for some cell types and QIAGEN recommends that the optimal cell sorting conditions be individually determined.

7. **Apply vessel to a magnetic separator for 5–10 min.**

Note: A clear supernatant indicates that the separation is complete.

Separation must be performed with the vessel held vertically so that the pellet forms on the side of the flask or tube. This ensures unselected cells do not contaminate the magnetic pellet.

8. **Once separation is complete, carefully remove the supernatant without disturbing the pellet.**

For optimal results, we recommend repeating step 7.

9. **Centrifuge the supernatant at 200–300 \times g for 5–10 min to pellet the cells, and resuspend the cells in fresh medium.**

Guidelines for the direct method of cell separation

Direct binding of primary antibodies to the BioMag particles may be useful for some applications. Washed BioMag particles (see step 1 of the indirect method for cell separation) should be incubated with the appropriate amount of primary antibody at 4°C for 20 min in 500 µl of sterile medium or buffer. Considerably larger amounts of primary antibody are required compared with the indirect method of cell separation. After separation (using a suitable magnetic separator), the BioMag particles/antibody complex should be washed 3 times using 1 ml sterile medium for each wash. Afterwards, the protocol for the indirect method of cell separation should be followed from step 5.

Ordering Information

Product	Contents	Cat. No.
BioMag Goat Anti-Mouse IgG (50 ml)	BioMag goat anti-mouse IgG secondary antibody suspension (1 mg/ml)	310004
BioMag Goat Anti-Mouse IgG (500 ml)	BioMag goat anti-mouse IgG secondary antibody suspension (1 mg/ml)	310007
Accessories		
Single-Tube Magnet	Magnet for separating magnetic particles in a 1.5 ml or 2 ml tube	36910
12-Tube Magnet	Magnet for separating magnetic particles in 12 x 1.5 ml or 2 ml tubes	36912
15 ml/50 ml Tube Magnet	Magnet for separating magnetic particles in 5 x 15 ml and 3 x 50 ml tubes	36935
Flask Magnet	Magnet for separating magnetic particles in a cell culture flask	36937

Patented or patent-pending and/or registered or registration-pending trademark of the QIAGEN Group: QIAGEN®, BioMag®. Ficoll is a registered trademark of Amersham Pharmacia Biotech.

Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

© 2002 QIAGEN, all rights reserved.

www.qiagen.com

QIAGEN:

Australia
Tel. 03-9489-3666
Fax 03-9489-3888

Canada
Tel. 800-572-9613
Fax 800-713-5951

France
Tel. 01-60-920-930
Fax 01-60-920-925

Germany
Tel. 02103-29-12400
Fax 02103-29-22000

Italy
Tel. 02-33430411
Fax 02-33430426

Japan
Tel. 03-5547-0811
Fax 03-5547-0818

Switzerland
Tel. 061-319-3031
Fax 061-319-3033

UK and Ireland
Tel. 01293-422-999
Fax 01293-422-922

USA
Tel. 800-426-8157
Fax 800-718-2056

Distributors:

Argentina Tecnolab S.A. [011] 4555 0010 **Austria/Slovenia** VWR International GmbH [01] 576 00 0 **Belgium/Luxemburg** Westburg b.v. 0800-1-9815 **Brazil** Uniscience do Brasil 011 3622 2320 **China** Gene Company Limited [852]2896-6283 **Cyprus** Scientronics Ltd [02] 357 22 765416 **Czech Republic** BIOCONSULT spol. s.r.o. [420] 2 417 29 792 **Denmark** VWR International ApS 43 86 87 88 **Egypt** Clinilab 52 57 212 **Finland** VWR International Oy [09] 804 551 **Greece** BioAnalytica S.A. [10]640 03 18 **India** Genetix [011]542 1714 or [011]515 9346 **Israel** Westburg (Israel) Ltd. 08 6650813/4 or 1-800 20 22 20 **Korea** ILS Laboratories, Inc. [02] 924-86 97 **Malaysia** RESEARCH BIOLABS SDN. BHD. [603]8070 3101 **Mexico** Quimica Valaner S.A. de C.V. [55] 55 25 57 25 **The Netherlands** Westburg b.v. [033]4950094 **New Zealand** Biolab Scientific Ltd. [09] 980 6700 or 0800 933 966 **Norway** VWR International AS 22 90 00 00 **Poland** Syngen Biotech Sp. z.o.o. [071] 351 41 06 or 0601 70 60 07 **Portugal** IZASA PORTUGAL LDA [21] 424 7312 **Singapore** Research Biolabs Pte Ltd 2731066 **Slovak Republic** BIOCONSULT Slovakia spol. s.r.o. [02] 5022 1336 **South Africa** Southern Cross Biotechnology (Pty) Ltd [021] 671 5166 **Spain** IZASA, S.A. [93] 902.20.30.90 **Sweden** VWR International AB [08] 621 34 00 **Taiwan** TAIGEN Bioscience Corporation [02] 2880 2913 **Thailand** Theera Trading Co. Ltd. [02] 412-5672 **In other countries contact:** QIAGEN, Germany

