

# Instructions for Q-Rex Gene Expression Plug-in version 2.0

## Adaptions for using the comparative $C_T$ method ( $2^{-\Delta\Delta C_q}$ )

In the Q-Rex Gene Expression update to version 2.0, the relative quantitation calculation is based on the more generalized mathematical model suggested by Pfaffl\*. This has the following advantages for the user:

- Normalization to more than one reference gene
- Adjustment of the real-time PCR efficiency

In case of 100% reaction efficiency for all tested targets and one reference gene, the relative quantity (RQ) results with version 2.0 are the same as with comparative  $C_T$  method ( $2^{-\Delta\Delta C_q}$ ) used in version 1.0. In fact, the  $2^{-\Delta\Delta C_q}$  method is simply a special case of the Pfaffl method, in which the efficiency of the target and reference gene is 2 (100%).

This guide provides step-by-step instructions for the setup with Q-Rex Gene Expression Plug-in version 2.0 to ensure the RQ values are the same as with version 1.0.

### Procedure

1. Open your finished experiment.
2. Make sure to define only one target as a reference gene (RG) in the sample layout step of the Step Marker. If you are working with SYBR Green, add and define new targets and assign them accordingly to the samples.

#### Probe

Define targets

Target	Acquisition	Type	
ref gene	Green	RG	<input type="checkbox"/>
target 1	Yellow	Test	<input type="checkbox"/>
target 2	Orange	Test	<input type="checkbox"/>

Add target

#### SYBR Green


Define targets

Target	Acquisition	Type	
Gene of i...	Green	Test	<input type="checkbox"/>
Referenc...	Green	RG	<input type="checkbox"/>


Add target

\* Pfaffl, M.W. (2001) A new mathematical model for relative quantification in real-time RT-PCR. Nucleic Acids Res. **29**, e45.

3. Add gene-expression analysis to the experiment by clicking the menu icon in the "Analysis" step of the Step Marker and selecting the gene expression analysis from the list.



4. Before you can get gene expression results, you need to define the efficiency and calibrator following the yellow highlighted fields. Check the "Use efficiency of 100% for all targets" box in the "Gene expression" panel.



5. Change to the "Sample result" view.

Notice that, in version 2.0 of the Q-Rex Gene Expression Plug-in, the results are grouped at acquisition and are ordered one after the other.


The RQ values are the same as in the previous version. However, the RQ min/max in the previous version was the standard deviation (SD) and in version 2.0 is the standard error (SE).

**Note:** Make sure that the technical replicates have exactly the same name, otherwise you will not see RQ results or the relative quantity plot.

## Exporting $\Delta\Delta Cq$ values

You can export  $\Delta\Delta Cq$  values with Q-Rex Software version 1.1.

1. To export  $\Delta\Delta Cq$  values, switch to the "Report/export" step in the Step Marker.
2. Click the "Export" tab.
3. Select the "Relative expression ( $2^{-\Delta\Delta Cq}$ )" export format.
4. Click "Export".
5. Enter a file name for the export and the location where it should be saved.  
The Q-Rex Software uses the export format as a default file name.
6. Click "OK" to save the export file.



---

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN® kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at [www.qiagen.com](http://www.qiagen.com) or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, Sample to Insight®, RNeasy® (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.  
08/2018 PROM-12987-001 © 2018 QIAGEN, all rights reserved.

Ordering [www.qiagen.com/shop](http://www.qiagen.com/shop) | Technical Support [support.qiagen.com](http://support.qiagen.com) | Website [www.qiagen.com](http://www.qiagen.com)