
april 2019.

QuantiFERON[®]-TB Gold Plus (QFT[®]-Plus) ELISA uputstvo za upotrebu

Verzija 1

IVD

Za korišćenje u in vitro dijagnostici

Interferon-gama (IFN- γ) test pune krvi koji meri reakcije na peptidne antigene ESAT-6 i CFP-10

REF

622120, 622822

QIAGEN GmbH, QIAGEN Strasse 1, 40724 Hilden,
Nemačka

R6 **MAT**

1083163SR

Sadržaj

Namena.....	5
Rezime i objašnjenje testa.....	5
Principi ispitivanja.....	7
Vreme potrebno za obavljanje ispitivanja.....	9
Komponente i skladištenje.....	10
Potreban materijal koji se ne isporučuje.....	12
Čuvanje i rukovanje uzorcima.....	13
Epruvete za prikupljanje uzoraka krvi.....	13
Komplet reagenasa.....	13
Rekonstituisani i neupotrebljeni reagensi.....	13
Upozorenja i mere opreza.....	14
Upozorenja.....	14
Mere opreza.....	15
Prikupljanje uzoraka i rukovanje njima.....	18
Uputstva za upotrebu.....	24
Faza 1 – Inkubacija krvi i prikupljanje plazme.....	24
Faza 2 – IFN- γ ELISA.....	25
Izračunavanje i tumačenje testa.....	30
Generisanje krive standarda.....	30
Kontrola kvaliteta testa.....	31
Tumačenje rezultata.....	31
Ograničenja.....	34

Karakteristike učinka	35
Kliničke studije	35
Karakteristike učinka ispitivanja	41
Tehničke informacije.....	46
Nejasni rezultati.....	46
Zgrušani uzorci plazme	46
Vodič za rešavanje problema	47
Reference.....	49
Simboli.....	58
Kontakt informacije	59
Skraćena procedura testiranja.....	60
Faza 1 – Inkubacija krvi	60
Faza 2 – IFN- γ ELISA.....	60
Važne promene	62
Istorija revizija priručnika	62

Namena

Ispitivanje QuantiFERON-TB Gold Plus (QFT-Plus) je in vitro dijagnostički test koji stimulira ćelije u hepariziranoj punoj krvi pomoću koktela peptida koji simuliraju proteine ESAT-6 i CFP-10. Otkrivanje proteina interferon- γ (IFN- γ) putem ELISA (Enzyme-Linked Immunosorbent Assay) ispitivanja služi za identifikaciju in vitro reakcija na ove peptidne antigene koji su povezani sa *Mycobacterium tuberculosis* infekcijom.

QFT-Plus je indirektan test za *M. tuberculosis* infekciju (uključujući oboljenje) i predviđen je za upotrebu u kombinaciji sa procenom rizika, radiografijom i drugim medicinskim i dijagnostičkim procenama.

Rezime i objašnjenje testa

Tuberkuloza je prenosiva bolest koju uzrokuje infekcija složenim *M. tuberculosis* (MTB) organizmima (*M. tuberculosis*, *M. bovis*, *M. africanum*) koja se najčešće prenosi na nove domaćine vazдушnim putem preko kapljičnih jezgara sa pacijenata obolelih od respiratorne tuberkuloze. Osoba koja je prvi put zaražena može da se razboli od tuberkuloze u roku od nekoliko nedelja ili par meseci, ali većina inficiranih osoba se ne razboli. Kod pojedinih osoba postoji latentna infekcija tuberkulozom (latent tuberculosis infection, LTBI), asimptomatsko stanje koje se ne prenosi na druge, a te osobe mogu da se razbole od tuberkuloze mesecima ili godinama kasnije. Glavna svrha dijagnostifikovanja LTBI stanja je razmatranje lečenja kojim se sprečava oboljevanje od tuberkuloze. Do nedavno, tuberkulinski kožni test (Tuberculin Skin Test, TST) predstavljao je jedini dostupan način za dijagnostifikovanje LTBI stanja. Kožna osetljivost na tuberkulin razvija se od 2 do 10 nedelja nakon inficiranja. Međutim, pojedine inficirane osobe, uključujući osobe sa širokim opsegom oboljenja koja umanjuju imunološke funkcije, ali i osobe bez ovih oboljenja, ne reaguju na tuberkulin. Suprotno tome, pojedine osobe kod kojih postoji mala verovatnoća da su zaraženi *M. tuberculosis* infekcijom pokazuju osetljivost na tuberkulin i imaju pozitivne TST rezultate

nakon vakcinacije bacilom Calmette-Guérin (BCG), infekciju mikobakterijom koja nije iz grupe *M. tuberculosis* ili druge neodređene faktore.

Mora se napraviti razlika između LTBI stanja i bolesti tuberkuloze, stanja koje se mora prijaviti i koje najčešće pogađa pluća i donje disajne puteve, iako sistemi drugih organa takođe mogu da budu zaraženi. Dijagnoza tuberkuloze se uspostavlja na osnovu istorijskih, fizičkih, radioloških, histoloških i mikrobakterioloških rezultata.

QFT-Plus je test za ćelijske imunske odgovore (cell-mediated immune, CMI) na peptidne antigene koji simuliraju mikobakterijske proteine. Ovi proteini, ESAT-6 i CFP-10 ne nalaze se ni u jednoj BCG vrsti, niti u većini netuberkuloznih mikobakterija, sa izuzetkom sledećih: *M. kansasii*, *M. szulgai* i *M. marinum* (1). Osobe inficirane složenim MTB organizmima najčešće imaju limfocite u krvi koji prepoznaju ove i druge mikobakterijske antigene. Ovaj proces prepoznavanja uključuje generisanje i sekreciju citokina, IFN- γ . Otkrivanje i naknadna kvantifikacija IFN- γ proteina čine osnovu ovog testa.

Antigeni koji se koriste u QFT-Plus su peptidni kokteli koji stimulišu proteine ESAT-6 i CFP-10. Veliki broj studija pokazao je da ovi peptidni antigeni stimulišu IFN- γ odgovore u T ćelijama kod osoba zaraženih bakterijom *M. tuberculosis*, ali generalno ne kod osoba koje nisu zaražene niti kod osoba koje su primile BCG vakcinu i koje nisu obolele ili kod kojih ne postoji rizik od LTBI infekcije (1–32). Međutim, medicinska lečenja ili stanja koja mogu da oslabe funkcije imunološkog sistema potencijalno mogu da smanje IFN- γ odgovore. Pacijenti sa određenim mikobakterijskim infekcijama takođe mogu da budu responsivni na ESAT-6 i CFP-10, jer se geni koji kodiraju ove proteine nalaze u *M. kansasii*, *M. szulgai* i *M. marinum* (1, 23). QFT-Plus je test za LTBI i korisno sredstvo za dijagnostifikovanje složene *M. tuberculosis* infekcije kod bolesnih pacijenata. Pozitivan rezultat podržava dijagnostifikovanje bolesti tuberkuloze, ali i infekcije koje izaziva druga mikobakterija (npr., *M. kansasii*) takođe mogu da dovedu do pozitivnih rezultata. Druge medicinske i dijagnostičke procene neophodne su da bi se potvrdila ili isključila tuberkuloza.

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Reference

1. Andersen, P. et al. (2000) Specific immune-based diagnosis of tuberculosis. *Lancet* 356, 1099.
2. Balcells, M.E. et al. (2008) A comparative study of two different methods for the detection of latent tuberculosis in HIV-positive individuals in Chile. *Int. J. Infect. Dis.* 12, 645.
3. Bartalesi, F. et al. (2009) QuantiFERON-TB Gold and TST are both useful for latent TB screening in autoimmune diseases. *Eur. Respir. J.* 33, 586.
4. Bocchino, M. et al. (2008) Performance of two commercial blood IFN-gamma release assays for the detection of *Mycobacterium tuberculosis* infection in patient candidates for anti-TNF-alpha treatment. *Eur. J. Clin. Microbiol. Infect. Dis.* 27,907.
5. Brock, I. et al. (2006) Latent tuberculosis in HIV positive, diagnosed by the *M. tuberculosis* specific interferon-gamma test. *Respir. Res.* 7, 56.
6. Chun, J.K. et al. (2008) The role of a whole blood interferon gamma assay for the detection of latent tuberculosis infection in bacille Calmette-Guerin vaccinated children. *Diagn. Microbiol. Infect. Dis.* 62, 389.
7. Connell, T.G. et al. (2008) A three-way comparison of tuberculin skin testing, QuantiFERON-TB gold and T-SPOT.TB in children. *PLoS ONE* 3, e2624. doi: 10.1371/journal.pone.0002624.
8. Detjen, A.K. et al. (2007) Interferon-gamma release assays improve the diagnosis of tuberculosis and nontuberculous mycobacterial disease in children in a country with a low incidence of tuberculosis. *Clin. Infect. Dis.* 45, 322.

-
9. Diel, R. et al. (2009) Comparative performance of tuberculin skin test, QuantiFERON-TB-Gold In-Tube assay, and T-Spot. TB test in contact investigations for tuberculosis. *Chest* 135, 1010.
 10. Diel, R. et al. (2008) Predictive value of a whole-blood IFN- γ assay for the development of active TB disease. *Am. J. Respir. Crit. Care Med.* 177, 1164.
 11. Diel, R. et al. (2006) Tuberculosis contact investigation with a new, specific blood test in a low-incidence population containing a high proportion of BCG-vaccinated persons. *Respir. Res.* 7, 77.
 12. Dogra, S. et al. (2007) Comparison of a whole blood interferon-gamma assay with tuberculin skin testing for the detection of tuberculosis infection in hospitalized children in rural India. *J. Infect.* 54, 267.
 13. Drobniowski, F. et al. (2007) Rates of latent tuberculosis in health care staff in Russia. *PLoS Med.* 4, e55.
 14. Gerogianni, I. et al. (2008) Whole-blood interferon-gamma assay for the diagnosis of tuberculosis infection in an unselected Greek population. *Respirology* 13, 270.
 15. Harada, N. et al. (2008) Comparison of the sensitivity and specificity of two whole blood interferon-gamma assays for *M. tuberculosis* infection. *J. Infect.* 56, 348.
 16. Higuchi, K. et al. (2009) Comparison of performance in two diagnostic methods for tuberculosis infection. *Med. Microbiol. Immunol.* 198, 33.
 17. Kang, Y.A. et al. (2005) Discrepancy between the tuberculin skin test and the whole-blood interferon gamma assay for the diagnosis of latent tuberculosis infection in an intermediate tuberculosis-burden country. *JAMA* 293, 2756.

-
18. Katiyar, S.K. et al. (2008) Use of the QuantiFERON-TB Gold In-Tube test to monitor treatment efficacy in active pulmonary tuberculosis. *Int. J. Tuberc. Lung Dis.* 12, 1146.
 19. Kipfer, B. et al. (2008) Tuberculosis in a Swiss army training camp: contact investigation using an Interferon gamma release assay. *Swiss. Med. Wkly.* 138, 267.
 20. Luetkemeyer, A. et al. (2007) Comparison of an interferon-gamma release assay with tuberculin skin testing in HIV-infected individuals. *Am. J. Respir. Crit. Care Med.* 175, 737.
 21. Mackensen, F. et al. (2008) QuantiFERON TB-Gold - A new test strengthening long-suspected tuberculous involvement in serpiginous-like choroiditis. *Am. J. Ophthalmol.* 146, 761.
 22. Manuel, O. et al. (2007) Comparison of Quantiferon-TB Gold with tuberculin skin test for detecting latent tuberculosis infection prior to liver transplantation. *Am. J. Transplant.* 7, 2797.
 23. Matulis, G. et al. (2007) Detection of latent tuberculosis in immunosuppressed patients with autoimmune diseases performance of a *Mycobacterium tuberculosis* antigen specific IFN-gamma assay. *Ann. Rheum. Dis.* 67, 84.
 24. Mirtskhulava, V. et al. (2008) Prevalence and risk factors for latent tuberculosis infection among health care workers in Georgia. *Int. J. Tuberc. Lung Dis.* 12, 513.
 25. Nakaoka, H. et al. (2006) Risk for tuberculosis among children. *Emerging Infect. Dis.* 12, 1383.
 26. Pai, M. et al. (2005) *Mycobacterium tuberculosis* infection in health care workers in rural India: comparison of a whole-blood, interferon-g assay with tuberculin skin testing. *JAMA* 293, 2746.

-
27. Ponce de Leon, D. et al. (2008) Comparison of an interferon-gamma assay with tuberculin skin testing for detection of tuberculosis (TB) infection in patients with rheumatoid arthritis in a TB-endemic population. *J Rheumatol.* 35, 776.
 28. Richeldi, L. et al. (2008) Prior tuberculin skin testing does not boost QuantiFERON-TB results in paediatric contacts. *Eur. Respir. J.* 32, 524.
 29. Rothel, J.S. and Andersen, P. (2005) Diagnosis of latent *Mycobacterium tuberculosis* infection: is the demise of the Mantoux test imminent? *Expert Rev. Anti Infect. Ther.* 3, 981.
 30. Schoepfer, A.M. et al. (2008) Comparison of interferon-gamma release assay versus tuberculin skin test for tuberculosis screening in inflammatory bowel disease. *Am. J. Gastroenterol.* 103, 2799.
 31. Silverman, M.S. et al. (2007) Use of an interferon-gamma based assay to assess bladder cancer patients treated with intravesical BCG and exposed to tuberculosis. *Clin. Biochem.* 40, 913.
 32. Stebler, A. et al. (2008) Whole-blood interferon-gamma release assay for baseline tuberculosis screening of healthcare workers at a Swiss university hospital. *Infect. Control Hosp. Epidemiol.* 29, 681.
 33. Turner, J. et al. (1996) Stimulation of human peripheral blood mononuclear cells with live *Mycobacterium bovis* BCG activates cytolytic CD8+ T cells in vitro. *Immunology* 87, 339.
 34. Brookes, R.H. et al. (2003) CD8+ T cell-mediated suppression of intracellular *Mycobacterium tuberculosis* growth in activated human macrophages. *Eur. J. Immunol.* 33, 3293.
 35. Stenger, S. et al. (1998) An antimicrobial activity of cytolytic T cells mediated by granulysin. *Science* 282, 121.

-
36. Lalvani, A. et al. (1998) Human cytolytic and interferon gamma-secreting CD8+ T lymphocytes specific for *Mycobacterium tuberculosis*. Proc. Natl. Acad. Sci. U.S.A. 95, 270.
 37. Lewinsohn, D.M. et al. (2001) Classically restricted human CD8+ T lymphocytes derived from *Mycobacterium tuberculosis*-infected cells: definition of antigenic specificity. J. Immunol. 166, 439.
 38. Lewinsohn, D.A. et al. (2007) Immunodominant tuberculosis CD8 antigens preferentially restricted by HLA-B. PLoS Pathol. 3, 1240.
 39. Day, C.L. et al. (2011) Functional capacity of *Mycobacterium tuberculosis*-specific T cell responses in humans is associated with mycobacterial load. J. Immunol. 187, 2222.
 40. Rozot, V. et al. (2013) *Mycobacterium tuberculosis*-specific CD8+ T cells are functionally and phenotypically different between latent infection and active disease. Eur. J. Immunol. 43, 1568.
 41. Nikolova, M. et al. (2013) Antigen-specific CD4- and CD8-positive signatures in different phases of *Mycobacterium tuberculosis* infection. Diagn. Microbiol. Infect. Dis. 75, 277.
 42. Chicchio, T. et al. (2014) Polyfunctional T-cells and effector memory phenotype are associated with active TB in HIV-infected patients. J. Infect. doi: 10.1016/j.jinf.2014.06.009. Epub.
 43. Ongaya, A. et al. (2013) Mycobacterium tuberculosis-specific CD8+ T cell recall in convalescing TB subjects with HIV co-infection. Tuberculosis 93, S60.
 44. Lanicioni, C. et al. (2012) CD8+ T cells provide an immunologic signature of tuberculosis in young children. Am. J. Respir. Crit. Care Med. 185, 206.

-
45. Long, G., Ji-Chun, M., Min, Jin-Long, L., Jin-Hui, T. (2014) Interferon- γ release assay for the diagnosis of latent *Mycobacterium tuberculosis* infection in children younger than 5 years: a meta-analysis. Clin. Pediatr. 53, 1255.
 46. Howley, M.M. et al. (2015) Evaluation of QuantiFERON-TB Gold In-Tube and tuberculin skin tests among immigrant children being screened for latent tuberculosis infection. Ped. Infect. Dis. 34, 35.
 47. Diel, R., Loddenkemper, R., Niemann, S., Meywald-Walter, K., and Nienhaus, A. (2011) Negative and positive predictive value of a whole-blood interferon- γ release assay for developing active tuberculosis. Am. J. Respir. Crit. Care Med. 183, 88.
 48. Machingadaize, S. et al. (2012) Predictive value of recent QuantiFERON conversion for tuberculosis disease in adolescents. Am. J. Respir. Crit. Care Med. 186, 1051.
 49. Riazi, S. et al. (2012) Rapid diagnosis of *Mycobacterium tuberculosis* infection in children using interferon-gamma release assays (IGRAs). Allergy Asthma Proc. 33, 217.
 50. Lighter-Fisher, J. and Surette, A-M. (2012) Performance of an interferon-gamma release assay to diagnose latent tuberculosis infection during pregnancy. Obstet. Gynecol. 119, 1088.
 51. Mathud, J.S. et al. (2014) Pregnancy differentially impacts performance of latent tuberculosis diagnostics in a high-burden setting. PLoS ONE 9, e92308.
 52. Hoffman, M. and Ravn, P. (2010) The use of interferon-gamma release assays in HIV-positive individuals. Eur. Infect. Dis. 4, 23.
 53. Cheallagh, C.N. et al. (2013) Interferon gamma release assays for the diagnosis of latent TB infection in HIV-infected individuals in a low TB burden country. PLoS ONE 8, e53330.

-
54. Ramos, J. M. et al. (2012) Contribution of interferon gamma release assays testing to the diagnosis of latent tuberculosis infection in HIV-infected patients: A comparison of QuantiFERON-TB gold in tube, T-SPOT.TB and tuberculin skin test. *BMC Infect. Dis.* 12, 169.
 55. Wolf, T. et al. (2013) Tuberculosis skin test, but not interferon- γ releasing assays is affected by BCG vaccination in HIV patients. *J. Infect.* 66, 376.
 56. Hsia, E.C. et al. (2012) Interferon- γ release assay versus tuberculin skin test prior to treatment with golimumab, a human anti-tumor necrosis factor antibody, in patients with rheumatoid arthritis, psoriatic arthritis, or ankylosing spondylitis. *Arthritis Rheum.* 64, 2068.
 57. Garcovich, S. et al. (2012) Clinical applicability of QuantiFERON-TB-Gold testing in psoriasis patients during long-term anti-TNF-alpha treatment: a prospective, observational study. *J. Eur. Acad. Dermatol. Ven.* 26, 1572.
 58. Kwakernaak, A.J. et al. (2011) A comparison of an interferon-gamma release assay and tuberculin skin test in refractory inflammatory disease patients screened for latent tuberculosis prior to the initiation of a first tumor necrosis factor α inhibitor. *Clin. Rheumatol.* 30, 505.
 59. Vinton, P. et al. (2009) Comparison of QuantiFERON-TB Gold In-Tube test and tuberculin skin test for identification of latent *Mycobacterium tuberculosis* infection in healthcare staff and association between positive test results and known risk factors for infection. *Infect. Control Hosp. Epidemiol.* 30, 215.
 60. de Perio, M.A., Tsevat, J., Roselle, G.A., Kralovic, S.M., and Eckman, M.H. (2009) Cost-effectiveness of interferon gamma release assays vs tuberculin skin tests in health care workers. *Arch. Intern. Med.* 169, 179.

-
61. Nienhaus, A. et al. (2008) Evaluation of the interferon- γ release assay in healthcare workers. *Int. Arch. Occup. Environ. Health* 81, 295.
 62. Nienhaus, A. et al. (2011) Systematic review of cost and cost-effectiveness of different TB-screening strategies. *BMC Health Serv. Res.* 11, 247.
 63. Centers for Disease Control and Prevention (2010) Updated guidelines for using interferon-gamma release assays to detect *Mycobacterium tuberculosis* infection — United States, 2010. *MMWR Recomm. Rep.* 59 (RR-5), 1.
 64. Dorman, S.E. et al. (2014) Interferon- γ release assays and tuberculin skin testing for diagnosis of latent tuberculosis infection in healthcare workers in the United States. *Am. J. Respir. Crit. Care Med.* 189, 77.
 65. Fong, K.S. et al. (2012) Challenges of interferon-gamma release assay conversions in serial testing of health care workers in a tuberculosis control program. *Chest* 142, 55.
 66. Thanassi, W. et al. (2012) Delineating a retesting zone using receiver operating characteristic analysis on serial QuantiFERON tuberculosis test results in US healthcare workers. *Pulm. Med.* doi: 10.1155/2012/291294. Epub.
 67. Behrman, A. et al. (2013) Protecting Health Care Workers from Tuberculosis, 2013: ACOEM Medical Center Occupational Health Section Task Force on Tuberculosis and Health Care Workers. *J. Occup. Environ. Med.* 55, 985.
 68. Nienhaus, A., Ringshausen, F.C., Costa, J.T, Schablon, A., and Tripodi, D. (2013) IFN- γ release assay versus tuberculin skin test for monitoring TB infection in healthcare workers. *Expert Rev. Anti Infect. Ther.* 11, 37.
 69. Arend, S.M. et al. (2007) Comparison of two interferon-gamma assays and tuberculin skin test for tracing TB contact. *Amer. J. Respir. Crit. Care Med.* 175, 618.

-
70. Mandalakas, A.M., Detjen, A.K., Hesselning, A.C., Benedetti, A., and Menzies, D. (2011) Interferon-gamma release assays and childhood tuberculosis: systematic review and meta-analysis. *Int. J. Tuberc. Lung Dis.* 15, 1018.
 71. Grinsdale, J.A., Ho, C.S., Banouvong, H., Kwamura, L.M. (2011) Programmatic impact of using QuantiFERON-TB Gold in routine contact investigation activities. *Int. J. Tuberc. Lung Dis.* 15, 1614.
 72. Walsh, M.C. et al. (2011) Sensitivity of interferon- γ release assays is not compromised in tuberculosis patients with diabetes. *Int. J. Tuberc. Lung Dis.* 15, 179.
 73. Faurholt-Jespén, D. et al. (2014) Diabetes is associated with lower tuberculosis antigen-specific interferon gamma release in Tanzanian tuberculosis patients and non-tuberculosis controls. *Scand. J. Infect. Dis.* 46, 384.
 74. Painter, J.A. et al. (2013) Tuberculosis screening by tuberculosis skin test or QuantiFERON-TB Gold In-Tube Assay among an immigrant population with a high prevalence of tuberculosis and BCG vaccination. *PLoS ONE* 8, e82727.
 75. Rogerson, T.E. et al. (2012) Tests for latent tuberculosis in people with ESRD: a systematic review. *Amer. J. Kidney Dis.* 61, 33.
 76. Pareek, M. et al. (2013) Community-based evaluation of immigrant tuberculosis screening using interferon γ release assays and tuberculin skin testing: observational study and economic analysis. *Thorax.* 68, 230.
 77. CDC, Tuberculosis — United States, 2018.
https://www.cdc.gov/mmwr/volumes/68/wr/mm6811a2.htm?s_cid=mm6811a2_w
Accessed 22 March 2019.

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© 2019 QIAGEN, sva prava zadržana.