

User-Developed Protocol:

Concentration of RNA in 96-well plates using the MinElute[®] 96 UF PCR Purification Kit

This procedure has been adapted by customers and is for concentration of RNA in 96-well plates using the MinElute 96 UF PCR Purification Kit. **The procedure has not been thoroughly tested and optimized by QIAGEN.**

IMPORTANT: Please read the *MinElute[®] 96 UF PCR Purification Handbook*, paying careful attention to the “Safety Information” and “Important Notes” sections, before beginning this procedure.

Equipment and reagents to be supplied by the user

- Multichannel pipet and pipet tips (pipet tips with aerosol barriers for preventing cross-contamination are recommended)
- Disposable gloves
- Vacuum manifold (e.g., QIAvac Multiwell, cat. no. 9014579)
- Vacuum pump capable of producing a vacuum of –800 to –900 mbar (e.g., QIAGEN Vacuum Pump)
- Vacuum Regulator (QIAGEN cat. no. 19530; for easy monitoring of vacuum pressures and easy releasing of vacuum).

Things to do before starting

- If using a microplate shaker for RNA elution in step 9, calibrate the shaker using the following steps:
 - Use a 96-well polystyrene microplate with 300 µl round-bottom wells.
 - Add 200 µl of a colored aqueous solution (e.g., bromophenol blue), to 2 wells.
 - Place the 96-well plate on the microplate shaker.
 - Set the speed to the lowest level and slowly increase the speed of the shaker. Ensure that the plate is fixed securely on top of the shaker.
 - The recommended shaking speed for elution is the maximum speed at which no liquid is splashed out of the wells.

Important points before starting

- MinElute 96 UF PCR Purification Plates are not certified as RNase free. However, no RNA degradation has been detected after concentration using this protocol.

Procedure

1. **Prepare the vacuum manifold according to the supplier's instructions.**
Place a waste tray inside the base of the manifold.
2. **Place the MinElute UF PCR Purification Plate on top of the vacuum manifold.**
3. **Pipet the RNA samples onto the MinElute 96 UF PCR Purification Plate.**
Note: Processing sample volumes larger than 150 µl may lead to increased processing time.
4. **Apply a vacuum and maintain at –800 mbar for 10 minutes or until the wells are completely dry. Switch off vacuum source.**
5. **Carefully remove the MinElute 96 UF PCR Purification Plate from the vacuum manifold.**
6. **Carefully tap the MinElute 96 UF PCR Purification Plate on a stack of clean absorbent paper to remove any liquid that might remain on the bottom of the plate.**
7. **Add 20 µl elution buffer (supplied with the original RNA isolation kit) to each well.**
8. **Elute RNA according to step 8a or 8b.**
- 8a. **Shake the MinElute 96 UF PCR Purification Plate on a microplate shaker for 2 min at the recommended speed (see “Things to do before starting”).**
Note: Ensure that the MinElute 96 UF PCR Purification Plate is fixed securely on top of the shaker.
- 8b. **Alternatively, purified RNA may be dissolved by pipetting samples up and down 20 times.**
9. **Recover the concentrated RNA by pipetting the eluate out of each well.**
For easier recovery of the eluates, the plate can be held at a slight angle.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

QIAGEN handbooks can be requested from QIAGEN Technical Service or your local QIAGEN distributor. Selected handbooks can be downloaded from www.qiagen.com/literature/handbooks/default.aspx. Material safety data sheets (MSDS) for any QIAGEN product can be downloaded from www.qiagen.com/ts/msds.asp.

Trademarks: QIAGEN[®], MinElute[®] (QIAGEN Group).
RN03 © 2006–2010 QIAGEN, all rights reserved.