

QIAcube | Protocol Sheet

General Information (April 2010)

Application	Cleanup
Kit	QIAquick® PCR Purification Kit
Sample material	Amplification reactions
Short protocol name	Large-volume samples
Version	1
Full protocol name	Purification of PCR products from 100–200 µl PCR samples
Editable parameters	Fill-up volume: 0–100 µl in increments of 5 µl; default 50 µl Elution volume: 30–100 µl in increments of 10 µl; default 50 µl
Required QIAcube® software versions	Firmware version FIW-50-001-J_FW_MB.hex and PLC program version FIW-50-002-G_PLC_MB.prs or higher; available from the QIAcube Web Portal

Shaker

Material	100–200 µl PCR samples; for volumes <200 µl, make sure to enter the volume of Buffer EB that the QIAcube should use to adjust the sample volumes to 200 µl before starting the protocol run
Vessel	2 ml safe-lock microcentrifuge tube*
Adapter	Shaker adapter for 2 ml microcentrifuge tubes (marked with "2")

* Sample Tubes RB, 2 ml (cat. no. 990381); see www.qiagen.com/MyQIAcube.

Disposable Tips


Disposable Filter-Tips, 1000 µl
Disposable Filter-Tips, 200 µl

Reagent Bottle Rack

Rack labeling strip	QIAquick
----------------------------	----------

Position	Reagent
1	Buffer PB
2	–
3	–
4	–
5	Buffer PE
6	Buffer EB

Rotor Adapter



Position	Labware	Lid position
1	QIAquick spin column	L1
2	–	–
3	1.5 ml collection tube†	L3

† Sarstedt, Micro tube 1.5 ml Safety Cap (see www.sarstedt.com).



QIAcube | Protocol Sheet

Microcentrifuge Tube Slots

	Position		
	A	B	C
Content			
Vessel			

Number of samples	Volume of reagent required for the indicated number of samples (μ l)		
	A	B	C
2			
3			
4			
5			
6			
7			
8			
9			
10			
12			

Comments

Before starting the protocol, be sure to enter the volume of Buffer EB that the QIAcube should use to make the sample volumes up to 200 μ l.

Press "Edit" and then select "Fill-up volume". For example, for sample volumes of 200 μ l enter "0 μ l" as the fill-up volume; for sample volumes of 150 μ l, enter "50 μ l" as the fill-up volume.

Fill-up volume range: 0–100 μ l; default 50 μ l.