

**User-developed  
protocol**

## **User-Developed Protocol:**

### **Whole genome amplification of 20 µg from genomic DNA using the REPLI-g<sup>®</sup> Midi Kit**

This procedure has been adapted by customers and is for whole genome amplification from genomic DNA using the REPLI-g Midi Kit with 20 µg yield. **The procedure has not been thoroughly tested and optimized by QIAGEN.**

**IMPORTANT:** Please consult the “Safety Information” and “Important Notes” sections in the *REPLI-g Mini/Midi Handbook* before beginning this procedure. For safety information on the additional chemicals mentioned in this protocol, please consult the appropriate material safety data sheets (MSDSs) available from the product supplier.

#### **Equipment and reagents to be supplied by user**

- Microcentrifuge tubes
- Microcentrifuge
- Water bath or heating block
- Vortexer
- Pipets and pipet tips
- Ice
- Nuclease-free water

#### **Important points before starting**

- This protocol is optimized for whole genome amplification from >10 ng genomic DNA template. The template DNA should be suspended in TE buffer. Smaller amounts (1–10 ng) of starting material can be used if the DNA is of sufficient quality.
- For best results, the template DNA should be >2 kb in length with some fragments >10 kb.
- REPLI-g Midi DNA Polymerase should be thawed on ice (see step 5). All other components can be thawed at room temperature (15–25°C).

## Things to do before starting

- Prepare Buffer DLB by adding 500 µl nuclease-free water to the tube; mix thoroughly and centrifuge briefly.  
**Note:** Reconstituted Buffer DLB can be stored for 6 months at –20°C. Buffer DLB is pH-labile. Avoid neutralization with CO<sub>2</sub>.
- All buffers and reagents should be vortexed before use to ensure thorough mixing.
- Set a water bath or heating block to 30°C.

## Procedure

1. **Prepare sufficient Buffer D1 (denaturation buffer) and Buffer N1 (neutralization buffer) for the total number of whole genome amplification reactions (see Tables 1 and 2).**  
**Note:** Buffer D1 and Buffer N1 should not be stored longer than 3 months.

**Table 1. Preparation of Buffer D1**

| Component                             | Volume       |
|---------------------------------------|--------------|
| Reconstituted Buffer DLB <sup>†</sup> | 9 µl         |
| Nuclease-free water                   | 32 µl        |
| <b>Total volume</b>                   | <b>41 µl</b> |

\* Volumes given are suitable for up to 15 reactions.

**Table 2. Preparation of Buffer N1**

| Component           | Volume       |
|---------------------|--------------|
| Stop solution       | 12 µl        |
| Nuclease-free water | 68 µl        |
| <b>Total volume</b> | <b>80 µl</b> |

\* Volumes given are suitable for up to 15 reactions.

2. **Place 2.5 µl template DNA into a microcentrifuge tube.**  
The amount of template DNA should be >10 ng.  
A DNA control reaction can be set up using 10 ng (1 µl) control genomic DNA (e.g., REPLI-g Human Control Kit, cat. no. 150090).
3. **Add 2.5 µl Buffer D1 to the DNA. Mix by vortexing and centrifuge briefly.**
4. **Incubate the samples at room temperature (15–25°C) for 3 min.**
5. **Add 5 µl Buffer N1 to the samples. Mix by vortexing and centrifuge briefly.**

- 6. Thaw REPLI-g Midi DNA Polymerase on ice. Thaw all other components at room temperature, vortex, and centrifuge briefly.**

The REPLI-g Midi Reaction Buffer may form a precipitate after thawing. The precipitate will dissolve by vortexing for 10 s.

- 7. Prepare a master mix on ice according to Table 3. Mix and centrifuge briefly.**

**Important:** Add the master mix components in the order listed in Table 3. After addition of water and REPLI-g Midi Reaction Buffer, briefly vortex and centrifuge the mixture before addition of REPLI-g Midi DNA Polymerase. The master mix should be kept on ice and used immediately upon addition of the REPLI-g Midi DNA Polymerase.

**Table 3. Preparation of Master Mix**

| <b>Component</b>             | <b>Volume/reaction</b>      |
|------------------------------|-----------------------------|
| REPLI-g Midi Reaction Buffer | 14.5 $\mu$ l                |
| REPLI-g Midi DNA Polymerase  | 0.5 $\mu$ l                 |
| <b>Total volume</b>          | <b>15 <math>\mu</math>l</b> |

- 8. Add 15  $\mu$ l of the master mix to 10  $\mu$ l denatured DNA (step 5).**

- 9. Incubate at 30°C for 8–16 h.**

Maximum DNA yield is achieved using an incubation time of 16 h.

After incubation at 30°C, heat the water bath or heating block up to 65°C if the same water bath or heating block will be used in step 11.

- 10. Inactivate REPLI-g Midi DNA Polymerase by heating the sample at 65°C for 3 min.**

- 11. Store amplified DNA at 4°C for short-term storage or –20°C for long-term storage.**

DNA amplified using the REPLI-g Midi Kit should be treated as genomic DNA with minimal freeze-thaw cycles. Storage of nucleic acids at low concentration over a long period of time may result in acid hydrolysis. We therefore recommend storage of nucleic acids at a concentration of at least 100 ng/ $\mu$ l.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at [www.qiagen.com](http://www.qiagen.com) or can be requested from QIAGEN Technical Services or your local distributor.

Material safety data sheets (MSDS) for any QIAGEN product can be downloaded from [www.qiagen.com/Support/MSDS.aspx](http://www.qiagen.com/Support/MSDS.aspx).

Trademarks: QIAGEN<sup>®</sup>, REPLI-g<sup>®</sup> (QIAGEN Group).

© 2007–2011 QIAGEN, all rights reserved.