

RNeasy[®] PowerPlant[®] Kit

The RNeasy PowerPlant Kit can be stored at room temperature (15–25°C) until the expiry date printed on the box label.

Further information

- Safety Data Sheets: www.qiagen.com/safety
- Technical assistance: support.qiagen.com

Notes before starting

- Solution MBL must be warmed at 55°C for 5–10 minutes to dissolve precipitates prior to use. Use Solution MBL while still warm.
 - Add the appropriate amount of β -mercaptoethanol (β ME) to Solution MBL to produce a final concentration of 10 μ l/ml of MBL/ β ME. For each prep, 600 μ l of MBL/ β ME will be needed. Alternatively, add 594 μ l of MBL and 6 μ l of β ME directly to each bead tube.
1. Place up to 50 mg of plant sample into a 2 ml PowerBead Tube (provided).
Note: We recommended cutting samples into smaller pieces before weighing/loading.
 2. Add 600 μ l of Solution MBL/ β ME to the PowerBead Tube. You may prepare Solution MBL/ β ME in advance or add 594 μ l of MBL and 6 μ l of β ME directly to each bead tube.
Note: If sample is high in phenolics and you are using the Phenolic Separation Solution, reduce Solution MBL/ β ME to 550 μ l and add 50 μ l of the Phenolic Separation Solution.
 3. For the highest yields of RNA, a high powered bead beater is recommended.
 - A. On the PowerLyzer[®] 24 Instrument, we recommend a starting setting of 1 cycle at 4200 rpm for 45 s for leaf tissue and seeds.
 - B. You may also homogenize using a Vortex Genie[®] 2 and Vortex Adapter (cat. no. 13000-V1-24) for soft leaf tissue only. Set the vortex on full speed for 10 min.
 4. Centrifuge at 13,000 $\times g$ for 2 min at room temperature.
 5. Transfer all the supernatant to a clean 2 ml collection tube (provided). Expect between 500 to 600 μ l of lysate.
 6. Add 150 μ l of Solution IRS and vortex briefly to mix. Incubate at 4°C for 5 min.

- Note:** For plant samples that still contain PCR inhibitors after RNA purification, try adding up to 200 μ l of Solution IRS.
7. Centrifuge the tubes at 13,000 \times g for 2 min. Avoiding the pellet, transfer the supernatant to a clean 2 ml collection tube. Transfer no more than 650 μ l at this step.
 8. Add 650 μ l of Solution PM3 and 650 μ l of Solution PM4. Vortex briefly to mix.
Note: To purify small RNAs, such as microRNAs and siRNAs, transfer the lysate to a larger tube to accommodate the higher volume (2.6 ml), and add an additional 650 μ l of 100% ethanol. You will need to supply the 100% ethanol.
 9. Load 650 μ l of supernatant onto an MB RNA Column and centrifuge at 13,000 \times g for 1 min. Discard the flow through and place the Spin Column back into the 2 ml collection tube. Repeat until all the supernatant has been loaded onto the Spin Column.
Note: A total of three loads for each sample processed are required (four loads if an additional volume of 100% ethanol is added for the miRNA and siRNA protocol).
 10. Shake to mix Solution PM5. Add 600 μ l of Solution PM5 to the MB RNA Column and centrifuge at 13,000 \times g for 1 min.
Optional: If the RNase-Free DNase Set (cat. no. 79254-50) was purchased separately, it should be incorporated after step 10.
 11. Discard the flow through, place the MB RNA Spin Column back into the 2 ml collection tube, and add 600 μ l of Solution PM4. Centrifuge at 13,000 \times g for 1 min.
 12. Discard the flow through and centrifuge again at 16,000 \times g for 2 min.
 13. Place the MB RNA Spin Filter into a clean 2 ml collection tube.
 14. Add 50-100 μ l of RNase-free water (provided) to the center of the white Spin Column membrane. Incubate at room temperature for 1 min.
Note: Eluting with 100 μ l of RNase-free water will maximize RNA yield. For more concentrated RNA, a **minimum** of 50 μ l of RNase-free water can be used.
 15. Centrifuge at 13,000 \times g for 1 min. Discard the MB RNA Spin Column. The RNA is now ready for downstream applications and may be stored at -65° C to -90° C.