

# HotStarTaq<sup>®</sup> Plus DNA Polymerase

HotStarTaq Plus DNA Polymerase (cat. nos. 203601, 203603, 203605, 203607 and 203609), including buffers and reagents, should be stored immediately upon receipt at  $-30$  to  $-15^{\circ}\text{C}$  in a constant-temperature freezer.

## Further information

- *HotStarTaq Plus PCR Handbook*: [www.qiagen.com/HB-0450](http://www.qiagen.com/HB-0450)
- Safety Data Sheets: [www.qiagen.com/safety](http://www.qiagen.com/safety)
- Technical assistance: [support.qiagen.com](http://support.qiagen.com)

## Notes before starting

- HotStarTaq Plus DNA Polymerase requires a heat-activation step of 5 min at  $95^{\circ}\text{C}$  (see step 5).
- The PCR Buffer provides a final concentration of 1.5 mM  $\text{MgCl}_2$  in the reaction mix, which will give satisfactory results in most cases. However, in some cases, reactions may be improved by increasing the final  $\text{Mg}^{2+}$  concentration. If a higher  $\text{Mg}^{2+}$  concentration is required, add the appropriate volume of 25 mM  $\text{MgCl}_2$  to the reaction mix as described in the *HotStarTaq Plus PCR Handbook*.
- If required, prepare a dNTP mix containing 10 mM of each dNTP. Store this mix in aliquots at  $-20^{\circ}\text{C}$ . High-quality, PCR-grade dNTP mix (10 mM) is available from QIAGEN (cat. no. 201900).
- HotStarTaq Plus DNA Polymerase is provided with Q-Solution<sup>®</sup>, which facilitates amplification of templates that have a high degree of secondary structure or that are GC-rich by modifying the melting behavior of DNA. When using Q-Solution for the first time for a particular primer-template pair, always perform parallel reactions with and without Q-Solution.

- HotStarTaq *Plus* DNA Polymerase is provided with CoralLoad® PCR Buffer, which contains a gel-loading reagent and gel-tracking dyes.
  - CoralLoad PCR Buffer must not be used in capillary sequencers.
  - It is not necessary to keep PCR tubes on ice as nonspecific DNA synthesis cannot occur at room temperature due to the inactive state of HotStarTaq *Plus* DNA Polymerase.
  - A No Template Control (NTC) should always be included.
1. Thaw 10x CoralLoad PCR Buffer or 10x PCR Buffer, dNTP mix, primer solutions and 25 mM MgCl<sub>2</sub> (if required). Mix the solutions thoroughly before use to avoid localized differences in salt concentration.
  2. Prepare a reaction mix according to Table 1. The reaction mix typically contains all the components needed for PCR except the template DNA. Prepare a volume of reaction mix 10% greater than that required for the total number of PCR assays to be performed.

**Table 1. Reaction setup using HotStarTaq *Plus* DNA Polymerase**

Component	Volume/reaction	Final concentration
<b>Reaction mix</b>		
10x PCR Buffer* or <b>Optional:</b> 10x CoralLoad PCR Buffer*	10 µl	1x
dNTP mix (10 mM of each)	2 µl	200 µM of each dNTP
Primer A	Variable	0.1–0.5 µM
Primer B	Variable	0.1–0.5 µM
HotStarTaq <i>Plus</i> DNA Polymerase	0.5 µl	2.5 units/reaction
<b>Optional:</b> 5x Q-Solution†	20 µl	1x
<b>Template DNA</b> (added at step 4)	Variable	≤1 µg/reaction
<b>Total reaction volume</b>	100 µl‡	

\* Contains 15 mM MgCl<sub>2</sub>.

† For templates with GC-rich regions or complex secondary structure.

‡ If using different reaction volumes, adjust the volume of each component accordingly.

3. Mix the reaction mix gently but thoroughly, for example, by pipetting up and down a few times. Dispense appropriate volumes into PCR tubes.

- Add template DNA ( $\leq 1 \mu\text{g}/100 \mu\text{l}$  reaction) to the individual PCR tubes containing the reaction mix. For RT-PCR, add an aliquot from the reverse transcriptase reaction. This should not exceed 10% of the final PCR volume.
- Program the thermal cycler according to the manufacturer's instructions.

**Note:** Each PCR program must start with an initial heat-activation step at  $95^\circ\text{C}$  for 5 min. Do not exceed the 5 min activation time. A typical PCR cycling program is outlined in Table 2. For maximum yield and specificity, temperatures and cycling times should be optimized for each new template target or primer pair.

**Table 2. Optimized cycling conditions**

Step	Time	Temperature	Comment
<b>Initial heat activation</b>	5 min	$95^\circ\text{C}$	Activates HotStarTaq <i>Plus</i> DNA Polymerase.
<b>3-step cycling:</b>			
Denaturation	0.5–1 min	$94^\circ\text{C}$	
Annealing	0.5–1 min	$50\text{--}68^\circ\text{C}$	Approximately $5^\circ\text{C}$ below $T_m$ of primers.
Extension	1 min	$72^\circ\text{C}$	For PCR products longer than 1 kb, use an extension time of approximately 1 min per kb DNA.
Number of cycles	25–35		
<b>Final extension</b>	10 min	$72^\circ\text{C}$	

- Place the PCR tubes in the thermal cycler and start the cycling program.
 

**Note:** After amplification, samples can be stored overnight at  $2\text{--}8^\circ\text{C}$ , or at  $-20^\circ\text{C}$  for longer storage.
- When using CoralLoad PCR Buffer, the PCR products can be directly loaded onto an agarose gel without prior addition of a PCR loading buffer and gel-tracking dyes. Refer to Table 3 to identify the gel-tracking dyes present in CoralLoad PCR Buffer according to migration distance in different percentage agarose gels.
 

**Note:** Due to the high viscosity of the solution, apply the solution slowly into the wells of the agarose gel.

**Table 3. Migration distance of gel-tracking dyes in CoralLoad PCR Buffer**

<b>% TAE (TBE) agarose gel</b>	<b>Red dye</b>	<b>Orange dye</b>
0.8	500 (270) bp	~80 (<10) bp
1.0	300 (220) bp	~40 (<10) bp
1.5	250 (120) bp	~20 (<10) bp
2.0	100 (110) bp	<10 (<10) bp
3.0	50 (100) bp	<10 (<10) bp



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN<sup>®</sup>, Sample to Insight<sup>®</sup>, CoralLoad<sup>®</sup>, HotStarTaq<sup>®</sup>, Q-Solution<sup>®</sup> (QIAGEN Group). 1101227 03/2016 HB-0655-002  
© 2016 QIAGEN, all rights reserved.