September 2020

QIAcube[®] Connect Safety Instructions and Quick Start Guide

For use with QIASphere Base

Contents

Introduction	4
Safety Information.....	5
Proper use	5
Electrical safety	8
Environment.....	9
Biological safety.....	10
Chemicals	11
Toxic fumes	11
Waste disposal	12
Mechanical hazards.....	12
Centrifuge	13
Heat hazard	14
Maintenance safety	15
Radiation safety	16
Translations of Warnings and Cautions	17
Warnings.....	17
Cautions.....	35
Symbols on the QIAcube Connect	42
Installation Procedures.....	47
Site requirements.....	47
Setting up the QIAcube Connect instrument	47
To set up the QIAcube Connect instrument	47

Connecting the QIAcube Connect via Wi-Fi	49
Connecting the QIAcube Connect via LAN cable	49
Connecting the QIASphere Base.....	50
Connecting the QIASphere Base.....	50
To define WLAN/Wi-Fi (WPAv2 only) setting.....	51
To define LAN (Static IP and DHCP) setting	51
Revision History	52

Introduction

The QIAcube® Connect is designed to perform fully-automated purification of nucleic acids and proteins in molecular biology applications. The system is intended for use by professional users trained in molecular biological techniques and the operation of the QIAcube Connect.

The QIAcube Connect is intended to be used only in combination with QIAGEN kits indicated for use with the QIAcube Connect for the applications described in the kit handbooks.

Safety Information

Before using the QIAcube Connect, it is essential that you read this guide carefully and pay attention to the safety information. The instructions and safety information in this guide must be followed to ensure safe operation of the instrument and to maintain the instrument in a safe condition.

<p>WARNING</p> 	<p>The term WARNING is used to inform you about situations that could result in personal injury to you or other persons.</p> <p>Details about these circumstances are given in a box like this.</p>
---	---

<p>CAUTION</p> 	<p>The term CAUTION is used to inform you about situations that could result in damage to the instruments or other equipment.</p> <p>Details about these circumstances are given in a box like this.</p>
---	--

The advice given in this guide is intended to supplement, not supersede, the normal safety requirements prevailing in the user's country.

Proper use

<p>WARNING</p> 	<p>Risk of personal injury and material damage (W1)</p> <p>Improper use of the QIAcube Connect may cause personal injuries or damage to the instrument. The QIAcube Connect must only be operated by qualified personnel who have been appropriately trained. Servicing of the QIAcube Connect must only be performed by a QIAGEN Field Service specialist.</p>
---	--

QIAGEN charges for repairs that are required due to incorrect maintenance.

WARNING 	<p>Risk of personal injury and material damage (W2)</p> <p>The QIAcube Connect is too heavy to be lifted by one person. To avoid personal injury or damage to the instrument, do not lift the instrument alone.</p>
---	--

WARNING 	<p>Risk of personal injury and material damage (W3)</p> <p>Do not attempt to move the QIAcube Connect during operation.</p>
---	--

In case of emergency, power OFF the QIAcube Connect at the power switch located in front of the instrument and unplug the power cord from the power outlet.

CAUTION 	<p>Damage to the instrument (C1)</p> <p>Avoid spilling water or chemicals onto the QIAcube Connect. Damage caused by water or chemical spillage will void your warranty.</p>
---	---

In case of emergency, power OFF the QIAcube Connect at the power switch located in front of the instrument and unplug the power cord from the power outlet.

CAUTION 	<p>Damage to the instrument (C2)</p> <p>Only use QIAGEN spin columns and QIAcube Connect specific consumables with the QIAcube Connect. Damage caused by use of other types of spin columns or chemistries will void your warranty.</p>
---	--

WARNING 	<p>Risk of personal injury and material damage (W4)</p> <p>Do not use damaged rotor adapters. The rotor adapters can only be used once. High g forces exerted in the centrifuge can cause damage to used rotor adapters.</p>
---	---

<p>CAUTION</p> 	<p>Risk of material damage (C3)</p> <p>Empty the tip disposal container prior to use to prevent a tip jam in the waste drawer. Failure to empty the waste container may block the robotic arm that could cause run failure or instrument damage.</p>
<p>WARNING</p> 	<p>Risk of personal injury and material damage (W5)</p> <p>To prevent plastic crash, load the tubes properly. After a plastic crash, sharp plastic particles could be inside the centrifuge. Be careful when handling items inside the centrifuge.</p>
<p>CAUTION</p> 	<p>Damage to the instrument (C4)</p> <p>Only use the correct volume of liquids. Exceeding the recommended volume of liquids may damage the centrifuge rotor or instrument.</p>
<p>WARNING</p> 	<p>Risk of fire or explosion (W6)</p> <p>When using ethanol or ethanol-based liquids on the QIAcube Connect, handle such liquids carefully and in accordance with the required safety regulations. If liquid has been spilled, wipe it off and leave the QIAcube Connect hood open to allow flammable vapors to disperse.</p>
<p>WARNING</p> 	<p>Risk of explosion (W7)</p> <p>The QIAcube Connect is intended for use with reagents and substances supplied with QIAGEN kits or other than outlined in respective Information for use. Use of other reagents and substances may lead to fire or explosion.</p>

If there are hazardous material inside the QIAcube Connect, you are responsible to perform the appropriate decontamination.

Note: Do not place items on top of the QIAcube Connect hoods.

<p>CAUTION</p> 	<p>Damage to the instrument (C5)</p> <p>Do not lean against the touchscreen when it is pulled out.</p>
---	---

Electrical safety

Note: Disconnect the line power cord from the power outlet before servicing.

<p>WARNING</p> 	<p>Electrical hazard (W8)</p> <p>Any interruption of the protective conductor (earth/ground lead) inside or outside the instrument or disconnection of the protective conductor terminal is likely to make the instrument dangerous. Intentional interruption is prohibited.</p> <p>Lethal voltages inside the instrument</p> <p>When the instrument is connected to line power, terminals may be live and opening covers or removing parts is likely to expose live parts.</p>
---	---

<p>WARNING</p> 	<p>Damage to electronics (W9)</p> <p>Before powering ON the instrument, make sure that the correct supply voltage is used. Use of incorrect supply voltage may damage the electronics. To check the recommended supply voltage, refer to the specifications indicated in the type plate of the instrument.</p>
---	---

<p>WARNING</p> 	<p>Risk of electric shock (W10)</p> <p>Do not open any panels on the QIAcube Connect.</p> <p>Risk of personal injury and material damage</p> <p>Only perform maintenance that is specifically described in this user manual.</p>
---	--

To ensure satisfactory and safe operation of the QIAcube Connect, follow these guidelines:

- The line power cord must be connected to a line power outlet that has a protective conductor (earth/ground).
- Do not adjust or replace internal parts of the instrument.
- Do not operate the instrument with any covers or parts removed.
- If liquid has spilled inside the instrument, power OFF the instrument, disconnect it from the power outlet, and contact QIAGEN Technical Services.

If the instrument becomes electrically unsafe, prevent other personnel from operating it, and contact QIAGEN Technical Services.

The instrument may be electrically unsafe when:

- It or the line power cord appears to be damaged.
- It has been stored under unfavorable conditions for a prolonged period.
- It has been subjected to severe transport stresses.
- Liquids come in contact directly with electrical components of the QIAcube Connect.

Environment

Operating Conditions

<p>WARNING</p> 	<p>Explosive atmosphere (W11)</p> <p>The QIAcube Connect is not designed for use in an explosive atmosphere.</p>
---	---

CAUTION 	Damage to the instrument (C6)
	Direct sunlight may bleach parts of the instrument and cause damage to plastic parts. The QIAcube Connect must be located out of direct sunlight.

Biological safety

Specimens and reagents containing materials from humans should be treated as potentially infectious. Use safe laboratory procedures as outlined in publications such as Biosafety in Microbiological and Biomedical Laboratories, HHS (www.cdc.gov/BMBL).

Samples may contain infectious agents. You should be aware of the health hazard presented by such agents and should use, store, and dispose of such samples according to the required safety regulations.

WARNING 	Samples containing infectious agents (W12)
	Some samples used with this instrument may contain infectious agents. Handle such samples with the greatest of care and in accordance with the required safety regulations. Always wear safety glasses, 2 pairs of gloves, and a lab coat. The responsible body (e.g., laboratory manager) must take the necessary precautions to ensure that the surrounding workplace is safe, and that the instrument operators are suitably trained and not exposed to hazardous levels of infectious agents as defined in the applicable Material Safety Data Sheets (MSDSs) or OSHA,* ACGIH, [†] or COSHH [‡] documents. Venting for fumes and disposal of wastes must be in accordance with all national, state, and local health and safety regulations and laws.

* OSHA: Occupational Safety and Health Administration (United States of America).

[†] ACGIH: American Conference of Government Industrial Hygienists (United States of America).

[‡] COSHH: Control of Substances Hazardous to Health (United Kingdom).

Chemicals

<p>WARNING</p> 	<p>Hazardous chemicals (W13)</p> <p>Some chemicals used with this instrument may be hazardous or may become hazardous after completion of the protocol run. Always wear safety glasses, gloves, and a lab coat. The responsible body (e.g., laboratory manager) must take the necessary precautions to ensure that the surrounding workplace is safe and that the instrument operators are not exposed to hazardous levels of toxic substances (chemical or biological) as defined in the applicable Material Safety Data Sheets (MSDSs) or OSHA,* ACGIH,[†] or COSHH[‡] documents. Venting for fumes and disposal of wastes must be in accordance with all national, state, and local health and safety regulations and laws.</p>
---	---

* OSHA: Occupational Safety and Health Administration (United States of America).

[†] ACGIH: American Conference of Government Industrial Hygienists (United States of America).

[‡] COSHH: Control of Substances Hazardous to Health (United Kingdom).

Toxic fumes

If working with volatile solvents or toxic substances, you must provide an efficient laboratory ventilation system to remove vapors that may be produced.

<p>WARNING</p> 	<p>Toxic fumes (W14)</p> <p>Do not use bleach to clean or disinfect the QIAcube Connect. Bleach in contact with salts from the buffers can produce toxic fumes.</p>
---	--

<p>WARNING</p> 	<p>Toxic fumes (W15)</p> <p>Do not use bleach to disinfect used labware. Bleach in contact with salts from the buffers used can produce toxic fumes.</p>
---	---

Waste disposal

Used labware, such as sample tubes, QIAGEN spin columns, filter-tips, buffer bottle and enzyme tubes, or rotor adapters, may contain hazardous chemicals or infectious agents from the purification process. These hazardous wastes must be collected and disposed of properly according to local safety regulations.

WARNING 	Hazardous chemicals and infectious agents (W16) The waste may contain toxic material and must be disposed of properly. Refer to your local safety regulations for proper disposal procedures.
---	---

Mechanical hazards

The hood of the QIAcube Connect must remain closed during operation of the instrument. Only open the hood when instructed to do so by the instruction for use.

While loading the worktable, always stand clear of the instrument. Do not lean on the worktable when the robotic arm of the instrument is moving to reach loading position with its lid open. Wait until the robotic arm completed its movements before you start to load or unload.

WARNING 	Moving parts (W17) Avoid contact with moving parts during operation of the QIAcube Connect. Under no circumstances should you place your hands under the robotic arm when it is lowering. Do not attempt to move any tip racks or tubes whilst the instrument is operating.
---	---

<p>WARNING</p> 	<p>Moving parts (W18)</p> <p>To avoid contact with moving parts during operation of the QIAcube Connect, the instrument must be operated with the hood closed. If the hood sensor or lock is not functioning correctly, contact QIAGEN Technical Services.</p>
---	--

Centrifuge

Make sure that the rotor and buckets are installed correctly. All buckets must be mounted before starting a protocol run, regardless of the number of samples to be processed. If the rotor or buckets show signs of mechanical damage or corrosion, do not use the QIAcube Connect; contact QIAGEN Technical Services.

<p>CAUTION</p> 	<p>Damage to the instrument (C7)</p> <p>The QIAcube Connect must not be used if the centrifuge lid is broken, or if the lid lock is damaged. Make sure that no loose material is inside the centrifuge during operation.</p> <p>Make sure that the rotor is installed correctly and that all buckets are properly mounted, regardless of the number of samples to be processed. Load the rotor only as instructed by the software.</p> <p>Only use rotors, buckets, and consumables designed for use with the QIAcube Connect. Damage caused by use of other consumables will void your warranty.</p> <p>We recommend replacing the centrifuge rotor and buckets after 20,000 cycles, which is equivalent to 9 years of usage with two runs per day for 220 days each year. For more information contact QIAGEN Technical Services.</p>
---	---

In case of breakdown, caused by power failure, the centrifuge lid can be manually opened to remove the samples.

WARNING 	Moving parts In case of breakdown caused by power failure, remove the power cord and wait 10 minutes before attempting to manually open the centrifuge lid.	(W19)
---	---	-------

CAUTION 	Damage to the instrument After a power failure, do not move the z-module (robotic arm) manually in front of the instrument. Damage may occur if the QIAcube Connect hood is closed and collides with the z-module.	(C8)
---	--	------

WARNING 	Risk of personal injury and material damage Raise the centrifuge lid carefully. The lid is heavy and may cause injury if it falls.	(W20)
---	--	-------

CAUTION 	Risk of overheating To ensure proper ventilation, maintain a minimum clearance of 10 cm at the sides and rear of the QIAcube Connect. Slits and openings that ensure the ventilation of the QIAcube Connect must not be covered.	(C9)
---	--	------

Heat hazard

The QIAcube Connect worktable contains a heated shaker.

WARNING 	Hot surface The shaker can reach temperatures of up to 70°C (158°F). Avoid touching it when it is hot. Carefully remove the samples after a run.	(W21)
---	--	-------

Maintenance safety

WARNING/ CAUTION 	Risk of personal injury and material damage (W22) Only perform maintenance that is specifically described in this user manual.
WARNING 	Risk of explosion (W23) When cleaning the QIAcube Connect with alcohol-based disinfectant, leave the QIAcube Connect hood open to allow flammable vapors to disperse. Only clean the QIAcube Connect when worktable components have cooled down.
WARNING 	Risk of fire (W24) Do not allow cleaning fluid or decontamination agents to come into contact with the electrical parts of the QIAcube Connect.
WARNING 	Risk of personal injury and material damage (W25) To prevent the rotor nuts from loosening during operation of the centrifuge, securely tighten the nuts using the rotor key supplied with QIAcube Connect.
WARNING 	UV irradiation hazard (W26) A mechanical lock ensures that the hood must be closed for operation of the UV LED. If the hood sensor or lock is not functioning correctly, contact QIAGEN Technical Services.
WARNING 	Risk of personal injury and material damage (W27) Make sure that lids from spin columns and 1.5 ml microcentrifuge tubes are in the correct position and pushed all the way down to the bottom of the slots on the sides of the rotor adapter. Incorrectly positioned lids can break off during centrifugation.

<p>WARNING</p> 	<p>Risk of personal injury and material damage (W28)</p> <p>Be sure the lid is completely removed from the spin column. Spin columns with partially removed lids may not be removed properly from the rotor, causing the protocol run to crash.</p>
---	--

<p>CAUTION</p> 	<p>Damage to the instrument(s) (C10)</p> <p>Do not use bleach, solvents, or reagents containing acids, alkalis, or abrasives to clean QIAcube Connect.</p>
---	---

<p>CAUTION</p> 	<p>Damage to the instrument(s) (C11)</p> <p>Do not use spray bottles containing alcohol or disinfectant to clean surfaces of the QIAcube Connect. Spray bottles should be used only to clean items that have been removed from the worktables.</p>
---	---

Radiation safety

<p>WARNING</p> 	<p>Risk of personal injury (W29)</p> <p>Do not expose your skin to UV-C light from the UV LED lamp.</p>
--	--

<p>WARNING</p> 	<p>Risk of personal injury (W30)</p> <p>Hazard Level 2 laser light: Do not stare into the light beam when using handheld barcode scanner.</p>
---	--

Translations of Warnings and Cautions

This subsection contains translations of the warnings and cautions used in this guide. Each warning or caution has a reference number in square brackets at the top right of its box.

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WARNING 	Risk of personal injury and material damage (W5) To prevent plastic crash, load the tubes properly. After a plastic crash, sharp plastic particles could be inside the centrifuge. Be careful when handling items inside the centrifuge.
---	--

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<p>WARNING</p> 	<p>Samples containing infectious agents (W12)</p> <p>Some samples used with this instrument may contain infectious agents. Handle such samples with the greatest of care and in accordance with the required safety regulations. Always wear safety glasses, 2 pairs of gloves, and a lab coat. The responsible body (e.g., laboratory manager) must take the necessary precautions to ensure that the surrounding workplace is safe, and that the instrument operators are suitably trained and not exposed to hazardous levels of infectious agents as defined in the applicable Material Safety Data Sheets (MSDSs) or OSHA,* ACGIH,† or COSHH‡ documents. Venting for fumes and disposal of wastes must be in accordance with all national, state, and local health and safety regulations and laws.</p>
---	---

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<p>WARNING</p> 	<p>Hazardous chemicals (W13)</p> <p>Some chemicals used with this instrument may be hazardous or may become hazardous after completion of the protocol run. Always wear safety glasses, gloves, and a lab coat. The responsible body (e.g., laboratory manager) must take the necessary precautions to ensure that the surrounding workplace is safe and that the instrument operators are not exposed to hazardous levels of toxic substances (chemical or biological) as defined in the applicable Material Safety Data Sheets (MSDSs) or OSHA,* ACGIH,† or COSHH‡ documents. Venting for fumes and disposal of wastes must be in accordance with all national, state, and local health and safety regulations and laws.</p>
---	---

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CAUTION 	Damage to the instrument (C7) <p>The QIAcube Connect must not be used if the centrifuge lid is broken, or if the lid lock is damaged. Make sure that no loose material is inside the centrifuge during operation.</p> <p>Make sure that the rotor is installed correctly and that all buckets are properly mounted, regardless of the number of samples to be processed. Load the rotor only as instructed by the software.</p> <p>Only use rotors, buckets, and consumables designed for use with the QIAcube Connect. Damage caused by use of other consumables will void your warranty.</p> <p>We recommend replacing the centrifuge rotor and buckets after 20,000 cycles, which is equivalent to 9 years of usage with two runs per day for 220 days each year. For more information contact QIAGEN Technical Services.</p>
---	--

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Installation Procedures

This section provides you with important information on how to install the QIAcube Connect, QIASphere Base and tablet including the required instrument setup and network settings. For detailed unpacking and installing information, refer to the respective user manual.

Note: Before using the QIASphere Base, it is essential that you read the supplier's user manual carefully and pay attention to the safety information. The safety information for QIASphere Base can be found at <https://www.eurotech.com/en/products/iot/multi-service-iot-edge-gateways/reliagate-10-12>.

Site requirements

The QIAcube Connect must be located out of direct sunlight, away from heat sources, and away from sources of vibration and electrical interference. Refer to user manual for the operating conditions (temperature and humidity) and detailed site, power and grounding requirements.

Setting up the QIAcube Connect instrument

This section describes important actions that must be performed before operating the QIAcube Connect.

To set up the QIAcube Connect instrument

1. Remove the power cord and the quick-start guide from the foam packing material on top of the QIAcube Connect.
2. Remove the USB flash drive, rotor key, rotor nut, Allen key, S2 shaker adapter, and shaker rack plugs from the waste drawer.
3. Carefully peel off the protective film from the QIAcube Connect hood. To remove the protective film, gently pull it towards you. Ensure to close the QIAcube Connect hood after removing the protective film.

4. To remove the foam protectors above the centrifuge, gently pull the top foam protector towards you. After removal of top foam protector, gently pull the middle foam protector towards you and find the barcode reader enclosed within middle foam. Gently remove the bottom foam protector above centrifuge.
5. To remove the foam protector around the robotic arm, gently pull the foam protector towards you. After removing the protector for the robotic arm, ensure to close the QIAcube Connect hood.
6. Before turning on the instrument, plug the Wi-Fi adapter into one of the USB ports behind the touchscreen.
7. Connect the QIAcube Connect to a power outlet using the provided power cord.
8. The centrifuge rotor and buckets are pre-installed in the QIAcube Connect. In case the centrifuge rotor and buckets have been manually removed (e.g., during maintenance), follow the procedures below to install the centrifuge rotor and buckets:

Note: The rotor can be mounted in only one orientation. The pin on the rotor shaft fits into a notch on the underside of the rotor, directly underneath rotor **Position 1**.

 - 8a. Line up **Position 1** of the rotor with the pin on the rotor shaft and carefully lower the rotor onto the shaft.
 - 8b. Install the rotor nut on top of the rotor and tighten using the provided rotor key. Ensure that the rotor is securely seated. Check that all buckets are properly suspended and can swing freely.

Note: To help prevent buckets from being loaded incorrectly, the side of the rotor bucket that must face towards the rotor shaft is marked with a **gray** line.

Mount all centrifuge buckets before starting a protocol run, even if fewer than 12 samples are to be processed.
9. **Optional:** Refer to the instrument user manual if you need to change the shaker adapter.
10. To power on the QIAcube Connect instrument, close the hood and press the power button. A beeper sounds and the startup screen appears.

-
11. To log in for the first time, enter **Admin** in both **User ID** and **Password** fields. After this login, you will have Admin rights with the option to set up other users. Refer to the instrument user manual for more details on user management.
 12. From the **Configuration** menu under the **System** tab, you can modify the **Set Date, Time** and **Name Settings** fields.
 13. Follow the instructions below to connect the QIAcube Connect instrument to a network. This can only be done by a user with an Admin role.

Connecting the QIAcube Connect via Wi-Fi

Note: The QIAcube Connect supports only WPA-PSK and WPA2-PSK modes. Additionally, the SSID of the Wi-Fi hotspot must be visible. A connection to a hidden SSID is not supported.

- 13a. From the **Configuration** menu under the **Wi-Fi** tab, you can set up the Wi-Fi network.
- 13b. Scan for available networks.
- 13c. From the list, select one of the available networks.
- 13d. Enter the Wi-Fi password and connect.

Connecting the QIAcube Connect via LAN cable

- 13e. From the **Configuration** menu under the **LAN** tab, enter the local network settings to set up the LAN network.

Note: You might need the help of an IT professional to enter detailed network settings.

Connecting the QIASphere Base

The QIASphere Base must be connected to the same network as the QIACube Connect instrument and tablet. The QIASphere Base does not need to be in proximity with the QIACube Connect instrument to be connected.

Note: The name **QIASphere Base** is abbreviated as **Q-Base** in the QIACube Connect instrument software and app.

Connecting the QIASphere Base

1. Power on the QIASphere Base by plugging the power cord into an appropriate power socket.

Note: If you use a LAN cable to connect to your network, plug the cable into ethernet port ETH1 prior to powering on the QIASphere Base.

Note: Bluetooth® is automatically turned off after 10 minutes.

2. Switch on the tablet and activate Bluetooth in the settings.
3. Launch the QIACube Connect App.
4. To proceed, read and accept the "Terms and Agreement".
5. Press the **Configuration** icon () to set up the network.

Note: For successful Bluetooth communication between QIASphere Base and tablet, the devices must be within a maximum distance of 10 meters.

6. After a few seconds, the QIASphere Base name will be listed under **Available Q-Base devices for configuration of network settings**. If the QIASphere Base cannot be found, refresh the Bluetooth search by pressing the **Configuration** icon again.
7. To pair QIASphere Base with the tablet via Bluetooth, select the QIASphere Base name from the list.
8. Enter the default login details. The **User Name** is Admin and the **Password** is Qiagen123.

9. Select **Login**. The tablet application prompts you to change the password.
10. Enter a new password of your choice according to the listed rules. Press **Reset**.
11. Follow the procedures below to define the LAN or Wi-Fi settings of the QIASphere Base.

Note: Do not click **Finish Configuration** before you enter the IP address of QIASphere Base into the QIAcube Connect instrument.

To define WLAN/Wi-Fi (WPAv2 only) setting

- 11a. Press the **WLAN** tab in the QIAcube Connect App to configure the Wi-Fi connection.
- 11b. Scan to search for available Wi-Fi networks.
- 11c. Select the appropriate Wi-Fi network and enter the password. The server, SSID and security type fields are automatically filled.
- 11d. Connect and proceed to **step 12**.

Note: QIASphere Base supports connecting to Wi-Fi networks with visible SSIDs only.

To define LAN (Static IP and DHCP) setting

- 11e. Select the **LAN** tab from the **Configuration** menu, enter the local network settings and set up the LAN Network. **Static IP Manual** is selected by default.

Note: When using DHCP-configured LAN, the network settings for QIASphere Base will be obtained automatically from the DHCP server. No further input is required.

Note: You might need the help of an IT professional to enter detailed network settings.

12. On the QIAcube Connect instrument screen, select the **Configuration** menu and then select the **Q-Base** tab.
13. Check the **Communication Enabled** box.
14. Enter the QIASphere Base IP address in the **Q-Base URL** field on the QIAcube Connect touchscreen. The IP address is shown in the QIAcube Connect App in the **Configuration** menu under the **Q-Base** tab. You can also scan the 2-D barcode.

Note: The test will fail if you press **Test** before pressing **Finish Configuration** on the tablet.

15. To continue, select **Apply**.
16. Press **Finish Configuration**.
Note: If you do not press **Finish Configuration**, a connection will be established after 10 minutes.
17. **Optional:** To test if the connection is successful, wait for 10 minutes and then press **Test**.
18. Reboot the QIAcube Connect instrument.

Revision History

Revision	Description
09/2020	Updated the steps in To set up the QIAcube Connect instrument section.

For further support with technical difficulties, please contact QIAGEN Technical Services at www.qiagen.com.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Notes

Notes

Trademarks: QIAGEN[®], Sample to Insight[®], QIAcube[®] (QIAGEN Group); Bluetooth[®] (Bluetooth SIG, Inc.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

HB-2619-004 09/2020 1117382 © 2020 QIAGEN, all rights reserved.