

# QIAseq™ Stranded Total RNA Library Kit

## Part 2: Strand-specific ligation

Immediately upon receipt, store the QIAseq Stranded Total RNA Kit (cat. nos. 180753 and 180755) at  $-30^{\circ}\text{C}$  to  $-15^{\circ}\text{C}$ . QIAseq Beads (cat. nos. 1107149, 1107460) should be stored at  $4^{\circ}\text{C}$  (**do not freeze**). If stored under these conditions, kits are stable until the date indicated on the QC label.

### Further information

- QIAseq Stranded Total RNA Library Handbook: [www.qiagen.com/HB-2465](http://www.qiagen.com/HB-2465)
- Safety Data Sheets: [www.qiagen.com/safety](http://www.qiagen.com/safety)
- Technical assistance: [support.qiagen.com](http://support.qiagen.com)

### Notes before starting

- To maximize output yields and minimize adapter dimer formation, use 2 ml tubes.
- **Important:** Ensure reactions are thoroughly mixed by pipetting up and down, prepared and incubated at recommended temperatures.
- Ligation Initiator is a very viscous solution. Pipet slowly and carefully to ensure the addition of the correct volume into the ligation reaction.

### Strand-specific ligation

1. Thaw and keep Ligation Initiator at room temperature until the ligation mix is prepared.
2. Vortex briefly and centrifuge the adapter plate. Remove the protective lid on the adapter plate, pierce the foil seal and dilute the QIAseq Adapters as described in Table 1 (Dilution of QIAseq adapters for total RNA input following mRNA enrichment)

or Table 2 (Dilution of QIAseq adapters for mRNA/depleted rRNA input).

- For each sample, add 2  $\mu$ l of diluted adapter. Ensure that only one adapter is used per sample, and be sure to track the barcodes. **Important:** To avoid cross-contamination, do not reuse diluted adapters.

**Table 1. Dilution of QIAseq Adapters for total RNA input following mRNA enrichment**

Input RNA (total RNA) amount	Volume of QIAseq Adapter	Volume of RNase-free water	Adapter dilution
100 ng	10 $\mu$ l	990 $\mu$ l*	1:100
500 ng	10 $\mu$ l	240 $\mu$ l*	1:25
1000 ng	10 $\mu$ l	115 $\mu$ l	1:12.5
5000 ng	10 $\mu$ l	40 $\mu$ l	1:5

**Table 2. Dilution of QIAseq Adapters for mRNA/depleted rRNA input**

Input RNA (total RNA) amount	Volume of QIAseq Adapter	Volume of RNase-free water	Adapter dilution
1 ng	10 $\mu$ l	1240 $\mu$ l*	1:100
10 ng	10 $\mu$ l	240 $\mu$ l*	1:25
50 ng	10 $\mu$ l	40 $\mu$ l	1:5
100 ng	10 $\mu$ l	15 $\mu$ l	1:2.5

\* Detailed dilution steps are described in the *QIAseq Stranded Total RNA Library Handbook*.

- Prepare ligations as described in Table 3 and incubate at 25°C for 10 min.

**Important:** Do not use a heated lid.

**Table 3. Strand-specific ligation reaction setup**

Component	Volume/reaction
Pre-mixed sample	50 $\mu$ l
Diluted adapter	2 $\mu$ l
Ultralow Input Ligation Buffer, 4X	25 $\mu$ l
Ultralow Input Ligase	5 $\mu$ l
Ligation Initiator	6.5 $\mu$ l
RNase-free water	11.5 $\mu$ l
<b>Total reaction volume</b>	<b>100 <math>\mu</math>l</b>

- 
5. After the ligation reaction, add 80  $\mu$ l QIAseq Beads. Mix well by pipetting up and down 10 times.
  6. Incubate for 5 min at room temperature.
  7. Place the tubes on a magnetic rack. After the solution has cleared (~10 min or longer), carefully remove and discard the supernatant.

**Important:** Do not discard the beads as they contain the DNA of interest.

8. With the tubes still on the magnetic stand, add 200  $\mu$ l of 80% ethanol. Rotate the tubes (2 to 3 times) to wash the beads. Carefully remove and discard the wash.
9. Repeat the ethanol wash.

**Important:** Completely remove all traces of the ethanol wash after this second wash. Remove the ethanol with a 200  $\mu$ l pipette first, and then use a 10  $\mu$ l pipette to remove any residual ethanol.

10. With the tubes (caps opened) still on the magnetic stand, air dry at room temperature for 5 to 10 min.

**Note:** Visually inspect that the pellet is completely dry, but avoid over-drying.

11. Remove the tubes from the magnetic stand, and elute the DNA from the beads by adding 92  $\mu$ l nuclease-free water. Mix well by pipetting.
  12. Return the tubes to the magnetic rack until the solution has cleared.
  13. Transfer 90  $\mu$ l supernatant to clean tubes.
  14. Add 108  $\mu$ l QIAseq Beads.
  15. Incubate for 5 min at room temperature.
  16. Place the tubes on a magnetic rack. After the solution has cleared (~10 min or longer), carefully remove and discard the supernatant.
- Important:** Do not discard the beads as they contain the DNA of interest.
17. With the tubes still on the magnetic stand, add 200  $\mu$ l of 80% ethanol. Rotate the tubes (2 to 3 times) to wash the beads. Carefully remove and discard the wash.
  18. Repeat the ethanol wash.
-

---

**Important:** Completely remove all traces of the ethanol wash after this second wash. Remove the ethanol with a 200 µl pipette first, and then use a 10 µl pipette to remove any residual ethanol.

19. With the beads (caps opened) still on the magnetic stand, air dry at room temperature for 5 to 10 min.

**Note:** Visually inspect that the pellet is completely dry, but avoid over-drying.

20. Remove the tubes from the magnetic stand, and elute the DNA from the beads by adding 25 µl nuclease-free water. Mix well by pipetting.

21. Return the tubes to the magnetic rack until the solution has cleared.

22. Transfer 23.5 µl to clean tubes/plate.

23. Proceed to CleanStart library amplification in part 3 of the protocol. Alternatively, the samples can be stored at  $-20^{\circ}\text{C}$  in a constant-temperature freezer overnight.



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN®, QIAseq™, Sample to Insight® (QIAGEN Group). 1107529 09/2017 HB-2451-001 © 2017 QIAGEN, all rights reserved.