

EpiTect[®] Fast DNA Bisulfite Kit – Part 2

See *Quick-Start Protocol: EpiTect Fast DNA Bisulfite Kit – Part 1* for instructions about kit storage and reagent preparation.

Further information

- *EpiTect Fast Bisulfite Conversion Handbook*: www.qiagen.com/HB-1211
- Safety Data Sheets: www.qiagen.com/safety
- Technical assistance: support.qiagen.com

Notes before starting

- Carrier RNA is not necessary if >100 ng DNA is used.
- Perform all centrifugation steps at room temperature (15–25°C) and maximum speed.
- Equilibrate samples and buffers to room temperature.

Cleanup of converted DNA

1. Upon completion of the bisulfite conversion (refer to Part 1), briefly centrifuge the PCR tubes. Transfer the reactions to clean 1.5 ml microcentrifuge tubes.
2. Add 310 µl freshly prepared Buffer BL to each sample (Table 1). For samples containing <100 ng DNA starting material, add dissolved carrier RNA to Buffer BL. Mix by vortexing and then centrifuge briefly.
3. Add 250 µl ethanol (96–100%) to each sample. Mix by pulse vortexing for 15 s and then centrifuge briefly to remove drops from inside the lid.
4. Place MinElute[®] DNA spin columns and collection tubes in a rack. Transfer the entire contents of each tube (step 3) to a corresponding spin column.
5. Centrifuge the spin columns for 1 min. Discard the flow-through and place the spin columns back into the collection tubes.
6. Add 500 µl Buffer BW to each spin column. Centrifuge for 1 min. Discard the flow-through and place the spin columns back into the collection tubes.



7. Add 500 μ l Buffer BD to each spin column, avoiding the transfer of any precipitates. Close the spin column lids, and incubate for 15 min at room temperature.
IMPORTANT: Minimize exposure of Buffer BD to air to prevent acidification.
8. Centrifuge the spin columns for 1 min. Discard the flow-through and place the spin columns back into the collection tubes.
9. Add 500 μ l Buffer BW to each spin column. Centrifuge for 1 min. Discard the flow-through and place the spin columns back into the collection tubes.
10. Repeat step 9 once.
11. Add 250 μ l ethanol (96–100%) to each spin column and centrifuge for 1 min.
12. Place the spin columns into new 2 ml collection tubes and centrifuge for 1 min to remove any residual liquid.
Optional: Incubate the spin columns with lids open in a heating block set to 60°C for 5 min to evaporate the liquid.
13. Place the spin columns into clean 1.5 ml microcentrifuge tubes (not provided). Add 15 μ l Buffer EB directly onto the center of each spin-column membrane and close the lids gently.
Note: As little as 10 μ l Buffer EB can be used for elution.
14. Incubate the spin columns at room temperature for 1 min. Centrifuge for 1 min to elute the DNA.
Note: Store purified DNA at 2–8°C for up to 24 h. For longer storage, we recommend storage at –20°C.

Table 1. Carrier RNA and Buffer BL volumes

Number of samples	1	4	8	16	24	48
Volume of Buffer BL*	350 μ l	1.4 ml	2.8 ml	5.6 ml	8.4 ml	16.8 ml
Volume of carrier RNA solution*	3.5 μ l	14 μ l	28 μ l	56 μ l	84 μ l	168 μ l

* Volumes give 10 μ g/ml carrier RNA in Buffer BL with a 10% surplus for pipetting inaccuracies.



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. Trademarks: QIAGEN®, Sample to Insight®, EpiTec®, MinElute® (QIAGEN Group). 1100634 02/2016 HB-1134-002 © 2016 QIAGEN, all rights reserved.