

**User-developed  
protocol**

## User-Developed Protocol:

### Isolation of plasmid DNA from yeast using the QIAprep<sup>®</sup> Spin Miniprep Kit

This procedure has been adapted from the QIAprep<sup>®</sup> Spin Miniprep Kit Protocol by Michael Jones, Chugai Institute for Molecular Medicine, Ibaraki, Japan. **It has not been thoroughly tested and optimized by QIAGEN.**

The procedure has been used successfully for isolation of plasmid DNA (6.6 kb low-copy plasmid pGAD 10, Clontech, containing inserts of average 1.6kb) from *Saccharomyces cerevisiae* strain PJ69-4A. The plasmid DNA obtained has been used successfully for PCR and transformation of *E. coli*.

Please be sure to read the *QIAprep Miniprep Handbook* and the detailed QIAprep Spin Miniprep Kit Protocol carefully before beginning this procedure, paying special attention to the important notes and the detailed QIAprep Spin Miniprep Kit Protocol using a microcentrifuge.

## Procedure

- 1. Inoculate a single colony into 2–5 ml of the appropriate selective media and grow the culture for 16–24 h at 30°C.**  
SD medium: 6,7g/liter Yeast Nitrogen Base (Difco 0919-15-3) supplemented with amino acids. To select for pGAD 10 do not add leucine.
- 2. Harvest the cells by centrifugation for 5 min at 5000 x g and resuspend cells in 250 µl Buffer P1 containing 0.1 mg/ml RNase A. Transfer the cell suspension to a 1.5 ml microcentrifuge tube.**
- 3. Add 50–100 µl of acid-washed glass beads (Sigma G-8772) and vortex for 5 min. Let stand to allow the beads to settle. Transfer supernatant to a fresh 1.5 ml microcentrifuge tube.**
- 4. Add 250 µl lysis buffer P2 to the tube and invert gently 4–6 times to mix. Incubate at room temperature for 5 min.**
- 5. Add 350 µl neutralization buffer N3 to the tube and invert immediately but gently 4–6 times.**
- 6. Centrifuge the lysate for 10 min at maximum speed in a tabletop microcentrifuge (13,000 rpm or  $\geq 10,000 \times g$ ). Meanwhile, place a QIAprep Spin Column in a 2 ml collection tube.**
- 7. Transfer the cleared lysate from step 6 to QIAprep Spin Column by decanting or pipetting.**
- 8. Centrifuge for 30–60 s (13,000 rpm or  $\geq 10,000 \times g$ ). Discard flow-through.**

**User-developed  
protocol**

9. **Wash QIAprep Spin Column by adding 0.75 ml of Buffer PE and centrifuging 30–60 s (13,000 rpm or  $\geq 10,000 \times g$ ).**
10. **Discard flow-through and centrifuge for an additional 1 min to remove residual wash buffer (13,000 rpm or  $\geq 10,000 \times g$ ).**

**IMPORTANT:** Residual wash buffer will not be completely removed unless the flow-through is discarded before this additional centrifugation. Residual ethanol from Buffer PE may inhibit subsequent enzymatic reactions.

11. **Place QIAprep Spin Column in a clean 1.5 ml microcentrifuge tube. To elute DNA, add 25  $\mu$ l of Buffer EB (10 mM Tris·Cl, pH 8.5) or water to the center of each QIAprep Spin Column, let stand for 1 min, and centrifuge for 1 min.**

Typical yield is up to 1  $\mu$ g. For subsequent PCR, use 0.1–1.0  $\mu$ l of the eluate. For subsequent transformation into *E. coli*, 2–3  $\mu$ l of eluate yields about 30 colonies.

QIAGEN handbooks can be requested from QIAGEN Technical Service or your local QIAGEN distributor.

Selected handbooks can be downloaded from [www.qiagen.com/literature/handbooks/default.asp](http://www.qiagen.com/literature/handbooks/default.asp).

Material safety data sheets (MSDS) for any QIAGEN product can be downloaded from [www.qiagen.com/ts/msds.asp](http://www.qiagen.com/ts/msds.asp).

Trademarks: QIAGEN<sup>®</sup>, QIAprep<sup>®</sup> (QIAGEN).

The PCR process is covered by US Patents 4,683,195 and 4,683,202 and foreign equivalents issued to Hoffmann-La Roche AG.

© 2001 QIAGEN, all rights reserved.