

## miRNeasy Serum/Plasma Kit

The miRNeasy Serum/Plasma Kit (cat. no. 217184) is shipped at ambient temperature. Store the RNeasy® MinElute® spin columns immediately at 2–8°C. QIAzol® Lysis Reagent can be stored at room temperature (15–25°C) or at 2–8°C. Store the remaining components dry at room temperature. All kit components are stable for at least 9 months under these conditions if not otherwise stated on label.

### Further information

- *miRNeasy Serum/Plasma Handbook*: [www.qiagen.com/HB-1002](http://www.qiagen.com/HB-1002)
- Safety Data Sheets: [www.qiagen.com/safety](http://www.qiagen.com/safety)
- Technical assistance: [support.qiagen.com](http://support.qiagen.com)

### Notes before starting

- This protocol is for purifying total RNA, including small RNAs, from small volumes (up to 200 µl) of serum, plasma, or other body fluids.
  - If necessary, redissolve any precipitate in Buffer RWT by warming.
  - Except for phase separation (step 6), all steps should be performed at room temperature (15–25°C). Work quickly.
  - Add ethanol (96–100%) to Buffer RWT and Buffer RPE concentrates before use (see bottle label for volume).
  - Before starting with step 1 for the first time, read the recommendations for preparing serum or plasma in the *miRNeasy Serum/Plasma Handbook*.
1. Prepare serum or plasma, or thaw frozen samples.
  2. Add 5 volumes QIAzol Lysis Reagent to the sample (e.g., for 200 µl sample, add 1 ml QIAzol Lysis Reagent). Mix by vortexing or pipetting up and down.
  3. Incubate the homogenate at room temperature (15–25°C) for 5 min.

4. Add chloroform of an equal volume to the starting sample and cap tube securely (e.g., for 200  $\mu$ l sample, add 200  $\mu$ l chloroform). Shake vigorously for 15 s.
5. Incubate at room temperature for 2–3 min.
6. Centrifuge for 15 min at 12,000  $\times g$  at 4°C.
7. Transfer the upper aqueous phase to a new collection tube (not supplied). Avoid transferring any interphase. Add 1.5 volumes of 100% ethanol (e.g., for 600  $\mu$ l aqueous phase, add 900  $\mu$ l ethanol). Mix thoroughly by pipetting.
8. Pipet up to 700  $\mu$ l sample, including any precipitate, into an RNeasy MinElute spin column in a 2 ml collection tube. Close the lid and centrifuge at  $\geq 8000 \times g$  for 15 s at room temperature. Discard the flow-through.
9. Repeat step 8 using the remainder of the sample.
10. Add 700  $\mu$ l Buffer RWT to the RNeasy MinElute spin column. Close the lid, and centrifuge for 15 s at  $\geq 8000 \times g$ . Discard the flow-through.
11. Pipet 500  $\mu$ l Buffer RPE onto the RNeasy MinElute spin column. Close the lid, and centrifuge for 15 s at  $\geq 8000 \times g$ . Discard the flow-through.
12. Add 500  $\mu$ l of 80% ethanol to the RNeasy MinElute spin column. Close the lid, and centrifuge for 2 min at  $\geq 8000 \times g$ . Discard the flow-through and the collection tube.
13. Place the RNeasy MinElute spin column in a new 2 ml collection tube (supplied). Open the lid of the spin column and centrifuge at full speed for 5 min to dry the membrane. Discard the flow-through and the collection tube.
14. Place the RNeasy MinElute spin column in a new 1.5 ml collection tube (supplied). Add 14  $\mu$ l RNase-free water directly to the center of the spin column membrane. Close the lid gently, and centrifuge for 1 min at full speed to elute the RNA.



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. Trademarks: QIAGEN<sup>®</sup>, Sample to Insight<sup>®</sup>, QiAzo<sup>®</sup>, MinElute<sup>®</sup>, RNeasy<sup>®</sup> (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

1127637 04/2022 HB-0982-003 © 2022 QIAGEN, all rights reserved.