

# RNeasy® Plus Mini Kit

The RNeasy Plus Mini Kit (cat. nos. 74134 and 74136) can be stored at room temperature (15–25°C) for at least 9 months if not otherwise stated on label.

## Further information

- *RNeasy Plus Mini Handbook*: [www.qiagen.com/HB-0405](http://www.qiagen.com/HB-0405)
- Safety Data Sheets: [www.qiagen.com/safety](http://www.qiagen.com/safety)
- Technical assistance: [support.qiagen.com](http://support.qiagen.com)

## Notes before starting

- If purifying RNA from cell lines rich in RNases, or tissue, add either 10 µl β-mercaptoethanol (β-ME), or 20 µl 2 M dithiothreitol (DTT), to 1 ml Buffer RLT Plus before use. Buffer RLT Plus containing DTT or β-ME can be stored at room temperature for up to 1 month.
- Add 4 volumes of ethanol (96–100%) to Buffer RPE for a working solution.
- Foaming can be reduced by adding Reagent DX (cat. no. 19088) at a final concentration of 0.5% (v/v) before disruption and homogenization.\*

\* This option not included in handbook; handbook to be updated.

1. **Cells:** Harvest a maximum of  $1 \times 10^7$  cells, either as a cell pellet, or lysed directly in the vessel. Add the appropriate volume of Buffer RLT Plus (see Table 1). Vortex for 30 s, or homogenize.

**Tissues:** Disrupt the tissue ( $\leq 30$  mg) and homogenize the lysate in the appropriate volume of Buffer RLT Plus (see Table 1). Centrifuge the lysate for 3 min at maximum speed.

Carefully remove the supernatant by pipetting and use it in step 2.

2. Transfer the homogenized lysate to a gDNA Eliminator spin column placed in a 2 ml collection tube (supplied).

3. Centrifuge for 30 s at  $\geq 8000 \times g$  ( $\geq 10,000$  rpm). Discard the column, and save the flow-through. Add 1 volume (usually 350  $\mu$ l or 600  $\mu$ l) of 70% ethanol to the flow-through, and mix well by pipetting. Do not centrifuge. Proceed immediately to step 4.
4. Transfer up to 700  $\mu$ l of the sample, including any precipitate, to an RNeasy spin column placed in a 2 ml collection tube (supplied). Close the lid, and centrifuge for 15 s at  $\geq 8000 \times g$ . Discard the flow-through.
5. Add 700  $\mu$ l Buffer RW1 to the RNeasy Mini spin column (in a 2 ml collection tube). Close the lid, and centrifuge for 15 s at  $\geq 8000 \times g$ . Discard the flow-through.
6. Add 500  $\mu$ l Buffer RPE to the RNeasy spin column. Close the lid, and centrifuge for 15 s at  $\geq 8000 \times g$ . Discard the flow-through.
7. Add 500  $\mu$ l Buffer RPE to the RNeasy spin column. Close the lid gently, and centrifuge for 2 min at  $\geq 8000 \times g$  ( $\geq 10,000$  rpm).

**Optional:** Place the RNeasy spin column in a new 2 ml collection tube (supplied). Centrifuge at full speed for 1 min to further dry the membrane.

8. Place the RNeasy spin column in a new 1.5 ml collection tube (supplied). Add 30–50  $\mu$ l RNase-free water directly to the spin column membrane. Close the lid, and centrifuge for 1 min at  $\geq 8000 \times g$  to elute the RNA.

**Optional:** Repeat elution with another volume of water or with RNA eluate.

**Table 1. Volumes of Buffer RLT Plus for sample disruption and homogenization**

Sample	Amount	Dish	Buffer RLT Plus	Disruption and homogenization
Pelleted cells	$<5 \times 10^6$	$<6$ cm	350 $\mu$ l	Add Buffer RLT Plus, vortex ( $\leq 1 \times 10^6$ cells); or use QIAshredder, TissueRuptor <sup>®</sup> , or needle and syringe
	$\leq 1 \times 10^7$	6–10 cm	600 $\mu$ l	
Animal tissues	$<20$ mg	–	350 $\mu$ l	TissueLyser LT; TissueLyser II; TissueRuptor, or mortar and pestle followed by QIAshredder or needle and syringe
	20–30 mg	–	600 $\mu$ l	

\* Use 600  $\mu$ l Buffer RLT Plus for tissues stabilized in RNAlater<sup>®</sup> Reagent, or for difficult-to-lyse tissues.



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. Trademarks: QIAGEN<sup>®</sup>, Sample to Insight<sup>®</sup>, RNeasy<sup>®</sup>, TissueRuptor<sup>®</sup> (QIAGEN Group). "RNAlater<sup>™</sup>" is a trademark of AMBION, Inc., Austin, Texas and is covered by various U.S. and foreign patents. 1101300 03/2016 HB-0574-002 © 2016 QIAGEN, all rights reserved.