September 2021

# EZ2<sup>®</sup> Connect and EZ2 Connect Fx User Manual

CE

REF

MAT

9003210, 9003220

QIAGEN GmbH  
QIAGEN Strasse 1, 40724 Hilden, GERMANY

1124549

# Contents

1	Introduction .....	7
1.1	About this user manual .....	7
1.2	General information.....	8
1.2.1	Technical assistance.....	8
1.2.2	Policy statement.....	8
1.3	Intended use of the EZ2.....	9
1.4	Requirements for EZ2 users.....	9
2	Safety Information .....	10
2.1	Proper use .....	11
2.2	Electrical safety .....	13
2.3	Operating conditions .....	15
2.4	Biological safety.....	16
2.5	Chemicals.....	17
2.6	Waste disposal .....	18
2.7	Mechanical hazards .....	18
2.8	Heat hazard .....	19
2.9	Radiation.....	19
2.10	Maintenance safety .....	20
2.11	Symbols on the EZ2 Connect instrument .....	22
3	General Description .....	24
3.1	Principle .....	24
3.2	External features of the EZ2 .....	25
3.2.1	Touchscreen.....	27
3.2.2	Hood .....	28
3.2.3	USB ports .....	29
3.2.4	RJ-45 Ethernet port and power cord socket .....	30
3.2.5	Ventilation openings .....	31
3.2.6	Bar code scanner .....	31
3.3	Internal features of the EZ2 .....	32

3.3.1	Pipettor head .....	33
3.3.2	Worktable .....	34
3.3.3	Magnet module .....	38
3.3.4	Camera .....	39
3.3.5	UV LED light .....	39
3.3.6	Internal light .....	40
4	Installation Procedures .....	41
4.1	Installation environment .....	41
4.1.1	Site requirements .....	41
4.1.2	Power requirements .....	42
4.1.3	Grounding requirements .....	43
4.2	Unpacking the EZ2 .....	44
4.3	Installing the EZ2 .....	46
4.3.1	Removal of the EZ2 accessories and shipping materials .....	46
4.3.2	Installation of the AC power cord .....	46
4.3.3	Installation of external barcode scanner (optional) .....	48
4.3.4	Installation of Wi-Fi adapter (optional when using QIAsphere) .....	48
4.3.5	Initial configuration of the EZ2 .....	48
4.3.6	Perform camera exposure calibration (for EZ2 Connect Fx) .....	50
5	Operating Procedures .....	52
5.1	General information .....	55
5.1.1	Entering text and numbers .....	58
5.2	Starting the EZ2 .....	60
5.3	Configuring the EZ2 .....	61
5.3.1	Setting basic system data .....	62
5.3.2	Changing Instrument settings .....	63
5.3.3	Managing users .....	64
5.3.4	Changing password .....	71
5.3.5	Managing user security .....	72
5.3.6	Installing new protocols .....	74
5.3.7	Updating software .....	74

5.3.8	Configuring network and QIAsphere Base connections .....	76
5.4	Setting up a protocol run .....	81
5.4.1	Selecting application type .....	81
5.4.2	Selecting a kit .....	82
5.4.3	Selecting a protocol .....	82
5.4.4	Defining parameters .....	83
5.4.5	Selecting sample positions .....	84
5.4.6	Entering sample IDs .....	86
5.4.7	Loading the EZ2 Connect Cartridge Rack .....	89
5.4.8	Loading the EZ2 Connect Tip Rack .....	90
5.5	Canceling the protocol run setup .....	91
5.6	Starting the protocol run and monitoring its progress .....	91
5.6.1	Load check .....	93
5.7	End of the protocol run .....	95
5.8	Aborting a protocol run .....	96
5.9	Sample recovery .....	98
5.10	Saving a run report .....	98
5.10.1	Run report contents .....	99
5.11	Setting up the worktable .....	99
5.11.1	Loading and unloading the cartridge rack .....	100
5.11.2	Loading and unloading the tip rack .....	101
5.12	Using the bar code scanner .....	103
5.13	Logging out and switching the instrument off .....	104
6	Maintenance Procedures .....	105
6.1	Cleaning agents .....	105
6.1.1	Disinfecting the EZ2 .....	107
6.1.2	Removing contamination .....	107
6.2	After run maintenance .....	108
6.3	Daily maintenance .....	111
6.4	Weekly maintenance .....	114
6.5	UV decontamination .....	117

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---

# 1 Introduction

Thank you for choosing the EZ2 Connect. We are confident it will become an integral part of your laboratory.

This user manual describes the EZ2 Connect.

Before using the EZ2 Connect, it is essential that you read this user manual carefully and pay attention to the safety information. The instructions and safety information in the user manual must be followed to ensure safe operation of the instrument and to maintain the instrument in a safe condition.

## 1.1 About this user manual

This user manual provides information about the EZ2 Connect and EZ2 Connect Fx (hereafter also referred to as EZ2) in the following sections:

- Introduction — contains the intended use and requirements for users
- General information — contains information regarding the intended use and requirements of the EZ2
- Safety Information — includes important information on any hazards related to the EZ2 and how to properly use the instrument
- General Description — an overview of the features of the EZ2
- Installation Procedures — instructions on how to set up the instrument before its first use
- Operating Procedures — includes instructions related to protocol runs
- Maintenance Procedures — contains information about cleaning and maintenance
- Troubleshooting — instructions on what to do in case of any problems with the EZ2
- Glossary — an alphabetical list of terms or words used in this user manual with explanations
- Technical Specification - Technical data

The appendices contain the following information:

- Appendix A - Legal Requirements for the EZ2
- Appendix B – EZ2 Accessories — an overview of the accessories available for the EZ2
- Appendix C - Ordering Information
- Document Revision History — the changes made to the user manual

---

## 1.2 General information

### 1.2.1 Technical assistance

At QIAGEN®, we pride ourselves on the quality and availability of our technical support. Our Technical Services Departments are staffed by experienced scientists with extensive practical and theoretical expertise in molecular biology and the use of QIAGEN products. If you have any questions or experience any difficulties regarding the EZ2 or QIAGEN products in general, do not hesitate to contact us.

QIAGEN customers are a major source of information regarding advanced or specialized uses of our products. This information is helpful to other scientists as well as to the researchers at QIAGEN. We therefore encourage you to contact us if you have any suggestions about product performance or new applications and techniques.

For technical assistance and more information, please see our Technical Support Center at [support.qiagen.com](https://support.qiagen.com) or call one of the QIAGEN Technical Service Departments or local distributors.

When contacting QIAGEN Technical Services about errors, please have the following information ready:

- EZ2 serial number, type, and version
- Error code (if applicable)
- Timepoint when the error occurred for the first time
- Frequency of error occurrence (i.e., intermittent or persistent error)
- Downloaded support package (see section 7.1.1)

### 1.2.2 Policy statement

It is the policy of QIAGEN to improve products as new techniques and components become available. QIAGEN reserves the right to change specifications at any time.

To produce useful and appropriate documentation, we appreciate your comments on this user manual. Please contact QIAGEN Technical Services.

### 1.3 Intended use of the EZ2

The EZ2 Connect and EZ2 Connect Fx systems are designed to perform automated isolation and purification of nucleic acids. The EZ2 is intended to be used only in combination with QIAGEN kits indicated for use with the EZ2 instrument for the applications described in the kit handbooks. The EZ2 system is intended for use by professional operators, such as technicians and physicians trained in molecular biological techniques and the operation of the EZ2 system.

### 1.4 Requirements for EZ2 users

The table below covers the general level of competence and training necessary for transportation, installation, use, maintenance, and servicing of the EZ2.

<b>Task</b>	<b>Personnel</b>	<b>Training and experience</b>
Delivery	No special requirements	No special requirements
Installation	Laboratory technicians or equivalent	Appropriately trained and experienced personnel familiar with the use of computers and automation in general
Routine use (running protocols)	Laboratory technicians or equivalent	Professional users, such as technicians or physicians, trained in molecular biology techniques
Routine maintenance	Laboratory technicians or equivalent	Professional users, such as technicians or physicians, trained in molecular biology techniques
Servicing and annual maintenance	QIAGEN Field Service specialists only	Specialists trained, certified, and authorized by QIAGEN

## 2 Safety Information

Before using the EZ2, it is essential that you read this user manual carefully and pay attention to the safety information. The instructions and safety information in the user manual must be followed to ensure safe operation of the instrument and to maintain the instrument in a safe condition.

Possible hazards that could harm the user or result in damage to the instrument are clearly stated at the appropriate places throughout this user manual.

If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

The following types of safety information appear in this user manual.

<p><b>WARNING</b></p> 	<p>The term <b>WARNING</b> is used to inform you about situations that could result in <b>personal injury</b> to you or others.</p> <p>Details about these circumstances are given in a box like this one.</p>
--	--

<p><b>CAUTION</b></p> 	<p>The term <b>CAUTION</b> is used to inform you about situations that could result in <b>damage to an instrument</b> or other equipment.</p> <p>Details about these circumstances are given in a box like this one.</p>
---	--

The advice given in this manual is intended to supplement, not supersede, the normal safety requirements prevailing in the user's country.

Please be aware that you may be required to consult your local regulations for reporting serious incidents that have occurred in relation to the device to the manufacturer and/or its authorized representative (only applicable for CE-marked devices with an authorized representative established in the EU) and the regulatory authority in which the user and/or the patient is established.

## 2.1 Proper use

<p><b>WARNING/ CAUTION</b></p> 	<p><b>Risk of personal injury and material damage</b></p> <p>Improper use of the EZ2 may cause personal injuries or damage to the instrument. The EZ2 must only be operated by qualified personnel who have been appropriately trained. Servicing of the EZ2 instrument must only be performed by a QIAGEN Field Service Specialist.</p>
--	--

<p><b>WARNING</b></p> 	<p><b>Risk of personal injury</b></p> <p>The EZ2 is too heavy to be lifted by one person. To avoid personal injury or damage to the instrument, do not lift the instrument alone. Use the handle attached to the box to lift the EZ2. After the EZ2 has been unpacked, two people must lift the instrument. Lift the instrument by placing your hands underneath the bottom of the instrument.</p>
---	--

<p><b>WARNING</b></p> 	<p><b>Risk of personal injury and material damage</b></p> <p>Do not attempt to move the EZ2 during operation.</p>
---	---

Perform the maintenance as described in the Maintenance Procedures section. QIAGEN charges for repairs that are required due to incorrect maintenance.

In case of emergency, power OFF the EZ2 at the power switch located in front of the instrument and unplug the power cord from the power outlet.

<p><b>CAUTION</b></p> 	<p><b>Damage to the instrument</b></p> <p>Avoid spilling water or chemicals onto the EZ2. Instrument damage caused by water or chemical spillage will void your warranty.</p>
---	---

<b>WARNING</b> 	<b>Risk of fire or explosion</b> When using ethanol or ethanol-based liquids on the EZ2, handle such liquids carefully and in accordance with the required safety regulations. If liquid has been spilled, wipe it off and leave the EZ2 hood open to allow flammable vapors to disperse.
---	--

<b>WARNING</b> 	<b>Risk of explosion</b> The EZ2 is intended for use with reagents and substances supplied with QIAGEN kits as outlined in respective Information for use. Use of other reagents and substances may lead to fire or explosion.
---	---

If hazardous material is spilled on or inside the EZ2, the user is responsible for carrying out appropriate decontamination.

**Note:** Do not place items on top of the EZ2.

<b>CAUTION</b> 	<b>Damage to the instrument</b> Ensure that the EZ2 is switched off before you manually move the mechanical components of the instrument.
---	--

<b>CAUTION</b> 	<b>Damage to the instrument</b> Do not lean against the instrument or touchscreen.
---	---

## 2.2 Electrical safety

**Note:** If operation of the instrument is interrupted in any way (for example, due to interruption of the power supply or a mechanical error), first switch off the EZ2 instrument, then disconnect the electrical cord from the power supply before attempting troubleshooting or service activity.

<p><b>WARNING</b></p> 	<p><b>Electrical hazard</b></p> <p>Any interruption of the protective conductor (earth/ground lead) inside or outside the instrument or disconnection of the protective conductor terminal is likely to make the instrument dangerous.</p> <p>Intentional interruption is prohibited.</p> <p><b>Lethal voltages inside the instrument</b></p> <p>When the instrument is connected to line power, terminals may be live and opening covers or removing parts is likely to expose live parts.</p>
---	---

<p><b>WARNING</b></p> 	<p><b>Damage to electronics</b></p> <p>Before powering ON the instrument, make sure that the correct supply voltage is used.</p> <p>Use of incorrect supply voltage may damage the electronics.</p> <p>To check the recommended supply voltage, refer to the specifications indicated in the type plate of the instrument.</p>
---	--

<p><b>WARNING</b></p> 	<p><b>Risk of electric shock</b></p> <p>Do not open any panels on the EZ2.</p> <p><b>Risk of personal injury and material damage</b></p> <p>Only perform maintenance that is specifically described in this user manual.</p>
---	--

To ensure satisfactory and safe operation of the EZ2, follow the advice below:

- The line power cord must be connected to a line power outlet that has a protective conductor (earth/ground).
- Place instrument in a location so that the power cord is accessible and can be connected/disconnected.
- Use only the power cord delivered by QIAGEN.
- Do not adjust or replace internal parts of the instrument.
- Do not operate the instrument with any covers or parts removed.
- If liquid has spilled inside the instrument, switch off the instrument, disconnect it from the power outlet, and contact QIAGEN Technical Services for advice before attempting troubleshooting or service activity.

If the instrument becomes electrically unsafe, prevent other personnel from operating it and contact QIAGEN Technical Services.

The instrument may be electrically unsafe when:

- The EZ2 or the line power cord appear to be damaged.
- The EZ2 has been stored in unfavorable conditions for a prolonged period.
- The EZ2 has been subjected to severe transport stresses.
- Liquids have come into direct contact with electrical components of the EZ2.
- The power cord has been exchanged with a non-official power cord.

<b>WARNING</b> 	<b>Electric hazard</b> Do not touch the EZ2 with wet hands.
---	--

<b>WARNING</b> 	<b>Electric hazard</b> Never install a fuse different from that specified in the user manual.
---	--

## 2.3 Operating conditions

Parameters such as temperature range and humidity range are described in the Technical Specification section.

<b>WARNING</b> 	<b>Explosive atmosphere</b> The EZ2 instrument is not designed for use in an explosive atmosphere.
---	---

<b>WARNING</b> 	<b>Risk of overheating</b> To ensure proper ventilation, maintain a minimum clearance of 10 cm at the sides and rear of the EZ2.  Slits and openings that ensure the ventilation of the instrument must not be covered.
---	--

<b>WARNING</b> 	<b>Risk of explosion</b> The EZ2 is intended for use with reagents and substances supplied with QIAGEN kits. Use of other reagents and substances may lead to fire or explosion.
---	---

<b>CAUTION</b> 	<b>Damage to the instrument</b> Direct sunlight may bleach parts of the instrument, cause damage to plastic parts or interfere with the correct function of the load check. The EZ2 must be located out of direct sunlight.
---	--

<p><b>CAUTION</b></p> 	<p><b>Damage to the instrument</b></p> <p>Do not use the EZ2 in the vicinity of sources of strong electromagnetic radiation (for example unshielded, deliberately operated high-frequency sources or mobile radio devices), because these can interfere with the proper operation.</p>
---	--

## 2.4 Biological safety

Specimens and reagents containing materials from humans should be treated as potentially infectious. Use safe laboratory procedures as outlined in publications such as Biosafety in Microbiological and Biomedical Laboratories, HHS (<https://www.cdc.gov/labs/pdf/CDC-BiosafetymicrobiologicalBiomedicalLaboratories-2009-P.pdf>). You should be aware of the health hazard presented by such agents and should use, store, and dispose of such samples according to the required safety regulations.

<p><b>WARNING</b></p> 	<p><b>Samples containing infectious agents</b></p> <p>Some samples used with the EZ2 instrument may contain infectious agents. Handle such samples with the greatest of care and in accordance with the required safety regulations.</p> <p>Always wear safety glasses, gloves, and a lab coat.</p> <p>The responsible body (for example, a laboratory manager) must take the necessary precautions to ensure that the surrounding workplace is safe, and that the instrument operators are suitably trained and not exposed to hazardous levels of infectious agents as defined in the applicable Material Safety Data Sheets (MSDSs) or the OSHA<sup>1,*</sup> ACGIH,<sup>†</sup> or COSHH<sup>‡</sup> documents.</p> <p>Venting for fumes and disposal of wastes must be in accordance with all national, state, and local health and safety regulations and laws.</p>
---	---

\* OSHA – Occupational Safety and Health Organization (United States of America)

† ACGIH – American Conference of Government Industrial Hygienists (United States of America)

‡ COSHH – Control of Substances Hazardous to Health (United Kingdom)

## 2.5 Chemicals

<p><b>WARNING</b></p> 	<p><b>Hazardous chemicals</b></p> <p>Some chemicals used with the EZ2 instrument may be hazardous or may become hazardous after completion of a purification.</p> <p>Always wear safety glasses, gloves, and a lab coat.</p> <p>The responsible body (for example, a laboratory manager) must take the necessary precautions to ensure that the surrounding workplace is safe, and that the instrument operators are not exposed to hazardous levels of toxic substances (chemical or biological) as defined in the applicable Material Safety Data Sheets (MSDSs) or the OSHA,* ACGIH,<sup>†</sup> or COSHH<sup>‡</sup> documents.</p> <p>Venting for fumes and disposal of wastes must be in accordance with all national, state, and local health and safety regulations and laws.</p>
---	---

\* OSHA: Occupational Safety and Health Administration (United States of America).

<sup>†</sup> ACGIH: American Conference of Government Industrial Hygienists (United States of America).

<sup>‡</sup> COSHH: Control of Substances Hazardous to Health (United Kingdom).

### Toxic Fumes

<p><b>WARNING</b></p> 	<p><b>Toxic fumes</b></p> <p>Do not use bleach to clean or disinfect the EZ2 instrument. Bleach in contact with salts from the buffers can produce toxic fumes.</p>
---	---

<p><b>WARNING</b></p> 	<p><b>Toxic fumes</b></p> <p>Do not use bleach to disinfect used labware. Bleach in contact with salts from the buffers can produce toxic fumes.</p>
---	--

**Note:** If you work with volatile solvents, toxic substances, etc., you must provide an efficient laboratory ventilation system to remove vapors that may be produced.

## 2.6 Waste disposal

Used consumables, such as reagent cartridges and disposable filter-tips, may contain hazardous chemicals or infectious agents from the purification process. Such waste must be collected and disposed of properly, according to local safety regulations.

For more information about how to dispose of the EZ2, see Appendix A: Waste Electrical and Electronic Equipment (WEEE).

<p><b>CAUTION</b></p> 	<p><b>Hazardous chemicals and infectious agents</b></p> <p>Waste may contain toxic or infectious material and must be disposed of properly. Refer to your local safety regulations for proper disposal procedures.</p>
---	--

## 2.7 Mechanical hazards

The hood of the EZ2 must remain closed during operation of the instrument. Only open the hood when instructed to do so by the instruction for use.

The worktable of the EZ2 instrument moves during operation of the instrument. While loading the worktable, always stand clear of the instrument. Do not lean on the worktable when the robotic arm of the instrument is moving to reach loading position with its lid open. Wait until the robotic arm completes its movements before you start to load or unload.

<p><b>WARNING</b></p> 	<p><b>Moving parts</b></p> <p>To avoid contact with moving parts during the operation of the EZ2, the instrument must be operated with the hood closed.</p> <p>If the hood sensor or lock is not functioning correctly, contact QIAGEN Technical Services.</p>
---	--

<p><b>WARNING</b></p> 	<p><b>Moving parts</b></p> <p>Avoid contact with moving parts during operation of the EZ2. Under no circumstances should hands be placed under the pipetting arm during movement. Do not attempt to remove any plasticware from the workdeck whilst the instrument is operating.</p>
---	--

## 2.8 Heat hazard

The EZ2 worktable contains a heating system.

<p><b>WARNING</b></p> 	<p><b>Hot surface</b></p> <p>The heating system can reach temperatures of up to 95°C (203°F). Avoid touching it when it is hot.</p>
---	---

## 2.9 Radiation

The EZ2 instrument has a UV LED lamp. The wavelength of the UV light produced by the UV LED lamp is 270 to 285 nm. This wavelength corresponds to ultraviolet light type C, which can be used for decontamination procedures. A mechanical lock ensures that the hood must be closed for operation of the UV LED. If the hood sensor or lock is not functioning correctly, contact QIAGEN Technical Services.

<p><b>WARNING</b></p> 	<p><b>UV radiation</b></p> <p>Avoid looking directly into UV light. Do not expose your skin to UV light.</p>
---	--

The EZ2 instrument has a 2D handheld barcode scanner to allow kit bar code and sample barcode scanning.

<p><b>WARNING</b></p> 	<p><b>Risk of personal injury</b></p> <p>Hazard Level 2 laser light: Do not stare into the light beam when using handheld barcode scanner.</p>
---	--

## 2.10 Maintenance safety

<b>WARNING/ CAUTION</b> 	<b>Risk of personal injury and material damage</b> Only perform maintenance that is specifically described in this user manual.
--	--

Perform the maintenance as described in the Maintenance Procedures section. QIAGEN charges for repairs that are required due to incorrect maintenance.

<b>WARNING/ CAUTION</b> 	<b>Risk of personal injury and material damage</b> Improper use of the EZ2 instrument may cause personal injuries or damage to the instrument.  The EZ2 must only be operated by suitably qualified personnel.  Servicing of the EZ2 instrument must only be performed by QIAGEN Field Service Specialists.
--	--

<b>WARNING</b> 	<b>Risk of fire</b> When cleaning the EZ2 instrument with alcohol-based disinfectant, leave the instrument door open to allow flammable vapors to disperse.  Only clean the EZ2 instrument with alcohol-based disinfectant when worktable components have cooled down.
---	---

<b>CAUTION</b> 	<b>Damage to the instrument</b> Do not use bleach, solvents, or reagents containing acids, alkalis, or abrasives to clean the EZ2 instrument.
---	--

<b>CAUTION</b> 	<b>Damage to the instrument</b> Do not use spray bottles containing alcohol or disinfectant to clean surfaces of the EZ2 instrument. Spray bottles should be used only to clean items that have been removed from the worktable and if permitted by local laboratory operating practices.
---	--

<b>CAUTION</b> 	<b>Damage to the instrument</b> After wiping the worktable with paper towels, make sure that no residual pieces of paper towel remain. Pieces of paper towel remaining on the worktable could lead to a worktable collision.
---	---

<b>WARNING/ CAUTION</b> 	<b>Risk of personal electric shock</b> Do not open any panels on the EZ2 instrument.  Only perform maintenance as described in this user manual.
--	---

## 2.11 Symbols on the EZ2 Connect instrument

The following symbols appear on the EZ2 Connect and EZ2 Connect Fx instruments.

Symbol	Location	Description
	Heating system — inside instrument	Heat hazard — the temperature of the heating system can reach up to 95°C.
	Near the tip rack	Biohazard — the tip rack may be contaminated with biohazardous material and must be handled with gloves.
	On the back of the instrument	UV radiation hazard — avoid looking directly into UV light. Do not expose your skin to UV light.
	On handheld barcode scanner	Hazard Level 2 laser light: Do not stare into the light beam when using handheld barcode scanner.
	Robotic arm — inside instrument	Crush hazard — the nozzle unit may crush your fingers or hand.
	Type plate on the back of the instrument	CE mark for Europe.
	Type plate on the back of the instrument	CSA mark for Canada and the USA.
	Type plate on the back of the instrument	RCM (former C-Tick) for Australia and New Zealand.
	Type plate on the back of the instrument	RoHS mark for China (the restriction of the use of certain hazardous substances in electrical and electronic equipment).
	Type plate on the back of the instrument	WEEE mark for Europe.
	Type plate on the back of the instrument	Legal manufacturer.

Symbol	Location	Description
	Type plate on the back of the instrument	Unique Device Identifier (UDI) as a 2D bar code in Data Matrix format.
	Type plate on the back of the instrument	Global Trade Item Number.
	Type plate on the back of the instrument	Serial number.
	Type plate on the back of the instrument	Consult instructions for use.
	Type plate on the back of the instrument	See warnings and precautions.

---

## 3 General Description

### 3.1 Principle

The EZ2 is available in different variants. This user manual focuses on the EZ2 Connect and EZ2 Connect Fx.

The EZ2 Connect performs fully automated nucleic acid purification from up to 24 samples in molecular biology applications using magnetic particle technology.

The EZ2 has preinstalled protocols that are used with QIAGEN kits to purify nucleic acids. The instrument's touchscreen display allows the user to easily select protocols. The intuitive software/user interface guides the user through the run setup process, which includes selecting variable parameters. Afterward, the operator loads labware and samples onto the EZ2 worktable, while following the instructions shown on the display. The order and contents of the labware are determined by each protocol. When the protocol run starts, the isolation of nucleic acids is carried out automatically in a modular workflow. The workflow includes sample lysis, nucleic acid binding to magnetic particles, contaminant removal by washing, and, finally, elution from magnetic beads. The execution of individual steps varies and is optimized for each protocol.

The aspiration and dispensation of samples and reagents and the separation of magnetic particles are performed by the 24-channel pipettor head. If required by the protocol, the temperature of the liquids is controlled by the heating system.

The EZ2 Connect Fx incorporates all the characteristics of the EZ2 Connect, as well as some additional features, which are particularly relevant for forensic workflows. These features include:

- An internal camera, which is used for load checks and bar code reading
- An external bar code reader, which is used for reading sample IDs and kit bar codes
- A sample recovery function for forensic protocols
- Extended user management
- User Interface with additional features
- Additional reporting functions
- UV schedule function

### 3.2 External features of the EZ2

Figure 1. The front of the EZ2 instrument.

- 1 Touchscreen
- 2 Hood
- 3 Power button
- 4 USB port

**Note:** Two additional USB ports are located on the rear of the touchscreen (not shown).

**Figure 2. The back of the EZ2 instrument.**

- 5 RJ-45 Ethernet port
- 6 Power cord socket — including instrument main fuse
- 7 Cooling air outlets
- 8 Instrument type plate

### 3.2.1 Touchscreen

The EZ2 has a 10.1-inch color touchscreen display with a resolution of 1280 x 800 pixels. The Graphical User Interface (GUI) is displayed on the touchscreen, allowing the user to operate the instrument, set up and start runs, perform maintenance procedures, monitor the instrument status, change settings, and view reports.

Figure 3. EZ2 touchscreen display with GUI.

### 3.2.2 Hood

The EZ2 hood protects the interior of the instrument from external contamination during protocol runs. Additionally, the hood protects operators from UV irradiation when the UV LED light is in use, that is during decontamination procedures.

The hood must be closed to enable a protocol run to start. The hood is locked at the beginning of a run and remains locked throughout the duration of the run. This protects users from moving parts on the worktable. The hood can be manually opened to gain access to the worktable when no protocol is running. During operation of the EZ2, the hood must remain closed and should only be opened when you are instructed to do so by the software.

<p><b>WARNING</b></p> 	<p><b>Moving parts</b></p> <p>To avoid contact with moving parts during the operation of the EZ2, the instrument must be operated with the hood closed.</p> <p>If the hood sensor or lock is not functioning properly, contact QIAGEN Technical Services.</p>
---	---


Figure 4. The EZ2 with the hood closed.

---

### 3.2.3 USB ports

The EZ2 has 3 USB ports. One is located next to the power button, on the front of the instrument. Two are located on the back of the touchscreen display.

The USB ports allow you to connect a USB drive to the EZ2. A USB drive that is connected to the instrument can be used to, for example, transfer report files. For more information on saving reports, refer to the instructions provided in the Saving a run report section.

In addition, you can use the USB drive to upload protocols, or update the software, if you have the relevant files on your USB drive. For more information on uploading protocols, refer to the Managing protocols section. For more information on software updates, refer to the [Updating software](#) section.

The bar code scanner is connected to the EZ2 using any one of the 3 available USB ports.

**Important:** Use only the USB flash drive provided by QIAGEN. Do not connect other USB flash drive devices to USB ports

**Important:** Do not remove the USB drive while downloading or transferring data or software to or from the instrument.

**Important:** Always turn off the EZ2 in order to plug or unplug the Wi-Fi USB device. Plug-and-play of the Wi-Fi USB device while the instrument is turned on is not supported.

### 3.2.4 RJ-45 Ethernet port and power cord socket

The RJ-45 Ethernet port is located on the back of the instrument. The port is used to connect the EZ2 to a local area network. See image below highlighted in white.

The power cord socket is located on the back of the EZ2 and is used for connecting the instrument to a power outlet via the supplied power cord. See image below highlighted in blue.

Figure 5. Location of RJ-45 port and power cord socket.

<p><b>WARNING</b></p> 	<p><b>Electrical hazard</b></p> <p>Any interruption of the protective conductor (earth/ground lead) inside or outside the instrument or disconnection of the protective conductor terminal is likely to make the instrument dangerous.</p> <p><b>Lethal voltages inside the instrument</b></p> <p>When the instrument is connected to line power, terminals may be live and opening covers or removing parts is likely to expose live parts.</p>
---	--

<p><b>WARNING</b></p> 	<p><b>Damage to electronics</b></p> <p>Before turning the instrument on, make sure that the correct supply voltage is used.</p> <p>Use of incorrect supply voltage may damage the electronics.</p> <p>To check the recommended supply voltage, refer to the specifications indicated on the type plate of the instrument.</p>
---	---

<p><b>WARNING</b></p> 	<p><b>Risk of electric shock</b></p> <p>Do not open any panels on the EZ2.</p> <p><b>Risk of personal injury and material damage</b></p> <p>Only perform maintenance that is specifically described in this user manual.</p>
---	--

### 3.2.5 Ventilation openings

EZ2 ventilation openings allow the internal components of the instrument to be cooled.

<p><b>CAUTION</b></p> 	<p><b>Risk of overheating</b></p> <p>To ensure proper ventilation, maintain a minimum clearance of 10 cm at the sides and rear of the EZ2.</p> <p>Slits and openings that ensure the ventilation of the instrument must not be covered.</p>
---	---

### 3.2.6 Bar code scanner

The bar code scanner delivered with the instrument can be connected to the EZ2 using any one of the 3 USB ports. The scanner is used for reading Q-Card and sample bar codes. For more information on how to use the bar code scanner, refer to the [Using the bar code scanner](#) section.

<b>WARNING</b> 	<b>Risk of personal injury</b> Hazard Level 2 laser light: Do not stare into the light beam when using handheld bar code scanner.
---	--

### 3.3 Internal features of the EZ2


Figure 6. Interior of the EZ2.

- 1 Pipettor head
- 2 Cartridge rack
- 3 Tip rack
- 4 Magnet module
- 5 Camera

Internal features not marked in the picture:

- Heating system
- UV LED lamp
- Internal light

### 3.3.1 Pipettor head

The pipettor head is mounted above the worktable and moves in the Z-direction (meaning up and down) to reach the sample and reagent tubes on the worktable. The worktable itself moves in the Y-direction (meaning front to back) so that the pipettor head is above the appropriate position in the cartridge or tip rack during each action that is performed by the instrument.

The pipettor head contains 24 high-precision syringe pumps that are connected to tip adapters which can be attached to filter-tips. The syringe pumps operate simultaneously and can aspirate or dispense small volumes of liquid (20–1000  $\mu$ l) through the attached filter-tips.

Another component of the pipettor head is the piercing unit that is located behind the tip adapters. The piercing unit is a row of 24 metal spikes that puncture the foil sealing the reagent cartridges. During operation, the piercing unit, controlled by the EZ2, opens all wells of the reagent cartridges in a dedicated order. The pipettor head then automatically picks up filter-tips from the tip rack and performs aspiration and dispensation operations at different locations on the worktable before ejecting the tips back into the tip rack at the end of the run.

Figure 7. EZ2 pipettor head.

<p><b>WARNING</b></p> 	<p><b>Moving parts</b></p> <p>To avoid contact with moving parts during the operation of the EZ2, the instrument must be operated with the hood closed.</p> <p>If the hood sensor or lock is not functioning correctly, contact QIAGEN Technical Services.</p>
---	--

<p><b>WARNING</b></p> 	<p><b>Moving parts</b></p> <p>Avoid contact with moving parts during operation of the EZ2. Under no circumstances should hands be placed under the pipetting arm during movement. Do not attempt to remove any plasticware from the workdeck whilst the instrument is operating.</p>
---	--

### 3.3.2 Worktable

The EZ2 worktable contains two movable EZ2 Connect Cartridge Rack that hold all the labware required for a protocol run, and the heating system, which controls the temperature of samples during a run.

#### EZ2 Connect Cartridge Rack


Figure 8. EZ2 Connect Cartridge Rack with two cartridges inserted.


**Figure 9. Cartridge rack inside the instrument.**

The EZ2 Connect Cartridge Rack consists of two separate parts. The left cartridge rack is used for cartridges in positions from 1 to 12. The right cartridge rack is used for cartridges in positions 13 to 24. On the worktable, the cartridge rack is located behind the tip rack. It holds up to 24 reagent cartridges.

For more information on how to load the EZ2 Connect Cartridge Rack, refer to the Loading the cartridge rack section.

Sealed reagent cartridges are prefilled and contain the reagents needed for a protocol run. Each cartridge consists of 10 sealed reagent wells and 2 empty heating positions. One heating position is a well and the other is a slot that can hold a tube.

**Figure 10. An EZ2 cartridge.**

## EZ2 Connect Tip Rack


Figure 11. EZ2 Connect Tip Rack.

The EZ2 Connect Tip Rack consists of two separate parts. The left part of the tip rack is used for labware in positions from 1 to 12. The right part of the tip rack is used for labware in positions 13 to 24.

The EZ2 Connect Tip Racks are located at the front of the worktable. Each consists of four rows and 12 columns.

Individual positions in the EZ2 Connect Tip Rack are marked by engravings. During run setup, the User Interface gives instruction to load specific positions of the EZ2 Connect Cartridge Rack with filter tips in tip holders, sample tubes, enzymes in tubes, or elution tubes.

Figure 12. Tip holders and filter tips.

**Important:** Only use QIAGEN recommended tubes for elution.

For information on how to load the EZ2 Connect Tip Rack, refer to the section Loading the tip rack.

## Heating system

The heating system is located under the rear part of cartridge rack. If needed, it heats up the wells in positions 11 and 12 of the cartridge.


Figure 13. Heatable positions in the cartridge rack.

<p><b>WARNING</b></p> 	<p><b>Hot surface</b></p> <p>The heating system can reach temperatures of up to 95°C (203°F). Avoid touching it when it is hot.</p>
---	---

## Tray

The tray is located underneath the tip and cartridge racks. Its role is to prevent contamination of the EZ2 that could be caused by liquids that are unintentionally spilled. The tray can be removed and cleaned as described in the Daily maintenance section.

**Figure 14. The EZ2 tray.**

### 3.3.3 Magnet module

The EZ2 magnet module consists of magnets that are used to capture particles that are present in the liquid aspirated into the filter-tips.

**Figure 15. The EZ2 magnet module.**

### 3.3.4 Camera


Figure 16. Camera module.

The EZ2 Connect Fx version has a built-in camera that is used for load checks and reading bar codes from cartridges.

Load checks are done before a protocol run starts. The camera checks if all the labware has been loaded in the correct positions, tubes are opened if required. The results of the load check are shown on the screen. For more information on the load check, refer to the Load check section.

The camera also reads the cartridge 2D bar codes. Information that is collected from the Q-Card is included in run reports.

### 3.3.5 UV LED light

The EZ2 is equipped with a UV LED light for decontamination. During the maintenance decontamination procedure, the UV LED moves over the worktable.

**Note:** The hood must be closed before starting a maintenance procedure and is automatically locked during the procedure.

<b>WARNING</b> 	<b>UV radiation</b> Do not expose your skin to UV light from the UV LED lamp.
---	--

<b>WARNING</b> 	<b>Moving parts</b> To avoid contact with moving parts during the operation of the EZ2, the instrument must be operated with the hood closed.  If the hood sensor or lock is not functioning correctly, contact QIAGEN Technical Services.
---	---

### 3.3.6 Internal light

The EZ2 has a built-in LED light. The internal light illuminates the worktable and informs about the current status of the run. There are two modes of the LED light:

- Blinking light – indicates that operator activity is required (for example, when an error occurs).
- Constant light – the default setting, used in all other situations.

## 4 Installation Procedures

This section provides instructions on installation environment requirements as well as unpacking, installing, and packing the EZ2.

### 4.1 Installation environment

The EZ2 is a plug-and-play instrument. Unpacking and installation procedures are easy to follow, but a person who is familiar with laboratory equipment should oversee the installation.

#### 4.1.1 Site requirements

The EZ2 must be located out of direct sunlight, away from heat sources and away from sources of vibration and electrical interference. Refer to Section Technical Specifications for the operating conditions (temperature and humidity). The site of installation should be free of excessive drafts, excessive moisture and excessive dust and should not be subject to large temperature fluctuations.

Use a level workbench that is large enough and strong enough to accommodate the EZ2. Refer to the Technical Specifications section for the weight and dimensions of the EZ2. Ensure that the workbench is dry, clean and vibration-proof and has additional space for accessories.

The EZ2 must be placed within approximately 1.5 m of a properly grounded (earthed) AC power outlet. The power line to the instrument should be voltage regulated and surge protected. Ensure that the EZ2 is positioned so that it is easy to access the power connector at the back of the instrument and the power switch on the front at all times, and that it is easy to turn the instrument off and disconnect it.

**Note:** It is recommended to plug the instrument directly into its own power outlet and to not share the power outlet with another lab equipment.

<p><b>WARNING</b></p> 	<p><b>Explosive atmosphere</b></p> <p>The EZ2 instrument is not designed for use in an explosive atmosphere.</p>
---	--

<p><b>CAUTION</b></p> 	<p><b>Risk of overheating</b></p> <p>To ensure proper ventilation, maintain a minimum clearance of 10 cm at the sides and rear of the EZ2.</p> <p>Slits and openings that ensure the ventilation of the instrument must not be covered.</p>
---	---

<p><b>WARNING</b></p> 	<p><b>Risk of personal injury and material damage</b></p> <p>The EZ2 is too heavy to be lifted by one person. To avoid personal injury or damage to the instrument, do not lift the instrument alone.</p>
---	---

<p><b>CAUTION</b></p> 	<p><b>Damage to the instrument</b></p> <p>Direct sunlight may bleach parts of the instrument and cause damage to plastic parts.</p> <p>The EZ2 must be located out of direct sunlight.</p>
---	--

<p><b>CAUTION</b></p> 	<p><b>Damage to the instrument</b></p> <p>Do not use the EZ2 in the vicinity of sources of strong electromagnetic radiation (for example unshielded, deliberately operated high-frequency sources or mobile radio devices), because these can interfere with the proper operation.</p>
---	--

#### 4.1.2 Power requirements

The EZ2 operates at: 100–240 V AC +/- 10%, 50/60 Hz, 1000 VA.

Make sure the voltage rating of the EZ2 is compatible with the AC voltage available at the installation site.

<p><b>WARNING</b></p> 	<p><b>Damage to electronics</b></p> <p>Before turning the instrument on, make sure that the correct supply voltage is used.</p> <p>Use of incorrect supply voltage may damage the electronics.</p> <p>To check the recommended supply voltage, refer to the specifications indicated on the type plate of the instrument.</p>
---	---

<p><b>WARNING</b></p> 	<p><b>Electrical hazards</b></p> <p>Any interruption of the protective conductor (earth/ground lead) inside or outside the instrument or disconnection of the protective conductor terminal is likely to make the instrument dangerous.</p> <p>Intentional interruption is prohibited.</p> <p><b>Lethal voltages inside the instrument</b></p> <p>When the instrument is connected to line power, terminals may be live and opening covers or removing parts is likely to expose live parts.</p>
---	--

4.1.3 Grounding requirements

To protect operating personnel, the National Electrical Manufacturers' Association (NEMA) recommends that the EZ2 be correctly grounded (earthed). The instrument is equipped with a 3-conductor AC power cord that, when connected to an appropriate AC power outlet, grounds (earths) the instrument. To preserve this protection feature, do not operate the instrument from an AC power outlet that has no ground (earth) connection.

<p><b>WARNING</b></p> 	<p><b>Electrical hazards</b></p> <p>Any interruption of the protective conductor (earth/ground lead) inside or outside the instrument or disconnection of the protective conductor terminal is likely to make the instrument dangerous.</p> <p>Intentional interruption is prohibited.</p> <p><b>Lethal voltages inside the instrument</b></p> <p>When the instrument is connected to line power, terminals may be live and opening covers or removing parts is likely to expose live parts.</p>
---	--

## 4.2 Unpacking the EZ2

<p><b>WARNING</b></p> 	<p><b>Risk of personal injury</b></p> <p>The EZ2 is too heavy to be lifted by one person. To avoid personal injury or damage to the instrument, do not lift the instrument alone.</p>
--	---

The following items are delivered:

- EZ2 instrument
- Quick-start guide and safety instructions
- Left and right EZ2 Connect Cartridge Rack
- Left and right EZ2 Connect Tip Rack
- Left and right EZ2 Connect Tip Rack - Flip Cap Tube (for EZ2 Connect Fx)
- Power cord set
- USB drive
- Silicone grease
- Certificate of manufacture
- Hand-held bar code scanner (for EZ2 Connect Fx)

### To unpack the EZ2:

1. Before unpacking the EZ2, move the package to the site of installation and check that the arrows on the package point upward. In addition, check whether the package is damaged. In case of damage, contact QIAGEN Technical Service.
2. Open the top of the transportation box and remove the top layer (PE foam).
3. Remove the accessories box together with the PE foam at its perimeter.
4. Remove the outer carton box by holding at the cut-out area and lifting up the outer carton box.
5. Remove the two pieces of protector from the unit.
6. Move the unit to the workbench or trolley from the packaging. When lifting the EZ2, slide your fingers under the side of the instrument and keep your back straight.

**Important:** Two persons are required to lift the EZ2.

**Important:** Do not hold the touchscreen display while unpacking or lifting the EZ2, this might damage the instrument.

7. Remove the transport lock for P-Axis by pulling the transport lock from bottom.

**Figure 17. Transport lock.**

8. Remove gel from the unit.
9. Remove the transport lock for Y-Axis by pushing the bottom of transport lock towards rear direction and pull it out from the rear. There is a total of two transport locks for Y-Axis.
10. Check if the packing list document is included after unpacking the EZ2.
11. Read the packing list to check that you have received all items. If anything is missing, contact QIAGEN Technical Services.
12. Check that the EZ2 is not damaged and that there are no loose parts. If anything is damaged, contact QIAGEN Technical Services. Make sure that the EZ2 has equilibrated to ambient temperature before powering on.

13. Retain the package in case you need to transport the EZ2 in the future. Refer to the Packing the EZ2 section for more details. Using the original package minimizes the possibility of damage during transportation of the EZ2.

## 4.3 Installing the EZ2

This section describes important actions that must be performed before operating the EZ2. These actions include:

- Removal of the EZ2 accessories and shipping material.
- Installation of the AC power cord
- Installation of external barcode scanner
- Installation of the Wi-Fi adapter
- Initial configuration
- Camera exposure calibration
- If an installation qualification/operational qualification (IQ/OQ) is required in your laboratory setup, this service can be ordered together with the instrument. For details, please contact QIAGEN Technical Service.

### 4.3.1 Removal of the EZ2 accessories and shipping materials

1. Remove the power cord, barcode scanner, and the quick-start guide from the foam packing material on top of the EZ2.
2. Remove the USB flash drive, cartridge and sample/tip racks
3. Ensure all shipping foams, transport locks and other packing material have been removed as described in the Unpacking the EZ2 section.

### 4.3.2 Installation of the AC power cord

1. Remove the power cord from the foam packing material on top of the EZ2.  
**Note:** Only use the power cord that is supplied with the EZ2.
2. Ensure that the power switch is set to OFF.

Figure 18. Location of the power button.

3. Check that the voltage rating on the label at the back of the EZ2 matches the voltage available at the installation site.
4. Plug the power cord into the instrument power-cord socket.
5. Plug the power cord into a grounded power outlet.
6. Do not switch on power to the instrument at this point. Power needs to be OFF for subsequent installation of USB devices in the following steps.

<p><b>WARNING</b></p> 	<p><b>Damage to electronics</b></p> <p>Before powering ON the instrument, make sure that the correct supply voltage is used.</p> <p>Use of incorrect supply voltage may damage the electronics.</p> <p>To check the recommended supply voltage, refer to the specifications indicated in the type plate of the instrument.</p>
---	--

<p><b>WARNING</b></p> 	<p><b>Electrical hazard</b></p> <p>Any interruption of the protective conductor (earth/ground lead) inside or outside the instrument or disconnection of the protective conductor terminal is likely to make the instrument dangerous. Intentional interruption is prohibited.</p> <p><b>Lethal voltages inside the instrument</b></p> <p>When the instrument is connected to line power, terminals may be live and opening covers or removing parts is likely to expose live parts.</p>
---	--

#### 4.3.3 Installation of external barcode scanner (optional)

1. Remove the barcode scanner from the box.
2. Before powering on the instrument connect the barcode scanner to any one of the 3 USB ports located either on the front of the instrument or on the rear of the touchscreen.

<b>WARNING</b> 	<b>Risk of personal injury</b> Hazard Level 2 laser light: Do not stare into the light beam when using handheld bar code scanner.
---	--

#### 4.3.4 Installation of Wi-Fi adapter (optional when using QIAsphere)

1. Remove the Wi-Fi adapter from the packaging
2. Before turning on the instrument, plug the Wi-Fi adapter into one of the 3 USB ports located on either the front of the instrument or on the rear of the touchscreen.

**Note:** It may be more convenient to use the USB ports on the rear of the touchscreen.

#### 4.3.5 Initial configuration of the EZ2

1. Follow steps 1 to 5 from the Starting the EZ2 section.
2. To log in for the first time, enter **Admin** in both **User ID** and **Password** fields then press **Log in**. After this login, you will have administrator rights with the option to set up other users. Refer to the Managing users section for further details.

**Note:** After first log in using the Admin account, the password should be changed according to the Password policy described in the Managing user security section and as per instruction in the Changing password section.

3. From the **Configuration** menu under the **System** tab, you can modify the Device Name, Date and Time fields. You also find information on Software version, Serial number, and Firmware version. For further details, see the [Setting basic system data](#) section.


Figure 19. The System configuration tab.

- You may adjust settings of the EZ2 Connect according to your preferences through the **Instrument Settings** tab under the **Configuration** menu. For further details see the Changing Instrument settings section.


Figure 20. Details on the Instrument setting tab.

#### 4.3.6 Perform camera exposure calibration (for EZ2 Connect Fx)

**Note:** Only Administrators can perform the camera exposure calibration.

<p><b>WARNING</b></p> 	<p><b>Moving parts</b></p> <p>To avoid contact with moving parts during the operation of the EZ2, the instrument must be operated with the hood closed.</p> <p>If the hood sensor or lock is not functioning correctly, contact QIAGEN Technical Services.</p>
---	--

1. Before first use in the final installation site, the camera exposure calibration must be performed.
2. From the **Maintenance** menu under the **Camera LED** tab follow the instructions on the user interface.


Figure 21. Camera exposure calibration screen.

3. During the calibration routine, the following will be displayed:


Figure 22. Camera exposure calibration is in progress message.

**Note:** For support or technical advice regarding this please contact QIAGEN Technical Service Department or see our Technical Support Center at <http://www.qiagen.com/service-and-support/technical-support>.

## 5 Operating Procedures

This section describes how to operate the EZ2 Instrument.

Before proceeding, it is recommended that you familiarize yourself with the features of the instrument as described in the External features of the EZ2 and Internal features of the EZ2 sections.

The EZ2 is intended to be used only in combination with QIAGEN kits indicated for use with the EZ2 instrument for the applications described in the kit handbooks.

The hood of the EZ2 must remain closed and will automatically lock during operation of the instrument. Only open the hood when instructed to do so by the instruction for use.

The worktable of the EZ2 instrument moves during operation of the instrument. Never open the EZ2 hood while the instrument is operating.

<b>WARNING</b> 	<b>Moving parts</b> To avoid contact with moving parts during the operation of the EZ2, the instrument must be operated with the hood closed.  If the hood sensor or lock is not functioning correctly, contact QIAGEN Technical Services.
---	---

<b>WARNING</b> 	<b>Moving parts</b> Avoid contact with moving parts during operation of the EZ2. Under no circumstances should hands be placed under the pipetting arm during movement. Do not attempt to remove any plasticware from the workdeck whilst the instrument is operating.
---	---

<b>WARNING</b> 	<b>Risk of personal injury and material damage</b> Do not attempt to move the EZ2 during operation.
---	--

<b>WARNING/ CAUTION</b> 	<b>Risk of personal injury and material damage</b> Improper use of the EZ2 may cause personal injuries or damage to the instrument. The EZ2 must only be operated by qualified personnel who have been appropriately trained. Servicing of the EZ2 must only be performed by QIAGEN Field Service specialists.
--	---

<b>CAUTION</b> 	<b>Damage to the instrument</b> Avoiding spilling water or chemicals onto the EZ2. Instrument damage caused by water or chemical spillage will void your warranty.
---	---

<b>WARNING</b> 	<b>Risk of fire or explosion</b> When using ethanol or ethanol-based liquids on the EZ2, handle such liquids carefully and in accordance with the required safety regulations. If liquid has been spilled, wipe it off and leave the EZ2 hood open to allow flammable vapors to disperse.
--	--

<b>WARNING</b> 	<b>Risk of explosion</b> The EZ2 is intended for use with reagents and substances supplied with QIAGEN kits as outlined in respective Information for use. Use of other reagents and substances may lead to fire or explosion.
---	---

<b>CAUTION</b> 	<b>Damage to the instrument</b> Ensure that the EZ2 is switched off before you manually move the mechanical components of the instrument.
---	--

<p><b>CAUTION</b></p> 	<p><b>Damage to the instrument</b></p> <p>Do not lean against the instrument or touchscreen.</p>
---	--

<p><b>WARNING</b></p> 	<p><b>Samples containing infectious agents</b></p> <p>Some samples used with the EZ2 instrument may contain infectious agents. Handle such samples with the greatest of care and in accordance with the required safety regulations.</p> <p>Always wear safety glasses, gloves, and a lab coat.</p> <p>The responsible body (for example, a laboratory manager) must take the necessary precautions to ensure that the surrounding workplace is safe, and that the instrument operators are suitably trained and not exposed to hazardous levels of infectious agents as defined in the applicable Material Safety Data Sheets (MSDSs) or the OSHA<sup>1,*</sup> ACGIH,<sup>†</sup> or COSHH<sup>‡</sup> documents.</p> <p>Venting for fumes and disposal of wastes must be in accordance with all national, state, and local health and safety regulations and laws.</p>
---	---

\* OSHA – Occupational Safety and Health Organization (United States of America)

† ACGIH – American Conference of Government Industrial Hygienists (United States of America)

‡ COSHH – Control of Substances Hazardous to Health (United Kingdom)

<p><b>CAUTION</b></p> 	<p><b>Hazardous materials and infectious agents</b></p> <p>The waste contains samples and reagents. This waste may contain toxic or infectious material and must be disposed of properly. Refer to your local safety regulations for proper disposal procedures.</p>
---	--

<b>WARNING</b> 	<b>Hot surface</b> The heating system can reach temperatures of up to 95°C (203°F). Avoid touching it when it is hot.
---	--

<b>WARNING</b> 	<b>UV radiation</b> Avoid looking directly into UV light. Do not expose your skin to UV light.
---	---

<b>WARNING</b> 	<b>Risk of personal injury</b> Hazard Level 2 laser light: Do not stare into the light beam when using handheld bar code scanner.
---	--

## 5.1 General information

The EZ2 is operated using a touchscreen display. Some of the actions that you can do using the Graphical User Interface (GUI) are:

- Guided run setups and maintenance
- Checking the status of the run and the instrument
- Generating and save run reports
- Changing instrument settings to customize your EZ2

Each screen of the GUI consists of three elements: the toolbar, the main content, and the footer.


Figure 23. The home screen (for EZ2 Connect Fx)


Figure 24. The home screen (for EZ2 Connect)

## Toolbar

The toolbar is used to access the main sections of the EZ2 software, to check the connection status of the LAN and Wi-Fi, and to log out of the application. The toolbar contains the same buttons throughout the whole application, but some buttons are disabled during protocol setup, protocol runs, and maintenance procedures.

Element	Description
 Setup	Tap this button to go to the <b>Home</b> screen, where you can start the setup process of a protocol run.
 Tools	Tap this button to go to the <b>Tools</b> section, where you can configure data exchange, and access maintenance procedures.
 Configuration	Tap this button to go to the <b>Configuration</b> section, where you can change settings, manage users, configure network, and QlAsphere Base connections, and change your password.
 LAN	Shows the status of the LAN connection.
 Wi-Fi	Shows the status of the Wi-Fi connection.
 Logout	Tap this button to log out.

Figure 25. Description of elements in the User Interface.

## Main content

The part of the screen where the main content of each view is displayed.

## Footer

The footer shows the current date and time, the software version, and the name of user that is currently logged in.

### 5.1.1 Entering text and numbers

An on-screen keyboard is used for entering text in editable fields of the EZ2 instrument software. To access the keyboard, tap the field that you want to edit. The keyboard appears.





Figure 26. The on-screen keyboard.

The default layout of the keyboard is lowercase QWERTY, with numbers from 1 to 0, commonly used special characters, a space bar, a **Shift**  key, a **Caps Lock**  key, and the **Special Characters**  key. To enter a character, tap the relevant letter, number or special character on the keyboard. To enter one uppercase letter, tap **Shift** . To enter multiple consecutive uppercase letters, tap **Caps Lock** . Tap **Caps Lock**  again to turn off the uppercase input mode. To see special characters, tap **Special Characters** . To go back to the alphabetic characters, tap **ABC** .


Figure 27. The on-screen keyboard in special character mode.

To remove one character to the left of the cursor, tap **Backspace** . To delete all characters from the field, tap **Clear All** .

Some fields have requirements or restrictions that have to be followed. If the entered text does not match the requirements of the field, an error message appears, and the input is not accepted. To proceed, modify the text, so that it follows the requirements.


Figure 28. Example of text field validation.

If certain characters are not allowed in a field, they are deactivated on the keyboard and they cannot be entered in the field.

## 5.2 Starting the EZ2

To start the EZ2, switch on the instrument and log in, then follow the steps below:

1. Close the hood of the instrument.
2. Press the power button. The startup screen appears, and the instrument is initialized. After the initialization is done, the **Login** screen appears.

Welcome to EZ2  
Log in to continue

Service login

User ID:

Password:

Log in

9/1/2021 9:26

**Figure 29. The login screen.**

3. Tap the **User ID** field and enter your user ID using the on-screen keyboard. For more information on using the on-screen keyboard, refer to the Entering text and numbers section.  
**Note:** If this is the first time the EZ2 is switched on, enter the default user ID, that is **Admin**.
4. Tap the **Password** field and enter your password using the on-screen keyboard.  
**Note:** If this is the first time the EZ2 is switched on, enter the default password, that is **Admin**.
5. Tap **Log in**. If the credentials you entered are correct, the **Home** screen appears. If the credentials you entered are incorrect, an error message is shown.  
**Note:** If the number of failed login attempts exceeds the limit set by your administrator, your account is locked. To activate your account, contact your administrator. If you are the only administrator and your account is locked, contact QIAGEN Technical Support.

### 5.3 Configuring the EZ2

Administrators of the EZ2 can adjust instrument settings, manage users, upload and delete protocols, update software and configure network connectivity.

**Note:** Users whose role is defined as **Operator** do not have access to settings and configuration features of the software.

### 5.3.1 Setting basic system data

To set the device name, date, time, date format and time zone, follow the steps below:

**Note:** Only Administrators can change system settings.

1. Tap the **Configuration** icon on the toolbar and tap the **System configuration** tab.


Figure 30. System configuration tab.

2. To set the device name, tap the **Device name** field, and enter a name using the on-screen keyboard. To save the name, tap **Accept**.

**Note:** The device name cannot be longer than 24 characters. The name cannot contain special characters or spaces.

3. To set the date, tap the **Date** field, and select the date using the date picker. To change the month, the year, or both, use the left and right arrows on either side of the month and year label. To select a specific date, tap the day on the calendar. To confirm your selection, tap **Accept**.


Figure 31. Select date dialog.

4. To set the time, tap the **Time** field and enter the time using the on-screen keyboard.
5. To set the date format, tap the **Date format** drop-down list and select one of the listed formats.
6. To use the 24-hour time format, select the **24-hour time format** checkbox. To use the 12-hour format, clear the **24-hour time format** checkbox.
7. To set the time zone automatically using the QIASphere Base connection, select the **Set time zone automatically** checkbox.
 

**Note:** This setting requires an active QIASphere Base connection. For more information on configuring QIASphere Base, refer to Configuring QIASphere Base. If QIASphere Base is not configured, the **Set time zone automatically** checkbox is disabled.
8. To save the settings, tap **Accept**.


### 5.3.2 Changing Instrument settings

You can adjust the settings of the EZ2 according to your preferences.

**Note:** Only Administrators can change Instrument settings.

To modify Instrument settings, follow the steps below:

1. Tap the **Configuration** icon on the toolbar and tap the **Instrument settings** tab.

**Figure 32. The Instrument settings tab.**

2. To adjust the volume of the audio, use the **Audio volume** slider. A sound is played each time you change the volume. To mute the device, tap **Audio** , or move the slider to the far-left position.
3. To adjust the brightness of the touchscreen, use the **Display brightness** slider. To set the lowest brightness, tap **Display brightness** , or move the slider to the far-left position.  
**Note:** The lowest brightness setting turns off the backlight of the display, which makes the content on the touchscreen barely visible.
4. To adjust the brightness of the internal light, use the **Internal light brightness** slider. To turn the internal light off, tap **Internal light brightness** , or move the slider to the far-left position.
5. Tap **Accept** to save the modified settings, or tap **Cancel** to restore previously saved settings.

### 5.3.3 Managing users

EZ2 user management allows you to create and edit user accounts with two different roles: administrator and operator.

**Note:** The EZ2 Connect Fx version has an additional advanced user role who can create shortcuts.


When you use the EZ2 for the first time, a default user (Admin) is pre-installed and configured.

**Note:** User management is only available to users with the administrator role.

## Adding a new user

**Note:** Only administrators can add new users.

1. Tap the **Configuration** icon on the toolbar and tap the **User management** tab. Existing users are shown in the table.

User ID	First Name	Last Name	Role	Status
Admin	Admin	Admin	Administrator	Active
Operator	Operator	Operator	Operator	Active

**Figure 33. The User Management tab.**

2. To add a new user, tap **New**. The **Create new account** dialog is shown.

### Create new account

Anonymous ID:	User ID:	
<input type="text" value="d6da443f"/>	<input type="text"/>	
First name:	Last name:	
<input type="text"/>	<input type="text"/>	
User role:	User status:	
<input type="text" value="Administrator"/>	<input type="text" value="Active"/>	
Password status: Not set		
<input type="button" value="Cancel"/>	<input type="button" value="Set password"/>	<input type="button" value="Accept"/>

Figure 34. Create new account dialog.

3. Tap the respective boxes to enter the **User ID**, **First name**, and **Last name** using the on-screen keyboard. The **Anonymous ID** is generated automatically and is used to identify users in reports and audit trails.

**Note:** Only administrators can view the full details of a user account. As a result, only administrators can identify users based on their anonymous IDs.

4. Select the **User role** and **User status** from their respective drop-down lists.
5. Tap **Set password**. The **Set user password** dialog appears. Enter the password in the **New password** field and again in the **Retype password** field. The password must meet the criteria displayed in the dialog. To modify the password policy, refer to the *Managing user security* section.

## Set user password

Password must consist of 8-40 characters and have at least one of:

- both upper-case and lower-case letters
- numerical digits
- special characters such as @, # or \$

New password:

Retype password:

Cancel

Set password

Figure 35. Set user password dialog.

6. Tap **Set password**. If the passwords match, the password status is set to active.

**Note:** Users can change their password later, see the Changing password section.

7. Tap **Accept**. The user is added.

**Note:** It is recommended to create at least one more account with administrator role in addition to the default Admin account. In case one administrator is locked due to wrong credentials, the other administrator can unlock the account.

### Editing an existing user account

**Note:** Only administrators can edit user accounts.

1. Tap the **Configuration** icon on the toolbar and tap the **User management** tab. Existing users are shown in the table.


Figure 36. The User Management tab.

2. Tap the table row that corresponds to the user whose profile you want to edit. The **Edit user account** dialog appears.

### Edit user account

Anonymous ID:

First name:

User role:

User ID:

Last name:

User status:

Password status: Active

Figure 37. Edit user account dialog.

3. To edit the **First name** or **Last name**, tap the respective boxes and modify the contents using the on-screen keyboard.

- To change the **User role** or **User status**, tap the respective drop-down lists and select an option.
- To change the password of the user, tap **Change password**. The **Change user password** dialog appears.

**Change user password**

Password must consist of 8-40 characters and have at least one of:

- both upper-case and lower-case letters
- numerical digits
- special characters such as @, # or \$

New password:

Retype password:

**Figure 38.** Change user password screen.

- Enter the password in the **New password** field and again in the **Retype password** field. The password must meet the criteria displayed in the dialog. To modify the password policy, refer to the Managing user security section.
- Tap **Change password**.
- To save the changes made to the user account, tap **Accept**.

**Note:** For cybersecurity reasons, an operator must change the password provided by the administrator upon first login.

### Deleting a user

**Note:** Only administrators can delete user accounts.

- Tap the **Configuration** icon on the toolbar and tap the **User management** tab. Existing users are shown in the table.


Figure 39. The User Management tab.

2. Tap the table row that corresponds to the user you want to delete.
3. Tap **Delete**. A confirmation dialog is shown.



Figure 40. Delete user confirmation dialog.

4. To confirm the account deletion, tap **Yes**.

## Deactivating a user

**Note:** Only administrators can deactivate user accounts.

1. Tap the **Configuration** icon on the toolbar and tap the **User management** tab. Existing users are shown in the table.

User ID	First Name	Last Name	Role	Status
Admin	Admin	Admin	Administrator	Active
Operator	Operator	Operator	Operator	Active

**Figure 41. The User Management tab.**


2. Tap the table row that corresponds to the user that you want to deactivate.
3. Tap **Edit**.
4. Tap the **User status** drop-down list and select **Inactive**.
5. Tap **Accept**.

### 5.3.4 Changing password

All active users can change their own passwords. Additionally, administrators can change the passwords of other users. For more information on how to change passwords of other users, refer to the Editing an existing user account section. **Note:** For cybersecurity reasons, an operator must change the password provided by the administrator upon first login.

To change your own password, follow the procedure below.

1. Tap the **Configuration** icon on the toolbar and tap the **Change password** tab.

The screenshot shows the 'Change password' tab in the QIAGEN configuration interface. At the top, there is a toolbar with icons for Setup, Tools, Configuration (highlighted), Service, LAN, Wi-Fi, and Logout. Below the toolbar is a navigation bar with tabs for User management, User security, Change password (selected), Settings, LAN, Wi-Fi, and Q-Base. The main content area displays password requirements: 'Password must consist of 8-40 characters and have at least one of: - both upper-case and lower-case letters - numerical digits - special characters such as @, # or \$'. There are three input fields labeled 'Current password:', 'New password:', and 'Retype password:'. At the bottom, there are 'Cancel' and 'Change password' buttons. The footer shows the date '2019/12/18 14:33', the software version '0.1.1-930-gb4101dd3', and the user 'Admin Admin'.

Figure 42. Change password tab.

2. Tap the **Current password** field and enter your current password using the on-screen keyboard.
3. Tap **Accept**.
4. Enter the password in the **New password** field and again in the **Retype password** field. The password must meet the criteria displayed on the screen.
5. Tap **Change password**.

### 5.3.5 Managing user security

Administrators can change the password policy for all user roles, as well as change the settings for password expiration and limit incorrect login attempts.

To modify settings related to user security, follow the steps below:

1. Tap the **Configuration** icon on the toolbar and tap the **User security** tab.

The screenshot shows the 'User security' configuration screen. At the top, there is a toolbar with icons for Setup, Tools, Configuration, Service, LAN, Wi-Fi, and Logout. Below the toolbar is a navigation bar with tabs for User management, User security, Change password, Settings, LAN, Wi-Fi, and Q-Base. The 'User security' tab is selected. The screen displays several configuration options: 'User role' is set to 'Administrator'; 'Password policy' is set to 'Strong policy'; 'Enable password expiration' is unchecked; 'Password expires every' is set to 90 days; 'Account locks after' is set to 3 incorrect login attempts. On the right side, there is a section for 'Password policies' with three options: 'No policy' (Any password up to 40 characters (could be empty)), 'Standard policy' (Password must consist of 8-40 characters), and 'Strong policy' (Password must consist of 8-40 characters and have at least one of: - both upper-case and lower-case letters, - numerical digits, - special characters such as @, # or \$). At the bottom, there are 'Cancel' and 'Accept' buttons. The footer shows the date and time (2019/12/18 14:44), software version (0.1.1-930-gb4101dd3), and the user name (Admin Admin).

Figure 43. User security tab.


2. To change the password policy, select the user role whose policy you want to change from the **User role** drop-down list.
3. Select an option from the **Password policy** drop-down list. There are three password policy options:
  - o No policy — there are no requirements for the password. The password length can be between 0 and 40 characters.
  - o Standard policy — the password length must be between 8 and 40 characters.
  - o Strong policy — the password length must be between 8 and 40 characters, and must contain either uppercase and lowercase letters, numerical digits or special characters.
4. To enable password expiration, select the **Enable password expiration** checkbox.
5. To specify the number of days after which user passwords expire, enter a value in the **Days** field.
6. To specify the number of incorrect login attempts after which a user account is locked, enter a value in the **Incorrect login attempts** field. It is recommended to set to at least 2 incorrect attempts. Otherwise, one typographical error will lock your account.
7. Tap **Accept** to save all changes, or tap **Cancel** to restore previously saved settings.

### 5.3.6 Installing new protocols

**Note:** Only administrators can install new protocol packages.

**Note:** Protocols tab is only available to users with administrator role.

1. Plug in a USB flash drive with a protocol package.
2. Tap **Protocols**.

**Figure 44. Protocols tab.**

3. Tap **Check for updates**.
4. Start the installation by choosing the protocol package and tap **Accept**.
5. After installation is finished, restart the device.

### 5.3.7 Updating software

**Note:** Only administrators can update the software.

**Note:** System Configuration is only available to users with administrator role.

1. Plug in a USB flash drive with the newest software version.
2. Tap the **Configuration** icon on the toolbar and tap the **System configuration** tab.


Figure 45. System Configuration tab.

3. Tap **Check for updates**.
4. Start the installation by choosing the software package and tap **Accept**.
5. After installation is finished, follow the restart instruction of the machine.

## 5.3.8 Configuring network and QIASphere Base connections

### Configuring a LAN connection

**Note:** LAN tab is only available to users with administrator role.

1. Tap the **Configuration** icon on the toolbar and tap the **LAN** tab.

The screenshot displays the QIAGEN configuration interface. At the top, there is a toolbar with icons for Setup, Maintenance, Configuration (highlighted), and Data. To the right are icons for Network and Logout. Below the toolbar is a horizontal menu with tabs for System configuration, Instrument settings, Protocols, User management, User security, LAN (selected), Wi-Fi, and QIASphere. The LAN configuration panel contains several input fields: IP address, Subnet mask, DNS server, Gateway, and Hardware address. A checkbox for 'DHCP enabled' is checked. At the bottom of the panel are 'Cancel' and 'Accept' buttons. The footer shows the date and time '8/25/2021 15:02' and the user 'Admin Admin'.

**Figure 46. The LAN tab.**

2. Enter **IP address**.
3. Enter **Subnet mask**.
4. Tap **Accept**.

## Configuring a Wi-Fi connection

**Note:** Only administrators can configure Wi-Fi settings.


**Note:** Make sure the Wi-Fi USB device was plugged-in before the machine was started.

The EZ2 can be connected to a Wi-Fi network using the Wi-Fi adapter inserted into one of the USB ports. The Wi-Fi adapter supports the 802.11b, 802.11g, and 802.11n Wi-Fi standards, and WEP, WPA-PSK, and WPA2-PSK encryption.

The SSID of the network must be visible. It is not possible to connect to a network with a hidden SSID.


To configure the Wi-Fi connection, follow the steps below:

1. Tap the **Configuration** icon on the toolbar and tap the Tap **Wi-Fi** tab.

**Figure 47. The Wi-Fi tab.**

2. To search for available networks, tap **Scan**. The signal strength, encryption type and name of each network is shown.

**Figure 48. List of available Wi-Fi networks.**

3. To select a network, tap the row that corresponds to the network you want to select.
4. To connect to the selected network, tap **Connect**.
5. Enter the network password using the on-screen keyboard and tap **Connect**. The instrument connects to the network, and the connection status is updated on the screen.

**Note:** If you failed to connect to the network several times despite entering the correct password, it is recommended to restart the instrument.

**Note:** In case of errors, refer to the Troubleshooting section for more information.


Figure 49. Successful connection to a Wi-Fi network.

- To disconnect from a network, tap the row that represents the network that the instrument is connected to. Then, tap **Disconnect**.
- To view additional information about the network, or to change its password, tap **Advanced**. To go back to the Wi-Fi tab, tap **Back**. To change the password, tap **Change password**. To disconnect from the network, tap **Disconnect**.


Figure 50. Network settings dialog.

## Configuring QIASphere Base connection

For information on how to configure QIASphere Base, refer to the QIASphere Base user manual.

For more details on the QIASphere Base network configuration and how to connect both devices with each other, refer to the QIASphere Base user manual which is available at [www.qiagen.com](http://www.qiagen.com).

Only users assigned the role Administrator can change the network configuration. It is recommended to consult your network administrator when configuring the network. For communication with QIASphere Base, the outbound TCP port 443 (https) is used; ping is supported.

1. Tap the **Configuration** icon on the toolbar and tap the **QIASphere** tab.


Figure 51. The QIASphere tab.

2. Enter the QIASphere Base URL.
3. Enter the QIASphere Base password.
4. Tick the **Enable QIASphere communication** box.
5. Tap **Accept**.
6. Tap **Test connection**.

**Note:** If you have any trouble regarding the QIASphere, please use the QIASphere Base user manual or contact the QIAGEN support.

## 5.4 Setting up a protocol run

To setup a protocol run, press the **Setup** tab on the **Home** screen.

Commonly used QIAGEN protocols are pre-installed on the EZ2 upon delivery. You can download all available QIAGEN protocols from the QIAGEN website.

The EZ2 software will guide you through the protocol run setup process. The steps of the protocol wizard vary, depending on the selected protocol and may differ from the figures included in this section.

**Important:** Before starting a protocol run, read the relevant QIAGEN kit handbook and follow the instructions on sample storage, handling and preparation provided in the handbook.

### 5.4.1 Selecting application type

To begin setting up a protocol run, tap an application in the **Applications** pane, a favorite protocol in the **Favorites** pane, or a recently used protocol in the **Recently Used** pane.


Figure 52. Select an application screen.

The software will automatically proceed to the next screen. Depending on the scanned Q-Card barcode, the software will automatically skip selections screens if the required information is provided by bar code scanning.

### 5.4.2 Selecting a kit

To select a kit that is to be used during the protocol run, tap a kit name in the **Select kit** pane.

**Note:** If there is only one available option, it is selected automatically.


Figure 53. Select kit step.

After you select a kit, the **Next** button is enabled. To proceed to the **Select protocol** step, tap **Next**.

To cancel the protocol, tap **Cancel**.

### 5.4.3 Selecting a protocol

To select a protocol that is compatible with the kit selected in the **Select kit** step, tap a protocol in the **Select protocol** pane.

**Note:** If there is only one available option, it is selected automatically.


Figure 54. Select protocol step.


After you select a protocol, the **Next** button is enabled. To proceed to the **Define parameters** step, tap **Next**. If there are no parameters to define, you will be taken to the **Select sample positions** step.

To return to the previous screen, tap **Back**.

To cancel the protocol, tap **Cancel**.

#### 5.4.4 Defining parameters

To set values for protocol parameters, tap the box next to each parameter in the **Define parameters** pane.

**Figure 55. Define parameters step.**

Additional information about the protocol is also displayed on the screen, e.g. references into the respective kit handbook with respect to storage, handling, and pretreatment (if required).

**Note:** Ensure you have read the respective information in the kit handbook.

To proceed to the **Select sample positions** step, tap **Next**.


To return to the previous screen, tap **Back**.

To cancel the protocol, tap **Cancel**.

#### 5.4.5 Selecting sample positions

To select the positions of your samples, tap the relevant rows on the worktable diagram or tap the corresponding row numbers underneath the diagram. The selected positions are highlighted. To select or deselect all positions, use the **Select all** toggle.

**Figure 56. Select sample positions step.**

After you select at least one sample position, the **Next** button is enabled. To proceed to the **Enter sample IDs** step, tap **Next**.

To return to the previous screen, tap **Back**.

To cancel the protocol, tap **Cancel**.

## 5.4.6 Entering sample IDs

The format of a sample ID is **YYYY-MM-DD\_hh-mm\_XX**, where the first 16 characters represent the current date and time, and the XX is the sample number. To generate the Sample IDs, you must press the **Generate missing sample IDs** button.

Enter sample IDs

Position	Sample ID	Note (optional)
1	2021-09-07_09-01_01	
2	2021-09-07_09-01_02	
3	2021-09-07_09-01_03	
4	2021-09-07_09-01_04	
5	2021-09-07_09-01_05	
6	2021-09-07_09-01_06	
7	2021-09-07_09-01_07	
8	2021-09-07_09-01_08	
9	2021-09-07_09-01_09	

Step 5 of 7

**Application:** DNA  
**Kit:** DNA FFPE  
**Protocol:** DNA FFPE  
**Elution volume:** 100µl  
**Rack type:** TipRack  
**Estimated run time:** 2 hrs 55 min 0 sec  
**Number of samples:** 24

**Pretreatment:** Pretreat the samples following the instructions in the kit handbook.


Cancel      Generate missing sample IDs      Back      Next

9/7/2021 9:02      Admin Admin

Figure 57. Enter sample IDs step.

## Modifying a sample ID

To modify a sample ID, tap it and use the on-screen keyboard to change the text.

The image shows a dialog box titled "Enter sample ID for sample at position 1". Inside the dialog, there is a text input field containing the text "Sample1". To the right of the text field are two icons: a blue 'X' for clearing the field and a grey back arrow for navigating back. Below the text field is a full QWERTY on-screen keyboard. At the bottom right of the dialog, there are two buttons: a grey "Cancel" button and a dark blue "Accept" button.

Figure 58. Changing the sample ID.

There is an 80-character limit for the **Sample ID** field. Tap **Accept** to save the changes, or tap **Cancel** to go back to the **Enter sample IDs** screen.

### Adding a note to a sample

Optionally, you can add a note to each sample. Tap the **Note (optional)** box next to the relevant sample ID, and use the on-screen keyboard to enter the note.

The screenshot shows a mobile application interface for entering a note. At the top, a grey header bar contains the text "Enter note for sample at position 1". Below this is a light blue text input field with a vertical cursor on the left and a clear button (a square with a left-pointing arrow and an 'x') on the right. Underneath the input field is a full QWERTY on-screen keyboard with light blue keys. The keyboard includes a home key on the bottom left, a spacebar, and a backslash key on the bottom right. At the bottom right of the screen, there are two buttons: a light blue "Cancel" button and a dark blue "Accept" button.

Figure 59. Adding a note to a sample.

There is an 80-character limit for the **Note (optional)** field. Tap **Accept** to save the changes, or tap **Cancel** to go back to the **Enter sample IDs** screen.

To proceed to the **Load the cartridge rack** step, tap **Next**.

To return to the previous screen, tap **Back**.

To cancel the protocol, tap **Cancel**.

## 5.4.7 Loading the EZ2 Connect Cartridge Rack

Remove the cartridge rack from the instrument and place it safely onto the bench. See the Loading and unloading the cartridge rack section.

The **Load the cartridge rack** step features instructions on how to load the cartridge rack. To highlight the well on the diagram and the corresponding row in the table, tap the well or the table row on the instrument touchscreen.

**Important:** Read the instructions carefully before loading the rack, and make sure to follow all directions including those in the respective kit handbook.

Well	Action
11	Load an open sample tube containing 4 discs of dried blood into the heating block.

Before you proceed to the instructions from the table above, invert the reagent cartridges 3 times to mix the magnetic particles. Place the cartridges in the same positions as the ones that you chose in the Select sample positions screen, then tap the cartridges to deposit the reagents at the bottom of their wells.

Step 6 of 7

12  
11  
10  
9  
8  
7  
6  
5  
4  
3  
2  
1

Cancel      Back      Next

8/31/2020 10:45      Software version: 0.2.9-93-g1d33f49c      Admin Admin

Figure 60. Load the cartridge rack step.

To proceed to the **Load the tip rack** step, tap **Next**.

To return to the previous screen, tap **Back**.

To cancel the protocol, tap **Cancel**.

### 5.4.8 Loading the EZ2 Connect Tip Rack

Remove the tip rack from the instrument and place it safely onto the bench. See the Loading and unloading the tip rack section.

The **Load the tip rack** step features instructions on how to load the tip rack. To highlight a position on the tip rack diagram and the corresponding row in the table, tap the position or the table row.

**Important:** Read the instructions carefully before loading the rack, and make sure to follow all directions including those in the respective kit handbook.

Row	Labware	Content
A	2.0 ml sample tube (ST)	Sample
C	Tip holder	1000 µl tip
D	1.5 ml elution tube (ET)	Empty

Load the labware in the same order as cartridges. Make sure to load the samples into the heating block and NOT into the tip rack!

Step 7 of 7

A B C D

Cancel Back Next

8/31/2020 10:47 Software version: 0.2.9-93-g1d33f49c Admin Admin

**Figure 61.** Load the tip rack step.

**Note:** In some applications, not all positions of the tip rack are to be used. This is indicated by the lack of this position on the left-hand side and a white circle on the right-hand side.

To proceed to the start of the protocol, tap **Next**. Before the run starts, you will see an overview of the selections you made during the run setup process.

To return to the previous screen, tap **Back**.

To cancel the protocol, tap **Cancel**.

## 5.5 Canceling the protocol run setup

You can cancel the protocol run setup process at any time. If you cancel the run setup, your progress is not saved, and the worktable does not move. If you have loaded anything onto the worktable, remove the labware.

To cancel the setup, tap **Cancel**. In the **Cancel run setup** dialog, tap **Yes** to confirm the cancellation, or tap **No** to go back to the run setup.


Figure 62. The Cancel run setup dialog.

## 5.6 Starting the protocol run and monitoring its progress

After successfully completing all the steps of the protocol run setup, you can start the run, and, optionally, add the protocol to your favorites. During the protocol run, you can monitor its progress. Ongoing steps are displayed on the screen.

To start the run and view its progress, follow the steps below:

1. Tap **Next** in the **Load the tip rack** step. The **Run setup selection overview** dialog is shown.


Figure 63. Run setup selection overview dialog.

2. If all the information in the overview is correct, tap **Create favorite** to save the protocol to your favorites, or tap **Start** to immediately proceed with the protocol run. To make modifications to any of the selections, tap **Return** to go back to the run setup.

**Note:** The **Estimated run time** does not include the time to complete the **Load Check**.

3. If the hood of the instrument is open, close it and tap **Proceed**.



Figure 64. Close the hood dialog.

4. **For E22 Connect Fx:** If your instrument has a built-in camera, a load check is performed. For more information about the load check, refer to the [Load check](#) section. The run will start after the check is completed. To stop the load check, tap **Abort**.

**Note:** Wait until the **Load Check** has successfully completed before leaving the instrument unattended.

5. Progress of the run and the elapsed run time are displayed on the **Protocol run in progress** screen.

**Figure 65.** Run progress screen.

### 5.6.1 Load check

**Note:** The load check feature is only available in EZ2 Connect Fx instruments.

The EZ2 Connect Fx comes with a built-in camera, which checks if you have loaded all of the cartridges and labware into the correct positions on the worktable. However, the content of the labware in the tip rack is not checked, so you must make sure that you have carefully followed the instructions related to the protocol you are running.

The load check starts after you tap **Start** in the **Run setup selection overview** dialog window. For more information on starting a run, refer to the Starting the protocol run and monitoring its progress section.

After the load check starts, the camera moves above the worktable and checks all positions on the cartridge rack and the tip rack, and the **Load check is in progress** screen is shown on the display. The positions that you chose in the **Select sample positions** screen are highlighted.


Figure 66. Load check is in progress screen.

### Successful load check

If all of the labware is loaded correctly, the load check ends successfully, and the run starts automatically.

### Failed load check

If the camera discovers one or more errors during the load check procedure, the **Load check failed** screen is shown. Incorrect labware placements are marked in red. To get more information about a particular load check error, tap one of the red positions. A dialog window with details about the error appears. To return to the loading instructions and start the load check procedure again, tap **Back**. The **Load the tip rack** screen is shown. Once you correctly load the worktable, tap **Next** in the **Load the tip rack** screen. The load check procedure starts again.

**Note:** In case of repeated load check failure, please recalibrate the camera (refer to [section 6.6](#)).

Load check failed


Figure 67. Load check failed screen.

## 5.7 End of the protocol run

When the protocol has successfully finished, the **Protocol run completed** screen is shown. Information about required cleanup/maintenance steps is also displayed. Follow the instructions to properly remove all labware from the instrument and clean the piercing unit, see the After run maintenance section. If after run maintenance is completed, tap the check box to transfer the maintenance status into the run report. For more information on how to save a run report, refer to the Saving a run report section.


Figure 68. Protocol run completed screen.

## 5.8 Aborting a protocol run

A protocol run can be stopped at any time. To stop the run, follow the steps below:

1. In the **Protocol run in progress** screen, tap **Abort run**. A confirmation dialog appears.
2. Tap **Yes** to stop the run, or tap **No** to go back to the **Protocol run in progress** screen.


Figure 69. Aborting protocol run dialog.

- When the run is aborted, the instrument finishes the current movement, and then attempts to dispense the contents of the pipettes into an empty tube, and to release tips into empty tip holders. Then, the worktable returns to its initial position. After this is done, a message is shown and the **Proceed to the summary** button is activated. Tap **Proceed to the summary**.

**Note:** If you tap **Abort** during a pause or when the machine is waiting to reach a specific temperature, the run is stopped immediately.


Figure 70. Protocol run interrupted dialog.

- Tap Finish to end the run and go back to the **Home** screen. A run report is generated. For more information on how to save a run report, refer to the Saving a run report section.



Figure 71. Protocol run aborted screen.

- Also, for aborted protocol runs, the after run maintenance needs to be performed. Please follow the instructions regarding the cleaning of the piercing unit from Daily maintenance section.

## 5.9 Sample recovery

In case of a run failure the EZ2 Connect Fx provides a recovery feature which allows localization and recovery of the sample.

Protocol run failed

Run failed at: 10:23  
Run duration: 7 sec

Finish the protocol run and perform the daily maintenance procedure.  
**State of the protocol after the recovery:**  
The sample is available at the following positions: Rack: A.  
The tip contents have been dispensed into the labware at the following position: Rack: D.  
Proceed with the following step from the Protocol Manual: 1.

Run steps:

✓ Preparation
! Bead preparation - Current step: Resuspending beads
Lyse
Binding
Wash
Elution

Finish

9/1/2021 10:24 Admin Admin

Figure 72. Protocol run failed screen.

## 5.10 Saving a run report

After a run ends successfully, fails, or is aborted, run report is generated in two formats: PDF and XML.

To save a run report, tap **Finish** in the **Protocol run completed**, **Protocol run failed**, or **Protocol run aborted** screens.

For more information on what the run report contains, refer to the Run report contents section.

### 5.10.1 Run report contents

An EZ2 run report is created by the software application after a run is completed, aborted or if it fails, after the user taps the **Finish** button on the screen that is shown after a run ends.

Each run report is saved in two formats: PDF and XML. Both formats include the same information, that is:

- The name of the user that was logged in when the run was started
- The serial number of the instrument and software version
- The duration of the run
- The time and date when the run was started and when it ended
- Protocol information:
  - Name
  - Version
  - Application
  - Parameters selected
- The kit name, material number, lot number, and expiration date
- The name of the run report file, which includes the date when the run ended and the serial number of the instrument
- The status of the run, which indicates if the run was completed successfully, if it failed, or if it was aborted
- The status of the cleanup procedure that is required after a run ends
- Information about the samples - their positions, names and any notes that were added by the user
- Validity of samples

### 5.11 Setting up the worktable

Setting up the EZ2 worktable consists of removing the tip rack and the cartridge rack from the instrument, loading the cartridges, tubes, tip holders and tips, and placing the racks back into the instrument. Some protocols require additional actions to be done before the start of a run, e.g. inversion of the cartridge to mix the magnetic beads. These actions are described in the Loading the cartridge rack and Loading the tip rack screens of the protocol run setup process, and in the kit handbooks.

### 5.11.1 Loading and unloading the cartridge rack

To load the cartridge rack, follow the steps below:

1. Follow the instructions of the respective kit handbook on how to treat the reagent cartridges before loading onto the cartridge rack.

**Note:** Remember to prepare the same number of reagent cartridges as the number of positions you chose in the **Select sample positions** step of the protocol run setup.

2. Tap the reagent cartridges until the reagents are deposited at the bottom of the wells.
3. Remove one or both sections (left or right) of the cartridge rack from the worktable, depending on which positions you chose in the **Select sample positions** step of the protocol run setup. To remove the cartridge rack sections, grab their handles and gently pull the rack up.
4. Slide the reagent cartridges, with the orientation depicted in the User Interface, into the cartridge rack in the direction of the arrow that is engraved on each cartridge rack section, until you feel resistance.
5. Once all the reagent cartridges are loaded, place each cartridge rack section onto the worktable.

At the end of a protocol run, you must remove the reagent cartridges from the cartridge rack. To do so, follow the steps below:

1. Remove one or both sections (left or right) of the cartridge rack from the worktable. To remove the cartridge rack sections, grab their handles and gently pull the rack up.
2. Slide out the reagent cartridges and dispose of them properly, according to local safety regulations.

<p><b>CAUTION</b></p> 	<p><b>Hazardous materials and infectious agents</b></p> <p>The waste contains samples and reagents. This waste may contain toxic or infectious material and must be disposed of properly. Refer to your local safety regulations for proper disposal procedures.</p>
---	--

3. Disinfect the cartridge rack and remove any spilled liquid or potential contamination that is present on the cartridge rack sections. For more information on disinfecting and removing contamination, refer to the Cleaning agents, Disinfecting the EZ2, and Removing contamination sections.
4. Place the cartridge rack sections back in the instrument.

### 5.11.2 Loading and unloading the tip rack

To load the tip rack, follow the steps below:

1. Remove one or both sections of the tip rack from the worktable, depending on which positions you chose in the **Select sample positions** step of the protocol run setup. To remove a tip rack section, grab both sides of the section and gently pull up.
2. Place the tips into their tip holders.
3. Load the tip holders with tips/labware into Row C.

**Note:** Information on whether the selected protocol requires one or two rows of tips is shown on the screen in the **Load the tip rack** step of the run setup process and provided in the kit handbook.

4. Load the labware into Row A and D (or C, if you are using the flip cap rack).

**Note:** Make sure that you follow any protocol-specific directions shown on the screen in the **Load the tip rack** step of the run setup process. You could be required to perform some additional actions. The instructions can also be found in the kit handbooks.

**Note:** Remove any caps from the labware. When using the flip cap rack, place the flip caps into the slots located below Row A or above Row C.

5. Once all the labware is loaded, place the tip rack sections onto the worktable.

At the end of a protocol run, you must remove the labware from the tip rack. To do so, follow the steps below:

**Note:** Remove the elution tube containing the eluant before removing the used labware from the tip rack.

1. Remove one or both sections of the tip rack from the worktable. To remove a tip rack section, grab both sides of the section and gently pull up.
2. Remove the labware from the tip rack and dispose of any unnecessary labware, according to local safety regulations.

**CAUTION****Hazardous materials and infectious agents**

The waste contains samples and reagents. This waste may contain toxic or infectious material and must be disposed of properly. Refer to your local safety regulations for proper disposal procedures.

3. Disinfect the tip rack and remove any spilled liquids or potential contamination that is present on the cartridge rack sections. For more information on disinfecting and removing contamination, refer to the Cleaning agents, Disinfecting the EZ2, and Removing contamination sections.
4. Place the tip rack sections back in the instrument.

## 5.12 Using the bar code scanner

**Important:** Basic version does not contain the 2D handheld barcode scanner, but it can be purchased separately. Usage of 2D handheld barcode scanner is optional

The 2D handheld barcode scanner can be used to scan the following barcodes during Protocol setup

- Kit Q-Cards
- Sample Bar Codes

Direct the handheld scanner at the required Q-Card/Barcode to scan the information into the Instrument when prompted by the User Interface


Figure 73. Scanning the information into the Instrument.

<p><b>WARNING</b></p> 	<p><b>Risk of personal injury</b></p> <p>Hazard Level 2 laser light: Do not stare into the light beam when using handheld bar code scanner.</p>
---	---

## 5.13 Logging out and switching the instrument off

To log out of the Software, tap the **Log out** button on the toolbar.

To switch the instrument off, press the power button.

**Note:** It is not recommended to switch the instrument off when a protocol run, or maintenance procedure are in progress

Figure 74. Location of the Log out button.

## 6 Maintenance Procedures

<b>WARNING/ CAUTION</b> 	<b>Risk of personal injury and material damage</b> Only perform maintenance that is specifically described in this user manual.
--	--

To ensure reliable operation of the EZ2, maintenance procedures must be performed. The procedures are shown in the table below. Each maintenance procedure must be carried out by appropriate personnel.

**Important:** The maintenance procedures as shown in the table below must be carried out to ensure reliable operation of the EZ2. Each maintenance procedure must be carried out by appropriate personnel.

Type of task(s)	Frequency	Personnel
After run maintenance	After each run.	Laboratory technicians or equivalent
Daily maintenance	At the end of each day, if at least one run was performed that day <b>Note:</b> Perform this procedure after you are done with the After run maintenance.	Laboratory technicians or equivalent
Weekly maintenance	Once a week. <b>Note:</b> Perform this procedure after you are done with the regular and daily procedures.	Laboratory technicians or equivalent
Annual maintenance and servicing	Annually or semiannually, depending on your requirements (for more information contact QIAGEN Technical Service.	QIAGEN Instrument Service Specialists only

Optionally, a UV decontamination procedure may be performed as required to reduce pathogen and nucleic acid decontamination.

### 6.1 Cleaning agents

The surfaces and removable parts of the EZ2 must be cleaned and disinfected with compatible detergents and disinfectants. Please follow the instructions provided by the manufacturer of such materials to safely clean the instrument.

**Note:** If you want to use different disinfectants than the recommended ones, make sure that they have the same composition.

If you are unsure about the suitability of the detergents or disinfectants to be used with the EZ2, do not use them.

General cleaning of the EZ2, with the exception of the hood, can be done using mild detergents, such as the Mikrozyd® AF sensitive, or 70% ethanol. 70% ethanol should only be used on the worktable. The hood should **ONLY** be cleaned with lint free tissue moistened with water.

<b>WARNING</b> 	<b>Toxic fumes</b> Do not use bleach to clean or disinfect the EZ2 instrument or used labware. Bleach in contact with salts from the buffers can produce toxic fumes.
---	--

<b>CAUTION</b> 	<b>Damage to the instrument</b> Do not use spray bottles containing alcohol or disinfectant to clean surfaces of the EZ2. Spray bottles should be used only to clean items that have been removed from the worktable and if permitted by local laboratory operating practices.
---	---

<b>WARNING</b> 	<b>Risk of fire</b> Do not allow cleaning fluid or decontamination agents to come into contact with the electrical parts of the EZ2.
---	---

<b>WARNING</b> 	<b>Risk of electric shock</b> Do not open any panels on the EZ2.  <b>Risk of personal injury and material damage</b> Only perform maintenance that is specifically described in this user manual.
---	---

<p><b>WARNING</b></p> 	<p><b>Risk of fire or explosion</b></p> <p>When using ethanol or ethanol-based liquids to clean the EZ2, handle such liquids carefully and in accordance with the required safety regulations. If any liquid spills, wipe it off and leave the EZ2 hood open to allow flammable vapors to disperse.</p>
---	---

### 6.1.1 Disinfecting the EZ2

Ethanol-based disinfectants can be used for disinfecting surfaces such as the worktable. An example of an ethanol-based disinfectant are the Mikrozyd® Liquid (Mikrozyd® Liquid consists of 25 g ethanol and 35 g 1-propanol per 100 g) or Mikrozyd AF wipes. These are available from Schülke & Mayr GmbH, see [www.schuelke.com](http://www.schuelke.com).

Disinfectants based on quaternary ammonium salt can be used for submerging worktable items. An example of such disinfectants are the Lysetol® AF or Gigasept® Instru AF (Schülke & Mayr GmbH, see [www.schuelke.com](http://www.schuelke.com)). These disinfectants consist of 14 g cocosporylene-diamine-guanidine diacetate, 35 g phenoxypropanols, and 2.5 g benzalkonium chloride per 100 g, with anticorrosion components, fragrance, and 15–30% nonionic surfactants.

**Note:** If you want to use different disinfectants than the recommended ones, make sure that they have the same composition.

**Note:** The hood should **ONLY** be cleaned with lint free tissue moistened with water.

<p><b>CAUTION</b></p> 	<p><b>Damage to the instrument</b></p> <p>Do not use spray bottles containing alcohol or disinfectant to clean surfaces of the EZ2. Spray bottles should be used only to clean items that have been removed from the worktable and if permitted by local laboratory operating practices.</p>
---	--

### 6.1.2 Removing contamination

The EZ2 could be contaminated during operation. To remove the contamination, use appropriate decontamination solutions.

In case of RNase contamination, the RNaseZap® RNase Decontamination Solution (Ambion, Inc., cat. no AM9780) can be used for cleaning surfaces and submerging worktable items. RNaseZap

can also be used to perform decontamination by spraying worktable items, if they have been removed from the instrument.

In case of nucleic acid contamination, DNA-ExitusPlus™ (AppliChem, cat. no. A7089,0100) can be used for cleaning surfaces and submerging worktable items. DNA-ExitusPlus can also be used to perform decontamination by spraying worktable items, if they have been removed from the instrument. Cleaning with DNA-ExitusPlus can leave a residue on surfaces so for this reason, after cleaning the items with DNA-ExitusPlus, it is required to clean the items with a wet cloth several times, or rinse them with running water, until the DNA-ExitusPlus is completely removed.

## 6.2 After run maintenance

After run maintenance is required after each run on the EZ2.

The EZ2 must only be operated by qualified personnel who have been appropriately trained.

Servicing of the EZ2 must only be performed by QIAGEN Field Service specialists.

<b>WARNING</b> 	<b>Moving parts</b> To avoid contact with moving parts during the operation of the EZ2, the instrument must be operated with the hood closed.  If the hood sensor or lock is not functioning correctly, contact QIAGEN Technical Services
---	--

<b>WARNING</b> 	<b>Moving parts</b> Avoid contact with moving parts during operation of the EZ2. Under no circumstances should hands be placed under the pipetting arm during movement. Do not attempt to remove any plasticware from the workdeck whilst the instrument is operating.
---	---

<b>WARNING/ CAUTION</b> 	<b>Risk of personal injury and material damage</b> Improper use of the EZ2 may cause personal injuries or damage to the instrument.
--	--

<b>WARNING</b> 	<b>Samples combining infectious agents</b> Some samples used with this instrument may contain infectious agents. Handle such samples with the greatest of care and in accordance with the required safety regulations.  Some chemicals used with the EZ2 instrument may be hazardous or may become hazardous after completion of a purification.  Always wear safety glasses, gloves and a lab coat.  Venting for fumes and disposal of waste must be in accordance with all national, state, and local health and safety regulations and laws.
---	--

If you are working with potentially infectious materials, such as human blood, serum or plasma, the EZ2 should be decontaminated after use (for more details, see the Disinfecting the EZ2, and Removing contamination sections).

After running a protocol, clean the piercing unit of the pipettor head.

1. Remove all sample preparation waste and discard it according to your local safety regulations.
2. Remove the eluates and store them appropriately
3. Close the hood


Figure 75. Protocol run completed screen.

4. To begin the procedure, tap **Move down**. The instrument lowers the piercing unit of the pipettor head.
5. Open the hood.
6. Wipe the piercing unit using a soft tissue moistened with 70% ethanol. The piercing unit is sharp. Two pairs of gloves are recommended.

Figure 76. Cleaning the EZ2 piercing unit.

7. Wipe the piercing unit with a soft tissue moistened with distilled water.
8. Close the hood.
9. To document the cleaning procedure in the run report, confirm that the maintenance has been completed by activating the checkbox.
10. On the touchscreen, tap **Finish**. The piercing unit returns to its home position.
11. Open the hood.
12. If contamination is visible on the worktable, clean it with 70% ethanol, and then with distilled water.
13. Wipe the surface of the instrument using a soft tissue moistened with 70% ethanol.

**Note:** The hood should **ONLY** be cleaned with lint free tissue moistened with water.

### 6.3 Daily maintenance

Daily maintenance is required after the last run of each day.

<p><b>WARNING</b></p> 	<p><b>Moving parts</b></p> <p>To avoid contact with moving parts during the operation of the EZ2, the instrument must be operated with the hood closed.</p> <p>If the hood sensor or lock is not functioning correctly, contact QIAGEN Technical Services</p>
---	---

<p><b>WARNING</b></p> 	<p><b>Moving parts</b></p> <p>Avoid contact with moving parts during operation of the EZ2. Under no circumstances should hands be placed under the pipetting arm during movement. Do not attempt to remove any plasticware from the workdeck whilst the instrument is operating.</p>
---	--

<b>WARNING/ CAUTION</b> 	<b>Risk of personal injury and material damage</b> Improper use of the EZ2 may cause personal injuries or damage to the instrument.
--	--

<b>WARNING</b> 	<b>Samples combining infectious agents</b> Some samples used with this instrument may contain infectious agents. Handle such samples with the greatest of care and in accordance with the required safety regulations.  Some chemicals used with the EZ2 instrument may be hazardous or may become hazardous after completion of a purification.  Always wear safety glasses, gloves and a lab coat.  Venting for fumes and disposal of waste must be in accordance with all national, state, and local health and safety regulations and laws.
---	--

**Important:** Before starting the daily maintenance procedure, complete the After run maintenance.

1. On the touchscreen, tap **Maintenance**.
2. Tap **Daily**. The date of the last finished procedure is shown on the screen.


Figure 77. Daily maintenance screen.

3. Clean the piercing unit (if this is not yet done during after run maintenance).
4. To prepare the piercing unit for cleaning tap **Move down**.
5. Open the hood.
6. Clean the piercing unit as recommended for After run maintenance.
7. Close the hood.
8. To return the piercing unit to home position, tap **Move up**.
9. Clean the worktable with 70% ethanol, and then with distilled water (if this is not yet done during after run maintenance).
10. To allow tray removal, tap **Move to back**.
11. Check if the tray is clean. If necessary, clean it with 70% ethanol, and then with distilled water.
12. Wipe the instrument racks with 70% ethanol, and then with distilled water.
13. Lightly wipe the O-rings of the tip adapters with a lint-free tissue to remove any residual liquid.

## 6.4 Weekly maintenance


<b>WARNING/ CAUTION</b> 	<b>Risk of personal injury and material damage</b> Improper use of the EZ2 may cause personal injuries or damage to the instrument.
--	--

<b>WARNING</b> 	<b>Samples combining infectious agents</b> Some samples used with this instrument may contain infectious agents. Handle such samples with the greatest of care and in accordance with the required safety regulations.  Some chemicals used with the EZ2 instrument may be hazardous or may become hazardous after completion of a purification.  Always wear safety glasses, gloves and a lab coat.  Venting for fumes and disposal of waste must be in accordance with all national, state, and local health and safety regulations and laws.
---	--

**Important:** Before starting the weekly maintenance procedure, complete the Daily maintenance.

To maintain good contact between tip adapters and filter-tips, and to prevent liquid leaking from the tips, a light application of grease should be applied to the O-rings of the tip adapters every week.

1. On the touchscreen, tap **Maintenance**.
2. Tap **Weekly maintenance**. The date of the last procedure, and the date of the next due date of the weekly maintenance are shown on the screen.

**Figure 78. Weekly maintenance screen.**

3. Put on gloves.
4. Apply a small amount of silicon grease into the end of a filter-tip.
5. Apply a small amount of silicon grease onto the surface of the O-rings, by just using a fingertip moistened with grease.
6. Place the filter-tip onto the pipettor head and rotate the filter-tip on the pipettor head to distribute the silicon grease evenly.
7. Make sure that the O-rings are only moistened with grease and that there are no clumps of grease visible.
8. Make sure that there is no grease present on parts other than the O-rings, especially on the bar on top and the opening of the pipettors.

**Note:** The filter-tips should sit flush against the upper white plastic bar if the O-rings are properly greased. There should not be a gap. Excess or insufficient grease can affect the performance of the EZ2.

**Note:** The opening in the nub of the pipettor should be checked after greasing to ensure no grease is present within the opening.


Figure 79. Placing a filter-tip onto the pipettor head.

## 6.5 UV decontamination

<b>WARNING/ CAUTION</b> 	<b>Risk of personal injury and material damage</b> Improper use of the EZ2 may cause personal injuries or damage to the instrument.
--	--

<b>WARNING</b> 	<b>Samples combining infectious agents</b> Some samples used with this instrument may contain infectious agents. Handle such samples with the greatest of care and in accordance with the required safety regulations.  Some chemicals used with the EZ2 instrument may be hazardous or may become hazardous after completion of a purification.  Always wear safety glasses, gloves and a lab coat.  Venting for fumes and disposal of waste must be in accordance with all national, state, and local health and safety regulations and laws.
---	--

<b>WARNING</b> 	<b>UV radiation</b> Avoid looking directly into UV light. Do not expose your skin to UV light.
---	---

**Important:** Before starting the UV decontamination, complete the After run maintenance.

1. On the touchscreen, tap **Maintenance**.
2. Tap **UV Run**. The date of the last procedure is shown on the screen.


Figure 80. Manual UV Run screen.

3. Select the number of decontamination cycles. The required decontamination time depends on the biological material processed on the device.
4. To begin the procedure, tap **Start**.

## 6.6 Servicing

Contact your local QIAGEN Technical Service or your local distributor for more information about flexible Service Support Agreements from QIAGEN.

<p><b>WARNING/ CAUTION</b></p> 	<p><b>Risk of personal injury and material damage</b></p> <p>Improper use of the EZ2 may cause personal injuries or damage to the instrument. The EZ2 must only be operated by qualified personnel who have been appropriately trained.</p> <p>Servicing of the EZ2 must only be performed by a QIAGEN Field Service specialist.</p>
--	--

## 7 Troubleshooting

This section contains information about what to do if an error occurs when using the EZ2 instrument.

### 7.1 Contacting QIAGEN Technical Services

Whenever encountering an EZ2 instrument error, be sure to have the following information at hand:

**Note:** Much of the information listed below can be found in the run report.

- Protocol name and version
- Software version
- Serial number of the instrument, this can be found on the type plate at the rear of the instrument or in each run report
- Sample input material and sample pre-treatment
- Detailed description of the error situation
- Support package

This information will help you and your QIAGEN Technical Service Specialist to deal most efficiently with your issue.

**Note:** Information about the latest software and protocol versions can be found at [www.qiagen.com](http://www.qiagen.com). In some cases, updates may be available for addressing specific problems.

#### 7.1.1 Creating a support package

The EZ2 can create a support package containing information about the device and the device status. This can help the product support to find issues.

1. Plug in a USB flash drive.
2. Tap **Data** on the toolbar.


Figure 81. Data button on the toolbar.

3. Tap **Create and download** in the Support Package section.

The screenshot shows the QIAGEN Data tab interface. At the top, there is a navigation bar with icons for Setup, Maintenance, Configuration, Data (highlighted), Network, and Logout. The main content area is divided into two sections. On the left, under the heading "Run reports available: 2", there are radio buttons for "All run reports" (selected) and "Recent reports". Below this is a "Number of reports:" field with the value "1" and a checked checkbox for "Delete reports after download". A "Download" button is at the bottom of this section. On the right, there are two panels: "Support Package" with a "Create and download" button, and "Audit trail" with a "Download" button. The footer shows the date and time "8/23/2021 16:10" and the user "Admin Admin".

**Figure 82. Data tab.**

4. The support package is now saved on the USB flash drive.

## 7.2 Operation

	Comments and suggestions
<b>Load check failed</b>	<p>Perform camera exposure calibration according to Section 4.3.6.</p> <p>Both racks must be loaded, even if just one rack is used. Racks have to be the same type.</p> <p>Check if there is no contamination (e.g. drops) on the labware.</p> <p><b>Note:</b> Avoid strong environmental light changes.</p> <p>Perform load check again.</p>
<b>Hood lock Error</b>	Check if the hood is firmly closed. Confirm by pressing down the hood with light force.
<b>Heater Fuse Error</b>	Restart the device (several times if necessary).
<b>Pipetting tips sticking to pipettor</b>	Use lab tissue to wipe the O-rings and make sure weekly maintenance was performed correctly.
<b>Pipette tip was crushed</b>	Tip may have been stuck to pipettor. See "pipetting tip sticking to pipettor".
<b>Imprecise pipetting</b>	If imprecise pipetting appears over several runs (elution volume): Make sure weekly maintenance was performed. Check if the tips sits firmly at the pipettor adapter.
<b>Leaking pipette tips</b>	Ensure weekly maintenance was performed. Check if the tips sit firmly at the pipettor adapter.
<b>Pipette tips not picked up by automatic pipettor</b>	Make sure that the tip rack is not damaged and is correctly positioned on the worktable.
<b>Cross-contamination</b>	Ensure maintenance was performed. Clean the piercer and worktable with 70 % ethanol. Start UV decontamination. Make sure samples and cartridge rack were handled properly.
<b>Update fails</b>	Restart the device and try to start the update again. Make sure you use the QIAGEN USB flash drive. USB flash drive should stay plugged in during the whole update procedure.
<b>USB device not detected</b>	Restart the device. Save the file(s) to the USB stick again. Check the USB stick on a PC to ensure it is functional. If the error persists, contact QIAGEN Technical Services.
<b>Connection problem</b>	Check if the ethernet cable is connected correctly. Check the Wi-Fi and LAN settings (Section 5.3.8).
<b>No W-LAN connection</b>	Check the Wi-Fi setting (Section 5.3.8). W-LAN adapter should be connected before you start the device. Reboot the device.
<b>Freeze during run</b>	Reboot the device. Restart the protocol run.
<b>Surface damaged</b>	Ensure that only the cleaning agents as assumed in Section 6.1 were used.
<b>Display does not turn on</b>	Do not touch the display with excessive force or use corrosive chemicals to clean the display surface. Contact QIAGEN Technical Services for repair.

## 8 Glossary

Term	Definition
<b>Bar code reader</b>	A handheld device that enables scanning of bar codes and conversion of them into data that is transmitted to the EZ2.
<b>EZ2 Connect Cartridge Rack</b>	A metal rack that accommodates reagent cartridges on the worktable.
<b>Connector panel</b>	The panel on the rear of the EZ2. It contains the power switch, the socket for the power cord, and the fuse box.
<b>Elution tube</b>	A polypropylene, screw-capped 1.5 ml tube for collecting purified nucleic acids. The recommended elution tubes are screw-capped, made of polypropylene, supplied by Sarstedt (cat. no. 72.692) and provided in EZ2 kits.
<b>Error code</b>	A number that represents an error of the EZ2.
<b>EZ2 kits</b>	Kits supplied by QIAGEN containing reagents, reagent cartridges and plasticware for use with EZ2 instruments.
<b>Filter-tip</b>	Labware that is picked up by a tip adapter during operation of the EZ2. Liquid is aspirated into and dispensed from a filter-tip. A filter-tip is also the location where separation of magnetic particles occurs.
<b>Heating system</b>	A component of the EZ2 that accommodates the heating positions of the reagent cartridges and heats samples.
<b>Hood</b>	The main door at the front of the EZ2. When open, it provides complete access to the worktable.
<b>O-ring</b>	A ring that is fitted at the bottom of a tip adapter. It is required for good contact between the tip adapter and a filter-tip.
<b>Pipettor head</b>	The component of the EZ2 that aspirates and dispenses liquid, and pierces cartridges via the piercing unit. The pipettor head moves up and down above the worktable and contains 24 syringe pumps, each of which is connected to a tip adapter.
<b>Protocol</b>	A set of instructions for the EZ2 that allows the instrument to automate a nucleic acid purification procedure.
<b>Reagent cartridge</b>	An item of labware that contains 10 wells and 2 heating positions. One heating position is a well, the other is a slot that can hold a tube. A reagent cartridge is prefilled with reagents and included in EZ2 kits.
<b>Report file</b>	A file generated by the EZ2 that contains system and run parameters.
<b>Sample tube</b>	A polypropylene, screw-capped 2 ml tube for holding a sample containing nucleic acids to be purified. Sample tubes are 2 ml in volume, screw-capped, made of polypropylene, supplied by Sarstedt (cat. no. 72.693), and provided in EZ2 kits.
<b>Tip adapter</b>	One of 24 metal probes installed on the pipettor head. During operation of the EZ2, the tip adapters pick up filter-tips from the worktable.
<b>Tip holder</b>	A polypropylene tube that holds a single filter-tip. Tip holders are loaded onto the tip rack.
<b>EZ2 Connect Tip Rack</b>	A metal rack that accommodates tip holders containing filter-tips on the worktable. The tip rack also accommodates sample tubes and elution tubes.
<b>Tray</b>	A metal tray that is located under the worktable. It collects any drops of liquid that may spill.
<b>UV LED lamp</b>	A light source of ultraviolet light for decontamination.
<b>Worktable</b>	The surface of the EZ2 that contains racks. The worktable is where samples, reagent cartridges and disposable labware are loaded. The worktable moves backward and forward to position samples and reagents under the pipettor head.

## 9 Technical Specifications

QIAGEN reserves the right to change specifications at any time.

### 9.1 Operating conditions

Power	100–240 V AC, 50/60 Hz, 1000 VA Mains supply voltage fluctuations are not to exceed 10% of nominal supply voltages.
Fuse	AC Inlet: T4A H 250 V Heater block (temperature fuse): 10A 250V 117C
Overvoltage category	II
Air temperature	18–30°C (64–86°F)
Relative humidity	10–75% RH
Altitude	Up to 2000 m (6500 ft.)
Place of operation	For indoor use only
Pollution level	2
Environmental class	3K21 (IEC 60721-3-3) 3M11 (IEC 60721-3-3)
Average noise level (over 8 hours)	Max. 70 dBA
IP Code (IEC 60529)	IP20

### 9.2 Transport conditions

Air temperature	–25–60°C (–13–140°F) in manufacturer's package <b>Note:</b> If the EZ2 is transported in temperatures below 0°C (32°F), it is recommended to wait 24 hours before switching on the instrument to allow it to reach the temperature conditions of the installation environment.
Relative humidity	5–85% RH
Environmental class	2K11 (IEC 60721-3-2) 2M4 (IEC 60721-3-2)

### 9.3 Storage conditions

Air temperature	5–40°C (41–104°F) in manufacturer's package
Relative humidity	5–85% RH

Environmental class	2K11 (IEC 60721-3-2) 2M4 (IEC 60721-3-2)
---------------------	---

## 9.4 Mechanical data and hardware features

Dimensions	Width: 720 mm (28.35 in.) Height: 575 mm (22.6 in.) Depth: 560 mm (22 in.)
Weight	70 kg (154.32 lb.)
Instrument features	Automated nucleic acid isolation using magnetic particles Desktop instrument Protocols stored on instrument Processes up to 24 samples in one run Aspirates and dispenses 24 samples or reagents simultaneously using a 24-channel pipettor head Separates magnetic particles using patented technology Controlled through a touchscreen Temperature control through a heating system
Data tracking	Bar code reader and on-screen keyboard enable data tracking of samples and consumables. System and run parameters are stored in a report file.
Pipettor head	Contains 24 high-precision syringe pumps, each containing a tip adapter that attaches to filter-tips. Syringe pumps are air-filled. Liquids containing salts, alcohol, solvents and/or magnetic particles can be aspirated and dispensed. Air gaps can be aspirated to prevent aspirated liquid from dripping. Filter-tips are picked up from the tip rack and ejected back into the tip rack. The pipettor head moves in the Z direction (up and down) above the worktable.
Heating system	Accommodates the heating positions of reagent cartridges and has a temperature range of between ambient temperature and 95°C (203°F). Heating block accuracy at 60°C is $\pm 2^\circ\text{C}$ .
Filter tips	Attach to the tip adapters of the pipettor head to allow liquid aspiration and dispensation. Capacity of 20–1000 $\mu\text{l}$ . The EZ2 accommodates up to 48 tip holders in 2 rows, each containing a filter-tip, in the tip rack on the worktable.

Labware	<p>Reagents are loaded onto the worktable in reagent cartridges. These cartridges are already prefilled by QIAGEN.</p> <p>Up to 24 reagent cartridges can be placed in the cartridge rack.</p> <p>Samples are loaded onto the worktable in 2 ml sample tubes.</p> <p>Steps that require heating occur on the heating system, which accommodates the heating positions of the reagent cartridges.</p> <p>Purified nucleic acids are collected in 1.5 ml elution tubes.</p>
UV LED lamp	UV LED wavelength: 253.7 nm
Capacity	Up to 24 samples per run
Display	10.1-inch color touchscreen. Display with resolution 1280 x 800 pixels.
Camera	<p>Monochrome camera. The USB interface provide power and communication.</p> <p>Sensor resolution is 0.34 MP.</p> <p>Width: 24 mm (0.94 in.)</p> <p>Height: 34 mm (1.34 in.)</p> <p>Depth: 39 mm (1.54 in.)</p>

---

# Appendix A

## Declaration of Conformity

### **Name and address of the legal manufacturer**

QIAGEN GmbH  
QIAGEN Strasse 1  
40724 Hilden  
Germany

An up-to-date declaration of conformity can be requested from QIAGEN Technical Services.

---

## Waste Electrical and Electronic Equipment (WEEE)

This section provides information about disposal of waste electrical and electronic equipment by users.

The crossed-out wheeled bin symbol (see below) indicates that this product must not be disposed of with other waste; it must be taken to an approved treatment facility or to a designated collection point for recycling, according to local laws and regulations.

The separate collection and recycling of waste electronic equipment at the time of disposal helps to conserve natural resources and ensures that the product is recycled in a manner that protects human health and the environment.

Recycling can be provided by QIAGEN upon request at additional cost. In the European Union, in accordance with the specific WEEE recycling requirements and where a replacement product is being supplied by QIAGEN, free recycling of its WEEE-marked electronic equipment is provided.

To recycle electronic equipment, contact your local QIAGEN sales office for the required return form. Once the form is submitted, you will be contacted by QIAGEN either to request follow-up information for scheduling collection of the electronic waste or to provide you with an individual quote.

---

## California Proposition 65 Warning

This product contains chemicals known to the State of California to cause cancer, birth defect or other reproductive harm.

## Liability Clause

QIAGEN shall be released from all obligations under its warranty in the event repairs or modifications are made by persons other than its own personnel, except in cases where the Company has given its written consent to perform such repairs or modifications.

All materials replaced under this warranty will be warranted only for the duration of the original warranty period, and in no case beyond the original expiration date of original warranty unless authorized in writing by an officer of the Company. Read-out devices, interfacing devices, and associated software will be warranted only for the period offered by the original manufacturer of these products. Representations and warranties made by any person, including representatives of QIAGEN, which are inconsistent or in conflict with the conditions in this warranty shall not be binding upon the Company unless produced in writing and approved by an officer of QIAGEN.

The EZ2 is equipped with an Ethernet port and a Wi-Fi USB device. The Purchaser of the EZ2 Connect is solely responsible for preventing any and all computer viruses, worms, trojans, malware, hacks, or any other type of cybersecurity breaches. QIAGEN assumes no liability for computer viruses, worms, trojans, malware, hacks, or any other type of cybersecurity breaches.

## Appendix B – EZ2 Accessories

### Ordering Information

<b>Product</b>	<b>Contents</b>	<b>Cat. no.</b>
EZ2 Connect	EZ2 Connect benchtop instrument to isolate nucleic acids from up to 24 samples in parallel, includes 1-year warranty on parts and labor	9003210
EZ2 Connect Fx	EZ2 Connect Fx benchtop instrument for isolation of DNA from casework and reference samples in forensic workflows, with internal camera for cartridge check and load checks, bar code reader, EZ2 Connect Tip Rack – Flip Cap Tube, and 1-year warranty on parts and labor.	9003220
<b>Accessories</b>		
EZ2 Connect Tip Rack	Tip Rack for EZ2 Connect, for use with screw-cap tubes	9027009
EZ2 Connect Tip Rack - Flip Cap Tube	Tip Rack for EZ2 Connect, for use with flip-cap tubes	9027010
EZ2 Connect Tip Rack - Large Volume	Tip Rack for EZ2 Connect, for use with large volume tubes	9027011
EZ2 Connect Cartridge Rack	Cartridge rack for EZ2 Connect, for use with prefilled kit cartridges	9027012
USB Flash Drive	USB flash drive provided by QIAGEN, for use with EZ2 USB ports	9026881
Barcode Reader	Barcode scanner to allow kit bar code and sample barcode scanning	9027101

For up-to-date licensing and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at [www.qiagen.com](http://www.qiagen.com) or can be requested from QIAGEN Technical Services or your local distributor.

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HB-2908-001 09/2021 © 2021 QIAGEN, all rights reserved.