

**User-developed  
protocol**

## User-Developed Protocol:

### Whole genome amplification from plasma and serum using the REPLI-g<sup>®</sup> Midi Kit

This procedure has been adapted by customers and is for whole genome amplification from plasma and serum using the REPLI-g Midi Kit. **The procedure has not been thoroughly tested and optimized by QIAGEN.**

**Note:** This protocol may be adapted for use with the REPLI-g Mini Kit, using the same reaction setup. In rare cases, potential inhibitors present in the starting material may have inhibitory effects on amplification when using the REPLI-g Mini Kit. In these cases, we recommend using the REPLI-g Midi Kit. Alternatively, upstream genomic DNA purification can be performed (e.g., using a QIAamp<sup>®</sup> Kit) with subsequent whole genome amplification of the purified DNA following the standard protocol in the *REPLI-g Mini/Midi Handbook*.

**IMPORTANT:** Please consult the “Safety Information” and “Important Notes” sections in the *REPLI-g Mini/Midi Handbook* before beginning this procedure. For safety information on the additional chemicals mentioned in this protocol, please consult the appropriate material safety data sheets (MSDSs) available from the product supplier.

#### Equipment and reagents to be supplied by user

- Microcentrifuge tubes
- Microcentrifuge
- Water bath or heating block
- Vortexer
- Pipets and pipet tips
- Ice
- Nuclease-free water
- TE buffer (10 mM Tris·Cl; 1 mM EDTA, pH 8.0)

#### Important points before starting

- REPLI-g Midi DNA Polymerase should be thawed on ice (see step 8). All other components can be thawed at room temperature.
- Buffer D2 should not be stored longer than 3 months.
- A DNA control reaction can be set up using 10 ng (1 µl) control genomic DNA (e.g., REPLI-g Human Control Kit, cat. no. 150090).

### Things to do before starting

- Prepare Buffer DLB by adding 500 µl nuclease-free water to the tube; mix thoroughly and centrifuge briefly.  
**Note:** Reconstituted Buffer DLB can be stored for 6 months at –20°C. Buffer DLB is pH-labile. Avoid neutralization with CO<sub>2</sub>.
- Set a water bath or heating block to 30°C.
- All buffers and reagents should be vortexed before use to ensure thorough mixing.

### Procedure

1. **Pipet a minimum of 50 µl plasma or serum sample into a microcentrifuge tube.**
2. **Centrifuge the sample at 6000 x g for 10 min and discard the supernatant.**
3. **Resuspend the cell pellet by adding 50 µl TE buffer and vortexing for 5 s. Centrifuge the pellet at 6000 x g for 10 min and discard the supernatant.**
4. **Resuspend the cell pellet by adding 10 µl TE buffer and vortexing for 5 s.**
5. **Prepare sufficient Buffer D2 (denaturation buffer) for the total number of whole genome amplification reactions (Table 1).**

**Note:** The total volume of Buffer D2 given in Table 1 is suitable for up to 6 reactions.

**Table 1. Preparation of Buffer D2**

<b>Component</b>	<b>Volume*</b>
DTT, 1 M	5 µl
Reconstituted Buffer DLB <sup>†</sup>	55 µl
<b>Total volume</b>	<b>60 µl</b>

\* Volumes given are suitable for up to 6 reactions. Excess Buffer D2 can be stored at -20°C for up to 3 months.

<sup>†</sup> Reconstitution of DLB is described in the “Things to do before starting” section.

6. **Add 10 µl Buffer D2 to each microcentrifuge tube containing cells obtained from plasma or serum. Mix by pipetting 3 times and place on ice for 10 min.**
7. **Add 10 µl Stop Solution to each microcentrifuge tube containing lysed cells and mix by pipetting 3 times.**  
**Note:** 10 µl lysed and neutralized cells are used in a 50 µl REPLI-g Midi reaction (Step10).
8. **Thaw REPLI-g Midi DNA Polymerase on ice. Thaw all other components at room temperature, vortex, and centrifuge briefly.**  
The REPLI-g Midi Reaction Buffer may form a precipitate after thawing. The precipitate will dissolve by vortexing for 10 s.

**9. Prepare a master mix on ice according to Table 2. Mix and centrifuge briefly.**

**IMPORTANT:** Add the master mix components in the order listed in Table 2. After addition of water and REPLI-g Midi Reaction Buffer, briefly vortex and spin down the mixture before addition of REPLI-g Midi DNA Polymerase. The master mix should be kept on ice and used immediately upon addition of the REPLI-g Midi DNA Polymerase.

**Table 2. Preparation of Master Mix**

<b>Component</b>	<b>Volume/reaction</b>
Nuclease-free water	10 $\mu$ l
REPLI-g Midi Reaction Buffer	29 $\mu$ l
REPLI-g Midi DNA Polymerase	1 $\mu$ l
<b>Total volume</b>	<b>40 <math>\mu</math>l</b>

**10. Add 10  $\mu$ l plasma or serum cell lysate (step 7) to 40  $\mu$ l master mix.**

**11. Incubate at 30°C for 8–16 h.**

Maximum DNA yield is achieved using an incubation time of 16 h. After incubation at 30°C, heat the water bath or heating block up to 65°C if the same water bath or heating block will be used in step 12.

**12. Inactivate REPLI-g Midi DNA Polymerase by heating the sample at 65°C for 3 min.**

**13. Store amplified DNA at 4°C for short-term storage or –20°C for long-term storage.**

DNA amplified using the REPLI-g Midi Kit should be treated as genomic DNA with minimal freeze-thaw cycles. Storage of nucleic acids at low concentration over a long period of time may result in acid hydrolysis. We therefore recommend storage of nucleic acids at a concentration of at least 100 ng/ $\mu$ l.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at [www.qiagen.com](http://www.qiagen.com) or can be requested from QIAGEN Technical Services or your local distributor.

Material safety data sheets (MSDS) for any QIAGEN product can be downloaded from [www.qiagen.com/Support/MSDS.aspx](http://www.qiagen.com/Support/MSDS.aspx).

Trademarks: QIAGEN®, QIAamp®, REPLI-g® (QIAGEN Group).

© 2005–2011 QIAGEN, all rights reserved .