

March 2015

BK Virus Probe ASR

Analyte Specific Reagent. Analytical and performance characteristics are not established.

Product	Contents	Cat. no.
BK Virus Probe ASR (24)	2 x 24 μ l BK Virus Probe ASR (1.5 μ M BK Virus Probe 1, 1.5 μ M BK Virus Probe 2) in buffered solution	1067392

Product description

The BK Virus Probe ASR contains a specific nucleic acid sequence for identification of BK virus DNA. The BK Virus Probe ASR is designed to detect a 274 bp region covering parts of the VP2 and VP3 genes in the BK virus genome in fluorescence channel Cycling Green (source 470 nm, detector 510 nm) of the real-time PCR cyclers. The DNA is synthetically derived.

Storage conditions

Upon receipt, store the BK Virus Probe ASR at -30°C to -15°C . Do not use after the expiration date printed on the outer box label. Repeated thawing and freezing (>2 x) should be avoided, as this may reduce reagent performance. If the reagents are to be used intermittently, they should be frozen in aliquots. Storage at $2-8^{\circ}\text{C}$ should not exceed a period of 5 hours.

Handling instructions

Store the probe away from any source of contaminating DNA, especially amplified DNA products. Use aerosol-barrier pipet tips for pipetting probe.



Quality control

This product was manufactured and released in accordance with the QIAGEN GmbH, Hilden, Quality Control and Quality Assurance procedures. The product has met QIAGEN GmbH, Hilden, quality requirements. Purity of the synthetic oligonucleotide is >80% as determined by high-performance liquid chromatography (HPLC).

Warnings and precautions

Dispose of in accordance with local regulations.

Safety information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient and compact PDF format at www.qiagen.com/safety where you can find, view, and print the SDS for each QIAGEN® kit and kit component.

Trademarks: QIAGEN® (QIAGEN Group).
THE PURCHASE OF THIS PRODUCT ALONE DOES NOT IMPLY ANY
LICENSE UNDER PATENTS OWNED BY ROCHE MOLECULAR SYSTEMS,
INC. OR F. HOFFMANN-LA ROCHE LTD.

