

## Application Note

# QIAxcel<sup>®</sup> system — mapping mutant gene loci in *Arabidopsis thaliana*

Nobuyoshi Mochizuki, Kyoto University, Graduate School of Science, Division of Biological Science, Kyoto, Japan

In this application note, we describe the assessment of simple sequence length polymorphism (SSLP) and cleaved amplified polymorphisms (CAPS) markers in the mapping of mutant gene loci and the homo/hetero examination of known mutant gene loci using the mutated CAPS or derived CAPS (dCAPS) markers. Both applications utilize the QIAxcel system and the QIAxcel DNA High Resolution Kit, with the OM500 method.

The QIAxcel system demonstrated superior separation capability, enabling the resolution of markers differing by only 4 base pairs in size. It is predicted that the QIAxcel system will be a powerful tool for increasing the speed of mapping and genotyping in the future.

## Introduction

There are several methods used in the identification of gene loci responsible for the mutation in *Arabidopsis thaliana*. Of these, mapping of the F<sub>2</sub> populations from crosses between different ecotypes has been a highly effective technique. However, as the entire genomic sequence became available, the increasing numbers of DNA markers made gene cloning much easier and faster using mapping and chromosomal walking. In addition, these DNA markers have often been used to detect polymorphisms of *Arabidopsis thaliana*. For polymorphism studies, we have routinely performed the following experiment types:

- Assessment of SSLP and CAPS markers in the mapping of mutant gene loci
- Homo/hetero examination of known mutant gene loci using the mutated CAPS or dCAPS markers
- Homo/hetero examination of T-DNA insertion in T-DNA insertion mutations

Design of PCR primers for the SSLP and dCAPS markers was extremely challenging, due to the sequence limitation, and the resulting PCR products were small (less than 200 bp). Consequently, the size differences among PCR products or restriction-enzyme-digested fragments are often just a few base pairs. Until recently, polyacrylamide gels or high concentration (2–4%) agarose gels were used to resolve such small size differences between DNA fragments. The introduction of the QIAxcel system has enabled the generation of rapid, reproducible results for these types of analysis.






Figure 2. Assessment of ecotypes based on the NGA707 marker.




Figure 4. Examination of mutations based on dCAPS.




Figure 3. Predicted sizes of restriction digested wild-type and mutant fragments.

## Conclusions

Although many DNA markers are recorded in the *Arabidopsis thaliana* database, markers with differences of only several bases pairs may be difficult to analyze. We anticipate that the use of the QIAxcel system, which demonstrates superior separation capability and simplicity, will enable the effective use of these markers and further increase the speed of mapping and genotyping.

## References

1. Bell, C.J., Ecker, J.R. (1994) Assignment of 30 microsatellite loci to the linkage map of *Arabidopsis*. *Genomics*, **1**, 137.
2. Konieczny, A., Ausubel, F.M. (1993) A procedure for mapping *Arabidopsis* mutations using co-dominant ecotype-specific PCR-based markers. *Plant J.* **4**, 403.
3. Neff, M.M., Turk, E., Kalishman, M. (2002) Web-based primer design for single nucleotide polymorphism analysis. *Trends Genet.* **18**, 6135.

## Ordering Information

Product	Contents	Cat. no.
QIAxcel Advanced system	Capillary electrophoresis device, including computer, and ScreenGel Software; 1-year warranty on parts and labor	9001941
QIAxcel DNA High Resolution Kit (1200)	QIAxcel DNA High Resolution Gel Cartridge, Buffers, Mineral Oil, QX Intensity Calibration Marker, 12-Tube Strips	929002

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at [www.qiagen.com](http://www.qiagen.com) or can be requested from QIAGEN Technical Services or your local distributor.

Visit [www.qiagen.com/mutant-mapping](http://www.qiagen.com/mutant-mapping) and find out how automated gel electrophoresis can benefit your lab!

Trademarks: QIAGEN®, QIAxcel®, Sample to Insight® (QIAGEN Group), (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

© 2016 QIAGEN, all rights reserved. PROM-3354-002

Ordering [www.qiagen.com/shop](http://www.qiagen.com/shop) | Technical Support [support.qiagen.com](mailto:support.qiagen.com) | Website [www.qiagen.com](http://www.qiagen.com)