
August 2021

QIAseq[®] 1-Step Amplicon Library Preparation Handbook

For the construction of Illumina[®]-compatible
libraries from multiplexed PCR amplicons

Contents

| | |
|---|----|
| Kit Contents | 4 |
| Shipping and Storage | 5 |
| Intended Use | 6 |
| Safety Information | 6 |
| Quality Control | 6 |
| Introduction | 7 |
| Principle and procedure | 7 |
| NGS adapter and index technologies | 9 |
| Starting materials | 11 |
| Equipment and Reagents to Be Supplied by User | 12 |
| Important Notes | 13 |
| Recommended library quantification method | 13 |
| Protocol: Library Construction Using the QIAseq 1-Step Amplicon Library Kit for Illumina Sequencing | 14 |
| Troubleshooting Guide | 20 |
| Appendix A: DNA Isolation and Quality Control | 22 |
| Appendix B: Amplicon Preparation and Quality Control | 24 |
| Appendix C: Optional A-Tailing Protocol | 26 |
| Appendix D: QIAseq Dual-Index Y-Adapters | 28 |
| Generation of sample sheets for Illumina instruments | 28 |
| Unique Dual-Index Y-Adapters | 28 |
| Combinatorial Dual-Index Y-Adapters | 51 |

| | |
|---|----|
| Appendix E: Adapter Indices for GeneRead Adapter I Sets A and B | 53 |
| Appendix F: Data Analysis | 55 |
| Ordering Information | 56 |
| Document Revision History | 60 |

Kit Contents

| | |
|---|---------------|
| QIAseq 1-Step Amplicon Library Kit | (12) |
| Catalog no. | 180412 |
| No. of reactions | 12 |
| 4x 1-Step Amplicon Buffer | 150 µl |
| 1-Step Amplicon Enzyme Mix | 24 µl |
| HiFi PCR Master Mix, 2x | 300 µl |
| Primer Mix Illumina Library Amp | 20 µl |
| RNase-free Water | 1.9 ml |
| Quick-Start Protocol | 1 |

| QIAseq 1-Step Amplicon Lib UDI/CDI Kit | UDI-A (96) | UDI-B (96) | UDI-C (96) | UDI-D (96) | CDI (96) |
|---|-------------------|-------------------|-------------------|-------------------|-----------------|
| Catalog no. | 180419 | 180420 | 180421 | 180425 | 180423 |
| No. of reactions | 96 | 96 | 96 | 96 | 96 |
| 4x 1-Step Amplicon Buffer | 1.2 ml | 1.2 ml | 1.2 ml | 1.2 ml | 1.2 ml |
| 1-Step Amplicon Enzyme Mix | 192 µl | 192 µl | 192 µl | 192 µl | 192 µl |
| HiFi PCR Master Mix, 2x | 1.25 ml | 1.25 ml | 1.25 ml | 1.25 ml | 1.25 ml |
| Primer Mix Illumina Library Amp | 150 µl | 150 µl | 150 µl | 150 µl | 150 µl |
| RNase-free Water | 5 x 1.9 ml | 5 x 1.9 ml | 5 x 1.9 ml | 5 x 1.9 ml | 5 x 1.9 ml |
| QIAseq UDI Y-Adapter Plate A/B/C/D (96) | 1 | 1 | 1 | 1 | – |
| QIAseq CDI Y-Adapter Plate (96) | – | – | – | – | 1 |
| QIAseq Y-Adapter Reference Card | 1 | 1 | 1 | 1 | 1 |
| Quick-Start Protocol | 1 | 1 | 1 | 1 | 1 |

| QIAseq CDI/UDI Y-Adapter Kit | CDI (24) | CDI (96) | UDI (24) | UDI A (96) | UDI B (96) | UDI C (96) | UDI D (96) |
|------------------------------|----------|----------|----------|------------|------------|------------|------------|
| Catalog no. | 180301 | 180303 | 180310 | 180312 | 180314 | 180316 | 180318 |
| Number of reactions | 24 | 96 | 24 | 96 | 96 | 96 | 96 |
| Adapter plate | 1 | 1 | 1 | 1 | 1 | 1 | 1 |
| Reference card | 1 | 1 | 1 | 1 | 1 | 1 | 1 |

Note: The QIAseq 1-Step Amplicon Lib CDI/UDI Kits (96) contain a QIAseq Y-Adapter Plate with either combinatorial dual-index adapters (CDI) or unique dual-index adapters (UDI).

To multiplex more than 96 libraries in a single sequencing run, combine kits with different UDI Y-adapter plates. For example, combining the QIAseq 1-Step Amplicon Lib UDI-A (or B or C or D) Kit (96) will allow the generation of 384 libraries with different sample indexes for 384-plex sequencing. For more information on QIAseq Y-Adapter Plates, please refer to Appendix D, page 28.

The QIAseq 1-Step Amplicon Library Kit (12) adapters are available in tube format. The GeneRead Adapter I Set A/B 12-plex (cat. 180985, 180986) can be ordered separately. The QIAseq 1-Step Amplicon Library Kit (12) is also fully compatible with all plate format QIAseq Y-Adapter Kits (24/96).

Shipping and Storage

The QIAseq 1-Step Amplicon Library Kit is shipped on dry ice and should be stored at -30 to -15°C upon arrival. When stored correctly, all reagents are stable for at least 6 months after delivery if not otherwise stated on the label.

QIAseq 1-Step Amplicon Lib CDI/UDI Kits (96) each contain an adapter plate shipped in a separate box. Store the adapter plate at -30 to -15°C upon arrival.

Intended Use

The QIAseq 1-Step Amplicon Library Kit is intended for molecular biology applications. This product is not intended for the diagnosis, prevention, or treatment of a disease.

All due care and attention should be exercised in the handling of the products. We recommend all users of QIAGEN® products to adhere to the NIH guidelines that have been developed for recombinant DNA experiments, or to other applicable guidelines.

Safety Information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient and compact PDF format at www.qiagen.com/safety where you can find, view and print the SDS for each QIAGEN kit and kit component.

Quality Control

In accordance with QIAGEN's ISO-certified Quality Management System, each lot of the QIAseq 1-Step Amplicon Library Kit is tested against predetermined specifications to ensure consistent product quality.

Introduction

While next-generation sequencing (NGS) has become a vital tool for life sciences and medical research, library preparation remains a key bottleneck in the NGS workflow. The QIAseq 1-Step Amplicon Library Kit is designed for the preparation of Illumina-compatible NGS libraries from multiplexed PCR or gene panel products, and employs an optimized 30-minute one-step library preparation protocol that reduces workflow duration, sample loss, and the potential for handling errors and cross-contamination. The kit accepts multiplexed PCR products from a variety of sources including the QIAGEN GeneRead® v2 panels, custom or lab-developed panels or multiplexed PCR assays, as well as other commercial panels. Optimized enzyme and buffer compositions ensure efficient library construction with a wide range of input amounts, and the entire protocol can optionally be performed at room temperature (15–25°C), enabling easy automation.

Principle and procedure

Purified amplicons from gene panels or multiplex PCR are converted to Illumina-compatible NGS libraries using a single enzymatic library construction step. During this reaction, amplicons are simultaneously prepared for ligation and barcoded adapters are ligated to both ends of the DNA inserts. The adapters contain sequences required for the PCR enrichment of the subsequent library, for flow-cell-binding during bridge amplification, and for sequencing primer binding sites for paired-end and multiplexed sequencing.




Figure 1. Scheme of optimized one-step amplicon library construction. Purified amplicons from multiplexed PCR or gene panels are converted to sequencing libraries by employing a 30-minute, one-tube library construction step. The libraries are purified with an easy and automatable size selection protocol, and the entire procedure can be performed at room temperature. For low-input applications, libraries can be amplified using the included HiFi PCR Master Mix.

Following library construction, excess adapters, adapter-dimers, and other reaction components are removed via precipitation onto Agencourt® AMPure® XP beads. This procedure is carried out at room temperature and can be easily automated on various liquid-handling platforms for high-throughput applications.

Following library purification, a high-fidelity library enrichment step can be performed to generate sufficient library from low amounts of starting material. This reaction relies on a high-fidelity DNA polymerase and optimized buffer conditions that ensure minimum GC bias and extremely low error rates.




Figure 2. Overview of a complete targeted resequencing NGS workflow including the QIAseq 1-Step Amplicon Library Kit. After DNA extraction and purification with an appropriate kit, target enrichment is performed with the QIAGEN Multiplex PCR Kit, GeneRead™ DNAseq Targeted Panels v2 or other gene panels or PCR products. NGS library construction is performed with the QIAseq 1-Step Amplicon Library Kit and data is analyzed with the CLC Genomics or Biomedical Workbench.

NGS adapter and index technologies

Sample multiplexing is one of the most important NGS tools for increasing throughput and reducing costs. It works by combining multiple samples to be processed together in a single sequencing run; as a consequence, sequencing reads need to be demultiplexed by reassigning each single read to its original source library. This is facilitated by the integration of index sequences into the individual adapter molecules.

QIAseq 1-Step Amplicon Lib CDI/UDI Kits (96) include a fully compatible indexing solution. We recommend using the QIAseq Dual-Index Y-Adapter Plates delivered with the kit. Each QIAseq 1-Step Amplicon Library Kit (96) includes one of the following:

- QIAseq Combinatorial Dual-Index (CDI) Y-Adapter Plate (96)
- QIAseq Unique Dual-Index (UDI) Y-Adapter Plate A, B, C, or D (96)

Combining QIAseq 1-Step Amplicon Lib UDI-A/B/C/D Kit (96) enables multiplexing of up to 384 samples per sequencing run. For more information on QIAseq Dual-Index Y-Adapters and index sequence motives, see Appendix D, page 28, and “Ordering Information”, page 56.

Adapter kits compatible with the QIAseq 1-Step Amplicon Library Kit (12) have to be purchased separately:

- GeneRead Adapter I Set A/B 12-plex (single-index adapters) (cat. no. 180985 or 180986)

GeneRead Single-Index adapter sets A and B are available in 12-plex format and can be combined with the QIAseq 1-Step Amplicon Library Kit (12). For more information on GeneRead Adapters and index sequence motives, see Appendix E, page 53, and “Ordering Information”, page 56.

CDI adapters use twelve i7 and eight i5 barcode motives that can be combined to form up to 96 combinatory dual indices. In contrast, QIAseq UDI Adapters use a fixed combination of 2 unique barcode motives per adapter molecule. Therefore, each single-index motive is only used exactly once on any UDI adapter plate.

Usage of UDI adapters effectively mitigates the risk of read misassignment due to index hopping. This is enabled by filtering misassigned reads during the demultiplexing of individual samples, thus generating highly accurate output data.

Compatible sequencing platforms

- Illumina NovaSeq™
- Illumina HiSeq®
- Illumina MiSeq®
- Illumina NextSeq®
- Illumina MiniSeq®

Starting materials

- PCR products generated with the GeneRead v2 DNAseq Targeted Panels
- PCR products generated with other custom or commercial gene panels
- PCR products generated with the QIAGEN Multiplex PCR Kit or other QIAGEN PCR reagents
- Multiplexed PCR amplicons generated with *Taq* or *Taq* derivatives
- Multiplexed PCR amplicons generated with Family B polymerases (see Appendix C)

Equipment and Reagents to Be Supplied by User

When working with chemicals, always wear a suitable lab coat, disposable gloves and protective goggles. For more information, consult the appropriate safety data sheets (SDSs), available from the product supplier.

For NGS library construction

- Only for 12-reaction kits: GeneRead Adapters (cat. no. 180985 or 180986; can be used with QIAseq 1-Step Amplicon Library Kit [12] for single-indexed libraries)
- Agencourt AMPure XP beads (cat. no. A63880 or A63881)
- 80% ethanol
- PCR tubes or plates
- Pipette tips and pipettes
- DNA LoBind tubes (from Axygen or Eppendorf)
- Vortexer
- Microcentrifuge
- Thermal cycler
- Magnetic racks for magnetic beads separation (e.g., Thermo Fisher Scientific/Life Technologies, DynaMag™-2 magnet, cat. no. 12321D)

For NGS library QC

- QIAxcel® Advanced Instrument (cat. no. 9001941) or similar capillary electrophoresis instrument
- GeneRead Library Quant Kit (cat. no. 180612)
- Real-time PCR machine for library quantification

Important Notes

PCR products should be free of contaminating primers, primer-dimers, or other PCR artifacts prior to library preparation.

Recommended library quantification method

QIAGEN's GeneRead Library Quant Kit (cat. no. 180612), which contains laboratory-verified forward and reverse primers together with a DNA standard, is highly recommended for accurate quantification of the prepared library.

We recommend purification of the PCR products prior to library preparation with Agencourt AMPure XP beads (see Appendix B).

Protocol: Library Construction Using the QIAseq 1-Step Amplicon Library Kit for Illumina Sequencing

This protocol describes library construction for sequencing on Illumina platforms.

Procedure

One-step end-repair and adapter ligation

1. Prepare a reaction mix for adapter ligation according to Table 1, adding the components to the PCR tube or plate containing purified PCR products.

Note: When using barcode adapters, open one adapter tube at a time and change gloves between pipetting the different barcode adapters to avoid cross-contamination.

Table 1. Reagents for adapter ligation of PCR product

| Component | Volume/reaction (µl) |
|---|----------------------|
| Purified PCR amplicons from gene panel or multiplex PCR product | Variable; 10–100 ng* |
| 4x 1-Step Amplicon Library Buffer | 12.5 |
| QIAseq CDI or UDI Y-Adapter Plate (96) <i>or</i> GeneRead 12-plex adapter (cat. nos. 180985 or 180986)† | 4 2 |
| 1-Step Amplicon Enzyme Mix | 2 |
| DNase-free water | Variable |
| Total | 50 |

* We recommend the quantification of PCR products via a microfluidics or capillary electrophoresis platform.

† If using adapters from another supplier, add the adapter to 1 µM final concentration or follow the supplier's directions.

2. Mix the components by pipetting up and down several times.

3. Program a thermal cycler to incubate at 25°C for 30 min. Optionally, the reaction can be incubated at room temperature.

IMPORTANT: Do not use a thermal cycler with a heated lid.

4. After the reaction is complete, place the reactions on ice and proceed with purification using Agencourt AMPure XP beads.

Cleanup of adapter-ligated DNA with Agencourt AMPure XP beads

5. Prepare 1.5 ml LoBind tubes for each ligation reaction and label tubes or prepare a 96-well plate (depending on availability of magnetic rack and individual preferences).
6. Transfer the 50 µl ligation reaction from step 4 to the prepared 1.5 ml LoBind tube or 96-well plate and add 50 µl nuclease-free water.

7. **For PCR products up to 350 bp in size:**

Add 40 µl (0.4x volume) Agencourt AMPure XP beads to each ligation reaction, mix well by pipetting, and then proceed to step 8.

For PCR products 350–1000 bp in size:

Add 60 µl (0.6x volume) Agencourt AMPure XP beads to each ligation reaction, mix well by pipetting, and then proceed to step 12 (**skip steps 8–11**).

8. Incubate the mixture for 5 min at room temperature. Pellet the beads on a magnetic stand (e.g., DynaMag) for 5 min.
9. Prepare a new LoBind tube for each ligation reaction or a new 96-well plate.
10. Carefully transfer 133 µl supernatant to the new tubes without disturbing the beads. This will leave behind about 7 µl supernatant. Discard the beads, which contain unwanted large DNA fragments. The large DNA fragments are generated by ligation of the adapter to nonspecific multiplex PCR products.
Note: Do not discard the supernatant.
11. Add 40 µl of resuspended Agencourt AMPure XP Beads slurry to the supernatant and mix well by pipetting.
12. Incubate the mixture for 5 min at room temperature.

13. Pulse-spin the tube or plate. Pellet the beads on a magnetic stand (e.g., DynaMag) for 5 min, then carefully remove and discard the supernatant. Be careful not to disturb the beads, because they contain the library.
Note: Do not discard the beads.
14. Wash the beads by adding 200 μ l of 80% ethanol. Pellet the beads on the magnetic stand and discard the supernatant. Repeat the wash once for a total of 2 ethanol washes.
15. Try to remove the residual ethanol as much as possible without disturbing the beads. Incubate the beads on the magnetic stand for 5–10 min or until the beads are dry. Over-drying of beads may result in lower DNA recovery. Remove from the magnetic stand.
16. Elute by resuspending in 26 μ l of nuclease-free water. Mix well by pipetting and incubate the tube or plate at room temperature for 2 min to elute the DNA from the beads.
17. Place the tube or plate back on the magnetic rack to pellet the beads. Incubate until the liquid is clear.
18. Use 23.5 μ l of the eluate in the library amplification procedure, or quantify and sequence directly if the amount of input amplicon was sufficient.

PCR amplification of purified library

19. Program a thermal cycler with a heated lid according to Table 2.

Table 2. Cycling conditions for the amplification of the DNA library

| Time | Temperature | Number of cycles |
|-----------------------------|-------------|------------------|
| Initial denaturation | | |
| 2 min | 98°C | 1 |
| Annealing | | |
| 20 s | 98°C | |
| 30 s | 60°C | 4–10* |
| 30 s | 72°C | |
| Final extension | | |
| 1 min | 72°C | 1 |
| ∞ | 4°C | Hold |

* **Note:** Cycle number depends on the amount and quality of input amplicon. In general, 4 cycles are sufficient for 20–500 ng of input PCR product and 10 cycles are sufficient for 1–20 ng of input PCR product. If input DNA is sufficient (>500 ng), library amplification can be omitted.

20. Mix the components in Table 3 in a 0.2 ml PCR tube or 96-well PCR plate.

Table 3. Reaction components for PCR amplification

| Component | Volume (µl) |
|----------------------------------|-------------|
| HiFi PCR Master Mix, 2x | 25 |
| Primer Mix (10 µM each) | 1.5 |
| Library DNA (from previous step) | 23.5 |
| Total | 50 |

21. Transfer the PCR plate to the thermal cycler and start the program.
22. Once PCR is complete, add 50 µl of resuspended Agencourt AMPure XP Beads to each reaction (50 µl) and pipet up and down thoroughly to mix the beads and PCR mix.
23. Incubate the mixture for 5 min at room temperature. Pellet the beads on a magnetic stand (e.g., DynaMag) and carefully discard the supernatant.
24. Wash the beads by adding 200 µl of 80% ethanol. Pellet the beads on the magnetic stand and discard the supernatant. Repeat the wash once for a total of two ethanol washes. Remove from the magnetic stand.
25. Incubate on the magnetic stand for 5–10 min or until the beads are dry. Over-drying of beads may result in lower DNA recovery. Remove from the magnetic stand.
26. Elute by resuspending in 30 µl nuclease-free water. Mix well by pipetting. Pellet the beads on the magnetic stand. Carefully transfer 28 µl supernatant to a clean LoBind 1.5 ml tube or PCR plate.
27. Assess the quality of the library using a capillary electrophoresis device or comparable method. Check for the expected size distribution of library fragments, for the absence of adapters, amplification primers, adapter-dimers or high-molecular-weight overamplification artifacts.

Note: The library should show a distribution reflecting the size of the input PCR amplicons plus 120 bp. The increase in library length reflects the addition of the sequencing adapters to the PCR amplicons.

Note: The median fragment size can be used in subsequent qPCR-based quantification methods to calculate library concentration.

28. Quantify the library using the GeneRead Library Quant Kit (cat. no. 180612 [not provided]) or a comparable method.
29. The purified library can be safely stored at -20°C in a DNA LoBind tube until ready to sequence.

Typical results

When high-quality, artifact-free PCR products are used; QIAseq 1-Step Amplicon libraries are typically free of adapter dimers, library amplification primers, excess adapters and high-molecular-weight amplification artifacts. When an appropriate cycle number for the input amount is chosen, library yield should be approximately 5 nM after purification, and the volume should be sufficient for quality control, library quantification and sequencing on most NGS platforms.

In the example experiment, as shown in Figure 3, reference DNA was amplified with the GeneRead DNaseq V2 Human Comprehensive Cancer Panel (06/2015 version). This panel consists of 4 separate PCR reactions, each including 1988 primer pairs. Amplicons are designed to capture 160 cancer-related genes comprising approximately 745 kb of the genome. After PCR, products were pooled and analyzed for yield and length distribution (Figure 3).




Figure 3. Electrophoresis trace of pooled multiplexed PCR amplicons. After amplification and pooling, a peak of expected size (160 bp), which represents the mean amplicon size of this set of PCR primers is seen. A second large peak centered at 35 bp is comprised of the electrophoresis marker, unextended primers and extended primer-dimers. During purification, these shorter products are removed, leaving only the desired PCR amplicons.

PCR products were purified with Agencourt AMPure XP Beads according to the *GeneRead DNaseq Targeted Panels V2 Handbook*, and 25 ng was used to generate libraries with the QIAseq 1-Step Amplicon Library Kit. Libraries were amplified with 4 cycles of PCR, and yielded a total of 28 μ l of 5 nM library after purification (Figure 4, 1:10 dilution shown). While sequencing quality, read length and the number of reads obtained will vary, depending on the sequencing platform, using a library quantification method such as the QIAGEN GeneRead DNaseq Library Quant Array enables more accurate clustering, optimizing data yield and quality.

Typical data generated with the QIAseq 1-Step Amplicon Library Kit contains minimal reads arising from adapter dimers, and when paired with a high-quality panel, often has >95% on-target reads, with even read distributions over target amplicons.




Figure 4. Electrophoresis trace of a purified library prepared using the QIAseq 1-Step Amplicon Library Kit. A characteristic size shift of 120 bp is observed due to the additional sequence added by the adapters. Typical libraries should be free of shorter peaks composed of excess primers, adapters or adapter-dimers and high-molecular-weight peaks, which can be introduced by the gene panel or through overamplification of the completed libraries.

Troubleshooting Guide

This troubleshooting guide may be helpful in solving any problems that may arise. For more information, see also the Frequently Asked Questions page at our Technical Support Center: www.qiagen.com/FAQ/FAQList.aspx. The scientists in QIAGEN Technical Services are always happy to answer any questions you may have about either the information and/or protocols in this handbook or sample and assay technologies (for contact information, visit www.qiagen.com).

Comments and suggestions

Low library yields

- | | |
|--|--|
| a) Suboptimal reaction conditions due to impurity with multiple PCR products | Make sure to use PCR purification methods that effectively remove impurities that could potentially inhibit QIAseq 1-Step Amplicon Library Kit enzymes. |
| b) Insufficient amount of starting DNA | Make sure at least 2 ng multiplex PCR products and correct cycle numbers for library amplification PCR are used for library prep. Consider increasing the amount of PCR products used or the number of cycles of library amplification |
| c) Over-drying of Agencourt AMPure XP beads | Over-drying Agencourt AMPure XP beads could decrease elution efficiency. Ensure that the beads are not dried for more than 10 min at room temperature. |

Unexpected signal peaks in capillary electrophoresis device traces

- | | |
|--|--|
| a) Presence of shorter peaks between 60 and 160 bp | These peaks represent library adapters (about 60 bp), adapter dimers (about 120 bp) or adapter-ligated PCR primer dimers (about 150–160 bp) that occur when there is no, or insufficient, adapter depletion after library preparation. As adapter dimers can form clusters on the flow cell and will be sequenced, this will reduce the capacity of the flow cell for the library fragments, even though a low ratio of adapter-dimers versus library will not be a problem. If necessary, repeat the protocol “Cleanup of adapter-ligated DNA with Agencourt AMPure XP beads”, page 15, to remove the residual adapters and adapter-dimers. |
| b) Presence of larger library fragments | If the fragment population shifts higher than expected after adapter ligation and PCR enrichment (e.g., more than the expected 120 bp shift), this can be due to insufficient depletion of the larger nonspecific fragments. If necessary, repeat the protocol “Cleanup of adapter-ligated DNA with Agencourt AMPure XP beads”, page 15, to remove the residual large fragments. |

Comments and suggestions

- c) Incorrect library fragment size after adapter ligation
- During library preparation, adapters of approximately 60 bp are ligated to both ends of the DNA library fragments. This should be reflected on a capillary electrophoresis device by a shift in size of all library fragments of 120 bp. If using library adapters from other suppliers, please refer to the size information given in the respective documentation. The absence of a clear size shift may indicate no, or only low, adapter ligation efficiency. Make sure to use the parameters and incubation times described in the handbook, as well as the correct amount of starting DNA.

Appendix A: DNA Isolation and Quality Control

High-quality DNA is essential for obtaining reliable sequencing results, and proper sample handling and DNA isolation procedures are critical to the success of the experiment. Residual traces of proteins, salts, extraction reagents or other contaminants can degrade DNA, interfere with DNA quantification or inhibit downstream PCR.

For the isolation and purification of high-quality genomic DNA, we recommend the following QIAGEN kits:

- QIAamp® DNA Mini Kit (cat. no. 51304) for the preparation of genomic DNA samples from fresh tissues
- GeneRead DNA FFPE Kit (cat. no. 180134) for the preparation of NGS-ready genomic DNA FFPE tissue samples
- MagAttract® HMW DNA Kit (cat. no. 67563) for isolation of high-molecular-weight genomic DNA
- QIAamp DNA Microbiome Kit (cat. no. 51704) for isolation of bacterial microbiome DNA from mixed samples
- QIAamp Circulating Nucleic Acid Kit (cat. no. 55114)

For accurate DNA quantification, we recommend the QIAxpert® (cat. no. 9002340).

Note: When using a UV-vis spectrophotometer to quantify DNA, ensure that samples have been treated to remove RNA, because the absorption spectra of RNA and DNA overlap significantly. For best results, DNA should be resuspended in DNase-free water or DNase-free 10 mM Tris buffer pH 8.0. Do not use DEPC-treated water.

For best results, all DNA samples should also demonstrate consistent quality according to the following criteria:

Concentration and purity determined by UV spectrophotometry

The concentration and purity of DNA should be determined with a UV-vis spectrophotometer. Prepare dilutions and measure absorbance in 10 mM Tris buffer pH 8.0, since the absorption spectra of nucleic acids is dependent on pH.

The $A_{260}:A_{280}$ ratio should be greater than 1.8.

The concentration determined by A_{260} should ideally be >2.5 $\mu\text{g}/\text{ml}$ DNA.

For accurate DNA quantification, we recommend the QIAxpert.

DNA integrity

For best results, genomic DNA should be greater than 2 kb in length, with many fragments greater than 10 kb. This can be checked by running a fraction of each DNA sample on a 1% agarose gel.

FFPE DNA

If FFPE DNA will be used, PCR conditions and cycle number for target enrichment may need to be optimized. This can be done with standard qPCR reagents, or for the GeneRead V2 DNAseq panels, the QIAGEN GeneRead DNA QuantiMIZE assay.

Appendix B: Amplicon Preparation and Quality Control

Amplicon preparation and quality control

The QIAseq 1-Step Amplicon Library Kit accepts amplicons generated from a wide variety of multiplex PCR reactions or gene panels. Regardless of the PCR system used, high-quality amplicons are key to the production of high-quality NGS libraries.

When conducting target enrichment via multiplex PCR, avoid overamplification into the plateau phase. When overamplified, amplicons produced in earlier cycles of PCR serve as both primer and template for further amplification, generating long, chimeric molecules. These are visible as high-molecular-weight products with a broad range of lengths and can be visible on the QIAxcel or Agilent® Bioanalyzer trace as a second peak much larger than the amplicons of interest, or as a high-molecular-weight smear sometimes extending into the wells on an agarose gel. These chimeric products can interfere with library preparation and sequencing, and if present, cycle number should be decreased appropriately.

Following PCR, amplicons should be purified with Agencourt AMPure XP Beads or a similar cleanup technology to remove excess primers, primer dimers and buffer components. As the QIAseq 1-Step Amplicon Library Kit accepts a wide range of DNA inputs, amplicon quantification after purification is not strictly necessary if a robust DNA extraction and PCR system is available, but is recommended. Amplicons can be quantified on a QIAxcel, agarose gel, Agilent Bioanalyzer, Qubit®, or NanoDrop®.

Low-diversity libraries

In contrast to a standard whole genome library, the base composition at the beginning of the reads in an amplicon-sequencing library reflects the base composition of the primers used to produce the PCR amplicons. In experiments involving hundreds or thousands of primer pairs,

the sequence diversity at the beginning of each read is high; however in low-multiplexing experiments, where only 8–12 primer pairs are used, sequence diversity is low.

Low sequence diversity can interfere with normal cluster calling, phasing and quality matrix calculation on Illumina platforms. We recommend users to follow all guidelines set forth by the instrument manufacturer for the sequencing of low-diversity libraries if lower-complexity amplicons are used. These may include the intentional underclustering of the flow cell or mixing amplicon libraries with the PhiX standard or barcoded high-diversity whole genome or RNAseq libraries.

Appendix C: Optional A-Tailing Protocol

While the QIAseq 1-step Amplicon Library Kit accepts PCR products from a wide range of sources, some care must be taken to confirm polymerase compatibility. The novel one-step reaction requires that PCR amplicons contain 3' A-overhangs for efficient ligation. *Taq* polymerase, the most commonly used thermostable DNA polymerase, and its derivatives, by default carry out this non-templated A-addition during the PCR reaction. *Taq* and *Taq* derivatives have attributes that make them amenable to multiplex PCR, and many commercial gene panels employ a *Taq*-based enzyme.

In contrast to *Taq*, other polymerases with strong 3'–5' exonuclease activities do not carry out this reaction. While these enzymes are not commonly used for multiplexed PCR, amplicons produced with such enzymes are still compatible with the QIAseq 1-Step Amplicon Library Kit, but require A-tailing prior to ligation. A suggested A-tailing protocol is given below.

Materials required for A-tailing

- MinElute® PCR Purification Kit (cat. no. 28004 or 28006), or Agencourt AMPure XP Beads (Beckman Coulter, cat. no. A63880)
- Klenow (3'→ 5' exo-) Low Concentration (Enzymatics, P7010-LC-L)
- 10X Blue Buffer (Enzymatics, B0110)
- 100 mM dATP Solution (Enzymatics, N2010-A-L)
- Pipette tips and pipettes

Procedure

1. Mix the components in Table 4 and add to the PCR tube or plate containing purified PCR products.

Table 4. Reaction components for A-tailing

| Component | Volume (µl) |
|--|--------------------|
| Purified PCR product | 20 |
| 10x Blue Buffer | 5 |
| 100 mM dATP | 1 |
| Klenow (3'→ 5' exo-) low concentration | 1 |
| DNase-free water | 23 |
| Total | 50 |

2. Incubate in a thermal cycler or heating block for 30 minutes at 37°C.
3. Purify with the MinElute PCR Purification Kit or Agencourt AMPure XP Beads as per the manufacturer's directions.

Appendix D: QIAseq Dual-Index Y-Adapters

Generation of sample sheets for Illumina instruments

Index sequences for QIAseq Unique and Combinatorial Dual-Index Y-Adapters are available for download at www.qiagen.com. Sequencing on the NextSeq, HiSeq X™, or HiSeq 3000/4000 system follow a different dual-indexing workflow than other Illumina systems. If you are manually creating sample sheets for these instruments, enter the reverse complement of the i5 index adapter sequence. If you are using Illumina Experiment Manager, BaseSpace, or Local Run Manager for run planning, the software will automatically reverse complement index sequences when necessary.

Ready-to-use sample sheets containing all QIAseq CDI and UDI Y-Adapter barcode sequences are available for MiSeq, NextSeq, MiniSeq, HiSeq, and NovaSeq instruments. These can be imported and edited using the Illumina Experiment Manager Software, Illumina Local Run Manager, or any text editor. Make sure to download the appropriate sample sheet for NextSeq, HiSeq X, or HiSeq 3000/4000 systems depending on whether you are using Local Run Manager or manually configuring the sequencing run.

Unique Dual-Index Y-Adapters

The layouts of the 24-plex and 96-plex (A/B/C/D) single-use UDI adapter plates are shown from Figure 5 to Figure 9. The index motives used in the QIAseq Unique Dual-Index Kits are listed in Table 5. To make sequencing preparation more convenient, you can download Illumina-compatible sample sheets for different sequencing instruments at www.qiagen.com.

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 |
|----------|---------|---------|---------|-------|-------|-------|-------|-------|-------|-------|-------|
| A | UDI 001 | UDI 009 | UDI 017 | empty | empty | empty | empty | empty | empty | empty | empty |
| B | UDI 002 | UDI 010 | UDI 018 | empty | empty | empty | empty | empty | empty | empty | empty |
| C | UDI 003 | UDI 011 | UDI 019 | empty | empty | empty | empty | empty | empty | empty | empty |
| D | UDI 004 | UDI 012 | UDI 020 | empty | empty | empty | empty | empty | empty | empty | empty |
| E | UDI 005 | UDI 013 | UDI 021 | empty | empty | empty | empty | empty | empty | empty | empty |
| F | UDI 006 | UDI 014 | UDI 022 | empty | empty | empty | empty | empty | empty | empty | empty |
| G | UDI 007 | UDI 015 | UDI 023 | empty | empty | empty | empty | empty | empty | empty | empty |
| H | UDI 008 | UDI 016 | UDI 024 | empty | empty | empty | empty | empty | empty | empty | empty |

Figure 5. QIAseq UDI Y-Adapter Plate (24) layout (UDI 1–24).

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|----------|---------|---------|---------|---------|---------|---------|---------|---------|---------|---------|---------|---------|
| A | UDI 001 | UDI 009 | UDI 017 | UDI 025 | UDI 033 | UDI 041 | UDI 049 | UDI 057 | UDI 065 | UDI 073 | UDI 081 | UDI 089 |
| B | UDI 002 | UDI 010 | UDI 018 | UDI 026 | UDI 034 | UDI 042 | UDI 050 | UDI 058 | UDI 066 | UDI 074 | UDI 082 | UDI 090 |
| C | UDI 003 | UDI 011 | UDI 019 | UDI 027 | UDI 035 | UDI 043 | UDI 051 | UDI 059 | UDI 067 | UDI 075 | UDI 083 | UDI 091 |
| D | UDI 004 | UDI 012 | UDI 020 | UDI 028 | UDI 036 | UDI 044 | UDI 052 | UDI 060 | UDI 068 | UDI 076 | UDI 084 | UDI 092 |
| E | UDI 005 | UDI 013 | UDI 021 | UDI 029 | UDI 037 | UDI 045 | UDI 053 | UDI 061 | UDI 069 | UDI 077 | UDI 085 | UDI 093 |
| F | UDI 006 | UDI 014 | UDI 022 | UDI 030 | UDI 038 | UDI 046 | UDI 054 | UDI 062 | UDI 070 | UDI 078 | UDI 086 | UDI 094 |
| G | UDI 007 | UDI 015 | UDI 023 | UDI 031 | UDI 039 | UDI 047 | UDI 055 | UDI 063 | UDI 071 | UDI 079 | UDI 087 | UDI 095 |
| H | UDI 008 | UDI 016 | UDI 024 | UDI 032 | UDI 040 | UDI 048 | UDI 056 | UDI 064 | UDI 072 | UDI 080 | UDI 088 | UDI 096 |

Figure 6. QIAseq UDI Y-Adapter Plate A (96) layout (UDI 1–96).

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|----------|------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|
| A | UDI 097 | UDI 105 | UDI 113 | UDI 121 | UDI 129 | UDI 137 | UDI 145 | UDI 153 | UDI 161 | UDI 169 | UDI 177 | UDI 185 |
| B | UDI 098 | UDI 106 | UDI 114 | UDI 122 | UDI 130 | UDI 138 | UDI 146 | UDI 154 | UDI 162 | UDI 170 | UDI 178 | UDI 186 |
| C | UDI 099 | UDI 107 | UDI 115 | UDI 123 | UDI 131 | UDI 139 | UDI 147 | UDI 155 | UDI 163 | UDI 171 | UDI 179 | UDI 187 |
| D | UDI 100 | UDI 108 | UDI 116 | UDI 124 | UDI 132 | UDI 140 | UDI 148 | UDI 156 | UDI 164 | UDI 172 | UDI 180 | UDI 188 |
| E | UDI 101 | UDI 109 | UDI 117 | UDI 125 | UDI 133 | UDI 141 | UDI 149 | UDI 157 | UDI 165 | UDI 173 | UDI 181 | UDI 189 |
| F | UDI 102 | UDI 110 | UDI 118 | UDI 126 | UDI 134 | UDI 142 | UDI 150 | UDI 158 | UDI 166 | UDI 174 | UDI 182 | UDI 190 |
| G | UDI 103 | UDI 111 | UDI 119 | UDI 127 | UDI 135 | UDI 143 | UDI 151 | UDI 159 | UDI 167 | UDI 175 | UDI 183 | UDI 191 |
| H | UDI 104 | UDI 112 | UDI 120 | UDI 128 | UDI 136 | UDI 144 | UDI 152 | UDI 160 | UDI 168 | UDI 176 | UDI 184 | UDI 192 |

Figure 7. QIAseq UDI Y-Adapter Plate B (96) layout (UDI 97–192).

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|----------|------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|
| A | UDI 193 | UDI 201 | UDI 209 | UDI 217 | UDI 225 | UDI 233 | UDI 241 | UDI 249 | UDI 257 | UDI 265 | UDI 273 | UDI 281 |
| B | UDI 194 | UDI 202 | UDI 210 | UDI 218 | UDI 226 | UDI 234 | UDI 242 | UDI 250 | UDI 258 | UDI 266 | UDI 274 | UDI 282 |
| C | UDI 195 | UDI 203 | UDI 211 | UDI 219 | UDI 227 | UDI 235 | UDI 243 | UDI 251 | UDI 259 | UDI 267 | UDI 275 | UDI 283 |
| D | UDI 196 | UDI 204 | UDI 212 | UDI 220 | UDI 228 | UDI 236 | UDI 244 | UDI 252 | UDI 260 | UDI 268 | UDI 276 | UDI 284 |
| E | UDI 197 | UDI 205 | UDI 213 | UDI 221 | UDI 229 | UDI 237 | UDI 245 | UDI 253 | UDI 261 | UDI 269 | UDI 277 | UDI 285 |
| F | UDI 198 | UDI 206 | UDI 214 | UDI 222 | UDI 230 | UDI 238 | UDI 246 | UDI 254 | UDI 262 | UDI 270 | UDI 278 | UDI 286 |
| G | UDI 199 | UDI 207 | UDI 215 | UDI 223 | UDI 231 | UDI 239 | UDI 247 | UDI 255 | UDI 263 | UDI 271 | UDI 279 | UDI 287 |
| H | UDI 200 | UDI 208 | UDI 216 | UDI 224 | UDI 232 | UDI 240 | UDI 248 | UDI 256 | UDI 264 | UDI 272 | UDI 280 | UDI 288 |

Figure 8. QIAseq UDI Y-Adapter Plate C (96) layout (UDI 193–288).

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|
| A | UDI 289 | UDI 297 | UDI 305 | UDI 313 | UDI 321 | UDI 329 | UDI 337 | UDI 345 | UDI 353 | UDI 361 | UDI 369 | UDI 377 |
| B | UDI 290 | UDI 298 | UDI 306 | UDI 314 | UDI 322 | UDI 330 | UDI 338 | UDI 346 | UDI 354 | UDI 362 | UDI 370 | UDI 378 |
| C | UDI 291 | UDI 299 | UDI 307 | UDI 315 | UDI 323 | UDI 331 | UDI 339 | UDI 347 | UDI 355 | UDI 363 | UDI 371 | UDI 379 |
| D | UDI 292 | UDI 300 | UDI 308 | UDI 316 | UDI 324 | UDI 332 | UDI 340 | UDI 348 | UDI 356 | UDI 364 | UDI 372 | UDI 380 |
| E | UDI 293 | UDI 301 | UDI 309 | UDI 317 | UDI 325 | UDI 333 | UDI 341 | UDI 349 | UDI 357 | UDI 365 | UDI 373 | UDI 381 |
| F | UDI 294 | UDI 302 | UDI 310 | UDI 318 | UDI 326 | UDI 334 | UDI 342 | UDI 350 | UDI 358 | UDI 366 | UDI 374 | UDI 382 |
| G | UDI 295 | UDI 303 | UDI 311 | UDI 319 | UDI 327 | UDI 335 | UDI 343 | UDI 351 | UDI 359 | UDI 367 | UDI 375 | UDI 383 |
| H | UDI 296 | UDI 304 | UDI 312 | UDI 320 | UDI 328 | UDI 336 | UDI 344 | UDI 352 | UDI 360 | UDI 368 | UDI 376 | UDI 384 |

Figure 9. QIAseq UDI Y-Adapter Plate D (96) layout (UDI 289–384).

Table 5. UDI motives used in the QIAseq UDI Y-Adapter Kits (24 and 96 A/B/C/D)

Unique Dual-Index adapters 1–24 are identical on the adapter plates of the QIAseq UDI Y-Adapter Kit (24) and QIAseq UDI Y-Adapter Kit A (96).

Note: Sequencing on the MiniSeq, NextSeq, HiSeqX, and HiSeq 3000/4000 systems follow a different dual-indexing workflow than other Illumina systems, which requires the reverse complement of the i5 index adapter sequence.

| Indices for entry on sample sheet | | | |
|-----------------------------------|---|---|--|
| Unique dual-index number | i5 bases for entry on sample sheet (MiSeq, HiSeq 2000/2500, NovaSeq 6000) | i5 bases for entry on sample sheet (iSeq100, MiniSeq, NextSeq, HiSeqX, HiSeq 3000/4000) | i7 bases for entry on sample sheet (all instruments) |
| UDI 001 | ATGGCCGACT | AGTCGGCCAT | TGAACGTTGT |
| UDI 002 | CGATGAGCAC | GTGCTCATCG | ACCAGACTTG |
| UDI 003 | GATAAGTCGA | TCGACTTATC | ACTGGCGAAC |
| UDI 004 | TCACGCCTTG | CAAGGCGTGA | GCGTTAGGCA |
| UDI 005 | AGGAACACAA | TTGTGTTCTT | TTATCGGCCT |
| UDI 006 | CTCAGTAGGC | GCCTACTGAG | GAGGTATAAG |
| UDI 007 | GAAGTGCCTG | CAGGCACCTC | TCAAGGATTC |
| UDI 008 | TCTCTCGCCT | AGGCGAGAGA | CGAACCGAGA |
| UDI 009 | AGGCACCTTC | GAAGGTGCCT | GAGCCAAGTT |
| UDI 010 | CTGTTGGTAA | TTACCAACAG | AAGGCCGTAG |
| UDI 011 | GCTGGTACCT | AGGTACCAGC | TTAGAGAAGC |
| UDI 012 | TAAGGAGCGG | CCGCTCCTTA | TCTAAGACCA |
| UDI 013 | AATCGCTCCA | TGGAGCGATT | TGTAACCACT |
| UDI 014 | CTCCTAATTG | CAATTAGGAG | CCGACACAAG |
| UDI 015 | GCCTCATAAT | ATTATGAGGC | CTCTGATGGC |
| UDI 016 | TGTATTGAGC | GCTCAATACA | CGGCCTGTTA |
| UDI 017 | AGCCATAACA | TGTTATGGCT | TGCATAGCTT |
| UDI 018 | CCACAAGTGG | CCACTTGTGG | AACCTTCTCG |
| UDI 019 | GTTATCACAC | GTGTGATAAC | AAGAGATCAC |
| UDI 020 | TACCGTTCTT | AAGAACGGTA | GCCTGAAGGA |

Table continues on next page

Table continued from previous page

| Indices for entry on sample sheet | | | |
|--|--|--|---|
| Unique dual-index number | i5 bases for entry on sample sheet (MiSeq, HiSeq 2000/2500, NovaSeq 6000) | i5 bases for entry on sample sheet (iSeq100, MiniSeq, NextSeq, HiSeqX, HiSeq 3000/4000) | i7 bases for entry on sample sheet (all instruments) |
| UDI 021 | AGGCCGTTAGG | CCTAACGCCT | ATTGTGCCTT |
| UDI 022 | CCGTAACGTC | GACGTACGG | TCCTCTACCG |
| UDI 023 | GTAATAGCCA | TGGCTATTAC | TACCATGAAC |
| UDI 024 | TAGCGCCGAT | ATCGGCGCTA | CATTGGCAGA |
| UDI 025 | CATTCTTGGGA | TCCAAGAATG | CACTGCTATT |
| UDI 026 | ATGCAAGGTT | AACCTGCAT | AATGGTAGGT |
| UDI 027 | CGCCAGACAA | TTGTCTGGCG | GATACCTATG |
| UDI 028 | GAAGGTTGGC | GCCAACCTTC | CACTAGGTAC |
| UDI 029 | TCGCATCACG | CGTGATGCGA | AGCTCGTCA |
| UDI 030 | CCGTCATGA | TCATGACCGG | TGTCAGTCTT |
| UDI 031 | ATCACAAGC | GCTTGTAAT | GATGAACAGT |
| UDI 032 | CAACCTGTAA | TTACAGGTTG | ACAATCGGCG |
| UDI 033 | GCCAGTCGTT | AACGACTGGC | GATTGAGTTC |
| UDI 034 | TGCCTTGTCG | CGACAAGGCA | GTAATGCCAA |
| UDI 035 | CTATCCGCTG | CAGCGGATAG | TCGTTGCGCT |
| UDI 036 | AATGCCGGA | TTCCGGCATT | AGGTGAGTAT |
| UDI 037 | CGGTTATCCG | CGGATAACCG | TCGATAATGG |
| UDI 038 | GCGGAAGAGT | ACTCTCCGC | GCGTCTCTTC |
| UDI 039 | TTGGTTAGTC | GACTAACCAA | GTCTCCTGCA |
| UDI 040 | TTCAGTGTGA | TCACACTGAA | GAGCTTCATT |

Table continues on next page

Table continued from previous page

| Indices for entry on sample sheet | | | |
|--|--|--|---|
| Unique dual-index number | i5 bases for entry on sample sheet (MiSeq, HiSeq 2000/2500, NovaSeq 6000) | i5 bases for entry on sample sheet (iSeq100, MiniSeq, NextSeq, HiSeqX, HiSeq 3000/4000) | i7 bases for entry on sample sheet (all instruments) |
| UDI 041 | AGAATTCTGG | CCAGAATTCT | AGGCCTACAT |
| UDI 042 | CATTGACTCT | AGAGTCAATG | TGTGGAACCG |
| UDI 043 | GCGGCTTCAA | TTGAAGCCGC | CGTATTAAGC |
| UDI 044 | TTATGGTCTC | GAGACCATAA | CCAGTGGTTA |
| UDI 045 | CGTAACCAGG | CCTGGTACG | GCGTTCGAGT |
| UDI 046 | AGCTCAGATA | TATCTGAGCT | CCTCCGGTT |
| UDI 047 | CCGGTGTTAC | GTAACACCGG | CACAAGACGG |
| UDI 048 | GACCTAACCT | AGGTTAGGTC | GCTTACACAC |
| UDI 049 | TTGTAGAAGG | CCTTCTACAA | AGGATGTCCA |
| UDI 050 | CCTAGCACTA | TAGTGCTAGG | CACCTTATGT |
| UDI 051 | ATCGTGTCT | AGAACACGAT | AAGCGGCTGT |
| UDI 052 | CCAACTTATC | GATAAGTTGG | TTCCTGTGAG |
| UDI 053 | GAAGCCAAGG | CCTGGCTTC | AGTACAGTTC |
| UDI 054 | TGGAGTCAA | TTGAACTCCA | TACAGCCTCA |
| UDI 055 | CTCAATCCT | AGGATTGAAG | GTTCTATTGG |
| UDI 056 | ATCTTGCGTG | CACGCAAGAT | ATATACCGGT |
| UDI 057 | CGTCTAAGGT | ACCTTAGACG | CCTCGGAATG |
| UDI 058 | GAGGTGAACA | TGTTACCTC | GTTCTGGAAC |
| UDI 059 | TCAGAACTAC | GTAGTCTGA | AGATTCACCA |
| UDI 060 | CGGATATTGA | TCAATATCCG | TCGGTCAGAT |

Table continues on next page

Table continued from previous page

| Indices for entry on sample sheet | | | |
|--|--|--|---|
| Unique dual-index number | i5 bases for entry on sample sheet (MiSeq, HiSeq 2000/2500, NovaSeq 6000) | i5 bases for entry on sample sheet (iSeq100, MiniSeq, NextSeq, HiSeqX, HiSeq 3000/4000) | i7 bases for entry on sample sheet (all instruments) |
| UDI 061 | AGGAGTAGAT | ATCTACTCCT | CACTCTCGCT |
| UDI 062 | CCGCCGAATA | TATTCGGCGG | GTTGGTCCAG |
| UDI 063 | GAGTCTATAC | GTATAGACTC | AGCTCGAAGC |
| UDI 064 | TTATTACCGG | CCGGTAATAA | AGAGGTTCTA |
| UDI 065 | CGCTCGTTAG | CTAACGAGCG | ATGACTCGAA |
| UDI 066 | AACAACGCTG | CAGCGTGTGT | GAACAATCCT |
| UDI 067 | CGCGGTATT | AATAGCCGCG | TGGCAAGGAG |
| UDI 068 | GCTCGACACA | TGTGTCGAGC | GAATATTGGC |
| UDI 069 | TTCTCCAAC | GTTGGAAGAA | CCGGAACCTA |
| UDI 070 | TTGGCGGTTG | CAACCGCCAA | ACTTGTTCGG |
| UDI 071 | AACAGGCAAT | ATTGCCTGTT | CAAGTCCAAT |
| UDI 072 | CAGAATGGCG | CGCCATTCTG | AACCGCAAGG |
| UDI 073 | GTTGAGATTC | GAATCTCAAC | ACGTTGACTC |
| UDI 074 | TGTGTGCGGA | TCCGCACACA | CCACTTAACA |
| UDI 075 | GTTCCGCGAA | TTCGCCGAAC | AGCAGTTCCT |
| UDI 076 | AGCTGTATTG | CAATACAGCT | TCGCCTTCGT |
| UDI 077 | CAGCGGATGA | TCATCCGCTG | TAGGACTGCG |
| UDI 078 | GTCCTTGGAT | ATCCAAGGAC | TCCGAGCGAA |
| UDI 079 | TCTAGATGCT | AGCATCTAGA | TTCGGTGTGT |
| UDI 080 | CGAGCCACAT | ATGTGGCTCG | ACAGGAGGAA |

Table continues on next page

Table continued from previous page

| Indices for entry on sample sheet | | | |
|--|--|--|---|
| Unique dual-index number | i5 bases for entry on sample sheet (MiSeq, HiSeq 2000/2500, NovaSeq 6000) | i5 bases for entry on sample sheet (iSeq100, MiniSeq, NextSeq, HiSeqX, HiSeq 3000/4000) | i7 bases for entry on sample sheet (all instruments) |
| UDI 081 | ATGGAATGGA | TCCATTCCAT | CCTCCATTAA |
| UDI 082 | CATTCTCAC | GTGAGGAATG | AGTCGCGGTT |
| UDI 083 | GCATAGGAAG | CTTCTATGC | CTCATCCAGG |
| UDI 084 | TGTTCTGTT | AACACGAACA | TGTGGTTGAA |
| UDI 085 | TAAGACCGTT | AACGGTCTTA | TTATGCGTGG |
| UDI 086 | ATGGTACCAG | CTGGTACCAT | GCGAATGTAT |
| UDI 087 | CCGACAGCTT | AAGCTGTCGG | GTC AAGCTCG |
| UDI 088 | GACGATATGA | TCATATCGTC | TAGAGTTGGA |
| UDI 089 | TTGTA CTCCA | TGGAGTACAA | CTGATGATCT |
| UDI 090 | GTGCACATAA | TTATGTGCAC | ACTAGGTGTT |
| UDI 091 | AGGACAAGTA | TACTTGTCCT | CTGTTAGCGG |
| UDI 092 | CCGATTCGAG | CTCGAATCGG | ATCGCACCAA |
| UDI 093 | GTAGGA ACTT | AAGTTCCTAC | CITACTGGT |
| UDI 094 | TACACTACGA | TCGTAGTGTA | CCTTAATGCG |
| UDI 095 | ATGACCTTGA | TCAAGGTCAT | TCTCGCTAG |
| UDI 096 | CTACGTGACG | CGTCACGTAG | TCTTCAGAGA |
| UDI 097 | AACAATCAGG | CCTGATTGTT | TACCGGTGGT |
| UDI 098 | CTGGTGTGCA | TGCACACCAG | AGGTGTTACG |
| UDI 099 | GCATATCCTT | AAGGATATGC | ACAGACCGAC |
| UDI 100 | TGTCCTGTAC | GTACAGGACA | CGAATACGTA |

Table continues on next page

Table continued from previous page

| Indices for entry on sample sheet | | | |
|--|--|--|---|
| Unique dual-index number | i5 bases for entry on sample sheet (MiSeq, HiSeq 2000/2500, NovaSeq 6000) | i5 bases for entry on sample sheet (iSeq100, MiniSeq, NextSeq, HiSeqX, HiSeq 3000/4000) | i7 bases for entry on sample sheet (all instruments) |
| UDI 101 | AGAACGTCGC | GCGACGTCT | TAGCATCGAT |
| UDI 102 | CACGGACTAG | CTAGTCCGTG | CCATGAGTCG |
| UDI 103 | GTTGAACACT | AGTGTCAAC | ACTAACATGC |
| UDI 104 | TCGCGTGGTA | TACCACGCGA | ACACTCTCTA |
| UDI 105 | AGCCACTATG | CATAGTGGCT | GCTCTTGCTT |
| UDI 106 | CCACCTACCA | TGGTAGGTGG | AATCTTGAGG |
| UDI 107 | GTTCCGGTGT | ACACCGGAAC | CTAACGGTC |
| UDI 108 | TAGGTCTGAC | GTCAGACCTA | TTGTGACCAA |
| UDI 109 | AGGAAGCATT | AATGCTTCT | TCACACACCT |
| UDI 110 | CCTTAGTTGG | CCAACAAAGG | CTGCAATTAG |
| UDI 111 | GTCCTATTCA | TGAATAGGAC | CTCCTTACTC |
| UDI 112 | TAAGATGGAC | GTCCATCTTA | GCAACGCAGA |
| UDI 113 | AGGCCATGGT | ACCATGGCCT | CCTTACCAAT |
| UDI 114 | CATTGGCCAA | TTGGCCAATG | TTAATCCTCG |
| UDI 115 | GCTATGAATC | GATTCATAGC | TTCCGAGTTC |
| UDI 116 | TTGGTCTCG | CGAGGACCAA | CTCGAGAGGA |
| UDI 117 | AGCGACATAC | GTATGTCGCT | TGTTGGCTGT |
| UDI 118 | CAAGTAGTCT | AGACTACTTG | CGTATCTGCG |
| UDI 119 | GTCAAGAAGA | TCTTCTGAC | CCATAGTATC |
| UDI 120 | TCCTGTTATG | CATAACAGGA | TGGACAGTAA |

Table continues on next page

Table continued from previous page

| Indices for entry on sample sheet | | | |
|--|--|--|---|
| Unique dual-index number | i5 bases for entry on sample sheet (MiSeq, HiSeq 2000/2500, NovaSeq 6000) | i5 bases for entry on sample sheet (iSeq100, MiniSeq, NextSeq, HiSeqX, HiSeq 3000/4000) | i7 bases for entry on sample sheet (all instruments) |
| UDI 121 | AAGTGC GATA | TATCGCACTT | GTACCTTGTT |
| UDI 122 | AGGCTACACG | CGTGTAGCCT | GAGTGCCTCT |
| UDI 123 | CTATATCGGC | GCCGATATAG | TAAGTAGCGG |
| UDI 124 | GCTAAGGTAA | TTACCTTAGC | CGTGGTGTTT |
| UDI 125 | TAACCTGGTT | AACCAGGTTA | CATTCCTGAA |
| UDI 126 | AGTTGGTCTA | TAGACCAACT | AAGATGCATG |
| UDI 127 | ATGCAGCTGG | CCAGTGCAT | CCTTGGAGCT |
| UDI 128 | CGTTGCCTTC | GAAGGCAACG | ACCGGAACAG |
| UDI 129 | GCGTGGAGAA | TTCTCCACGC | GAATGGAAGC |
| UDI 130 | TACGCCTCCT | AGGAGGCGTA | GTTCTCCATA |
| UDI 131 | AATTCGGTAG | CTACCGAATT | GTCACTATGT |
| UDI 132 | ATTGTGCAAC | GTTGACAAT | TGGTAGAACT |
| UDI 133 | CAACCTTGCG | CGCAAGGTTG | ACGCCTATGG |
| UDI 134 | GCACTGCGTA | TACGCAGTGC | AATCCGTAC |
| UDI 135 | TGCTAGTAGT | ACTACTAGCA | GTTGAGGCTA |
| UDI 136 | AAGTCACGGA | TCCGTGACTT | TATCAACTGG |
| UDI 137 | AGCGATTGAA | TTCAATCGCT | AAGAGGAGAT |
| UDI 138 | CTACCTCTCT | AGAGAGGTAG | GTCTTCTCGG |
| UDI 139 | GACAACGTGC | GACAGTTGTC | GAAGCCACTC |
| UDI 140 | TCCATTGCGG | CCGCAATGGA | GTAGGACACA |

Table continues on next page

Table continued from previous page

| Indices for entry on sample sheet | | | |
|--|--|--|---|
| Unique dual-index number | i5 bases for entry on sample sheet (MiSeq, HiSeq 2000/2500, NovaSeq 6000) | i5 bases for entry on sample sheet (iSeq100, MiniSeq, NextSeq, HiSeqX, HiSeq 3000/4000) | i7 bases for entry on sample sheet (all instruments) |
| UDI 141 | AGCCTCGCAA | TTGCGAGGCT | CTCCTCGTAT |
| UDI 142 | AATACAGGCT | AGCCTGTATT | CCACATGATT |
| UDI 143 | CGGACCGTTA | TAACGGTCCG | AGACGGTTGG |
| UDI 144 | GCGCTTATGC | GCATAAGCGC | CTAGGTTGAC |
| UDI 145 | TTAACACGAG | CTCGTGTTAA | AAGCGTACCA |
| UDI 146 | CGCCTCTAGA | TCTAGAGGCG | TCATGTTGGT |
| UDI 147 | AATCGACCTT | AAGGTCGATT | TTGGAATGGT |
| UDI 148 | CCGCAATAAC | GTTATTGCGG | GTGTATGTTG |
| UDI 149 | GTTCCAACGA | TCGTTGGAAC | TCCTGTCAAC |
| UDI 150 | TGTTAGACCG | CGGTCTAACA | TAATCAGGCA |
| UDI 151 | AACCTCATAG | CTATGAGGTT | GTAGTGGATT |
| UDI 152 | ATGAATCCAC | GTGGATTCAT | AATTGCGCAT |
| UDI 153 | CGGCTTAATT | AATTAAGCCG | GACAATAACG |
| UDI 154 | GAGTGCAGG | CCTGCAACTC | ACAGTTAAGC |
| UDI 155 | TCCACGAACA | TGTTCTGTTGA | AGCCACACTA |
| UDI 156 | TGACGGAGGA | TCCTCCGTC | CAATCGTCTT |
| UDI 157 | AATGAGTACG | CGTACTCATT | AGGAGCTTGT |
| UDI 158 | CGTCTCCGA | TCGGAAGACG | TTGAGCGGAG |
| UDI 159 | GACAGAGATT | AATCTCTGTC | AGTAGCTCTC |
| UDI 160 | TTACGCTAAC | GTTAGCGTAA | CACGCTGTCA |

Table continues on next page

Table continued from previous page

| Indices for entry on sample sheet | | | |
|--|--|--|---|
| Unique dual-index number | i5 bases for entry on sample sheet (MiSeq, HiSeq 2000/2500, NovaSeq 6000) | i5 bases for entry on sample sheet (iSeq100, MiniSeq, NextSeq, HiSeqX, HiSeq 3000/4000) | i7 bases for entry on sample sheet (all instruments) |
| UDI 161 | CTCCTCGAAG | CTTCGAGGAG | AAGACCTCTT |
| UDI 162 | ATACCGCAGA | TCTGCGGTAT | GACCTCTTCT |
| UDI 163 | CCTATCTGAT | ATCAGATAGG | TACTTCCTTG |
| UDI 164 | GATCGGTTAC | GTAACCGATC | TGCGATACGC |
| UDI 165 | TGGTGAGGTG | CACCTACCA | GCAGGCTTAA |
| UDI 166 | AACCGGCGTA | TACGCCGGTT | TAAGCTTGTG |
| UDI 167 | AATACCGATC | GATCGGTATT | ATGGTCCGCT |
| UDI 168 | CGATACTCAA | TTGAGTATCG | ATGTCAGAAG |
| UDI 169 | GTAAGGCGGT | ACCGCCTTAC | GACGAAGGTC |
| UDI 170 | TTC AAGTCTG | CGACCTTGAA | ATCACCGTGA |
| UDI 171 | TATCCGAGTA | TACTCGGATA | GCTACAGTGT |
| UDI 172 | AGCGCGCTTA | TAAGCGCGCT | CGTCGAATAT |
| UDI 173 | CCGGAGACAT | ATGTCTCCGG | CAACCATCGG |
| UDI 174 | GAGATAACTG | CAGTTATCTC | CGGTCCATTC |
| UDI 175 | TTGTAAGCGC | GCGCTTACAA | AGAAGAGCCA |
| UDI 176 | CAAGAGGAGG | CCTCTCTTG | CTATGCAATG |
| UDI 177 | AACCTTAGGA | TCCTAAGGTT | CACTGAACCG |
| UDI 178 | CTGGCAACTC | GAGTTGCCAG | TACTGTGTGA |
| UDI 179 | GAACCTGTG | CAACAAGTTC | GCATTCTGTT |
| UDI 180 | TGTGCAAGAT | ATCTTGACACA | CTCCGCTAAG |

Table continues on next page

Table continued from previous page

| Indices for entry on sample sheet | | | |
|--|--|--|---|
| Unique dual-index number | i5 bases for entry on sample sheet (MiSeq, HiSeq 2000/2500, NovaSeq 6000) | i5 bases for entry on sample sheet (iSeq100, MiniSeq, NextSeq, HiSeqX, HiSeq 3000/4000) | i7 bases for entry on sample sheet (all instruments) |
| UDI 181 | AATCGAGAGA | TCTCTCGATT | TCGCTTGAGA |
| UDI 182 | AGCGTGTCAG | CTGACACGCT | AACTAGCCTT |
| UDI 183 | CTTGGTGATT | AATCACCAAG | TTGCTCAGG |
| UDI 184 | GAAGCAGCAA | TTGCTGCTTC | CTCTACAACA |
| UDI 185 | TTCCGTCGAC | GTCGACGGAA | TGAGTGTGTT |
| UDI 186 | CGAGATGCCA | TGGCATCTCG | TAGTTAGTCG |
| UDI 187 | AAGTTCGTGC | GCACGAACTT | GCCTGATCCT |
| UDI 188 | CGTCCATAAG | CTTATGGACG | CGAGTACAGG |
| UDI 189 | TTGTGGCATA | TATGCCACAA | GCCTAGATTA |
| UDI 190 | AGATCGGAAT | ATTCCGATCT | TCGGCACTGT |
| UDI 191 | CATTCTACTG | CAGTAGAATG | CCGTGCAAGA |
| UDI 192 | ATCGCCGTAG | CTACGGCGAT | CTGGCTGGTT |
| UDI 193 | ATCCTTACAC | GTGTAAGGAT | CGTTAGGATT |
| UDI 194 | CGCAAGGACT | AGTCCTTGCG | TTCCATTACG |
| UDI 195 | GCTGGCGTTA | TAACGCCAGC | TAGCGGTAAC |
| UDI 196 | TACTTAGAGG | CCTCTAAGTA | GTAGCCAGGA |
| UDI 197 | ATGGCGATGC | GCATCGCCAT | AGGATACTCT |
| UDI 198 | CATTGGTGCG | CGCACCAATG | TATCTCCAG |
| UDI 199 | GCGAGATATA | TATATCTCGC | TAAGTCGTTT |
| UDI 200 | TGACTGCTAT | ATAGCAGTCA | TCCGGATTGA |

Table continues on next page

Table continued from previous page

| Indices for entry on sample sheet | | | |
|--|--|--|---|
| Unique dual-index number | i5 bases for entry on sample sheet (MiSeq, HiSeq 2000/2500, NovaSeq 6000) | i5 bases for entry on sample sheet (iSeq100, MiniSeq, NextSeq, HiSeqX, HiSeq 3000/4000) | i7 bases for entry on sample sheet (all instruments) |
| UDI 201 | AACGTCCGCT | AGCGGACGTT | ACGTCTTGTT |
| UDI 202 | CGCACATGTC | GACATGTGCG | ATGAAGTGCG |
| UDI 203 | GCACACCTGA | TCAGGTGTGC | CGATCACTGC |
| UDI 204 | TTGTCCAGAG | CTCTGGACAA | CCTATCGGAA |
| UDI 205 | AGCCTTCTCG | CAGGAAGGCT | CAGAGAGCTT |
| UDI 206 | CCTTAGCCCA | TGGCGTAAGG | GCAACTTGCG |
| UDI 207 | GAATACGTAC | GTACGTATTC | TATGGAGGAC |
| UDI 208 | TTGGCACCGT | ACGGTGCCAA | TGAGATCAGA |
| UDI 209 | ATTAGGTGGC | GCCACCTAAT | TCAGCCTATT |
| UDI 210 | CGATCAAGAA | TTCTTGATCG | GTTGTGAGCG |
| UDI 211 | GCTGTCTTCT | AGAAGACAGC | TCAGTAACAC |
| UDI 212 | TACATGTCTG | CAGACATGTA | AAGGCTCAGA |
| UDI 213 | AACCAGTTGA | TCAACTGGTT | GTGTGGTGGT |
| UDI 214 | CCGTAAGCT | AGCTTACCGG | CCGAGCTTAG |
| UDI 215 | GTTCGAATAG | CTATTGGAAC | ATCACGCTTC |
| UDI 216 | TGTCAGGCTC | GAGCCTGACA | TAGCTATGCA |
| UDI 217 | CAACAGTGTT | AACACTGTTG | TGTTCTCAT |
| UDI 218 | AAGAGAGGAA | TTCTTCTT | CATACCTTCT |
| UDI 219 | CGGTGTAGC | GCTACAACCG | GCCTTCAATG |
| UDI 220 | GCCTGAAGTG | CAC TTCAGGC | CTTGACCAGC |

Table continues on next page

Table continued from previous page

| Indices for entry on sample sheet | | | |
|--|--|--|---|
| Unique dual-index number | i5 bases for entry on sample sheet (MiSeq, HiSeq 2000/2500, NovaSeq 6000) | i5 bases for entry on sample sheet (iSeq100, MiniSeq, NextSeq, HiSeqX, HiSeq 3000/4000) | i7 bases for entry on sample sheet (all instruments) |
| UDI 221 | TTACGACACT | AGTGTGCTAA | CTACACACAA |
| UDI 222 | CGCCTAGATC | GATCTAGGCG | TAGGCTGAAT |
| UDI 223 | AATCTGGATG | CATCCAGATT | TCGGAGTCTT |
| UDI 224 | CGACGGTACA | TGTACCGTCG | AACATCGCGG |
| UDI 225 | GTAGTATTGC | GCAATACTAC | GTTGTCTTAC |
| UDI 226 | TCCAGCGGAT | ATCCGCTGGA | GTGGCAACTA |
| UDI 227 | CAACCACCTC | GAGGTGGTTG | GAGCAGGCAT |
| UDI 228 | AGCTTAGGCG | CGCCTAAGCT | AACGGCACCT |
| UDI 229 | CCGTTCTCTT | AAGGAACCGG | AGTAACCTTG |
| UDI 230 | GACATTGAAC | GTTCAATGTC | TCTCATAAGC |
| UDI 231 | TTAGAGGCGA | TCGCCTCTAA | TGCTTGCCAA |
| UDI 232 | CAAGCCGAAC | GTTGCGCTTG | CGGTTCTGT |
| UDI 233 | AGGAGAACGG | CCGTTCTCTT | CCAAGTAGAT |
| UDI 234 | CCTGTTAGAC | GTCTAACAGG | AAGGTTGGCG |
| UDI 235 | GTTCTACGTT | AACGTAGAAC | TGCTCTGGTC |
| UDI 236 | TAAGTCCACA | TGTGGACTTA | ACTGTAACGA |
| UDI 237 | CAAGAACCAT | ATGGTTCTTG | GATTCCAGGT |
| UDI 238 | AGTTGATGAC | GTCATCAACT | TTCACCAGAT |
| UDI 239 | CCTACTCTTG | CAAGAGTAGG | ACTTCCAAGG |
| UDI 240 | GAACAATCCA | TGGATTGTTT | CCGAATATTC |

Table continues on next page

Table continued from previous page

| Indices for entry on sample sheet | | | |
|--|--|--|---|
| Unique dual-index number | i5 bases for entry on sample sheet (MiSeq, HiSeq 2000/2500, NovaSeq 6000) | i5 bases for entry on sample sheet (iSeq100, MiniSeq, NextSeq, HiSeqX, HiSeq 3000/4000) | i7 bases for entry on sample sheet (all instruments) |
| UDI 241 | TTCTGTTGGT | ACCAACAGAA | CTCTATCCA |
| UDI 242 | CATCGTCAGG | CCTGACGATG | TCACAGCGGT |
| UDI 243 | ATGCATGAAG | CTTCATGCAT | CCTCTGTCGT |
| UDI 244 | CGTGAATCGC | GCGATTCACG | TCTGTTCTCG |
| UDI 245 | GAGCAGCCTT | AAGGCTGCTC | GATACTTCAC |
| UDI 246 | TCGATTACCA | TGGAATCGA | AGTGCTGATA |
| UDI 247 | CAGTCCAATT | AATTGGACTG | ATCCTTCGGT |
| UDI 248 | AGAGGCTTGG | CCAAGCCTCT | GACAACGATT |
| UDI 249 | CAGGCTCTCA | TGAGAGCCTG | GAACCGGTAG |
| UDI 250 | GTTGCCTCTC | GAGAGCGAAC | AGCAATGAGC |
| UDI 251 | TCGGACTAAT | ATTAGTCCGA | CAAGACTCCA |
| UDI 252 | CGAGATCTTC | GAAGATCTCG | ACCGTGTAGG |
| UDI 253 | ATAACCGGAC | GTCCGGTTAT | AGGCACAGGT |
| UDI 254 | CGTGTAGTTA | TAACTACACG | CGACAGATCG |
| UDI 255 | GAACATAGGT | ACCTATGTTT | ACGCGACAAC |
| UDI 256 | TCTAACATCG | CGATGTTAGA | ACTTGCCTTA |
| UDI 257 | AACGGTGGCA | TGCCACCGTT | CACCACTCAT |
| UDI 258 | AGGACGGTGT | ACACCGTCTT | CTTCGTAACT |
| UDI 259 | CTGTGACCTG | CAGGTCACAG | CAGTATTCGG |
| UDI 260 | GCTGTAACAA | TTGTTACAGC | CAGTCTGGAC |

Table continues on next page

Table continued from previous page

| Indices for entry on sample sheet | | | |
|--|--|--|---|
| Unique dual-index number | i5 bases for entry on sample sheet (MiSeq, HiSeq 2000/2500, NovaSeq 6000) | i5 bases for entry on sample sheet (iSeq100, MiniSeq, NextSeq, HiSeqX, HiSeq 3000/4000) | i7 bases for entry on sample sheet (all instruments) |
| UDI 261 | TACGGACGTC | GACGTCCGTA | TACCGTTCTA |
| UDI 262 | CCTAAGGAGC | GCTCCTTAGG | GTGTCCACAG |
| UDI 263 | ATAAGGCCAG | CTGGCCTTAT | TTACGACTGT |
| UDI 264 | CTCATCTGTA | TACAGATGAG | GACGCGAATG |
| UDI 265 | GAAGGCATCT | AGATGCCTTC | CAACGTACGC |
| UDI 266 | TCTCTACTGC | GCAGTAGAGA | AGCTCAGGAA |
| UDI 267 | AACCGAACAA | TTGTTCGGTT | GATAGGCGGT |
| UDI 268 | ATCTCGCCAC | GTGGCGAGAT | AGTAGGAAGT |
| UDI 269 | CCATGCAACG | CGTTGCATGG | CATGTTGTAG |
| UDI 270 | GAATGGTGTA | TACACCATTC | CACATTCTTC |
| UDI 271 | TATATGCCGT | ACGGCATATA | GCAGCTCGTA |
| UDI 272 | CTCGATAGAT | ATCTATCGAG | GTTCAGACGG |
| UDI 273 | AACACAAGAG | CTCTTGTT | TCCTGGAAGT |
| UDI 274 | CGCAATCGGT | ACCGATTGCG | GCATTGTTAG |
| UDI 275 | GTTGCGTAGA | TCTACGCAAC | GACCTACAGC |
| UDI 276 | TAGAGTGATC | GATCACTCTA | CACCGACGTA |
| UDI 277 | AAGACGCAGC | GCTGCGTCTT | CTCTCACCTT |
| UDI 278 | AACTTCTCGA | TCGAGAAGTT | CTCGTTCAAT |
| UDI 279 | CGCAACTGAG | CTCAGTTGCG | TGGTGGCAAG |
| UDI 280 | GCTCCGCAAT | ATTGCGGAGC | GATTGCTTGA |

Table continues on next page

Table continued from previous page

| Indices for entry on sample sheet | | | |
|--|--|--|---|
| Unique dual-index number | i5 bases for entry on sample sheet (MiSeq, HiSeq 2000/2500, NovaSeq 6000) | i5 bases for entry on sample sheet (iSeq100, MiniSeq, NextSeq, HiSeqX, HiSeq 3000/4000) | i7 bases for entry on sample sheet (all instruments) |
| UDI 281 | GTAACCTCCG | CGGAAGTTAC | CCGTTAAGGT |
| UDI 282 | CTCACGACTA | TAGTCGTGAG | TGCTGAGAGG |
| UDI 283 | AACCAACGGC | GCCGTGGTT | TTGTCACTTG |
| UDI 284 | CCTGCCTGTA | TACAGGCAGG | GCTGTTATGT |
| UDI 285 | TACGCTGCAG | CTGCAGCGTA | GCAGCAGTTG |
| UDI 286 | AATGTTGCGA | TCGCAACATT | GCAGATCAAT |
| UDI 287 | CGACGTTCTG | CAGAACGTCG | TGGTTCACGG |
| UDI 288 | AATAGGACAC | GTGTCCTATT | TCGACCGCAT |
| UDI 289 | ATGTGCCTCA | TGAGGCACAT | TAACCTAGGT |
| UDI 290 | CGACTCCGTT | AACGGAGTCG | AACTCATGCG |
| UDI 291 | GCTGTTGTGG | CCACAACAGC | CCGGATGAAC |
| UDI 292 | TACCAATCAC | GTGATTGGTA | CGTTGCCGTA |
| UDI 293 | ATGCTTACG | CGTAAGACAT | GCTCTACGGT |
| UDI 294 | CGCAACAATA | TATTGTTGCG | TGCATTGGCG |
| UDI 295 | GAACGAAGAC | GTCTTCGTTT | CGATTGTGAC |
| UDI 296 | TCGAGGACGT | ACGTCCTCGA | GACTGCACTA |
| UDI 297 | ATTATGAGCG | CGCTCATAAT | GTTAACTGCT |
| UDI 298 | CGCGTTATAA | TTATAACGCG | TCGGACCTTG |
| UDI 299 | GCGTGTCATGT | ACATGCACGC | TGCAGCAAGC |
| UDI 300 | TAAGCGGCTC | GAGCCGCTTA | CACATGCGAA |

Table continues on next page

Table continued from previous page

| Indices for entry on sample sheet | | | |
|--|--|--|---|
| Unique dual-index number | i5 bases for entry on sample sheet (MiSeq, HiSeq 2000/2500, NovaSeq 6000) | i5 bases for entry on sample sheet (iSeq100, MiniSeq, NextSeq, HiSeqX, HiSeq 3000/4000) | i7 bases for entry on sample sheet (all instruments) |
| UDI 301 | AACATGGAGA | TCTCCATGTT | CAGACGTAAT |
| UDI 302 | CCGAGTCTCT | AGAGACTCGG | ATTCGGTACG |
| UDI 303 | GTA CTCTAC | GTAGAAGTAC | TTAGCACGGC |
| UDI 304 | TGTTACATG | CATGTGAACA | GAGGATAGTA |
| UDI 305 | AAGTAACGC | GCGTTACCTT | AACTGTGGTT |
| UDI 306 | CCGCCTACT | AGTAAGGCGG | ATTACCTCGG |
| UDI 307 | GTTGAGGCAG | CTGCCTCAAC | CGTGTATAC |
| UDI 308 | TGCGGACCTA | TAGGTCGCCA | CTTGCTCACA |
| UDI 309 | AGAAGCGACA | TGTCGCTTCT | CAACACCTGT |
| UDI 310 | CAGGATAATC | GATTATCCTG | CAATTGCTCG |
| UDI 311 | GCTCCTACAG | CTGTAGGAGC | CATAGACAAC |
| UDI 312 | TCAACAGGT | ACCTGTTGAA | TTGGTGCTA |
| UDI 313 | CCTCGTCCAT | ATGGACGAGG | TATGTCCTGT |
| UDI 314 | AGCGTTGGTT | AACCAACGCT | GCCAAATCGT |
| UDI 315 | CATTGAACA | TGTTGAATG | TAGGCGATCG |
| UDI 316 | GCTTACCGAC | GTCGGTAAGC | ATGAGTGTAC |
| UDI 317 | TTAGCTTAGG | CCTAAGCTAA | CCGAAGGATA |
| UDI 318 | CCGACACACA | TGTGTGTCGG | AGTCCACTGT |
| UDI 319 | ATTCGCTGAT | ATCAGCGAAT | GCGGCTAATT |
| UDI 320 | CCAAGAGGCA | TGCCTCTGG | TCTAACTCAG |

Table continues on next page

Table continued from previous page

| Indices for entry on sample sheet | | | |
|--|--|--|---|
| Unique dual-index number | i5 bases for entry on sample sheet (MiSeq, HiSeq 2000/2500, NovaSeq 6000) | i5 bases for entry on sample sheet (iSeq100, MiniSeq, NextSeq, HiSeqX, HiSeq 3000/4000) | i7 bases for entry on sample sheet (all instruments) |
| UDI 321 | GACGCAGTTC | GAAC TGCGTC | CAAGCTGAGC |
| UDI 322 | TGGAACTCGG | CCGAGTCCA | CCAGAGCACA |
| UDI 323 | CCACACCAAT | ATTGGTGTGG | TGTACAAGGT |
| UDI 324 | AGTTCTCGGC | GCCGAGAACT | TAGAATGCCT |
| UDI 325 | CTTGACGACG | CGTCGTCAAG | TGCTTACTG |
| UDI 326 | GAGGTCGCTA | TAGCGACCTC | ATGACTAAGC |
| UDI 327 | TCAGTAGCAT | ATGCTACTGA | ATGTAGGCAA |
| UDI 328 | CTAACGTGGA | TCCACGTTAG | GCGAAGAGGT |
| UDI 329 | ATGCCAACCG | CGGTTGGCAT | CGGTGGTCT |
| UDI 330 | CGGTCGATTC | GAATCGACCG | CTGTCGTGG |
| UDI 331 | GAAGTACAGT | ACTGTACTTC | TGATCGACAC |
| UDI 332 | TCTGCAGTAA | TTACTGCAGA | CCACCAGCTA |
| UDI 333 | CTATCCTAGC | GCTAGGATAG | CACGGTTCGT |
| UDI 334 | AACACTCCTT | AAGGAGTGTT | AGTGAGAGCT |
| UDI 335 | CCGAACCTAA | TTAGGTTCCGG | TTGCATGCGG |
| UDI 336 | GTCTAGTCGC | GCGACTAGAC | TATACGTGTC |
| UDI 337 | TGGATGTACG | CGTACATCCA | TGACGCGTTA |
| UDI 338 | CTACCAGCGT | ACGCTGGTAG | TACAGAACGT |
| UDI 339 | AAGGATTACAG | CTGAATCCTT | CTTGTCAGGT |
| UDI 340 | CGAGGTGTGT | ACACACCTCG | ATCCACAGCG |

Table continues on next page

Table continued from previous page

| Indices for entry on sample sheet | | | |
|--|--|--|---|
| Unique dual-index number | i5 bases for entry on sample sheet (MiSeq, HiSeq 2000/2500, NovaSeq 6000) | i5 bases for entry on sample sheet (iSeq100, MiniSeq, NextSeq, HiSeqX, HiSeq 3000/4000) | i7 bases for entry on sample sheet (all instruments) |
| UDI 341 | GTAGACGCTC | GAGCGTCTAC | CCTATCCATC |
| UDI 342 | TCGTCCGTCA | TGACGGACGA | ACCGCGAGTA |
| UDI 343 | CCGTGATAGG | CCTATCACGG | AAGTTCTGGT |
| UDI 344 | AGGATGACCT | AGGTCATCCT | ACAGGTATCG |
| UDI 345 | CCTCGAGTAC | GTAICTGAGG | ATGACGGATT |
| UDI 346 | GTCACTGAGG | CCTCAGTGAC | GTCTGAGTAG |
| UDI 347 | TACGGTTAGA | TCTAACCGTA | TGCCAGATGT |
| UDI 348 | CAACGAGAAT | ATTCTCGTTG | GCTAAGCATT |
| UDI 349 | AATACACCGG | CCGGTGTATT | ACAGCATGGT |
| UDI 350 | CCGATCCATC | GATGGATCGG | ATAGAGACCG |
| UDI 351 | GAATCTCGCT | AGCGAGATTC | ATATCGCGTA |
| UDI 352 | TGACCGGCAA | TTGCCGGTCA | TTAAGGAGGT |
| UDI 353 | CATGATAGCA | TGCTATCATG | CTGTGCGACT |
| UDI 354 | AACAGCTTCG | CGAAGCTGTT | TCCGTATGCT |
| UDI 355 | CTAGTGCTTA | TAAGCACTAG | CCATCGATGT |
| UDI 356 | TGTGATACGT | ACGTATCACA | GTGAGCCGTT |
| UDI 357 | ATGAGCGTAT | ATACGCTCAT | TGCCGTTAAT |
| UDI 358 | CTAGATATGG | CCATATCTAG | CGGATGTGGT |
| UDI 359 | CGCTATGCTG | CAGCATAGCG | TCGCGTGTG |
| UDI 360 | TACTACGTGA | TCACGTAGTA | CCGCGATCAT |

Table continues on next page

Table continued from previous page

| Indices for entry on sample sheet | | | |
|--|--|--|---|
| Unique dual-index number | i5 bases for entry on sample sheet (MiSeq, HiSeq 2000/2500, NovaSeq 6000) | i5 bases for entry on sample sheet (iSeq100, MiniSeq, NextSeq, HiSeqX, HiSeq 3000/4000) | i7 bases for entry on sample sheet (all instruments) |
| UDI 361 | ATGTGGAGGT | ACCTCCACAT | CGCGTTATCG |
| UDI 362 | CCATGGCTCA | TGAGCCATGG | GTAGCCTCCT |
| UDI 363 | CCAATCACGC | GCGTGATTGG | ACTAGACACT |
| UDI 364 | TTAGATCCAG | CTGGATCTAA | CGATTCTGTG |
| UDI 365 | AGGAATATCG | CGATATTCTT | GAAGAGATGT |
| UDI 366 | CCTCCTATGT | ACATAGGAGG | AGATCCGACG |
| UDI 367 | TAGAGACACG | CGTGTCTCTA | CCAGGACATT |
| UDI 368 | CCAGCTCAGT | ACTGAGCTGG | ACGTGGCATT |
| UDI 369 | ATGGCTCATA | TATGAGCCAT | AAGCAGGACG |
| UDI 370 | CGGAGTGAAG | CTTCACTCCG | ACGAGTCGGT |
| UDI 371 | TACCTATGGT | ACCATAGGTA | AGTGTACGCG |
| UDI 372 | ATGAGACAGT | ACTGTCTCAT | ACCGACCATT |
| UDI 373 | CTAAGAGTTG | CAACTCTTAG | TTGCTAACGT |
| UDI 374 | TAACCGTATG | CATACGGTTA | CTTGATACTG |
| UDI 375 | AGAGTCCATG | CATGGACTCT | CTGATAAAGT |
| UDI 376 | CTAGACCGCA | TGCGGTCTAG | ATAGCTTACG |
| UDI 377 | TATGGCTTGT | ACAAGCCATA | GTCCATGAGT |
| UDI 378 | CGTTGTCTCT | AGGAACAACG | ACTCCAGTCG |
| UDI 379 | CCGACATTAG | CTAATGTCGG | TCTCAGCACG |
| UDI 380 | TGTGAAGGCA | TGCCTTCACA | ATCGTGATGT |
| UDI 381 | AGCATCGTCT | AGACGATGCT | ACGCAATCCG |
| UDI 382 | CCGACTAGGA | TCCTAGTCGG | GAGATCGGCT |
| UDI 383 | AACATTACCG | CGGTAATGTT | CTACGTCTCG |
| UDI 384 | CCTAATTCGT | ACGAATTAGG | CTCAGGCTGT |

Combinatorial Dual-Index Y-Adapters

The layout of the 24-plex and 96-plex single-use CDI adapter plate is shown in Figure 10 and Figure 11. The index motives used in the QIAseq Combinatorial Dual-Index Kits are listed in Table 6. To make sequencing preparation more convenient, you can download Illumina-compatible sample sheets for different sequencing instruments at www.qiagen.com.

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|
| A | 501/ 701 | 501/ 702 | 501/ 703 | 501/ 704 | 501/ 705 | 501/ 706 | 501/ 707 | 501/ 708 | 501/ 709 | 501/ 710 | 501/ 711 | 501/ 712 |
| B | 502/ 701 | 502/ 702 | 502/ 703 | 502/ 704 | 502/ 705 | 502/ 706 | 502/ 707 | 502/ 708 | 502/ 709 | 502/ 710 | 502/ 711 | 502/ 712 |
| C | 503/ 701 | 503/ 702 | 503/ 703 | 503/ 704 | 503/ 705 | 503/ 706 | 503/ 707 | 503/ 708 | 503/ 709 | 503/ 710 | 503/ 711 | 503/ 712 |
| D | 504/ 701 | 504/ 702 | 504/ 703 | 504/ 704 | 504/ 705 | 504/ 706 | 504/ 707 | 504/ 708 | 504/ 709 | 504/ 710 | 504/ 711 | 504/ 712 |
| E | 505/ 701 | 505/ 702 | 505/ 703 | 505/ 704 | 505/ 705 | 505/ 706 | 505/ 707 | 505/ 708 | 505/ 709 | 505/ 710 | 505/ 711 | 505/ 712 |
| F | 506/ 701 | 506/ 702 | 506/ 703 | 506/ 704 | 506/ 705 | 506/ 706 | 506/ 707 | 506/ 708 | 506/ 709 | 506/ 710 | 506/ 711 | 506/ 712 |
| G | 507/ 701 | 507/ 702 | 507/ 703 | 507/ 704 | 507/ 705 | 507/ 706 | 507/ 707 | 507/ 708 | 507/ 709 | 507/ 710 | 507/ 711 | 507/ 712 |
| H | 508/ 701 | 508/ 702 | 508/ 703 | 508/ 704 | 508/ 705 | 508/ 706 | 508/ 707 | 508/ 708 | 508/ 709 | 508/ 710 | 508/ 711 | 508/ 712 |

Figure 10. QIAseq CDI Y-Adapter Plate (96) layout (CDI 1–96).

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|----------|---------|---------|---------|-------|-------|-------|-------|-------|-------|-------|-------|-------|
| A | 501/701 | 501/702 | 501/703 | empty | empty | empty | empty | empty | empty | empty | empty | empty |
| B | 502/701 | 502/702 | 502/703 | empty | empty | empty | empty | empty | empty | empty | empty | empty |
| C | 503/701 | 503/702 | 503/703 | empty | empty | empty | empty | empty | empty | empty | empty | empty |
| D | 504/701 | 504/702 | 504/703 | empty | empty | empty | empty | empty | empty | empty | empty | empty |
| E | 505/701 | 505/702 | 505/703 | empty | empty | empty | empty | empty | empty | empty | empty | empty |
| F | 506/701 | 506/702 | 506/703 | empty | empty | empty | empty | empty | empty | empty | empty | empty |
| G | 507/701 | 507/702 | 507/703 | empty | empty | empty | empty | empty | empty | empty | empty | empty |
| H | 508/701 | 508/702 | 508/703 | empty | empty | empty | empty | empty | empty | empty | empty | empty |

Figure 11. QIAseq CDI Y-Adapter Plate (24) layout (CDI 1–24).

Table 6. CDI motives used in the QIAseq CDI Y-Adapter Kits (24 and 96)

Note: Sequencing on the MiniSeq, NextSeq, HiSeqX, and HiSeq 3000/4000 systems follow a different dual-indexing workflow than other Illumina systems, which requires the reverse complement of the i5 index adapter sequence.

| Indices for entry on sample sheet | | | | |
|-----------------------------------|---|---|-------------------|------------------------------------|
| D50X barcode name | i5 bases for entry on sample sheet (MiSeq, HiSeq 2000/2500, NovaSeq 6000) | i5 bases for entry on sample sheet (iSeq100, MiniSeq, NextSeq, HiSeqX, HiSeq 3000/4000) | D50X barcode name | i7 bases for entry on sample sheet |
| D501 | TATAGCCT | AGGCTATA | D701 | ATTACTCG |
| D502 | ATAGAGGC | GCCTCTAT | D702 | TCCGGAGA |
| D503 | CCTATCCT | AGGATAGG | D703 | CGCTCATT |
| D504 | GGCTCTGA | TCAGAGCC | D704 | GAGATTCC |
| D505 | AGGCGAAG | CTTCGCCT | D705 | ATTCAGAA |
| D506 | TAATCTTA | TAAGATTA | D706 | GAATTCGT |
| D507 | CAGGACGT | ACGTCCTG | D707 | CTGAAGCT |
| D508 | GTACTGAC | GTCAGTAC | D708 | TAATGCGC |
| | | | D709 | CGGCTATG |
| | | | D710 | TCCGCGAA |
| | | | D711 | TCTCGCGC |
| | | | D712 | AGCGATAG |

Appendix E: Adapter Indices for GeneRead Adapter I Sets A and B

The index sequences used in the GeneRead Adapter I Set A 12-plex are listed in Table 7. Indices 1–12 correspond to the respective Illumina adapter indices.

Table 7. Adapter indices

| Adapter name | Indices |
|-----------------------|----------------|
| Adapter Bc1 Illumina | ATCACG |
| Adapter Bc2 Illumina | CGATGT |
| Adapter Bc3 Illumina | TTAGGC |
| Adapter Bc4 Illumina | TGACCA |
| Adapter Bc5 Illumina | ACAGTG |
| Adapter Bc6 Illumina | GCCAAT |
| Adapter Bc7 Illumina | CAGATC |
| Adapter Bc8 Illumina | ACTTGA |
| Adapter Bc9 Illumina | GATCAG |
| Adapter Bc10 Illumina | TAGCTT |
| Adapter Bc11 Illumina | GGCTAC |
| Adapter Bc12 Illumina | CTTGTA |

The index sequences used in the GeneRead Adapter I Set B 12-plex are listed in Table 8. Indices 13–16, 18–23, 25, and 27 correspond to their respective Illumina adapter indices.

Table 8. Adapter indices

| Adapter name | Indices |
|-----------------------|----------------|
| Adapter Bc13 Illumina | AGTCAA |
| Adapter Bc14 Illumina | AGTTCC |
| Adapter Bc15 Illumina | ATGTCA |
| Adapter Bc16 Illumina | CCGTCC |
| Adapter Bc18 Illumina | GTCCGC |
| Adapter Bc19 Illumina | GTGAAA |
| Adapter Bc20 Illumina | GTGGCC |
| Adapter Bc21 Illumina | GTTTCG |
| Adapter Bc22 Illumina | CGTACG |
| Adapter Bc23 Illumina | GAGTGG |
| Adapter Bc25 Illumina | ACTGAT |
| Adapter Bc27 Illumina | ATTCCT |

Appendix F: Data Analysis

After sequencing, multiplex data can be analyzed using QIAGEN's CLC Genomics or Biomedical Workbench or the cloud-based GeneRead DNaseq Sequence Variant Analysis Software. Our data analysis software will perform quality control, read trimming (removing primer sequences), mapping to a reference genome, and variant identification. Please refer to the corresponding documentation for data analysis.

Ordering Information

| Product | Contents | Cat. no. |
|---|---|-----------------|
| QIAseq 1-Step Amplicon Library Kit (12) | Reagents for Illumina Amplicon Seq library preparation, and library amplification | 180412 |
| QIAseq 1-Step Amplicon Lib UDI-A Kit (96) | Reagents for Illumina Amplicon Seq library preparation and library amplification, including 96-plex Illumina Adapters | 180419 |
| QIAseq 1-Step Amplicon Lib UDI-B Kit (96) | Reagents for Illumina Amplicon Seq library preparation and library amplification, including 96-plex Illumina Adapters | 180420 |
| QIAseq 1-Step Amplicon Lib UDI-C Kit (96) | Reagents for Illumina Amplicon Seq library preparation and library amplification, including 96-plex Illumina Adapters | 180421 |
| QIAseq 1-Step Amplicon Lib UDI-D Kit (96) | Reagents for Illumina Amplicon Seq library preparation and library amplification, including 96-plex Illumina Adapters | 180425 |
| QIAseq 1-Step Amplicon Lib CDI Kit (96) | Reagents for Illumina Amplicon Seq library preparation and library amplification, including 96-plex Illumina Adapters | 180423 |
| GeneRead Adapter I Set A 12-plex (144) | 12 barcoded adapters for ligation to DNA library, for use with Illumina instruments | 180985 |
| GeneRead Adapter I Set B 12-plex (144) | 12 barcoded adapters for ligation to DNA library, for use with Illumina instruments | 180986 |
| Related products | | |
| GeneRead DNaseq Library Quant Array | Reagents for NGS sample library quantification following targeted enrichment | 180601 |

| Product | Contents | Cat. no. |
|--|--|----------|
| GeneRead qPCR SYBR® Green Mastermix | Master mix for use with the GeneRead Library Quant Arrays and Kit | Varies* |
| QIAquick® PCR Purification Kit (50) | QIAquick Spin Columns, buffers, collection tubes (2 ml) for purification of PCR products <150 bp | 28104 |
| QIAseq Y-Adapter Kits for Illumina | | |
| QIAseq CDI Y-Adapter Kit (24) | Combinatorial Dual-Index Adapters for Illumina | 180301 |
| QIAseq CDI Y-Adapter Kit (96) | Combinatorial Dual-Index Adapters for Illumina | 180303 |
| QIAseq UDI Y-Adapter Kit (24) | Unique Dual-Index Adapters for Illumina (1–24) | 180310 |
| QIAseq UDI Y-Adapter Kit A (96) | Unique Dual-Index Adapters for Illumina (1–96) | 180312 |
| QIAseq UDI Y-Adapter Kit B (96) | Unique Dual-Index Adapters for Illumina (97–192) | 180314 |
| QIAseq UDI Y-Adapter Kit C (96) | Unique Dual-Index Adapters for Illumina (193–288) | 180316 |
| QIAseq UDI Y-Adapter Kit D (96) | Unique Dual-Index Adapters for Illumina (289–384) | 180318 |
| QIAGEN panels for target enrichment | | |
| GeneRead Custom DNaseq Gene Panels | Pools containing primer sets for targeted enrichment of a customized panel of genes or genomic regions | 181902 |
| GeneRead DNaseq Panel PCR Kit V2 (12) | PCR chemistry for use with the GeneRead DNaseq Panel V2 System | 181940 |
| GeneRead DNaseq Panel PCR Kit V2 (96) | PCR chemistry for use with the GeneRead DNaseq Panel V2 System | 181942 |

* See www.qiagen.com.

| Product | Contents | Cat. no. |
|---|---|-----------------|
| GeneRead DNA QuantiMIZE Array Kit | qPCR arrays for optimizing the amount of input DNA and PCR cycling conditions for targeted enrichment of FFPE DNA | 180642 |
| QIAGEN PCR reagents for target enrichment | | |
| QIAGEN Multiplex PCR Kit (100) | For 100 x 50 µl multiplex PCR reactions: 2x QIAGEN Multiplex PCR Master Mix (providing a final concentration of 3 mM MgCl ₂ , 3 x 0.85 ml), 5x Q-Solution® (1 x 2.0 ml), RNase-free water (2 x 1.7 ml) | 206143 |
| QIAGEN Multiplex PCR <i>Plus</i> Kit (100) | For 100 x 50 µl multiplex PCR reactions: 2x Multiplex PCR Master Mix (3 x 0.85 ml), 5x Q-Solution (1 x 2 ml), RNase-free water (2 x 1.9 ml), 10x CoralLoad® Dye (1 x 1.2 ml) | 206152 |
| REPLI-g® whole genome amplification for sensitive applications | | |
| REPLI-g Single Cell Kit (96) | Polymerase, buffers, and reagents for whole genome amplification from limited input materials or single cells | 150345 |
| REPLI-g FFPE Kit (100) | Polymerase, buffers, and reagents for whole genome amplification from FFPE samples | 150245 |
| QIAGEN kits for genomic DNA isolation and purification | | |
| QIAamp DNA Mini Kit (50) | For 50 DNA preps: 50 QIAamp Mini Spin Columns, QIAGEN Proteinase K, collection tubes (2 ml), reagents, and buffers | 51304 |
| GeneRead DNA FFPE Kit | QIAamp MinElute columns, Proteinase K, UNG, collection tubes (2 ml), buffers, deparaffinization solution, RNase A | 180134 |

| Product | Contents | Cat. no. |
|--------------------------------|---|----------|
| MagAttract HMW DNA Kit (48) | For 48 DNA preps: MagAttract Suspension G, Buffer ATL, Buffer AL, Buffer MB, Buffer MW1, Buffer PE, Proteinase K, RNase A, Buffer AE, Nuclease-free Water | 67563 |
| QIAamp DNA Microbiome Kit (50) | For 50 DNA preps: 50 QIAamp UCP Mini Columns, 50 Pathogen Lysis Tubes L, buffers, reagents, collection tubes (2 ml) | 51704 |

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Document Revision History

| Date | Changes |
|---------|---|
| 08/2021 | Removed the QIAseq 1-Step Amplicon Library Kit (96) (cat. no. 180415) product for discontinuation. Added information about QIAseq Y-adapter kits and NGS adapter and index technologies. Removed information exclusive to formerly recommended adapters. Updated information related to Illumina instruments. |

Notes

Limited License Agreement for the QIAseq 1-Step Amplicon Library Kit

Use of this product signifies the agreement of any purchaser or user of the product to the following terms:

1. The product may be used solely in accordance with the protocols provided with the product and this handbook and for use with components contained in the kit only. QIAGEN grants no license under any of its intellectual property to use or incorporate the enclosed components of this kit with any components not included within this kit except as described in the protocols provided with the product, this handbook, and additional protocols available at www.qiagen.com. Some of these additional protocols have been provided by QIAGEN users for QIAGEN users. These protocols have not been thoroughly tested or optimized by QIAGEN. QIAGEN neither guarantees them nor warrants that they do not infringe the rights of third-parties.
2. Other than expressly stated licenses, QIAGEN makes no warranty that this kit and/or its use(s) do not infringe the rights of third-parties.
3. This kit and its components are licensed for one-time use and may not be reused, refurbished, or resold.
4. QIAGEN specifically disclaims any other licenses, expressed or implied other than those expressly stated.
5. The purchaser and user of the kit agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. QIAGEN may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to the kit and/or its components.

For updated license terms, see www.qiagen.com.

Oligonucleotide sequences © 2007–2013 Illumina, Inc. All rights reserved.

Trademarks: QIAGEN®, Sample to Insight®, QIAamp®, QIAquick®, QIAseq®, QIAxcel®, QIAxper®, CoralLoad®, GeneRead™, QIAGEN GeneRead®, MagAttract®, MinElute®, Q-Solution®, REPLI-g®, (QIAGEN Group); Agencourt®, AMPure® (Beckman Coulter, Inc.); Agilent® (Agilent Technologies, Inc.); HiSeq®, HiSeq X™, Illumina®, MiniSeq®, MiSeq®, NextSeq®, NovaSeq™ (Illumina, Inc.); NanoDrop®, Qubit®, SYBR® (Thermo Fisher Scientific or its subsidiaries). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

08/2021 HB-2028-002 © 2021 QIAGEN.com, all rights reserved.

Ordering www.qiagen.com/shop | Technical Support support.qiagen.com | Website www.qiagen.com