

EZ1&2™ Virus Mini Kit v2.0

For use with EZ1 instrument

For usage of EZ1&2 Virus Mini Kit v2.0 with EZ2 instruments, please refer to the handbook (www.qiagen.com/HB-2978) and quick-start protocol (www.qiagen.com/HB-2977).

The EZ1&2 Virus Mini Kit v2.0 (cat. no. 955134) can be stored at room temperature (15–25°C) until the indicated expiration date. Do not freeze the reagent cartridges.

Further information

- *EZ1&2 Virus Handbook*: www.qiagen.com/HB-0113
- Safety Data Sheets: www.qiagen.com/safety
- Technical assistance: support.qiagen.com

Notes before starting

- The EZ1&2 Virus Mini Kit v2.0 requires use of the EZ1® Advanced XL Virus Card v2.0 with the EZ1 Advanced XL, or use of the EZ1 Advanced Virus Card v2.0 with the EZ1 Advanced, or use of the EZ1 Virus Card v2.0 with the BioRobot® EZ1.
- The buffer in well 1 of the reagent cartridge may form a precipitate upon storage. If necessary, redissolve by warming at 37°C and then place at room temperature.
- EZ1® instruments should only be switched on after an EZ1 Card is completely inserted; otherwise, essential instrument data could be lost. EZ1 Cards should not be exchanged while the instrument is switched on.
- Before loading reagent cartridges into the EZ1 instrument, invert the cartridges three times to mix the magnetic particles then tap to deposit the reagents at the bottom of their wells. Check that the magnetic particles are completely resuspended.

1. If necessary, pretreat samples according to the appropriate procedure, as described in the *EZ1&2 Virus Handbook*. See Table 1 for information on which sample types require pretreatment as well as accessories required for pretreatment.
2. Carry out nucleic acid purification using an EZ1 instrument with the appropriate EZ1 Card, according to the *EZ1&2 Virus Handbook*.

Note: Specific guidance regarding the amount of starting material for different sample types is provided in the *EZ1&2 Virus Handbook*.

Table 1. Sample materials covered by the EZ1&2 Virus Mini Kit v2.0

Sample type	Pretreatment needed	Pretreatment accessories	Ordering information
Serum, plasma, cerebrospinal fluid, transport media	No	–	–
Urine, whole blood, dried swabs*	Yes*	Buffer ATL	Cat. no. 939011
Stool*	Yes*	Buffer ASL	Cat. no. 19082
Respiratory samples*	Yes*	Sputasol	www.thermofisher.com

* Information about specific pretreatment procedures is provided in the *EZ1&2 Virus Handbook*.



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN®, Sample to Insight®, BioRobot®, EZ1®, EZ1&2™ (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

1125908 11/2021 HB-0787-003 © 2021 QIAGEN, all rights reserved.