

A NEW IDH1/2 PCR ASSAY FOR ONE-STEP DETECTION OF 12 IDH1 AND IDH2 MUTATIONS IN GLIOMA

H. Girardi ¹, F. Monville ¹, C. Poggionovo ¹, A. Cateau ¹, M. Giry ², J. Voss ³, R. Jenkins ³, K. Mokhtari ², H. Peyro-Saint-Paul ¹, C. Giannini ³
¹ Qiagen Marseille, France - ² Pitié -Salpêtrière, Paris, France - ³ Mayo Clinic Rochester, MN, USA

BACKGROUND

Isocitrate dehydrogenase (IDH) mutational status is a strong diagnostic and prognostic marker in glioma which will probably be introduced in the next WHO classification system.

Current IDH mutations screening is performed with an IHC assay specific for IDH1 R132H, the most common mutation. Sequencing is recommended as a second-step test for IHC-negative or –equivocal cases. However, sequencing is not readily accessible in all centers, and its use generally leads to additional delay in providing a comprehensive IDH1/2 mutational status assessment. Moreover, both IHC and Sanger sequencing lack of standardization can lead to inter-laboratory variability.

Beyond IDH1 R132H, 11 other IDH mutations, 6 IDH1 and 5 IDH2, have been reported so far in literature, in large cohorts of gliomas.

Recent data indicate that mutant specific IDH1-inhibitors impair growth of mutant-IDH1 gliomas in mice.

A real-time PCR assay was designed to detect the 12 IDH1/2 mutations in one single step in FFPE samples, and identify the most frequent ones

OBJECTIVES

- Establish the analytical performance of the new IDH1/2 PCR assay
- Validate the IDH1/2 PCR assay performance on FFPE glioma clinical samples by comparing PCR IDH mutational status to IHC and Sanger sequencing.

MATERIAL & METHODS

IDH1/2 one-step qPCR assay:

PCR Clamping was used for the qualitative detection of IDH1 R132H and 11 additional IDH1/2 mutations. ARMS PCR technology was combined to selectively identify the most frequent IDH1 (R132H / R132C) and IDH2 (R172K) mutations (Table 1).

Evaluation of analytical sensitivity:

Limit of Detection (LOD) (min. % mutant DNA detected in a WT background) was determined following CLSI/NCCLS EP17-A guidelines.
 5 low positive samples (2.5-10-15 and 20%) obtained by mixing IDH mutant plasmid DNA with glioma IDH1/2 WT DNA were tested per mutation (n= 30 to 110 measurements per mutation and mutation percentage)

Validation of the IDH1/2 PCR assay:

171 FFPE glioma samples series: 121 samples retrospectively collected in a reverse chronological order from 2 academic centers (n1=102; n2=19) and 50 additional commercial samples. No specific selection criteria beyond tumor characteristics assessed by local pathologists.
 Samples selection: < 10 yrs; ≥ 50mm² tissue area with ≥ 40% tumoral cells.
 DNA extracted from 10 µm FFPE sections using the QIAamp DNA FFPE Tissue Kit (Qiagen)
 qPCRs performed on 25 ng DNA acc. to *therascreen* IDH1/2 RGQ PCR Kit protocol (Qiagen) and run on the Rotor-Gene Q 5plex HRM instrument (Qiagen)
 IHC was performed locally using the IDH1 R132H monoclonal antibody Clone H09 (Dianova)
 Bidirectional Sequencing (central) using recommended primers for IDH1 codon 132 and IDH2 codon 172⁽³⁾ and newly designed primers for IDH1 codon 100.

22 Synthetic samples (30% and 40% Mutant DNA in WT DNA) were additionally tested for the 11 rare mutations. Synthetic samples were processed similarly to clinical samples

RESULTS

IDH1/2 ASSAY CONTENT


Gene / Codon	Mutation	Base change
IDH1 / R132	R132H*	395 G>A
	R132C*	394 C>T
	R132S	394 C>A
	R132G	394 C>G
	R132L	394 G>T
R132V	394_395 CG>GT	
IDH1 / R100	R100Q	299 G>A
	R172K*	515 G>A
IDH2 / R172	R172M	515 G>T
	R172W	514 A>T
	R172S	516 G>T
	R172G	514 A>G

Table 1 - IDH1/2 mutations detected and identified* with the IDH1/2 PCR Assay


Detection of 12 mutations (6 within IDH1 codon R132, 5 within the homologous codon 172 of IDH2, and one within IDH1 codon 100)

Identification of 3 major IDH1/2 mutations (IDH1 R132H / R132C and IDH2 R172K)

ANALYTICAL SENSITIVITY



CLINICAL VALIDATION COHORT



- 147 samples met the selection criteria - The histological subtype distribution reflected typical clinical distribution with > 40% GBMs.
- IDH1/2 PCR success was 100% on samples collected in academic centers

CLINICAL PERFORMANCE

Table 2 - Agreement between IHC and IDH1/2 PCR for detection of IDH1 R132H mutation.

IDH1/2 PCR	IHC		TOTAL
	R132H POS	R132H NEG	
R132H	60	0	60
NON-R132H	1	72	73
TOTAL	61	72	133

Table 3 - Agreement between Sequencing and IDH1/2 PCR for IDH1/2 mutation status assessment.

IDH1/2 PCR	SANGER SEQUENCING		TOTAL
	IDH1/2 MUT +	IDH1/2 MUT -	
IDH1/2 MUT +	67	5	72
IDH1/2 MUT -	0	61	61
TOTAL	67	66	133

Table 4 - Number and types of IDH mutations detected by IHC, Sequencing or IDH1/2 PCR (n=133)

	IDH1 R132H	IDH1 R132C	IDH1 R132 OTHER	IDH1 R100	IDH2 R172K	IDH2 R172 OTHER	WT	% Mutated Cases
IHC	61	0	0	0	0	0	72	46 %
SEQUENCING	58	2	6*	0	1	0	66	50 %
IDH1/2 PCR	60	3	7	0	1	1	61	54 %

* 3 R132S; 2 R132G; 1 R132L

- Overall concordance between IHC and PCR for IDH1 R132H detection was 99% (Table 2) - The only PCR/IHC discordant case was a sample of the commercial series.
- Comparison to the reference Sanger Sequencing with a typical LOD of 15-20 % showed a high overall concordance of 96% - PCR detected 5 additional mutated cases (2 IDH1 R132H, 1 IDH1 R132C, 1 IDH1 R132, 1 IDH2 172) not detected by Sequencing (Table 3)
- Positive agreement between PCR and IHC was 98.4 % [91.3;99.7] and 100% between PCR and Sequencing [94.6;100] thus exceeding the target (PPA ≥95% with lower confidence limit ≥90%)
- All but one IHC-positive cases (n=61) were concordantly identified as IDH1 R132H by PCR whereas 3 (incl. the PCR-neg case) were not detected by Sequencing (Table 2-4)
- Out of the IHC-negative cases (n=72), 9 rare mutations were additionally detected by Sequencing (12 %) and 12 (10 IDH1, 2 IDH2) by PCR (17 %) (Table 4)

SYNTHETIC SAMPLES

- The IDH1/2 assay correctly detected the 11 IDH1/2 rare mutations at the two tested mutation-rates (30% and 40%).

SUMMARY AND CONCLUSION

The newly developed IDH1/2 PCR assay showed a high technical success.

The PCR assay analytical sensitivity (different for each mutation) was high - LOD was < 5% for all but one (rare IDH2) mutations, making this assay more sensitive than published references for sequencing techniques.

Positive concordance with IHC (R132H) and Sequencing was high (98% & 100% resp.)

The PCR assay detected 12 rare IDH1/2 mutations (9%)

Additional work is ongoing to confirm the rare mutations identified with the kit with more sensitive sequencing techniques

The IDH1/2 PCR assay can reliably be performed on FFPE glioma samples of less than 10 yrs of age - Analysis of a small series of samples > 10 yrs (n=17) also indicated a high technical success rate of IDH1/2 PCR (94%) and a high concordance between PCR with IHC (15/16) and Sequencing (14/16).

This new highly standardized IDH1/2 assay is able to detect the major IDH1R132H mutation and 11 rare IDH1/2 mutations in one step and should facilitate the implementation of IDH testing in pathological laboratories to optimize diagnosis of glioma in routine clinical setting

REFERENCES

- Parsons DW et al. 2008; Science 321(5897):1807-12
- Yan H et al. 2009 N Engl J Med 360(8):765-73.
- Preusser M et al. 2011; Clin Neuropath. 30:217-230
- Hartmann C, et al. Acta Neuropathol. 2009
- Pusch S et al. 2011; Neuropathol Appl Neurobiol. 37(4):428-30
- Rohle D et al. 2013; Science. 340(6132):626-30
- van den Bent MJ et al. 2013; J Neurooncol. 112(2):173-8