

ipsogen[®] RT Kit Handbook



Version 1

IVD

In vitro diagnostics



REF

679923



QIAGEN GmbH, QIAGEN Strasse 1, 40724 Hilden, GERMANY

R3

MAT

1072504EN



QIAGEN Sample and Assay Technologies

QIAGEN is the leading provider of innovative sample and assay technologies, enabling the isolation and detection of contents of any biological sample. Our advanced, high-quality products and services ensure success from sample to result.

QIAGEN sets standards in:

- Purification of DNA, RNA, and proteins
- Nucleic acid and protein assays
- microRNA research and RNAi
- Automation of sample and assay technologies

Our mission is to enable you to achieve outstanding success and breakthroughs. For more information, visit www.qiagen.com.

Contents

Intended Use	4
Summary and Explanation	4
Principle of the Procedure	4
Materials Provided	5
Kit contents	5
Materials Required but Not Provided	6
Warnings and Precautions	7
General precautions	7
Reagent Storage and Handling	8
Specimen Handling and Storage	8
Procedure	9
Sample RNA preparation	9
Protocol	
■ Reverse transcription	9
Troubleshooting	11
Quality Control	11
Limitations	11
Performance Characteristics	11
References	11
Symbols	12
Contact Information	12
Ordering Information	13

Intended Use

The *ipsogen* RT Kit allows reverse transcription of total RNA to be used in molecular diagnostic testing.

Note: Performance characteristics for the *ipsogen* RT Kit have only been established with BCR-ABL Mbc and ABL transcripts. The user is responsible for establishing appropriate performance characteristics for other target transcripts.

Summary and Explanation

Reverse transcription of RNA is required when quantifying RNA (e.g., by RT-PCR or real-time RT-PCR) or cloning a sequence of RNA. Reverse transcriptases are used *in vitro* for first-strand cDNA synthesis with RNA as the starting template. The efficiency of the reaction is highly dependent on the quality and quantity of the starting RNA template. It is important to have intact RNA as starting template. Even trace amounts of contaminating RNases in the RNA sample can cause RNA cleavage, resulting in shortened cDNA products. Chemical impurities, such as protein, poly-anions (e.g., heparin), salts, EDTA, ethanol, phenol, and other solvents, can affect the activity and processivity of the reverse transcriptase.

To ensure reproducible and efficient reverse transcription, it is important to determine the quality and quantity of the starting RNA. For best results, we recommend starting with RNA purified using silica-gel–membrane technology. For example, the QIAGEN® RNeasy® Mini Kit (cat. no. 74104), or RNeasy Midi Kit (cat. no. 75144) can be used to isolate RNA from a variety of starting materials and provide high-quality RNA highly suited for use in reverse-transcription and RT-PCR applications.

Principle of the Procedure

Reverse transcriptase is a multifunctional enzyme with 3 distinct enzymatic activities: an RNA-dependent DNA polymerase, a hybrid-dependent exoribonuclease (RNase H), and a DNA-dependent DNA polymerase. *In vivo*, the combination of these 3 activities allows transcription of the single-stranded RNA genome into double-stranded DNA for retroviral infection. The RNA-dependent DNA-polymerase activity (reverse transcription) transcribes cDNA from an RNA template. This activity allows synthesis of cDNA for cloning, PCR, and RNA sequencing.

Materials Provided

Kit contents

<i>ipsogen</i> RT Kit	(33)
Catalog no.	679923
Number of reactions	33
Reverse Transcriptase	36 μ l
5x RT Buffer for reverse transcription	180 μ l
dNTP Mix*	72 μ l
Random Primer†	190 μ l
RNase Inhibitor	18 μ l
DTT‡	45 μ l
<i>ipsogen</i> RT Kit Handbook (English)	1

* Deoxynucleotides 10 mM each.

† Random nonamer oligonucleotide.

‡ Dithiothreitol.

Materials Required but Not Provided

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, consult the appropriate safety data sheets (SDSs), available from the product supplier.

Consumables

- Nuclease-free aerosol-resistant sterile PCR pipet tips with hydrophobic filters
- 0.5 ml or 0.2 ml RNase- and DNase-free PCR tubes
- Ice

Reagents

- Nuclease-free PCR grade water
- Reagents for 1.2% formaldehyde agarose gel electrophoresis

Equipment

- Microliter pipets* dedicated for PCR (1–10 μ l; 10–100 μ l; 100–1000 μ l)
- Benchtop centrifuge* with rotor for 0.2 ml/0.5 ml reaction tubes (capable of attaining 10,000 rpm)
- Spectrophotometer,* or Agilent® BioAnalyzer®,* for RNA quantitation
- Equipment* for pulsed-field gel electrophoresis
- Thermomixer, heated orbital incubator, heating block, or water bath* (reverse transcription step)

* Ensure that instruments have been checked and calibrated according to the manufacturer's recommendations.

Warnings and Precautions

For in vitro diagnostic use

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient and compact PDF format at www.qiagen.com/safety where you can find, view, and print the SDS for each QIAGEN kit and kit component.

Discard sample and assay waste according to your local safety regulations.

General precautions

qPCR tests require good laboratory practices, including equipment maintenance, that are dedicated to molecular biology and compliant with applicable regulations and relevant standards.

This kit is intended for in vitro diagnostic use. Reagents and instructions supplied in this kit have been validated for optimal performance. Further dilution of the reagents or alteration of incubation times and temperatures may result in erroneous or discordant data. All reagents are formulated specifically for use with this kit. For optimal performance of the procedure, no substitutions should be made.

Determining transcript levels using qPCR requires both the reverse transcription of the mRNA and the amplification of the generated cDNA by PCR. Therefore, the entire assay procedure must be performed under RNase-/DNase-free conditions.

Use extreme caution to prevent:

- RNase/DNase contamination, which might cause degradation of the template mRNA and the generated cDNA
- mRNA or PCR carryover contamination resulting in false positive signal

We therefore recommend the following.

- Use nuclease-free labware (e.g., pipets, pipet tips, reaction vials) and wear gloves when performing the assay.
- Use fresh aerosol-resistant pipet tips for all pipetting steps to avoid cross-contamination of the samples and reagents.
- Prepare pre-PCR master mix with dedicated material (pipets, tips, etc.) in a dedicated area where no DNA matrixes (cDNA, DNA, plasmid) are introduced. Add template in a separate zone (preferably in a separate room) with specific material (pipets, tips, etc.).

Reagent Storage and Handling

The kit is shipped on dry ice and must be stored at -30°C to -15°C upon receipt.

- Gently mix and centrifuge the tubes before opening.
- Store all kit components in original containers.

These storage conditions apply to both opened and unopened components. Components stored under conditions other than those stated on the labels may not perform properly and may adversely affect the assay results.

Expiration dates for each reagent are indicated on the individual component labels. Under correct storage conditions, the product will maintain performance until the expiration date printed on the label.

Specimen Handling and Storage

Whole blood samples should be anti-coagulated with potassium EDTA and stored at $2-8^{\circ}\text{C}$ for no more than 5 days before RNA extraction.

Procedure

Sample RNA preparation

RNA extraction must be performed with a validated procedure (QIAGEN RNeasy Mini Kit, cat. no. 74104 or RNeasy Midi Kit, cat. no. 75144; or Life Technologies TRIzol[®], cat. nos. 15596-026 and 15596-018).

The performance of an assay is dependent on the concentration and quality of input RNA. We therefore recommend qualifying the purified RNA, prior to downstream analysis, by agarose* gel electrophoresis, Agilent BioAnalyzer, or spectrophotometry.[†]

Protocol: Reverse transcription

Things to do before starting

- Thaw all necessary components and place them on ice.
- Mix contents of tubes well (do not vortex) and spin briefly (10 s at 10,000 rpm) to collect the liquid at the bottom of the tube.
- Adjust RNA samples to 0.1 $\mu\text{g}/\mu\text{l}$ with nuclease-free water.

Note: A no template control (NTC) generated during reverse transcription using nuclease-free water as template can be used to check reverse transcription quality.

Procedure

1. Incubate 1 μg of each RNA sample to be tested (10 μl) for 5 min at 65°C.
2. Immediately cool on ice for 5 min.
3. Centrifuge briefly (10 s at 10,000 rpm) to collect the liquid at the bottom of the tube. Keep on ice.
4. Prepare the reverse transcription premix on ice, and keep on ice (see Table 1).

* When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles.

[†] Optical density measured at 260 and 280 nm: OD of 1.0 at 260 nm is equivalent to approximately 40 $\mu\text{g}/\text{ml}$ single-stranded RNA. An A_{260}/A_{280} ratio between 1.8 and 2.1 is indicative of highly purified RNA.

Table 1. Preparation of reverse transcription premix

Premix component	Volume per sample (μl)*	Final concentration
5x Reverse transcriptase buffer	5.0	1x
dNTP (10 mM each)	2.0	0.8 mM
Random nonamer (100 μ M)	5.25	21 μ M
RNase Inhibitor (40 U/ μ l)	0.5	0.8 U/ μ l
Reverse transcriptase (200 U/ μ l)	1.0	8 U/ μ l
DTT	1.25	–
RT premix volume per sample	15	

* Prepare $n + 1$, where n is the number of RNA samples

- 5. Mix with care (do not vortex), centrifuge briefly, and add 15 μ l of the premix to each RNA sample (for 40 ng/ μ l) and to the water control (NTC).**
- 6. Mix each tube with care (do not vortex) and centrifuge briefly.**
- 7. Run the reverse transcription program (see Table 2) using a thermomixer, heated orbital incubator, heating block, or water bath.**

Table 2. Settings for reverse transcription

Reverse transcription 1	25°C for 10 min
Reverse transcription 2	50°C for 60 min
Inactivation	85°C for 5 min
Cooling	4°C for 5 min

- 8. Centrifuge briefly (10 s at 10,000 rpm) to collect the cDNA at the bottom of the tube.**
- 9. Keep on ice, or store at –20°C, until qPCR is performed.**

Troubleshooting

For information on troubleshooting for this kit, see the Frequently Asked Questions page at our Technical Support Center: www.qiagen.com/FAQ/FAQList.aspx. The scientists in QIAGEN Technical Services are always happy to answer any questions you may have about either the information and protocol in this handbook or sample and assay technologies (for contact information, see “Contact Information”, page 12).

Quality Control

This kit is manufactured according to ISO 13485 standard. Certificates of Analysis are available upon request at www.qiagen.com/support/.

Limitations

The users must be trained and familiar with this technology prior to the use of this device. This kit should be used following the instructions given in this manual, in combination with a validated instrument mentioned in “Materials Required but Not Provided”, page 6.

Any diagnostic results generated must be interpreted in conjunction with other clinical or laboratory findings. It is the user’s responsibility to validate system performance for any procedures used in their laboratory which are not covered by the QIAGEN performance studies.

Attention should be paid to expiration dates printed on the box and labels of all components. Do not use expired components.

Performance Characteristics

Performance characteristics for the *ipsogen* RT Kit have only been established with BCR-ABL Mbc and ABL transcripts. The user is responsible for establishing appropriate performance characteristics for other target transcripts.

For further information regarding performance characteristics of the *ipsogen* RT Kit, see the *ipsogen BCR-ABL1 Mbc IS-MMR DX Kit Handbook*, available at www.qiagen.com.

References

QIAGEN maintains a large, up-to-date online database of scientific publications utilizing QIAGEN products. Comprehensive search options allow you to find the articles you need, either by a simple keyword search or by specifying the application, research area, title, etc.

For a complete list of references, visit the QIAGEN Reference Database online at www.qiagen.com/RefDB/search.asp or contact QIAGEN Technical Services or your local distributor.

Symbols

The following symbols may appear on the packaging and labeling:

 Σ <N>	Contains reagents sufficient for <N> reactions
	Use by
	In vitro diagnostic medical device
	Catalog number
	Lot number
	Material number
	Global Trade Item Number
	Temperature limitation
	Manufacturer
	Consult instructions for use

Contact Information

For technical assistance and more information, please see our Technical Support Center at www.qiagen.com/Support, call 00800-22-44-6000, or contact one of the QIAGEN Technical Service Departments or local distributors (see back cover or visit www.qiagen.com).

Ordering Information

Product	Contents	Cat. no.
<i>ipsogen</i> RT Kit (33)	For 33 reactions: Reverse transcriptase, 5x RT buffer, dNTP mix, Random primer, RNase Inhibitor, DTT	679923
Rotor-Gene® Q MDx — for IVD-validated real-time PCR analysis in clinical applications		
Rotor-Gene Q MDx 5plex HRM Platform	Real-time PCR cyclers and High Resolution Melt analyzer with 5 channels (green, yellow, orange, red, crimson) plus HRM channel, laptop computer, software, accessories, 1-year warranty on parts and labor, installation and training not included	9002032
Rotor-Gene Q MDx 5plex HRM System	Real-time PCR cyclers and High Resolution Melt analyzer with 5 channels (green, yellow, orange, red, crimson) plus HRM channel, laptop computer, software, accessories, 1-year warranty on parts and labor, installation and training	9002033

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

This page intentionally left blank

This product is intended for in vitro diagnostic use. *ipsogen* products may not be resold, modified for resale or used to manufacture commercial products without written approval of QIAGEN.

Information in this document is subject to change without notice. QIAGEN assumes no responsibility for any errors that may appear in this document. This document is believed to be complete and accurate at the time of publication. In no event shall QIAGEN be liable for incidental, special, multiple, or consequential damages in connection with, or arising from the use of this document.

ipsogen products are warranted to meet their stated specifications. QIAGEN's sole obligation and the customer's sole remedy are limited to replacement of products free of charge in the event products fail to perform as warranted.

This product contains SuperScript® III Reverse Transcriptase that is the subject of one or more issued US patents or pending US patent applications, and corresponding non-US equivalents, owned by Life Technologies Corporation, and is sold under an agreement between Life Technologies Corporation and *ipsogen*. The purchase price of this product includes limited, nontransferable rights under the foregoing patents to use this amount of the product to practice the claims in said patents solely for activities of the purchaser in measuring BCR-ABL p210 transcripts. No other rights are conveyed, including no right to use this product for applications in forensics. Further information about acquiring rights under patents owned by Life Technologies Corporation may be obtained by contacting the Licensing Department, Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, CA 92008. (760) 603-7200. E-mail: Outlicensing@lifetech.com.

Trademarks: QIAGEN®, *ipsogen*®, RNeasy®, Rotor-Gene® (QIAGEN Group); SuperScript® (Life Technologies Corporation); Agilent®, Bioanalyzer® (Agilent Technologies, Inc.); TRIzol® (Molecular Research Center, Inc.).

Limited License Agreement

Use of this product signifies the agreement of any purchaser or user of the *ipsogen* RT Kit to the following terms:

1. The *ipsogen* RT-Dx Kit may be used solely in accordance with the *ipsogen RT Kit Handbook* and for use with components contained in the Kit only. QIAGEN grants no license under any of its intellectual property to use or incorporate the enclosed components of this Kit with any components not included within this Kit except as described in the *ipsogen RT Kit Handbook* and additional protocols available at www.qiagen.com.
2. Other than expressly stated licenses, QIAGEN makes no warranty that this Kit and/or its use(s) do not infringe the rights of third-parties.
3. This Kit and its components are licensed for one-time use and may not be reused, refurbished, or resold.
4. QIAGEN specifically disclaims any other licenses, expressed or implied other than those expressly stated.
5. The purchaser and user of the Kit agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. QIAGEN may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to the Kit and/or its components.

For updated license terms, see www.qiagen.com.

© 2015 QIAGEN, all rights reserved.

www.qiagen.com

Australia ■ techservice-au@qiagen.com

Austria ■ techservice-at@qiagen.com

Belgium ■ techservice-bnl@qiagen.com

Brazil ■ suportetecnico.brasil@qiagen.com

Canada ■ techservice-ca@qiagen.com

China ■ techservice-cn@qiagen.com

Denmark ■ techservice-nordic@qiagen.com

Finland ■ techservice-nordic@qiagen.com

France ■ techservice-fr@qiagen.com

Germany ■ techservice-de@qiagen.com

Hong Kong ■ techservice-hk@qiagen.com

India ■ techservice-india@qiagen.com

Ireland ■ techservice-uk@qiagen.com

Italy ■ techservice-it@qiagen.com

Japan ■ techservice-jp@qiagen.com

Korea (South) ■ techservice-kr@qiagen.com

Luxembourg ■ techservice-bnl@qiagen.com

Mexico ■ techservice-mx@qiagen.com

The Netherlands ■ techservice-bnl@qiagen.com

Norway ■ techservice-nordic@qiagen.com

Singapore ■ techservice-sg@qiagen.com

Sweden ■ techservice-nordic@qiagen.com

Switzerland ■ techservice-ch@qiagen.com

UK ■ techservice-uk@qiagen.com

USA ■ techservice-us@qiagen.com

