

# Quick-Start Protocol

## miRNeasy Mini Kit

March 2016

The miRNeasy Mini Kit (cat. no. 217004) can be stored dry at room temperature (15–25°C) for at least 9 months if not otherwise stated on label. QIAzol® Lysis Reagent can be stored at room temperature or at 2–8°C.

### Further information

- *miRNeasy Mini Handbook*: [www.qiagen.com/HB-1277](http://www.qiagen.com/HB-1277)
- Safety Data Sheets: [www.qiagen.com/safety](http://www.qiagen.com/safety)
- Technical assistance: [support.qiagen.com](http://support.qiagen.com)

### Notes before starting

- This protocol is for purifying total RNA, including small RNAs, from animal cells and tissue.
  - If necessary, redissolve any precipitate in Buffer RWT by warming.
  - Except for phase separation (step 5), all steps should be performed at room temperature (15–25°C). Work quickly.
  - Add ethanol (96–100%) to Buffer RWT and Buffer RPE concentrates before use (see bottle label for volume).
  - Before starting with step 1 for the first time, select disruption and homogenization methods according to recommendations in the *miRNeasy Mini Kit Handbook*.
1. Add 700 µl QIAzol Lysis Reagent to the sample and disrupt and homogenize using an appropriate method.
  2. Incubate the homogenate at room temperature (15–25°C) for 5 min.
  3. Add 140 µl chloroform and cap tube securely. Shake vigorously for 15 s.
  4. Incubate at room temperature for 2–3 min.
  5. Centrifuge for 15 min at 12,000 x g at 4°C.

6. Transfer the upper aqueous phase to a new collection tube. Avoid transferring any interphase. Add 1.5 volumes (usually 525  $\mu$ l) of 100% ethanol, and mix thoroughly by pipetting.
7. Pipet up to 700  $\mu$ l sample, including any precipitate, into an RNeasy<sup>®</sup> Mini column in a 2 ml collection tube. Close the lid and centrifuge at  $\geq 8000 \times g$  for 15 s at room temperature. Discard the flow-through.
8. Repeat step 7 using the remainder of the sample.
9. **Optional:** Perform DNase digest according to instructions in Appendix B of the handbook (not required for detecting **mature miRNA** using the miScript PCR system).
10. Add 700  $\mu$ l Buffer RWT to the RNeasy Mini column. Close the lid, and centrifuge for 15 s at  $\geq 8000 \times g$ . Discard the flow-through.  
This step is optional if working with cultured cells.
11. Pipet 500  $\mu$ l Buffer RPE onto the RNeasy Mini column. Close the lid, and centrifuge for 15 s at  $\geq 8000 \times g$ . Discard the flow-through.
12. Add 500  $\mu$ l Buffer RPE to the RNeasy Mini column. Close the lid, and centrifuge for 2 min at  $\geq 8000 \times g$ .
13. **Optional:** Place the RNeasy Mini column into a new 2 ml collection tube. Centrifuge at full speed for 1 min to further dry the membrane.
14. Transfer the RNeasy Mini column to a new 1.5 ml collection tube. Pipet 30–50  $\mu$ l RNase-free water directly onto the RNeasy Mini column membrane. Close the lid, and centrifuge for 1 min at  $\geq 8000 \times g$  to elute.
15. If expected RNA yield is  $>30 \mu$ g, repeat step 14 using an additional 30–50  $\mu$ l RNase-free water or using the eluate from step 14 (if high RNA concentration is required). Reuse the collection tube from step 14.



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. Trademarks: QIAGEN<sup>®</sup>, Sample to Insight<sup>®</sup>, QIAzol<sup>®</sup>, RNeasy<sup>®</sup> (QIAGEN Group). 1101293 03/2016 HB-0604-002 © 2016 QIAGEN, all rights reserved.