

QIAseq™ UPX 3' Targeted RNA Panel

Part 1: Cell lysis, reverse transcription

Upon receipt, store the QIAseq UPX 3' Targeted RNA Panel at –30 to –15°C. The QIAseq Beads and Cell Index (ID) RT Plate should be stored at 4°C.

Further information

- *QIAseq UPX 3' Targeted RNA Library Handbook*: www.qiagen.com/HB-2480
- Safety Data Sheets: www.qiagen.com/safety
- Technical assistance: support.qiagen.com

Notes before starting

- If using the multi-use 96-well Cell ID RT Plate, consult the handbook.
- The recommended starting amount is 1–100 cells.
- For single cells, a minimum of 24 cells must be multiplexed per sample index.
- If purified RNA is used, consult the handbook.
- Ensure reactions are thoroughly mixed (12 times), prepared and incubated at recommended temperatures. Do not vortex.

Cell lysis

1. Prepare the Cell Lysis Premix on ice as described in Table 1.

Table 1 Preparation of Cell Lysis Premix

Component	Volume 8 wells	Volume 24 wells	Volume 96 wells	Volume 384 wells
Cell Lysis Buffer	10 μ l	30 μ l	120 μ l	480 μ l
RNase Inhibitor	2.5 μ l	7.5 μ l	30 μ l	120 μ l
Nuclease-free Water	17.5 μ l	52.5 μ l	210 μ l	840 μ l
Total volume	30 μl	90 μl	360 μl	1440 μl

2. Aliquot 3 μ l of the premix into the required wells of the single-use Cell ID RT Plate.
3. Capture cells into the plate containing the Cell Lysis Premix.
4. Incubate for 15 min on ice.
5. Freeze at -80°C or proceed immediately with "Reverse transcription".

Reverse transcription

6. Prepare the RT Premix on ice as described in Table 2.

Table 2 Preparation of RT Premix for single-use UPX 3' Cell ID RT Plates

Component	Volume 8 wells	Volume 24 wells	Volume 96 wells	Volume 384 wells
3' Tar RT Buffer	10 μ l	30 μ l	120 μ l	480 μ l
EZ Reverse Transcriptase	2.5 μ l	7.5 μ l	30 μ l	120 μ l
Nuclease-free Water	5 μ l	15 μ l	60 μ l	240 μ l
ERCC diluted 1.25 x 10 ⁶ -fold* or Nuclease-free Water	2.5 μ l	7.5 μ l	30 μ l	120 μ l
Total volume	20 μl	60 μl	240 μl	960 μl

*For single-cell analysis, ERCC can be added as a control. Prepare the ERCC by diluting the stock 1.25x10⁶-fold using 0.1X TE Buffer.

- Aliquot 2 μ l of the RT Premix into the wells.
- Incubate as follows: 5 min at 25°C, 1 h at 42°C, 5 min at 95°C, and hold at 4°C.
- Combine the synthesized cDNA from the different wells into one, 2 ml microcentrifuge tube. Up to 96 wells can be combined in one tube.
- Add 1.3X (volume) QIAseq Beads to the combined cDNA synthesis reactions as described in Table 3. Mix well by pipetting up and down 12 times.

Table 3 Addition of QIAseq Beads for cDNA cleanup

Number of wells combined	QIAseq Bead volume
8	52 μ l
24	156 μ l
96*	624 μ l

* When working with 384 wells, perform the cleanup as 4 sets of 96 wells. Combine supernatants at step 18.

- Incubate for 5 min at room temperature.
- Place the tubes on a magnetic rack for 5 min. After the solution has cleared, with the beads still on the magnetic stand, carefully remove and discard the supernatant.
- With the beads still on the magnetic stand, add 200 μ l 80% ethanol. Rotate the tube (3 times) to wash the beads. Carefully remove and discard the wash.

14. Repeat the ethanol wash. Completely remove all traces of the ethanol wash.

15. With the tube (caps opened) still on the magnetic stand, air dry for 10 min.

Note: Visually inspect that the pellet is completely dry.

16. Remove the beads from the magnetic stand and elute the DNA from the beads by adding 52 μ l nuclease-free water. Mix well by pipetting.

17. Return the tube/plate to the magnetic rack until the solution has cleared.

18. Transfer 50 μ l of the supernatant to clean tubes/plate.

Important: With 384 wells, mixture can be combined across all four sets to consolidate all 384 indices into one tube.

19. Add 1.3X volume QIAseq Beads to the supernatant. Mix well by pipetting up and down 12 times.

Note: QIAseq Bead volume will be 65 μ l if supernatant is 50 μ l.

20. Incubate for 5 min at room temperature.

21. Place the tubes on a magnetic rack for 5 min. After the solution has cleared, with the beads still on the magnetic stand, carefully remove and discard the supernatant.

22. With the beads still on the magnetic stand, add 200 μ l 80% ethanol. Rotate the tube (3 times) to wash the beads. Carefully remove and discard the wash.

23. Repeat the ethanol wash. Completely remove all traces of the ethanol wash.

24. With the tube (caps opened) still on the magnetic stand, air dry for 10 min.

Note: Visually inspect that the pellet is completely dry.

25. Remove the beads from the magnetic stand, and elute the DNA from the beads by adding 12 μ l nuclease-free water. Mix well by pipetting.

26. Return the tube/plate to the magnetic rack until the solution has cleared.

27. Transfer 10 μ l of the supernatant to clean tubes/plate.

Note: From this point of the protocol, the procedures assume that all cDNA wells (8, 24, 96 or 384) have been combined into a single tube.

28. Proceed with "Single primer extension" in Quick-Start Protocol Part 2. Alternatively, the samples can be stored at -30 to -15°C in a constant-temperature freezer.





Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN®, Sample to Insight®, QIAseq™ (QIAGEN Group). 1113984 05/2018 HB-2473-003 © 2018 QIAGEN, all rights reserved.