

User-Developed Protocol:

Isolation of plasmid DNA from *Bacillus subtilis* using the QIAGEN® Plasmid Midi Kit

This procedure has been adapted by customers from the QIAGEN® Plasmid Midi Kit Protocol. **It has not been thoroughly tested and optimized by QIAGEN.**

The procedure has been used successfully for isolation of high- and low-copy-number plasmids from various *Bacillus subtilis* strains. Yield of plasmid DNA was typically 10–20 µg plasmid DNA from 100 ml culture.

Please be sure to read the *QIAGEN Plasmid Purification Handbook* and the detailed QIAGEN Plasmid Midi Kit Protocol carefully before beginning this procedure.

Procedure

- 1. Pick a single colony from a selective plate and inoculate a starter culture of 10 ml LB medium containing the appropriate antibiotic. Grow overnight at 37°C with vigorous shaking (240 rpm).**
- 2. Dilute the miniculture 1:50 to 1:100 into 100 ml selective LB medium. Grow at 37°C for 3–4 hours with vigorous shaking (~240 rpm).**
The culture should reach an A_{600} of 0.8–1.2 units/ml.
- 3. Harvest the cells by centrifugation at 3000 x g for 15 min at 4°C.**
- 4. Resuspend the bacterial pellet in 4 ml Buffer P1 containing 5 mg/ml lysozyme.**
Ensure that RNase A (100 µg/ml) has been added to Buffer P1.
- 5. Incubate at 37°C for 30 min.**
- 6. Add 4 ml Buffer P2, mix gently but thoroughly by inverting 4–6 times, and incubate at room temperature for 5 min.**
Check Buffer P2 before use for SDS precipitation due to low storage temperatures. If necessary, dissolve the SDS by warming to 37°C.
- 7. Add 4 ml chilled Buffer P3, mix immediately but gently by inverting 4–6 times, and incubate on ice for 15 min.**
- 8. Centrifuge at $\geq 20,000$ x g for 30 min at 4°C. Remove supernatant containing plasmid DNA promptly.**
- 9. Centrifuge again at $\geq 20,000$ x g for 15 min at 4°C. Remove supernatant containing plasmid DNA promptly.**
- 10. Equilibrate a QIAGEN-tip 100 by applying 4 ml Buffer QBT, and allow the column to empty by gravity flow.**

User-developed protocol

11. Apply the supernatant from step 9 to the QIAGEN-tip and allow it to enter the resin by gravity flow.
12. Wash the QIAGEN-tip with 2 x 10 ml Buffer QC.
13. Elute DNA with 5 ml Buffer QF.
14. Precipitate DNA by adding 3.5 ml room-temperature isopropanol to the eluted DNA. Mix and centrifuge immediately at $\geq 15,000 \times g$ for 30 min at 4°C. Carefully decant the supernatant.
15. Wash the DNA pellet with 2 ml of room-temperature 70% ethanol and centrifuge at $\geq 15,000 \times g$ for 10 min. Carefully decant the supernatant without disturbing the pellet.
16. Air-dry the pellet for 5–10 min, and redissolve the DNA in a suitable volume of buffer (e.g., TE, pH 8.0, or 10 mM Tris·Cl, pH 8.5).

QIAGEN handbooks can be requested from QIAGEN Technical Service or your local QIAGEN distributor. Selected handbooks can be downloaded from www.qiagen.com/literature/handbooks/default.asp. Material safety data sheets (MSDS) for any QIAGEN product can be downloaded from www.qiagen.com/ts/msds.asp.

Trademark: QIAGEN® (QIAGEN).

© 2001 QIAGEN, all rights reserved.