

Quick-Start Protocol

Investigator[®] 24plex GO! Kit, Part 3

Protocol for buccal swabs

All components of the Investigator 24plex GO! Kit should be stored at –30 to –15°C. Avoid repeated thawing and freezing. The primer mix, allelic ladder, and DNA Size Standard must be stored protected from light. DNA samples and post-PCR reagents (allelic ladder and DNA size standard) should be stored separately from the PCR reagents. Under these conditions, the components are stable until the expiration date indicated on the kit.

Further information

- *Investigator 24plex GO! Handbook*: www.qiagen.com/HB-1913
- Safety Data Sheets: www.qiagen.com/safety
- Technical assistance: support.qiagen.com

Notes before starting

- Set up all reaction mixtures in an area separate from that used for DNA isolation and PCR product analysis (post-PCR).
- Use disposable tips with hydrophobic filters to minimize cross-contamination risks.
- Before opening the tubes, thaw PCR components, vortex, and then centrifuge briefly to collect the contents at the bottom of the tubes.
- For the swab lysis protocol, Investigator STR GO! Lysis Buffer (QIAGEN, cat. no. 386516) must be ordered separately.

Procedure

1. Place the swab into a 2 ml microcentrifuge tube.
2. Carefully, cut, break off, or eject the end part of the swab.

Note: Prepare a blank swab as negative control.

3. Add 500 μ l STR GO! Lysis Buffer to the sample.
4. Incubate at 95°C for 5 min shaking at 1200 rpm in a thermomixer.

Optional: Incubate at room temperature (15–25°C) for 5 min shaking at 1200 rpm in a thermomixer.

5. Prepare a master mix according to Table 1.

As some loss of reagents can occur during transfers, prepare the mix with additional reactions included. Also include positive and negative control reactions. The master mix contains all of the components needed for PCR except the template (sample) DNA.

Table 1. Master mix setup

Component	Volume per reaction
Fast Reaction Mix 2.0	7.5 μ l
Primer Mix	12.5 μ l
Total volume	20.0 μ l

6. Vortex the reaction mix thoroughly, and dispense 20 μ l into PCR tubes or the wells of a PCR plate.
7. Mix the swab lysate thoroughly and transfer 2 μ l swab lysate to each reaction.
8. Prepare the positive and negative controls.

Positive control: Use 1 μ l Control DNA.

Negative control: Use a blank swab lysate.

9. Program the thermal cycler according to the manufacturer's instructions, using the conditions given in Table 2.

Note: If using the GeneAmp 9700 thermal cycler with an Aluminum block, use "Std Mode", or with a Silver block or Gold-plated Silver block, use "Max Mode". Do not use "9600 Emulation Mode".

Table 2a. Standard cycling conditions

Temperature	Time	Number of cycles
98°C*	30 s	
64°C	40 s	3 cycles
72°C	5 s	
96°C	10 s	
61°C	40 s	24 cycles
72°C	5 s	
68°C	5 min	–
60°C	5 min	–
10°C	∞	–

* Hot-start to activate DNA polymerase.

Table 2b. Optional cycling conditions

Temperature	Time	Number of cycles
98°C*	30 s	
64°C	40 s	3 cycles
72°C	5 s	
96°C	10 s	
61°C	40 s	24 cycles
72°C	5 s	
68°C	2 min	–
60°C	2 min	–
10°C	∞	–

* Hot-start to activate DNA polymerase.

Table 2b details previously published cycling conditions, which may continue to be used if incomplete adenylation is not visible within the electropherograms.

10. After the cycling protocol is completed, store samples at –30 to –15°C protected from light, or proceed directly with electrophoresis.

Document Revision History

Date	Changes
04/2021	This is the initial release of the third part of the Investigator 24plex GO! Kit Quick-Start Protocol, which is divided into 3 parts for printing purposes. The first and second parts of the Quick-Start Protocol are HB-1907-006 and HB-2889-001, respectively.



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN®, Sample to Insight®, Investigator® (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

1124318 04/2021 HB-2895-001 © 2021 QIAGEN, all rights reserved.