

QIAamp[®] UCP DNA Micro Kit

Store QIAamp UCP MinElute[®] columns immediately upon receipt at 2–8°C. Store the remaining kit components dry at room temperature (15–25°C). All kit components are stable for at least 12 months under these conditions.

Further information

- *QIAamp UCP DNA Micro Kit Handbook*: www.qiagen.com/HB-2174
- Safety Data Sheets: www.qiagen.com/safety
- Technical assistance: support.qiagen.com

Notes before starting

- Perform all centrifugation steps at room temperature (15–25°C).
- Set a thermal mixer to 56°C for use in step 4.
- Dilute Buffers AUW1 and AUW2 according to the labels on the bottles.
- If Buffer AUL or Buffer AUT contains precipitates, dissolve by heating to 70°C with gentle agitation.

Procedure

1. Pipet 1–100 µl sample into a 1.5 ml microcentrifuge tube (not provided). Add buffer AUT to a final volume of 100 µl.
2. Add 10 µl Proteinase K.
3. Add 100 µl of Buffer AUL to the sample. Close the cap and mix by pulse-vortexing for 15 s.
4. Incubate the sample at 56°C for 10 min. Briefly centrifuge the 1.5 ml tube to remove drops from the inside of the lid.

5. Add 50 μ l ethanol to the lysate. Close the cap, and mix thoroughly by pulse-vortexing for 15 s. Incubate for 3 min at room temperature. Briefly centrifuge the 1.5 ml tube to remove drops from inside the lid.
6. Carefully apply the mixture from step 5 to the QIAamp UCP MinElute spin column without wetting the rim. Close the cap and centrifuge at 6000 \times g for 1 min. Place the spin column in a clean 2 ml collection tube (provided), and discard the tube containing the filtrate. Close each spin column to avoid aerosol formation during centrifugation.
Note: Flow-through containing Buffer AUL or AUW1 is not compatible with bleach.
7. Carefully open the spin column and add 500 μ l Buffer AUW1 without wetting the rim. Close the cap and centrifuge at 6000 \times g for 1 min. Place the spin column in a clean 2 ml collection tube (provided), and discard the collection tube containing the filtrate.
8. Carefully open the spin column and add 500 μ l Buffer AUW2 without wetting the rim. Close the cap and centrifuge at 6000 \times g for 1 min. Place the spin column in a clean 2 ml collection tube (provided), and discard the collection tube containing the filtrate.
9. To completely dry the membrane, centrifuge the spin column at full speed (20,000 \times g) for 3 min. Place the spin column in a clean 1.5 ml elution tube (not provided) and discard the collection tube.
10. Carefully apply 20–100 μ l of Buffer AUE to the center of the spin column membrane. Close the lid and incubate at room temperature for 1 min.
Important: Ensure that the elution buffer is equilibrated to room temperature.
11. Centrifuge at full speed (20,000 \times g) for 1 min to elute the DNA.
12. Repeat steps 10 and 11 to increase yield. Re-use the first eluate if higher concentrations are required.



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN[®], QIAamp[®], MinElute[®], Sample to Insight[®] (QIAGEN Group). © 2016 QIAGEN, all rights reserved. 1103789 07/2016 HB2171-001