

Efficient Identification and Characterization of Pathogenic Bacteria in Veterinary Diagnostics

Urška Zajc, Jana Avberšek, Mateja Pate, and Matjaž Ocepek

University of Ljubljana, Veterinary Faculty, Institute of Microbiology and Parasitology,
Laboratory for Molecular Bacteriology, Ljubljana, Slovenia

Abstract: Isolates of Shiga toxin-producing *Escherichia coli* (STEC), *Campylobacter* sp., *Clostridium perfringens*, *Clostridium difficile*, methicillin-resistant *Staphylococcus aureus* (MRSA), and *Mycobacterium avium* were successfully identified using the QIAxcel® System for separation and analysis of either single or multiplex PCR products. Furthermore, the associated ScreenGel Software enabled automated analysis and identification of pathogens, thereby minimizing manual intervention and yielding reproducible and highly accurate results. Results demonstrate the utility of PCR product analysis using the QIAxcel as a method for the efficient and reliable identification of pathogenic bacteria significant in both veterinary and human medicine.

Introduction

Rapid and accurate identification of pathogens in veterinary medicine is cornerstone for a successful battle against animal infectious diseases, which can also endanger human populations. The use of antibiotics has led to the development of an increasing number of antibiotic-resistant bacteria that may cause significant health problems in humans and animals. The latter are a significant source of resistant strains. Thus, the definition of an approach to rapidly and efficiently identify pathogens is of great clinical and commercial interest in order to take appropriate preventative measures to the benefit of animal and public health.

Material and methods

Isolates of Shiga toxin-producing *Escherichia coli* (STEC), *Campylobacter* sp., *Clostridium perfringens*, *Clostridium difficile*, methicillin-resistant *Staphylococcus aureus* (MRSA) and *Mycobacterium avium* were collected in the scope of routine and research work performed at the home research institute. DNA for PCR assays was extracted from cultures with a simple lysis: the colonies were suspended in PCR-grade water, boiled at 100°C for 15 minutes and centrifuged at 12,000 rpm for 2 min. The supernatant was used for PCR without further purification. Some of the pathogens, such as STEC and thermotolerant campylobacters, were detected also directly

from enrichment media. In those cases, DNA was extracted using a commercial DNA extraction kit adapted to food samples.

For each pathogen, references related to individual PCR assays are cited in the respective sections. The subsequent electrophoretic analysis of PCR products was identical for all bacterial isolates. Using the QIAxcel capillary electrophoresis system, the amplicons were separated with the QX DNA Screening Kit employing the AM420 method with the following electrophoresis parameters: alignment marker injection at 4kV for 20 s, sample injection ►

Contents

Introduction	1
Materials and methods	1
Results and discussion	2
<i>Campylobacter</i> sp.	2
<i>Clostridium difficile</i>	3
<i>Clostridium perfringens</i>	3
Shiga toxin-producing <i>Escherichia coli</i>	4
<i>Mycobacterium avium</i>	5
Methicillin-resistant <i>Staphylococcus aureus</i>	6
Conclusion	6
References	7

at 5 kV for 10 s, and separation at 5 kV for 420 s. The alignment markers 15 bp–1 kb or 15 bp–3 kb and the DNA size markers 50–800 bp or 100 bp–2.5 kb were run simultaneously.

The selected genes listed in Table 1 were targeted in pathogen-specific PCR assays and subsequently identified by the QIAxcel ScreenGel Software. The established PCR protocols followed by an optimized analysis procedure on the QIAxcel Advanced System yielded clear results, as seen in Figures 1–6.

Table 1. Overview of genes targeted in specific PCR assays.

Bacterium	Gene	Amplicon size (bp)
<i>Escherichia coli</i> (STEC)	<i>stx1</i>	180
	<i>stx2</i>	255
	<i>eae</i>	384
	<i>stx2f</i>	428
	<i>hlyA</i>	534
<i>Campylobacter coli</i>	<i>glyA</i>	126
<i>Campylobacter upsaliensis</i>	<i>glyA</i>	204
<i>Campylobacter lari</i>	<i>glyA</i>	251
<i>Campylobacter jejuni</i>	<i>hipO</i>	323
<i>Campylobacter fetus</i> subsp. <i>fetus</i>	<i>sapB2</i>	435
<i>Campylobacter</i> sp.	23S rRNA	650
<i>Clostridium perfringens</i>	<i>etx</i>	655
	<i>cpb2</i>	567
	<i>iA</i>	446
	<i>cpa</i>	324
	<i>cpe</i>	233
	<i>cpb</i>	190
<i>Clostridium perfringens</i>	<i>netB</i>	384
<i>Clostridium difficile</i>	<i>tcdA</i>	369
	<i>tpi</i>	230
	<i>tcdB</i>	160
	<i>tcdA</i> deleted	110
Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA)	<i>nuc</i>	255
	<i>mecA</i>	527
	16S rRNA	886
<i>Mycobacterium avium</i> subsp. <i>avium</i>	IS901	1108
<i>Mycobacterium avium</i> subsp. <i>avium</i> / <i>Mycobacterium avium</i> subsp. <i>hominissuis</i>	IS901/ FR300	1742/300
<i>Mycobacterium avium</i>	IS1245	427

Results and discussion

Campylobacter sp.

Thermotolerant campylobacters are the causative agents of intestinal campylobacteriosis, the leading bacterial zoonosis in Europe and worldwide (1). *Campylobacter jejuni* and *C. coli* are the two most important species among thermotolerant campylobacters and responsible for the vast majority of cases of this food-borne illness. Another two species, *C. lari* and *C. upsaliensis*, seem to be isolated only sporadically. *C. jejuni* and *C. coli* are commensals in the intestinal tract of a wide range of birds and mammals, including domestic animals used for food production (especially poultry), and are therefore widespread in the environment. Contaminated poultry meat is considered a major source of infection for humans. Consumption and handling of undercooked contaminated poultry meat was identified as an important risk factor (2). While the usual manifestation is an acute self-limiting gastroenteritis, severe infections might also occur.

C. fetus subsp. *fetus* is usually found in the intestinal tract of cattle and sheep, and may cause sporadic abortion in both these animals (3). The pathogen is also involved in human infections.

Members of the genus *Campylobacter* have fastidious growth requirements, making conventional detection and identification procedures problematic. As a consequence, rapid and reliable detection procedures are essential. A colony multiplex PCR has been developed and optimized to simultaneously identify the 23S rRNA from *Campylobacter* sp., the *hipO* gene (hippuricase) from *C. jejuni*, the *glyA* gene (serine hydroxymethyltransferase) from *C. coli*, *C. lari*, and *C. upsaliensis*, and the *sapB2* gene (surface layer protein) from *C. fetus* subsp. *fetus*. The multiplex PCR protocol is capable of detecting the type strains and clinical isolates from all 5 species with a high degree of specificity (4). In our hands, this PCR assay proved to be a useful tool for differentiating campylobacters isolated from animal clinical samples, food, and environmental samples (Figure 1). Accurate and simple to perform, the assay offers an effective alternative to traditional biochemical methods for the identification and differentiation of *C. jejuni*, *C. coli*, *C. lari*, *C. upsaliensis*, and *C. fetus* subsp. *fetus*.




Figure 1. PCR product analysis of thermotolerant campylobacters on QIAxcel. Lanes 1 and 5: *C. lari*, lane 2: *C. coli*, lanes 3 and 4: *C. jejuni*, lane 6: *C. upsaliensis*, lane 7: DNA size marker 50–800 bp. Amplicon of the 23S rRNA gene amplified as an internal control is visible on all lanes with the exception of lane 5.

Clostridium difficile

This anaerobic sporegenic Gram-positive bacterium is known as an important cause of hospital-acquired diarrhea in humans. Recent reports have indicated an increased frequency and severity of community-acquired *C. difficile*-associated disease. Three toxins, toxin A (TcdA), toxin B (TcdB) and binary toxin, are known to be primarily responsible for the pathogenicity of this organism. Apart from toxigenic strains, non-toxigenic strains lacking all toxin genes (*tcdA*, *tcdB*, and *cdtA/B*) are known. These strains do not induce any of the symptoms specific for *C. difficile*-associated diseases. Toxigenic strains of *C. difficile* have also been recognized as a cause of disease in different animal species. On the other hand, a high percentage of subclinical animal carriers are reported. The overlap of animal strains with those from some human cases and isolation of the organism from food and the environment suggests that animals may be a reservoir and *C. difficile* a zoonosis (5).

Since toxigenic culture is a method of choice for the detection of *C. difficile* (providing isolates suitable for epidemiological analysis), a PCR or real-time PCR-toxigenic culture scheme may be proposed as a diagnostic approach for animal *C. difficile* intestinal infections/subclinical colonization, offering combined species identification and toxin type characterization of isolate (6).

In our laboratory, a conventional multiplex PCR described by Leme et al. (7) is used for the identification of isolates and toxin type characterization. Differentiation between the three main toxin types (A+B+, A-B+, and A-B-) of *C. difficile* can be performed with this multiplex PCR assay (Figure 2).




Figure 2. PCR product analysis of *Clostridium difficile* strains on QIAxcel. The genes identified are *tcdB* and *tpi* in lane 1 and *tcdB*, *tpi* and *tcdA* in lane 2. Lane 3 contains the negative control and lane 4 displays the DNA size marker 50–800 bp.

Clostridium perfringens

This Gram-positive, rod-shaped, anaerobic, spore-forming, non-motile bacterium is widely distributed in the environment and in the gut of animals and humans. *C. perfringens* is one of the most common causes of food poisoning. *C. perfringens* strains are classified into 5 types (A, B, C, D, and E) based on the production of 4 major toxins and on their fatal effect in mouse lethality tests. Type A is mediated by alpha-toxin encoded by the *cpa* gene and causes disease in both humans and animals. In humans, it is related to food poisoning and strains often produce enterotoxin encoded by the *cpe* gene. *C. perfringens* types B, C, D, and E are associated with animal diseases. Types B, C, and D are mediated by beta-toxin (encoded by *cpb*) and epsilon-toxin (encoded by *etx*). Beta-toxin in combination with epsilon-toxin gives rise to type B, but when expressed ►

individually each gives rise to types C and D, respectively (8). The iota-toxin is encoded by the *iA* gene and associated with type E. In addition, a novel toxin, beta2-toxin, has been found to be encoded by *cpb2*, but its association with a disease in animals remains unclear (9, 10). In poultry, *C. perfringens* type A can cause necrotic enteritis by expressing the virulence factor *netB*, which encodes a toxin (11). Outbreaks are most common in broiler chickens aged 2–6 weeks and have devastating economic effects on producers because the disease causes high mortality and considerable loss of productivity. Therefore, it is of commercial interest to rapidly identify *C. perfringens* infection in the poultry industry (12).

The results of a PCR assay used to genotype *C. perfringens* isolates are presented in Figure 3. The assay identifies genes *cpa*, *cpb*, *etx*, *iA*, *cpb2* and *cpe* encoding the alpha-, beta-, epsilon-, iota-, and beta2-toxins, and enterotoxin (13, 14). A previously described PCR protocol (11) was used to detect the *netB* gene.




Figure 3. PCR product analysis of *Clostridium perfringens* strains on QIAxcel. The genes identified are *cpb2* and *cpa* in lane 1, *etx*, *cpa*, and *cpb* in lane 2, *cpb2*, *iA*, and *cpa* in lane 3, and *cpa* and *cpe* in lane 4. Lane 5 contains the negative control and lane 6 displays the DNA size marker 50–800 bp.

Shiga toxin-producing *Escherichia coli* (STEC)

A group of pathogenic *E. coli*, STEC (also called verotoxin-producing) can cause bloody diarrhea and more severe diseases in humans, such as hemolytic-uremic syndrome (HUS). The infection occurs after consumption of contaminated food (including ground beef, leafy greens, and sprouts) and water.

STEC possess unique virulence factors that facilitate effective colonization of the human gastrointestinal tract and subsequent release of Shiga toxin and are often associated with foodborne outbreaks. The major virulence genes of STEC are *stx1*, *stx2*, and *eae*, which encode toxins Stx1, Stx2 and intimin, respectively (15, 16). Toxin Stx1 is almost identical to the toxin produced by *Shigella dysenteriae*, though the second is immunologically distinct (17). There is an increasing demand for improved diagnostic procedures for the detection of STEC not only in fecal samples (the organism is present in the gastrointestinal tract of healthy ruminants) and food (e.g., meat and dairy products), but also in drinking water. The organism propagates in the animal gut and is shed in feces, which represents the major source of food contamination (18). A large number of serogroups are Shiga toxin producers; however, human infections are associated with a minor number of them (19).

In our laboratory, detection and genetic characterization of STEC from enrichment cultures is performed by multiplex PCR targeting *stx1*, *stx2*, *eae*, *hlyA* and the *stx2f* subtype of *stx2* gene (20, 21; Figure 4).




Figure 4. PCR product analysis of STEC strains on QIAxcel. The genes identified are *hly* and *eae* in lane 1, *hly*, *eae*, and *stx2* in lane 2, *stx2* in lane 4, *hly*, *eae*, *stx2*, and *stx1* in lane 5, and *stx2f* and *eae* in lane 6. The sample in lane 3 was negative. Lane 7 contains the negative control and lane 8 displays the DNA size marker 50–800 bp.

Mycobacterium avium

This bacterial species comprises organisms that range from ubiquitous mycobacteria causing opportunistic infections in a variety of hosts to obligate pathogens of birds and ruminants. On the basis of biochemical and DNA analysis, pathogenicity, and host-range characteristics, this species is currently divided into three subspecies: *M. avium* subsp. *avium*, *M. avium* subsp. *paratuberculosis*, and *M. avium* subsp. *silvaticum* (22) while a fourth, *M. avium* subsp. *hominissuis* (23), has not been formally proposed yet (24; www.bacterio.net).

M. avium subsp. *avium* (MAA) is the causative agent of avian tuberculosis. The bacterium may infect many animal species but birds (especially poultry) are particularly susceptible to infection which often leads to fatal outcome. In farm animals, particularly pigs and cattle, MAA causes avian tuberculosis with tuberculous lesions mostly localized in the head and mesenteric lymph nodes (25). *M. avium* subsp. *hominissuis* (MAH) was proposed to distinguish organisms found in humans and pigs from those isolated from birds. Predominantly found in the environment, MAH isolates are weakly virulent for birds but are frequently encountered in tuberculous lesions in different animals (especially pigs, rarely in cattle, deer, and other animals) and in humans (23). From the veterinary laboratory and disease investigation standpoint, *M. avium* subsp. *paratuberculosis* (MAP) which causes chronic enteritis in ruminants (paratuberculosis) is often referred to separately from other subspecies of *M. avium*.

The genome of *M. avium* contains several insertion sequences, short DNA fragments which are capable of transposition. These insertion sequences can be species-specific, strain-specific, or specific for a certain group of strains and represent a suitable tool in epidemiological research as genetic markers (26). The insertion sequence IS901 is present in all MAA isolates that are virulent for mice (27) and in all field MAA isolates that are fully virulent for birds (25), while IS1245 is found in the majority of *M. avium* isolates (28, 29) except in MAP. IS901 is used routinely in our laboratory to differentiate between MAA (IS1245+, IS901+) and MAH (IS1245+, IS901-), following the protocol of Kunze et al. (27). IS901 and its flanking regions are targeted, generating PCR products of 1108 bp or 1742 bp in case of MAA isolates (depending on amplification of IS901 itself or its flanking regions, respectively) and a PCR product of 300 bp in MAH isolates from amplifying flanking regions (Figure 5).

However, a protocol for detection of IS1245 (28) is usually applied to first identify isolates to the species level (Figure 6).




Figure 5. PCR product analysis of *M. avium* subsp. *avium* and *M. avium* subsp. *hominissuis* isolates on QIAxcel. Lanes 1–3: *M. avium* subsp. *hominissuis* (IS901-), lanes 4–5: *M. avium* subsp. *avium* (IS901+), lane 6: negative control, lane 7: DNA size marker 100 bp–2.5 kb.




Figure 6. PCR product analysis of *M. avium* isolates on QIAxcel. Lanes 1–3 show positive samples (IS1245+) while lanes 4 and 5 show negative samples. Lane 6 is a positive control, lane 7 contains the negative control, and lane 8 displays the DNA size marker 50–800 bp.

Methicillin-resistant *Staphylococcus aureus* (MRSA)

First reported in 1961, MRSA is the most common multidrug-resistant pathogen causing nosocomial infections in Europe (30). Hospitals and other healthcare facilities have been the traditional locations for MRSA infections for many decades. Since the 1990s, an increasing number of community-acquired MRSA infections have been documented in some countries. In recent years, studies have concentrated on livestock and the term livestock-associated MRSA is often used to describe these isolates (31). This pathogen has been isolated from both healthy and diseased animals including dogs, cats, horses, cows, pigs, and poultry.

Until recently, MRSA was generally considered a human problem, with most efforts devoted to controlling MRSA within the human population (32). However, reports of a porcine-associated MRSA strain causing infections in humans in close contact with pigs (33) have highlighted the risk that animals, particularly food-producing animals, might pose to public health. In addition, the recent discovery of a novel bovine-associated MRSA strain infecting cattle and humans (34) has heightened concerns about livestock acting as a potential reservoir of resistant pathogens. Rapid detection and accurate identification of MRSA are therefore crucial for the effective management of infections caused by this organism.

In our laboratory, a multiplex PCR assay using primer sets previously described (35, 36) is implemented to identify MRSA isolates. The assay targets *nuc*, which encodes the thermostable nuclease specific for *S. aureus*, *mecA*, which encodes the low-affinity penicillin-binding protein PBP 2A, and 16S rRNA genes. Detailed instructions to perform this assay are available in the protocol established by the DTU National Food Institute of the Technical University of Denmark (37). The results of the PCR product analysis on QIAxcel are shown in Figure 7.




Figure 7. QIAxcel analysis of *Staphylococcus aureus* isolates. Lanes 1 and 2: *S. aureus* (*nuc*+), lane 3: MRSA (*nuc*+, *mecA*+), lane 4: negative control, lane 5: DNA size marker 100 bp–2.5 kb. Amplicon of the 16S rRNA gene amplified as an internal control is visible in lanes 1–3.

Conclusion

With a growing need for fast and accurate identification of pathogenic bacteria, the QIAxcel System for analysis of PCR products can make an important contribution to improving veterinary diagnostics.

References

1. World Health Organization (2011) World health statistics 2011. Geneva: WHO Press.
2. Corry, J.E. and Atabay, H.I. (2001) Poultry as a source of *Campylobacter* and related organisms. Symp. Ser. Soc. Appl. Microbiol, **30**, 96S.
3. Skirrow, M.B. (1994) Diseases due to *Campylobacter*, *Helicobacter* and related bacteria. J. Comp. Pathol. **111**, 113.
4. Wang, G. et al. (2002) Colony multiplex PCR assay for identification and differentiation of *Campylobacter jejuni*, *C. coli*, *C. lari*, *C. upsaliensis*, and *C. fetus* subsp. *fetus*. J. Clin. Microbiol. **40**, 4744.
5. Keessen, E.C., Gaastra, W., and Lipman, L.J. (2011) *Clostridium difficile* infection in humans and animals, differences and similarities. Vet. Microbiol. **153**, 205.
6. Avberšek, J., Zajc, U., Mičunović, J., and Očepek, M. (2013) Improved detection of *Clostridium difficile* in animals by using enrichment culture followed by LightCycler real-time PCR. Vet. Microbiol. **164**, 93.
7. Lemeë, L. et al. (2004) Multiplex PCR targeting *tpi* (triose phosphate isomerase), *tcdA* (Toxin A), and *tcdB* (Toxin B) genes for toxigenic culture of *Clostridium difficile*. J. Clin. Microbiol. **42**, 5710.
8. Quinn, P.J. et al. (2011) Veterinary Microbiology and Microbial Disease. 2nd ed. West Sussex: Wiley-Blackwell.
9. Gibert, M., Jolivet-Reynaud, C., and Popoff, M.R. (1997) Beta2 toxin, a novel toxin produced by *Clostridium perfringens*, Gene **203**, 65.
10. Farzan, A. et al. (2013) An investigation into the association between *cpb2*-encoding *Clostridium perfringens* type A and diarrhea in neonatal piglets. Can. J. Vet. Res. **77**, 45.
11. Keyburn, A.I. et al. (2008) *NetB*, a new toxin that is associated with avian necrotic enteritis caused by *Clostridium perfringens*, PLoS Pathog. **4**, e26. doi:10.1371/journal.ppat.0040026.
12. Cooper, K.K., Songer, J.G., and Uzal, F.A. (2013) Diagnosing clostridial enteric disease in poultry. J. Vet. Diagn. Invest. **25**, 314.
13. Meer, R.R. and Songer, J.G. (1997) Multiplex polymerase chain reaction assay for genotyping *Clostridium perfringens*. Am. J. Vet. Res. **58**, 702.
14. Garmory, H.S., Chanter, N., French, N.P., Bueschel, D., Songer, J.G., and Titball, R.W. (2000) Occurrence of *Clostridium perfringens* beta2-toxin amongst animals, determined using genotyping and subtyping PCR assays. Epidemiol. Infect. **124**, 61.
15. Farrokh, C. et al. (2013) Review of Shiga-toxin-producing *Escherichia coli* (STEC) and their significance in dairy production, Int. J. Food Microbiol. **162**, 190.
16. Bai, J., Paddock, Z.D., Shi, X., Li, S., An, B., and Nagaraja, T.G. (2012) Applicability of a multiplex PCR to detect the seven major Shiga toxin-producing *Escherichia coli* based on genes that code for serogroup-specific O-antigens and major virulence factors in cattle feces. Foodborne Pathog. Dis. **9**, 541.
17. Scheutz, F. et al. (2012) Multicenter evaluation of a sequence-based protocol for subtyping Shiga toxins and standardizing Stx nomenclature. J. Clin. Microbiol. **50**, 2951.
18. Ferens, W.A. and Hovde, C.J. (2011) *Escherichia coli* O157:H7: animal reservoir and sources of human infection. Foodborne Pathog. Dis. **8**, 465.
19. European Food Safety Authority. (2010) The Community summary report on trends and sources of zoonoses, zoonotic agents and food-borne outbreaks in the European Union in 2008. The EFSA Journal **8**, 1496.
20. Paton, A.W. and Paton, J.C. (1998) Detection and characterization of Shiga toxin-producing *Escherichia coli* by using multiplex PCR assays for *stx1*, *stx2*, *eaeA*, enterohemorrhagic *E. coli* *hlyA*, *rfbO111*, and *rfbO157*. J. Clin. Microbiol. **36**, 598.
21. Schmidt et al. (2000). A new Shiga toxin 2 variant (Stx2f) from *Escherichia coli* isolated from pigeons. Appl. Environ. Microbiol. **66**, 1205.
22. Thorel, M.F., Krichevsky, M., and Lévy-Frébault, V.V. (1990) Numerical taxonomy of mycobactin-dependent mycobacteria, emended description of *Mycobacterium avium*, and description of *Mycobacterium avium* subsp. *avium* subsp. nov., *Mycobacterium avium* subsp. *paratuberculosis* subsp. nov., and *Mycobacterium avium* subsp. *silvaticum* subsp. nov. Int. J. Syst. Bacteriol. **40**, 254.
23. Mijis, W., et al. (2002) Molecular evidence to support a proposal to reserve the designation *Mycobacterium avium* subsp. *avium* to birdtype isolates and *M. avium* subsp. *hominissuis* for the human/porcine type of *M. avium*. Int. J. Syst. Evol. Microbiol. **52**, 1505.
24. Turenne, C.Y., Wallace, R. Jr, and Behr, M.A. (2007) *Mycobacterium avium* in the postgenomic area. Clin. Microbiol. Rev. **20**, 205.
25. Pavlik, I., Svastova, P., Bartl, J., Dvorska, L., and Rychlik, I. (2000) Relationship between IS901 and the *Mycobacterium avium* complex strains isolated from birds, animals, humans and the environment and virulence for poultry. Clin. Diagn. Lab. Immunol. **7**, 212.
26. Grange, J.M. (1996) Mycobacteria and human disease. 2nd ed. London: Arnold.
27. Kunze, Z.M., Portaels, F., and McFadden, J.J. (1992) Biologically distinct subtypes of *Mycobacterium avium* differ in possession of insertion sequence IS901. J. Clin. Microbiol. **30**, 2366.
28. Guerrero, C., Bernasconi, C., Burki, D., Bodmer, T., and Telenti, A. (1995) A novel insertion element from *Mycobacterium avium*, IS1245, is a specific target for analysis of strain relatedness. J. Clin. Microbiol. **33**, 304.
29. Roiz, M.P., Palenque, E., Guerrero, C., and Garcia, M.J. (1995) Use of restriction fragment length polymorphism as a genetic marker for typing *Mycobacterium avium* strains. J. Clin. Microbiol. **33**, 1389.
30. Köck, R. et al. (2010) Methicillin-resistant *Staphylococcus aureus* (MRSA): burden of disease and control challenges in Europe. Euro. Surveill. **15**, 19688.
31. Köck, R., Mellmann, A., Schaumburg, F., Friedrich, A.W., Kipp, F., and Becker, K. (2011) The epidemiology of methicillin-resistant *Staphylococcus aureus* (MRSA) in Germany. Dtsch. Arztebl. Int. **108**, 761.
32. Garcia-Alvarez, L., Dawson, S., Cookson, B., and Hawkey, P. (2012) Working across the veterinary and human health sectors. J. Antimicrob. Chemother. **67**, i37.
33. Huijsdens, X.W. et al. (2006). Community-acquired MRSA and pig-farming. Ann. Clin. Microbiol. Antimicrob. **5**, 26.
34. Garcia-Alvarez L. et al. (2011) Methicillin-resistant *Staphylococcus aureus* with a novel *mecA* homologue in human and bovine populations in the UK and Denmark: a descriptive study. Lancet Infect. Dis. **11**, 595.
35. Salisburry, S. (1997) Genotypic detection of methicillin-resistant staphylococci by multiplex PCR assay. LabMedica Int. **1997**, 16.
36. Poulsen, A.B., Skov, R., and Pallesen, L.V. (2003) Detection of methicillin resistance in coagulase-negative staphylococci and in staphylococci directly from simulated blood cultures using the EVIGENE MRSA Detection Kit. J. Antimicrob. Chemother. **51**, 419.
37. National Food Institute, Technical University of Denmark (2008). Multiplex PCR for the detection of *mecA* gene. www.crl-ar.eu/data/images/meca-pcr_protocol%2006.02.08.pdf.

Ordering Information

Product	Contents	Cat. no.
QIAXcel Advanced Instrument	Capillary electrophoresis device: includes computer, QIAXcel ScreenGel software, and 1-year warranty on parts and labor	9001941
QIAXcel DNA High Resolution Kit (1200)	QIAXcel DNA High Resolution Gel Cartridge, Buffers, Mineral Oil, QX Intensity Calibration Marker, 12-Tube Strips	929002
QIAXcel DNA Screening Kit (2400)	QIAXcel DNA Screening Gel Cartridge, Buffers, Mineral Oil, QX Intensity Calibration Marker, 12-Tube Strips	929004
QX Alignment Marker 15 bp/3 kb (1.5 ml)	Alignment marker with 15 bp and 3 kb fragments	929522
QX Alignment Marker 15 bp/1 kb (1.5 ml)	Alignment marker with 15 bp and 1 kb fragments	929521
QX DNA Size Marker 50–800 bp (50 µl)	DNA size marker with fragments of 50, 100, 150, 200, 250, 300, 400, 500, 600, 700, and 800 bp; concentration 100 ng/µl	929561
QX DNA Size Marker 100 bp–2.5 kb (50 µl)	DNA size marker with fragments of 100, 200, 300, 400, 500, 600, 700, 800, 1000, 1200, 1500, 2000, and 2500 bp; concentration 100 ng/µl	959559

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Visit www.qiagen.com/QIAXcel and find out how automated gel electrophoresis can benefit your lab!

Trademarks: QIAGEN®, QIAXcel® (QIAGEN Group).
1083558 06/2014 © 2014 QIAGEN, all rights reserved.

www.qiagen.com

Australia ■ 1-800-243-800
Austria ■ 0800-281011
Belgium ■ 0800-79612
Brazil ■ 0800-557779
China ■ 800-988-0325
Denmark ■ 80-885945
Finland ■ 0800-914416

France ■ 01-60-920-930
Germany ■ 02103-29-12000
Hong Kong ■ 800 933 965
India ■ 1-800-102-4114
Ireland ■ 1800 555 049
Italy ■ 800-787980
Japan ■ 03-6890-7300

Korea (South) ■ 080-000-7145
Luxembourg ■ 8002 2076
Malaysia ■ 603-7981-5510
Mexico ■ 01-800-7742-436
The Netherlands ■ 0800-0229592
Norway ■ 800-18859
Singapore ■ 1800-742-4368

Spain ■ 91-630-7050
Sweden ■ 020-790282
Switzerland ■ 055-254-22-11
Taiwan ■ 0800-665-957
UK ■ 0808-234-3665

