

# TissueLyser LT Adapter, 12-Tube

For disruption of up to 12 samples on the TissueLyser LT

## Kit Contents

TissueLyser LT Adapter, 12-Tube	
Catalog no.	69980
Base	1
Top	1
Insert with 12 Aluminum Tube Holders	1
Product Sheet	1

## Storage

The TissueLyser LT Adapter should be stored at ambient temperature. Prior to sample disruption, the insert of the TissueLyser LT Adapter (including aluminum tube holders) can be placed on dry ice. Do not place the adapter in liquid nitrogen.

## Product Description

The TissueLyser LT Adapter allows disruption of up to 12 biological samples on the TissueLyser LT. It consists of the following components:

- A base which attaches to the piston of the TissueLyser LT
- An insert with 12 aluminum tube holders, each of which holds a 2 ml tube containing sample and bead
- A top which screws on top of the base

The TissueLyser LT moves the adapter up and down very rapidly to cause sample disruption.

## Procedure

For instructions on installing and cleaning the TissueLyser LT Adapter, please refer to the "Installation Procedures" and "Maintenance Procedures" chapters of the *TissueLyser LT User Manual*.



## Ordering Information

Product	Contents	Cat. no.
TissueLyser LT Adapter, 12-Tube	Adapter for disruption of up to 12 samples on the TissueLyser LT	69980
<b>Related products</b>		
TissueLyser LT	Bead mill for low- to medium-throughput sample disruption	85600
Stainless Steel Beads, 5 mm (200)	Stainless steel beads, suitable for use with TissueLyser systems	69989
Tungsten Carbide Beads, 3 mm (200)	Tungsten carbide beads, suitable for use with TissueLyser systems	69997
Sample Tubes RB (2 ml)	1000 x 2 ml sample tubes, suitable for use with the TissueLyser LT	990381

The TissueLyser LT is intended for molecular biology applications. This product is neither intended for the diagnosis, prevention, or treatment of a disease, nor has it been validated for such use either alone or in combination with other products.

Trademarks: QIAGEN® (QIAGEN Group).

1056714 05/2009 © 2009 QIAGEN, all rights reserved.

[www.qiagen.com](http://www.qiagen.com)

Canada = 800-572-9613

Ireland = 1800 555 049

Norway = 800-18859

China = 0086-21-3865-3865

Italy = 800-787980

Singapore = 65-67775366

Denmark = 80-885945

Japan = 03-6890-7300

Spain = 91-630-7050

Australia = 1-800-243-800

Finland = 0800-914416

Korea (South) = 1544 7145

Sweden = 020-790282

Austria = 0800/281010

France = 01-60-920-930

Luxembourg = 8002 2076

Switzerland = 055-254-22-11

Belgium = 0800-79612

Germany = 02103-29-12000

Mexico = 01-800-7742-639

UK = 01293-422-911

Brazil = 0800-557779

Hong Kong = 800 933 965

The Netherlands = 0800 0229592

USA = 800-426-8157

