

## QX Intensity Calibration Marker

### For calibration of signal intensity

QX Intensity Calibration Marker is for calibration of signal intensity for each new gel cartridge in automated analysis of DNA fragments and RNA and is intended for use with all QIAxcel® Kits. The volume of QX Intensity Calibration Marker provided is sufficient for the number of sample runs specified in the kit handbook.

### Storage conditions

QX Intensity Calibration Marker can be stored at 2–8°C for up to 6 months.

### Quality control

In accordance with QIAGEN's ISO-certified Quality Management System, each lot of QX Intensity Calibration Marker is tested against predetermined specifications to ensure consistent product quality.

### Safety information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient and compact PDF format at [www.qiagen.com/safety](http://www.qiagen.com/safety) where you can find, view, and print the SDS for each QIAGEN kit and kit component.

### Product use limitations

QX Intensity Calibration Marker is intended to be used only in combination with QIAxcel Kits for applications described in the respective QIAxcel Kit handbooks.

### Procedure

For complete instructions for use, please refer to the *QIAxcel DNA Handbook*, which can be downloaded from the resources tab at [www.qiagen.com/p/QIAxcel](http://www.qiagen.com/p/QIAxcel).



## Ordering Information

Product	Contents	Cat. no.
QX Intensity Calibration Marker (1 ml)	1 ml QX Intensity Calibration Marker	929500

QX Intensity Calibration Marker is intended for molecular biology applications. This product is not intended for the diagnosis, prevention, or treatment of a disease.

Trademarks: QIAGEN®, QIAxcel® (QIAGEN Group).  
1089526 11/2014 © 2014 QIAGEN, all rights reserved.

[www.qiagen.com](http://www.qiagen.com)  
Australia ■ 1-800-243-800  
Austria ■ 0800-281011  
Belgium ■ 0800-79612  
Brazil ■ 0800-557779  
Canada ■ 800-572-9613  
China ■ 800-988-0325

Denmark ■ 80-885945  
Finland ■ 0800-914416  
France ■ 01-60-920-930  
Germany ■ 02103-29-12000  
Hong Kong ■ 800 933 965  
India ■ 1-800-102-4114  
Ireland ■ 1800 555 049

Italy ■ 800-787980  
Japan ■ 03-6890-7300  
Korea (South) ■ 080-000-7145  
Luxembourg ■ 8002 2076  
Mexico ■ 01-800-7742-436  
The Netherlands ■ 0800-0229592  
Norway ■ 800-18859

Singapore ■ 1800-742-4368  
Spain ■ 91-630-7050  
Sweden ■ 020-790282  
Switzerland ■ 055-254-22-11  
Taiwan ■ 0080-665-1947  
UK ■ 0808-2343665  
USA ■ 800-426-8157

