

ISO 18385 Forensic DNA Grade Certificate

The following QIAGEN products are manufactured in compliance to ISO18385:2016:

Investigator STR GO! Lysis Buffer (200), cat. no. 386516

Investigator STR GO! Punch Buffer (200), cat. no. 386526

Investigator STR GO! Punch Buffer (1000), cat. no. 386528

Our manufacturing facilities and production processes are ISO18385-compliant and meet the following requirements:

- Quality management system
- Staff exclusion database service
- Risk management of human DNA contamination by performing risk assessments for potential human DNA contamination and implementation of risk mitigations
- Validated cleaning and environmental monitoring procedures
- Batch-release testing with validated test methods

