


pQE-40 Vector

Positions of elements in bases

Vector size	4031
Start of numbering at <i>Xho</i> I (CTCGAG)	1-6
T5 promoter/lac operator element	7-87
T5 transcription start	61
6xHis-tag coding sequence	127-144
Dehydrofolate reductase coding sequence	157-714
Multiple cloning site	721-762
Lambda <i>t_o</i> transcriptional termination region	778-872
<i>rrnB</i> T1 transcriptional termination region	1634-1732
ColE1 origin of replication	2238
β-lactamase coding sequence	3826-2966



pQE-40

