


pQE-70 Vector

Positions of elements in bases

Vector size (bp)	3426
Start of numbering at <i>Xho</i> I (CTCGAG)	1–6
T5 promoter/lac operator element	7–87
T5 transcription start	61
6xHis-tag coding sequence	133–150
Multiple cloning site	113–132
Lambda <i>t</i> ₀ transcriptional termination region	173–267
<i>rnnB</i> T1 transcriptional termination region	1029–1127
ColE1 origin of replication	1603
β-lactamase coding sequence	3221–2361



pQE-70

