

QX DNA Alignment Markers

For calibration of migration time variation

QX DNA Alignment Markers are for calibration of migration time variation across all channels and are intended for use with QIAxcel® Kits for automated analysis of DNA fragments using QIAxcel or QIAxcel Advanced systems.

Storage conditions

QX DNA Alignment Markers can be stored at 2–8°C for up to 1 year.

Quality control

In accordance with QIAGEN's ISO-certified Quality Management System, each lot of QX DNA Alignment Marker is tested against predetermined specifications to ensure consistent product quality.

Safety information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient and compact PDF format at www.qiagen.com/safety where you can find, view, and print the SDS for each QIAGEN kit and kit component.

Product use limitations

QX DNA Alignment Markers are intended to be used only in combination with QIAxcel Kits for applications described in the respective QIAxcel Kit handbooks.

Procedure

For complete instructions for use, please refer to the *QIAxcel DNA Handbook*, which can be downloaded from the resources tab at www.qiagen.com/p/QIAxcel.



Ordering Information

Product	Contents	Cat. no.
QX DNA Alignment Marker 15 bp/500 bp (1.5 ml)	Alignment marker with 15 bp and 500 bp fragments	929520
QX DNA Alignment Marker 15 bp/600 bp (1.5 ml)	Alignment marker with 15 bp and 600 bp fragments	929530
QX DNA Alignment Marker 15 bp/1 kb (1.5 ml)	Alignment marker with 15 bp and 1 kb fragments	929521
QX DNA Alignment Marker 15 bp/3 kb (1.5 ml)	Alignment marker with 15 bp and 3 kb fragments	929522
QX DNA Alignment Marker 15 bp/5 kb (1.5 ml)	Alignment marker with 15 bp and 5 kb fragments	929524
QX DNA Alignment Marker 15 bp/10 kb (1.5 ml)	Alignment marker with 15 bp and 10 kb fragments	929523
QX DNA Alignment Marker 50 bp/1 kb (1.5 ml)	Alignment marker with 50 bp and 1 kb fragments	929526
QX DNA Alignment Marker 50 bp/5 kb (1.5 ml)	Alignment marker with 50 bp and 5 kb fragments	929529

QX DNA Alignment Markers are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

Trademarks: QIAGEN®, QIAxcel® (QIAGEN Group).
1081893 09/2014 © 2014 QIAGEN, all rights reserved.

www.qiagen.com
Australia ■ 1-800-243-800
Austria ■ 0800-281011
Belgium ■ 0800-79612
Brazil ■ 0800-557779
Canada ■ 800-572-9613
China ■ 800-988-0325

Denmark ■ 80-885945
Finland ■ 0800-914416
France ■ 01-60-920-930
Germany ■ 02103-29-12000
Hong Kong ■ 800 933 965
India ■ 1-800-102-4114
Ireland ■ 1800 555 049

Italy ■ 800-787980
Japan ■ 03-6890-7300
Korea (South) ■ 080-000-7145
Luxembourg ■ 8002 2076
Mexico ■ 01-800-7742-436
The Netherlands ■ 0800-0229592
Norway ■ 800-18859

Singapore ■ 1800-742-4368
Spain ■ 91-630-7050
Sweden ■ 020-790282
Switzerland ■ 055-254-22-11
Taiwan ■ 0080-665-1947
UK ■ 0808-2343665
USA ■ 800-426-8157

