

# Take the next step towards efficiency

Automate your sample preparation  
with the QIAcube<sup>®</sup> HT



QIAcube HT: automated mid- to high-throughput nucleic acid purification in 96-well format

- Silica-based membrane technology with dedicated purification kits ensures high-quality of purified nucleic acids
- Convenient, flexible and easy-to-use instrument with small footprint
- Innovative design features increase safety and minimize cross-contamination
- Efficient worktable decontamination with UV light
- Optimized use of plasticware with a smart and economical tip reuse strategy
- Flexible throughput of 24–96 samples in a single run, with sample volumes of up to 200  $\mu$ l
- Intuitive software with pre-installed protocols

## Simple and reliable automated processing for cost and time saving

		QIAamp 96 DNA QIAcube HT Kit	DNeasy <i>mericon</i> 96 QIAcube HT Kit	QIAamp 96 Virus QIAcube HT Kit	RNeasy 96 QIAcube HT Kit	New! miRNeasy 96 Advanced QIAcube HT Kit	QIAamp 96 PowerFecal QIAcube HT Kit	New! DNeasy 96 PowerSoil Pro QIAcube HT Kit
Nucleic acid	DNA	•	•	•			•	•
	RNA			•	•	•		
	miRNA				•	•		
Sample type	Cells and tissue	•		•	•			
	Whole blood	•		•				
	Buffy coat	•						
	Serum and plasma			•		•		
	Food		•					
	Stool			•			•	•
	Soil							•
	Plant tissue*	•	•		•			

\*Available as custom protocol

## Ordering Information

Product	Kit description	Cat. no.
QIAcube HT System	Robotic workstation with UV lamp, HEPA filter, laptop, QIAcube HT operating software, startup pack, installation and training, 1-year warranty on parts and labor	9001793
QIAamp 96 DNA QIAcube HT Kit (5)	For 480 preps: QIAamp 96 plates, QIAGEN Proteinase K, buffers	51331
DNeasy <i>mericon</i> 96 QIAcube HT Kit (5)	For 480 preps: DNeasy 96 plates, QIAGEN Proteinase K, buffers	69571
QIAamp 96 Virus QIAcube HT Kit (5)	For 480 preps: QIAamp 96 plates, QIAGEN Proteinase K, carrier RNA, buffers	57731
RNeasy 96 QIAcube HT Kit (5)	For 480 preps: RNeasy 96 plates, RNase-free water, buffers	74171
miRNeasy 96 Advanced QIAcube HT Kit (5)	For 480 preps: RNeasy 96 plates, RNase-free water, buffers	217261
QIAamp 96 PowerFecal QIAcube HT Kit (5)	For 480 preps: QIAamp 96 plates, QIAGEN Proteinase K, buffers	51531
DNeasy 96 PowerSoil Pro QIAcube HT Kit (5)	For 480 preps: QIAamp 96 plates, buffers*	47021
QIAcube HT Plasticware	Plasticware for 480 typical preps on QIAcube HT	950067

\*For mechanical disruption, PowerBead Pro Plates or PowerBead Pro Tubes must be ordered separately.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at [www.qiagen.com](http://www.qiagen.com) or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, Sample to Insight®, QIAamp®, QIAcube®, DNeasy®, *mericon*®, PowerFecal®, PowerSoil®, RNeasy® (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

1118446 08/2019 PROM-14446-002 © 2019 QIAGEN, all rights reserved.

Ordering [www.qiagen.com/shop](http://www.qiagen.com/shop) | Technical Support [support.qiagen.com](mailto:support.qiagen.com) | Website [www.qiagen.com](http://www.qiagen.com)