

March 2015

Important Note

Dear *artus* CMV RGQ MDx Kit user,

We would like to inform you about important updates to the information in *artus*® CMV RGQ MDx Kit Instructions for Use (Handbook).

The *artus* CMV RGQ MDx Kit (cat. no. 4503245) can now be used with the new Rotor-Gene® Q Software version 2.3 on Rotor-Gene Q MDx instruments. This new software is compatible with both Windows® XP and Windows 7 operating systems. If your Rotor-Gene Q MDx instrument was delivered with a previous software version, you can easily upgrade by downloading Rotor-Gene Q Software version 2.3 from the [Rotor-Gene Q MDx \(US\)](#) product page on www.qiagen.com. You can find the new software in the **Product Resources** section under the **Operating Software** tab.

For Rotor-Gene Q Software version 2.3, a new assay package for automated quantitation of human cytomegalovirus (CMV) DNA has been released, Rotor-Gene Q *artus* CMV Assay Package 3.1.1 (cat. no. 9023822). This assay package is available to order on CD.

Note: Rotor-Gene Q *artus* CMV Assay Package 1.2.7 (cat. no. 9022860) is not compatible with Rotor-Gene Q Software version 2.3.

Upgrading to Rotor-Gene Q Software version 2.3 is recommended but not mandatory as QIAGEN® will continue to support Rotor-Gene Q Software version 2.1 and Rotor-Gene Q *artus* CMV Assay Package 1.2.7.

If you have any questions, please contact your local QIAGEN Technical Service Department.

Best regards,

Your QIAGEN team

Trademarks: QIAGEN®, *artus*®, Rotor-Gene® (QIAGEN Group); Windows® (Microsoft Corporation).
1093689 03/2015 © 2015 QIAGEN, all rights reserved.